# Biodesign®

## TISSUE GRAFT

FP0005-04N

**LOT** Batch code

**REF** Catalogue number

Consult instructions for use

Contains biological material of animal origin

Country of manufacture

Date of manufacture

Do not resterilize

Do not reuse

Do not use if package is damaged and consult instructions for use

Keep dry

Manufacturer

Medical device

Double sterile barrier system

Sterilized using ethylene oxide

Unique device identifier

Use-by date

**MANUFACTURER**  
COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.

www.cookbiotech.com  
© COOK BIOTECH 2024

March 2024

ENGLISH

## BIODESIGN® TISSUE GRAFT

### DEVICE DESCRIPTION:

The Cook® Biodesign® Tissue Graft is a dried single or multi-layered small intestinal submucosa (SIS) sheet. The graft is used to reinforce soft tissue and can be cut to size to accommodate the patient's anatomy. The Tissue Graft is provided sterile for single use only.

### INTENDED USE

The Cook® Biodesign® Tissue Graft is intended for implantation to reinforce soft tissue. The graft is supplied sterile in peel-open packages and is intended for one-time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**MR** Magnetic resonance safe

**1-LAYER TISSUE GRAFT** This symbol means the following: 1-Layer Tissue Graft

**4-LAYER TISSUE GRAFT** This symbol means the following: 4-Layer Tissue Graft

This product is intended for use by trained medical professionals.

### CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine materials.

### PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- The single-layer graft should not be used in applications requiring high strength.
- The graft may not have sufficient strength to support stresses encountered in some ventral hernias or large-area, body-wall repairs.
- Ensure that the graft is rehydrated prior to suturing or stapling.
- Patients undergoing radiation therapy may not experience normal wound healing.
- **The graft's performance has not been evaluated with suture spacing greater than 2 mm.**
- **Ensure that all layers of the graft are secured when suturing or stapling.**

### POTENTIAL COMPLICATIONS

Complications that can occur with the use of any surgical graft material may include, but are not limited to:

- Abscess
- Adhesion
- Allergic Reaction
- Bleeding
- Delayed or failed incorporation
- Erosion
- Extrusion
- Fever
- Fistula Formation
- Induration
- Infection
- Inflammation
- Migration
- Pain
- Recurrence
- Seroma formation
- Tissue trauma

If any of the following conditions occur and cannot be resolved, device removal should be considered:

- Acute or chronic inflammation (initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction
- Infection
- Seroma formation

### STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

### STERILIZATION

This graft has been sterilized with ethylene oxide.

### INSTRUCTIONS FOR USE

#### Required Materials

- A sterile basin
- Hydration fluid: room temperature sterile saline or sterile lactated Ringer's solution.

**NOTE: Handle the graft using aseptic technique. Minimize contact with latex gloves.**

1. Remove the packaging containing the graft from the envelope.
2. Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.

3. Open the inner pouch and remove the graft.
4. Trim the graft to fit the patient's anatomy, if necessary, providing a small allowance for overlap with the adjacent tissue.
5. Place the graft into the sterile basin and hydrate the graft until the desired handling characteristics are achieved. A hydration time of greater than one minute is not required. Multiple grafts may be hydrated simultaneously in the same basin.
6. Prepare the graft site using standard surgical techniques.
7. Suture or staple the graft into place, avoiding excess tension. **NOTE:** Surgical experience indicates that suturing or stapling graft sheets with close tissue approximation produces better outcomes.
8. Complete the surgical procedure.
9. Discard any unused portions of the graft according to institutional guidelines for medical waste.

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