# Biodesign®

## DURAL GRAFT

FP0030-01P

**LOT** Batch code

**REF** Catalogue number

**i** Consult instructions for use

**BIO** Contains biological material of animal origin

**US** Country of manufacture

**STERILIZE** Do not resterilize

**⌘** Do not re-use

**⊠** Do not use if package is damaged and consult instructions for use

Keep dry

Manufacturer

**MD** Medical device

Double sterile barrier system

**STERILE EO** Sterilized using ethylene oxide

**UDI** Unique device identifier

Use-by date

**MANUFACTURER**  
COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.

www.cookmedical.com  
© COOK BIOTECH 2024

January 2024

ENGLISH

## BIODESIGN® DURAL GRAFT

### DESCRIPTION

The Biodesign® Dural Graft is a nonporous sheet of extracellular collagen matrix derived from porcine small intestinal submucosa. The graft is non-pyrogenic and has sufficiently low endotoxin levels to make it suitable for applications where it will contact cerebrospinal fluid.

### INTENDED USE

The graft is intended for use as a dura substitute for repairing dura mater. This graft is supplied sterile in peel-open packages and is intended for one-time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**MR** This symbol means the following: Magnetic Resonance Safe

**DURAL GRAFT** This symbol means the following: Dural Graft

This product is intended for use by trained medical professionals.

### CONTRAINDICATIONS

The graft should not be implanted into areas with known or suspected infection or in necrotic tissue. The graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

### PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to failure and/or transmission of disease.
- **Do not resterilize.** Do not subject the graft to steam, gas (ethylene oxide), or radiation sterilization as these may damage the graft.
- Graft performance when placed in tissue having compromised circulation has not been established.
- Discard all open and unused portions of the graft.
- The graft is sterile and non-pyrogenic if the package is dry, unopened, and undamaged.
- Discard the graft if damage or contamination is observed, or if the graft is past its expiration date.
- Prior to touching the graft, rinse surgical gloves thoroughly to remove any glove powder.
- Ensure that the graft is hydrated before suturing.
- Ensure that all layers of the graft are secured when suturing.

### POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials.

- Adhesion
- Allergic reaction
- Calcification
- CSF leak
- Inflammation
- Infection

### STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

### STERILIZATION

The graft has been sterilized with ethylene oxide.

### CLINICAL EVALUATIONS

The Biodesign Dural Graft has been evaluated in a prospective multicenter clinical study.<sup>1</sup>

### INSTRUCTIONS FOR USE

#### Required Materials

- A sterile basin
- Sterile forceps
- Hydration fluid: sterile, room temperature, saline or lactated Ringer's solution

**NOTE: Always handle the graft using aseptic technique. Minimize contact with latex gloves.**

#### Procedure

1. Choose the graft size based on the size of the dura mater defect, ensuring that the graft covers the entire defect and extends slightly beyond the wound margins.
2. Remove the packaging containing the graft from the box.
3. Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
4. Open the inner pouch and remove the graft. Place the graft into the sterile basin on the sterile field. Multiple grafts may be hydrated simultaneously in the same basin.
5. Add enough hydration fluid in the basin to completely cover the graft(s).
6. Hydrate the graft in the fluid until the desired handling characteristics are achieved. A hydration time of greater than two minutes is not required.
7. If necessary, trim the graft to fit the patient's defect, providing a small allowance for overlap of the wound margins.

**NOTE: The graft may be trimmed prior to hydration. Be sure to hydrate the graft prior to suturing it in place. See steps 5-6.**

8. Transfer the graft to the defect site and suture into place, avoiding excess tension.  
**NOTE:** Suturing the graft with close tissue approximation produces better outcomes.
9. Complete the surgical procedure.
10. Discard unused portions of the graft according to institutional guidelines for disposal of medical waste.

### REFERENCES

1. Bejani GK, Zabramski J; Durasis Study Group. Safety and efficacy of the porcine small intestinal submucosa dural substitute: results of a prospective multicenter study and literature review. *J Neurosurg* 2007;106:1028-1033.

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