**COOK**

MEDICAL

**CE**  
2797

**EN**  
3

**Biodesign®  
Dural Graft**

Instructions for Use

**CS**  
4

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**DA**  
5

**Biodesign®  
Dura-implantat**

Brugsanvisning

**DE**  
7

**Biodesign®  
Duraimplantat**

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**EL**  
8

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
# BIODESIGN® DURAL GRAFT

### DESCRIPTION

Cook® Biodesign® Dural Graft is a nonporous sheet of extracellular collagen matrix derived from porcine small intestinal submucosa. The graft is non-pyrogenic and has sufficiently low endotoxin levels to make it suitable for applications where it will contact cerebrospinal fluid. Biodesign® Dural Graft is comprised of an extracellular matrix (ECM) and is fully remodeled during the healing process. Animal studies have demonstrated that complete resorption of the graft, when used as indicated, is evident by 60 days.<sup>1,2</sup>

### INTENDED USE

The graft is intended for use as a dura substitute for repairing dura mater. This graft is supplied sterile in peel-open packages and is intended for one-time use.

Biodesign® Dural Graft is designed and labeled specifically for the above indication. Use of this product for indications other than specified above is not recommended and may result in increased complications.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**MR** This symbol means the following: Magnetic Resonance Safe

**DURAL GRAFT** This symbol means the following: Dural Graft

This product is intended for use by trained medical professionals.

### CONTRAINDICATIONS

This graft should not be implanted into areas with known or suspected infection or in necrotic tissue. This graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

### PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Do not subject this graft to steam, gas (ethylene oxide), or radiation sterilization as these may damage the graft.
- Graft performance when placed in tissue having compromised circulation has not been established.
- Discard all open and unused portions of the graft sheets.
- Graft is sterile and non-pyrogenic if the package is dry, unopened, and undamaged.
- Discard graft if damage or contamination is observed, or if the graft is past its expiration date.
- Prior to touching the graft, wash surgical gloves thoroughly to remove the glove powder.
- Ensure that graft is rehydrated before suturing.
- Ensure that all layers of the graft are secured when suturing.

### POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials.

- Infection, adhesion, CSF leak, and calcification
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation.)
- Allergic reaction

### STORAGE

This graft should be stored in a clean, dry location at controlled room temperature. Do not freeze this graft.

### STERILIZATION

This graft has been sterilized with ethylene oxide.

### CLINICAL EVALUATIONS

The Biodesign® Dural Graft has been evaluated in a prospective multicenter clinical study.<sup>3</sup>

### INSTRUCTIONS FOR USE

#### Required Materials

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Rehydration fluid: at least 100 mL of sterile, room temperature, saline or sterile, room temperature lactated Ringer's solution for each graft sheet

**NOTE:** Always handle the graft sheets using aseptic technique. Minimize contact with latex gloves.

1. Using aseptic technique, remove the graft inner pouch from its outer bag, and place the inner pouch in the sterile field.
2. Open the inner pouch carefully, and aseptically remove the graft sheet with the sterile forceps.
3. Place the graft sheet into the sterile dish in the sterile field. (Multiple graft sheets may be rehydrated simultaneously in the same dish.)
4. Add to the dish at least 100 mL of the rehydration fluid for each graft sheet.
5. Allow graft sheets to rehydrate for at least two minutes.
6. Using aseptic technique, trim the graft sheet to fit the site, providing a small allowance for overlap. **NOTE:** An alternative method is to cut the graft sheet to size prior to rehydration. If this method is selected, be sure to rehydrate the graft sheet prior to suturing it into place. See step 5.
7. Using aseptic technique, transfer the graft sheet to the graft site and suture into place, avoiding excess tension. **NOTE:** Surgical experience indicates that suturing graft sheets with close tissue approximation produces better outcomes.
8. Complete the standard surgical procedure.
9. Discard any unused portions of the graft sheet following standard technique for disposal of medical waste.

## REFERENCES

1. Cobb MA, Badyalak SF, Janas W, Boop FA. Histology after dural grafting with small intestinal submucosa. *Surg Neurol.* 1996;46(4):389-393; discussion 393-394.
2. Cobb MA, Badyalak SF, Janas W, Simmons-Byrd A, Boop FA. Porcine small intestinal submucosa as a dural substitute. *Surg Neurol.* 1999;51(1):99-104.
3. Bejjani GK, Zabramski J, Durasis Study Group. Safety and efficacy of the porcine small intestinal submucosa dural substitute: results of a prospective multicenter study and literature review. *J Neurosurg.* 2007;106(6):1028-1033.

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
#### REFERENSER

1. Cobb MA, Badylak SF, Janas W, Boop FA. Histology after dural grafting with small intestinal submucosa. *Surg Neurol.* 1996;46(4):389-393; discussion 393-394.
2. Cobb MA, Badylak SF, Janas W, Simmons-Byrd A, Boop FA. Porcine small intestinal submucosa as a dural substitute. *Surg Neurol.* 1999;51(1):99-104.
3. Bejjani GK, Zabramski J, Durasis Study Group. Safety and efficacy of the porcine small intestinal submucosa dural substitute: results of a prospective multicenter study and literature review. *J Neurosurg.* 2007;106(6):1028-1033.

**Keep dry**

**Manufacturer**

**Use-by date**

**Caution**

**Do not re-use**

**Magnetic Resonance Safe**

**Date of manufacture**

**STERILE**

**EO**

**Sterilized using ethylene oxide**

**MANUFACTURER**

COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.

**EC REP**

**EC REPRESENTATIVE**  
COOK BIOTECH EUROPE APS  
Sandet, DK-4632  
Bjaeverskov, DENMARK  
Tel: +45 56868686

[www.cookmedical.com](http://www.cookmedical.com)  
© COOK BIOTECH 2024

January 2024