Biodesign®

HIATAL HERNIA GRAFT

FP0032-03M

Caution

Date of manufacture

Do not re-use

Keep dry

Manufacturer

STERILE EO

Sterilized using ethylene oxide

Temperature limit

Use-by date

MANUFACTURER
COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.

www.cookbiotech.com
© COOK BIOTECH, INC. 2022

February 2022

BIODESIGN® HIATAL HERNIA GRAFT

DESCRIPTION

The Biodesign® Hiatal Hernia Graft is a multi-layered sheet of extracellular collagen matrix derived from porcine small intestinal submucosa. The Biodesign® Hiatal Hernia Graft is implanted by a trained medical professional to support weakened soft tissue during surgery to repair hiatal hernias.

INTENDED USE

The Biodesign® Hiatal Hernia Graft is used for implantation to reinforce soft tissue where weakness exists, including repair of hiatal hernias. This graft is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This symbol means the following: Magnetic Resonance Safe.

HIATAL HERNIA GRAFT This symbol means the following: Hiatal Hernia Graft.

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

The graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to graft failure and/or transmission of disease.
- Do not resterilize.** Discard all open and unused portions.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that the graft is hydrated prior to cutting or suturing. Insufficient hydration may cause device damage during laparoscopic delivery.
- Ensure that all layers of the graft are secured when suturing.
- Ensure that the graft does not touch the esophagus.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials.

- Adhesion
- Allergic reaction
- Fistula formation
- Graft detachment/migration
- Hematoma
- Infection
- Inflammation
- Recurrence of tissue defect
- Seroma formation

The following complications are associated with use of graft materials in hiatal hernia repair.

- Dysphagia
- Esophageal erosion
- Perforation
- Stricture formation

STORAGE

Store the graft in a clean, dry location at room temperature.

STERILIZATION

The graft has been sterilized with ethylene oxide.

SUGGESTED INSTRUCTIONS FOR USING THE HIATAL HERNIA GRAFT

Required Materials

- Sterile basin
- Sterile smooth forceps
- Hydration fluid: room temperature sterile saline or sterile lactated Ringer's solution

NOTE: Handle the graft using aseptic technique, minimizing contact with latex gloves.

- Remove the packaging containing the graft from the envelope.
- Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
- Place the graft into the sterile basin on the sterile field. (Multiple grafts may be hydrated simultaneously in the same basin).
- Place enough hydration fluid in the basin to completely cover the graft(s).
- Hydrate the graft in the fluid until the desired handling characteristics are achieved. A hydration time of greater than one minute is not required.
- Prepare the graft site using standard surgical techniques.

NOTE: Suturing the graft with close tissue approximation produces better outcomes. Maximum graft-tissue approximation is only achieved with complete crural closure.

- If necessary, trim the graft to the patient's anatomy, providing a small allowance for overlap with adjacent tissues.


NOTE: The graft may be trimmed prior to hydration to fit the patient's anatomy. Be sure to hydrate the graft prior to suturing it into place. See step 5.

NOTE: For laparoscopic placement, hydrate the graft prior to loading it into the delivery port. See step 5.


- Suture the graft into place, avoiding excess tension.

NOTE: Clinical literature suggests the use of tacks or staples at the diaphragm may be associated with cardiovascular complications.¹

NOTE: The 7X10-U graft has a larger hole in the corner to identify the graft orientation. The larger hole should be placed at the patient's upper left.

NOTE: The 7X10-AP graft should be oriented such that the wide arm reinforces the posterior space while the small arm reinforces the anterior space.

- Complete the standard surgical procedure.

- Discard any unused portions of the graft according to institutional guidelines for disposal of medical waste.

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