Biodesign®

HIATAL HERNIA GRAFT

FP0032-03N

	Batch code		Do not use if package is damaged and consult instructions for use
	Catalogue number		Keep dry
	Consult instructions for use		Manufacturer
	Contains biological material of animal origin		Medical device
	Country of manufacture		Double sterile barrier system
	Date of manufacture		Sterilized using ethylene oxide
	Do not re-sterilize		Unique device identifier
	Do not re-use		Use-by date

MANUFACTURER
COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.

BIODESIGN® HIATAL HERNIA GRAFT

DESCRIPTION

The Biodesign® Hiatal Hernia Graft is a multi-layered sheet of extracellular collagen matrix derived from porcine small intestinal submucosa. The Biodesign Hiatal Hernia Graft is implanted by a trained medical professional to support weakened soft tissue during surgery to repair hiatal hernias.

INTENDED USE

The Biodesign Hiatal Hernia Graft is used for implantation to reinforce soft tissue where weakness exists, including repair of hiatal hernias. This graft is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

MR This symbol means the following: Magnetic Resonance Safe.

HIATAL HERNIA GRAFT This symbol means the following: Hiatal Hernia Graft.

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

The graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine material.

PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to graft failure and/or transmission of disease.
- Do not resterilize.** Discard all open and unused portions.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that the graft is hydrated prior to cutting or suturing. Insufficient hydration may cause device damage during laparoscopic delivery.
- Ensure that all layers of the graft are secured when suturing.
- Ensure that the graft does not touch the esophagus.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials.

- Adhesion
- Allergic reaction
- Fistula formation
- Graft detachment/migration
- Hematoma
- Infection
- Inflammation
- Recurrence of tissue defect
- Seroma formation

The following complications are associated with use of graft materials in hiatal hernia repair.

- Dysphagia
- Esophageal erosion
- Perforation
- Stricture formation

STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

STERILIZATION

The graft has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

Required Materials

- Sterile basin
- Sterile smooth forceps
- Hydration fluid: room temperature sterile saline or lactated Ringer's solution

NOTE: Handle the graft using aseptic technique, minimizing contact with latex gloves.

Procedure

- Remove the packaging containing the graft from the envelope.
- Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
- Place the graft into the sterile basin on the sterile field. Multiple grafts may be hydrated simultaneously in the same basin.
- Place enough hydration fluid in the basin to completely cover the graft(s).

- Hydrate the graft in the fluid until the desired handling characteristics are achieved. A hydration time of greater than one minute is not required.
- Prepare the graft site using standard surgical techniques.

NOTE: Suturing the graft with close tissue approximation produces better outcomes. Maximum graft-tissue approximation is only achieved with complete crural closure.

- If necessary, trim the graft to the patient's anatomy, providing a small allowance for overlap with adjacent tissues.


NOTE: The graft may be trimmed prior to hydration to fit the patient's anatomy. Be sure to hydrate the graft prior to suturing it into place. See step 5.

NOTE: For laparoscopic placement, hydrate the graft prior to loading it into the delivery port. See step 5.


- Suture the graft into place, avoiding excess tension.

NOTE: Clinical literature suggests the use of tacks or staples at the diaphragm may be associated with cardiovascular complications.¹

NOTE: The 7X10-U graft has a larger hole in the corner to identify the graft orientation. The larger hole should be placed at the patient's upper left.

NOTE: The 7X10-AP graft should be oriented such that the wide arm reinforces the posterior space while the small arm reinforces the anterior space.

- Complete the standard surgical procedure.

- Discard any unused portions of the graft according to institutional guidelines for disposal of medical waste.

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