Manufacturer
Fabricante
Fabricant
Fabricante

Temperature limit
Limite de temperatura
Limite de température
Limite de temperatura

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BIODESIGN® HERNIA GRAFT

INTENDED USE

The Cook® Biodesign® Hernia Graft is intended for implantation to reinforce soft tissues where weakness exists. Indications for use include the repair of a hernia or body wall defect. The graft is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

HERNIA GRAFT This symbol means the following: Hernia Graft

This graft is intended for use by trained medical professionals.

CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to graft failure and/or transmission of disease.
- Do not resterilize.** Discard all open and unused portions of the graft.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that the graft is hydrated prior to cutting, suturing, stapling, tacking or loading the graft laparoscopically.
- Ensure that all layers of the graft are secured when suturing, stapling, or tacking.
- Place the graft in maximum possible contact with healthy, well-vascularized tissue to encourage cell ingrowth and tissue remodeling.
- Suturing, stapling, or tacking more than one graft together may decrease graft performance.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the graft.
- Extended hydration or excessive handling could lead to partial delamination of superficial layers of the graft.
- Care should be taken to avoid damage to the graft when loading laparoscopically. Loading the graft through a 10 mm or larger port is recommended.
- If wound is left open, keep graft moist to prevent dryness.

POTENTIAL COMPLICATIONS

Possible adverse reactions with the use of any prosthesis may include, but are not limited to:

- adhesion
- inflammation
- bowel erosion
- premature degradation
- fistula formation
- recurrence of tissue defect
- hematoma
- seroma formation
- infection

Complications, such as delayed wound infection, premature degradation, hernia recurrence, bowel erosion, and the need for re-operation, should be reasonably expected in patients who are critically ill or who have severely contaminated abdomens.

STORAGE

Store the graft in a clean, dry location at room temperature.

STERILIZATION

This graft has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

Required Materials

- Sterile forceps
- Hydration fluid: room temperature, sterile saline or sterile lactated Ringer's solution.

NOTE: Always handle the graft using aseptic technique, minimizing contact with latex gloves.

- Remove the packaging containing the graft from the envelope.
- Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
- Hydrate the graft in room temperature, sterile saline or sterile lactated Ringer's solution until the desired handling characteristics are achieved. A hydration time of greater than one minute is not required.
- Prepare the hernia repair site using standard surgical techniques.

GRAFT PLACEMENT

- Preperitoneal/retractorctus (Rives and Stoppa technique) or intraperitoneal placement (underlay) is recommended for open procedures.
- Intraperitoneal placement (underlay) is recommended for a laparoscopic approach.
- Trim the graft to fit the site, providing an allowance for overlap. Fundamental surgical principles suggest recurrence can be minimized if the graft overlaps the surrounding tissue by at least 4-5 cm in all directions.¹

NOTE: If the graft is cut too small for the defect, excess tension may be placed on the suture line. This can result in hernia recurrence or development of a defect in the adjacent tissues.

- Suture, staple, or tack the graft into place, avoiding excess tension.

NOTE: Excessive handling may lead to partial delamination of a superficial layer. This can be managed by securing the detached portion to adjacent tissue with sutures or staples.

NOTE: Surgical experience indicates that suturing or stapling the graft with close tissue approximation produces better outcomes. Leave less than or equal to 3 cm between sutures with bite depth of 1 cm. Use permanent or long-term absorbable sutures or tacks according to surgeon preference.

NOTE: Interrupted sutures can provide additional security against recurrence of tissue defect in the event of suture failure.

NOTE: The Hernia Graft has been designed with an open matrix structure to optimize tissue incorporation and angiogenesis.

NOTE: Minimizing the bridging of defects may improve outcomes.

- Complete the surgical procedure.
- Discard any unused portions of the graft according to institutional guidelines for medical waste.
- Place closed suction drains for 2-6 weeks. Remove when output is less than 20 mL/24 hours for at least two (2) consecutive days or until drain is dry.

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