# Biodesign<sup>®</sup>

## HERNIA GRAFT

FP0036-02Q

	<b>Batch code</b>		<b>Do not use if package is damaged and consult instructions for use</b>
	<b>Catalogue number</b>		<b>Keep dry</b>
	<b>Consult instructions for use</b>		<b>Manufacturer</b>
	<b>Contains biological material of animal origin</b>		<b>Medical device</b>
	<b>Country of manufacture</b>		<b>Double sterile barrier system</b>
	<b>Date of manufacture</b>		<b>Sterilized using ethylene oxide</b>
	<b>Do not re-sterilize</b>		<b>Unique device identifier</b>
	<b>Do not re-use</b>		<b>Use-by-date</b>

**MANUFACTURER**  
COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.

www.cookbiotech.com  
© COOK BIOTECH 2024

January 2024

## BIODESIGN<sup>®</sup> HERNIA GRAFT

### INTENDED USE

The Biodesign<sup>®</sup> Hernia Graft is intended for implantation to reinforce soft tissues where weakness exists. Indications for use include the repair of a hernia or body wall defect. The graft is supplied sterile and is intended for one-time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**MR** This symbol means the following: Magnetic Resonance Safe.

**HERNIA GRAFT** This symbol means the following: Hernia Graft

This graft is intended for use by trained medical professionals.

### CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

### PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to graft failure and/or transmission of disease.
- Do not re-sterilize.** Discard all open and unused portions of the graft.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that the graft is hydrated prior to cutting, suturing, stapling, tacking or loading the graft laparoscopically.
- Ensure that all layers of the graft are secured when suturing, stapling, or tacking.
- Place the graft in maximum possible contact with healthy, well-vascularized tissue to encourage cell ingrowth and tissue remodeling.
- Suturing, stapling, or tacking more than one graft together may decrease graft performance.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the graft.
- Extended hydration or excessive handling could lead to partial delamination of superficial layers of the graft.
- Care should be taken to avoid damage to the graft when loading laparoscopically. Loading the graft through a 10 mm or larger port is recommended.
- If wound is left open, keep graft moist to prevent dryness.

### POTENTIAL COMPLICATIONS

Possible adverse reactions with the use of any prosthesis may include, but are not limited to:

- adhesion
- inflammation
- bowel erosion
- premature degradation
- fibula formation
- recurrence of tissue defect
- hematoma
- seroma formation
- infection

Complications, such as delayed wound infection, premature degradation, hernia recurrence, bowel erosion, and the need for re-operation, should be reasonably expected in patients who are critically ill or who have severely contaminated abdomens.

### STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

### STERILIZATION

This graft has been sterilized with ethylene oxide.

### INSTRUCTIONS FOR USE

#### Required Materials

- Sterile forceps
- Hydration fluid: room temperature, sterile saline or sterile lactated Ringer's solution.

**NOTE: Always handle the graft using aseptic technique, minimizing contact with latex gloves.**

#### Procedure

- Remove the packaging containing the graft from the envelope.
- Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
- Hydrate the graft in room temperature, sterile saline or sterile lactated Ringer's solution until the desired handling characteristics are achieved. A hydration time of greater than one minute is not required.
- Prepare the hernia repair site using standard surgical techniques.

#### GRAFT PLACEMENT

- Preperitoneal/retrorectus (Rives and Stoppa technique) or intraperitoneal placement (underlay) is recommended for open procedures.
- Intraperitoneal placement (underlay) is recommended for a laparoscopic approach.

- Trim the graft to fit the site, providing an allowance for overlap. Fundamental surgical principles suggest recurrence can be minimized if the graft overlaps the surrounding tissue by at least 4-5 cm in all directions.!

**NOTE:** If the graft is cut too small for the defect, excess tension may be placed on the suture line. This can result in hernia recurrence or development of a defect in the adjacent tissues.

- Suture, staple, or tack the graft into place, avoiding excess tension.

**NOTE:** Excessive handling may lead to partial delamination of a superficial layer. This can be managed by securing the detached portion to adjacent tissue with sutures or staples.

**NOTE:** Surgical experience indicates that suturing or stapling the graft with close tissue approximation produces better outcomes. Leave less than or equal to 3 cm between sutures with bite depth of 1 cm. Use permanent or long-term absorbable sutures or tacks according to surgeon preference.

**NOTE:** Interrupted sutures can provide additional security against recurrence of tissue defect in the event of suture failure.

**NOTE:** The Hernia Graft has been designed with an open matrix structure to optimize tissue incorporation and angiogenesis.

**NOTE:** Minimizing the bridging of defects may improve outcomes.

- Complete the surgical procedure.
- Discard any unused portions of the graft according to institutional guidelines for medical waste.
- Place closed suction drains for 2-6 weeks. Remove when output is less than 20 mL/24 hours for at least two (2) consecutive days or until drain is dry.

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