COOK[®]
BIOTECH

EN
3

**Biodesign[®]
Anal Fistula Plug**

Instructions for Use

ES
5

**Tapón para fístulas anales
Biodesign[®]**

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FR
8

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
BIODESIGN® ANAL FISTULA PLUG

INTENDED USE

The Biodesign® Anal Fistula Plug is for implantation to reinforce soft tissue where a rolled configuration is required, for repair of anorectal fistulas. The plug is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

MR Magnetic resonance safe

ANAL FISTULA PLUG This symbol means the following: Anal Fistula Plug

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

- The plug is derived from a porcine source and should not be used in patients sensitive to porcine materials.
- **Not for vascular use.**

PRECAUTIONS

- The plug is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- The plug is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the plug if mishandling has caused possible damage or contamination, or if the plug is past its expiration date.
- Do not implant the plug in a grossly infected or abscessed fistula tract.
- **A seton should be used until there is no evidence of acute inflammation, purulence, or excessive drainage. Allow the tract to mature and stabilize for six to eight weeks before placing the plug.**
- Ensure that the plug is hydrated prior to placement, cutting, or suturing.
- Expect some drainage for two to four weeks. Drainage can continue for up to 12 weeks after the procedure as the plug is incorporated and the fistula tract closes.

POTENTIAL COMPLICATIONS

Complications that can occur with the plug include:

- | | | |
|-----------------------------------|----------------------|-------------|
| • Abscess | • Extrusion | • Migration |
| • Allergic reaction | • Fistula recurrence | • Redness |
| • Bleeding | • Induration | • Sepsis |
| • Delayed or failed incorporation | • Infection | • Swelling |
| | • Inflammation | |

STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

STERILIZATION

This plug has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

Required Materials

- Sterile basin
- Hydration fluid: room temperature sterile saline or sterile lactated Ringer's solution.
- Resorbable suture, such as: 0 chromic, 2-0 or 0 coated polyglycolic acid suture (coated PGA).
- 10 cc syringe
- Brush for tract debridement
- A 5 Fr catheter
- Hydrogen peroxide or saline

NOTE: Always handle the plug using aseptic technique, minimizing contact with latex gloves.

PREPARATION

1. Remove the packaging containing the plug from the box.
2. Using aseptic technique, open the outer pouch and pass the inner pouch containing the plug onto the sterile field.
3. Open the inner pouch and remove the tray containing the plug. Open the tray and place the plug into a sterile basin.
4. Add enough hydration fluid to the basin to fully submerge the plug. Allow the plug to hydrate until the desired handling characteristics are achieved. Do not hydrate the plug for more than one minute.
5. Remove the hydrated plug from the basin using forceps.

6. Attach a suitable resorbable suture (approximately 30 cm in length) around the tail (narrow end) of the plug to use as a tether to pull the plug into the fistula tract.
7. Prepare the patient using standard surgical techniques appropriate for anal fistula repair.

NOTE: The recommended practice for preoperative bowel preparation in elective colorectal surgery includes mechanical bowel cleansing through the use of enemas and cathartic agents, and administration of prophylactic antimicrobial agents, oral, or intravenous. Insufficient cleansing or inadequate antibacterial prophylaxis can predispose the patient to infections.¹⁻³

PROCEDURE

1. Perform under local, regional, or general anesthesia.
2. If a seton is not in place, identify the internal (rectal) fistula opening.

NOTE: Failure to locate the internal opening may lead to persistence of the fistula. If the internal opening cannot be reliably identified, an alternative method of treatment should be considered.
3. Insert a fistula probe, seton, or suitable instrument through the fistula tract, entering through the external opening and exiting via the internal opening.
4. Attach a brush to the instrument at the internal opening.
5. If a seton is in place, cut the seton and attach a brush or similar instrument to the seton at the internal (rectal) opening.
6. Pull the brush into the tract and gently debride to remove non-vascularized tissue. Do not enlarge the tract.


NOTE: Failure to debride the tract may increase infection risk.

NOTE: Blood should be noted on the bristles of the brush and at both the internal and external fistula openings after adequate debridement.
7. Flush the tract thoroughly with hydrogen peroxide or sterile saline using the 10 cc syringe and 5 Fr catheter.
8. Use the tether to pull the plug into the internal opening and through the fistula tract. Ensure the plug is drawn completely into the tract.


IMPORTANT: The internal opening is the high-pressure zone of the fistula. The thicker end of the plug must be securely in the internal opening to prevent ingress of fecal debris. The higher pressures within the rectum and anal canal assist in maintaining the plug in the fistula tract.
9. Once the plug is properly positioned, trim away and discard any remaining portion of the plug that is not implanted within the fistula tract.

IMPORTANT: Trim the plug at the level of the bowel wall at the internal opening to minimize contact with bowel contents.
10. Using a suitable resorbable suture, secure the internal end of the plug to the adjacent tissue, obtaining adequate bites of sphincter muscle.

NOTE: Inadequate fixation of the plug to the deep tissue layers may result in early fistula recurrence.

NOTE: No part of the plug should be visible at the internal opening.

NOTE: If the internal opening is dimpled or recessed, consider limited mobilization of the mucosal edges prior to suture placement to ensure adequate coverage of the plug.

11. Trim any excess external plug flush with the skin. No fixation is required. The external opening may be enlarged slightly to facilitate drainage.

NOTE: Do not close the external opening when suturing the plug to tissue. This is to allow continued drainage of the fistula tract. Complete obstruction of the external opening may result in accumulation of fluid, infection, or abscess.
12. Place a sterile dressing over the external opening.

NOTE: Expect continued drainage of the fistula tract for up to 12 weeks after plug placement.

POST-OPERATIVE CARE

To provide the best environment for tissue integration into the plug, patient activity should be minimized. Provide patients with a list of post-procedure care recommendations. The following patient guidelines should be considered.

1. Patients should avoid any strenuous physical activity beyond a gentle walk for at least six weeks.
2. Patients should avoid any heavy lifting over 10 lbs (5 kg) for at least six weeks.
3. Patients should abstain from sexual intercourse and other forms of vaginal/rectal insertion (e.g. tampons) for at least six weeks.
4. Patients should bathe standing up to keep the area clean.
5. Patients should observe a liquid diet for the initial two days post-procedure.
6. Patients should begin a high fiber diet after two days post-procedure.
7. Patients should use a stool softener for at least two weeks post-procedure.
8. Patients should take over-the-counter pain medication as needed.
9. Patients should avoid the use of steroid creams.
10. Patients should expect drainage for two to four weeks. Drainage can continue for up to twelve weeks after the procedure as the plug incorporates and the fistula tract closes.

REFERENCES

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
3. Berrios-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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RÉFÉRENCES

- Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
- Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
- Berrios-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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REFERÊNCIAS

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
3. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

Batch code

Catalogue number

Consult instructions for use

Contains biological material
of animal origin

Country of manufacture

Date of manufacture

Do not resterilize

Do not reuse

Do not use if package is damaged
and consult instructions for use

Keep dry

Manufacturer

Medical device

Double sterile barrier system

Sterile

Sterilized using ethylene oxide

Unique device identifier

Use-by date

MANUFACTURER

COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.