
COOK®
MEDICAL

CE 0123

EN 3	Biodesign® Peyronie's Repair Graft Instructions for Use
CS 4	Biodesign® Štěp pro opravu Peyronieovy choroby Pokyny k použití
DA 6	Biodesign® Peyronie's ophelende transplantat Brugsanvisning
DE 7	Biodesign® Peyronie-Reparaturimplantat Gebrauchsanweisung
EL 9	Biodesign® Μόσχευμα αποκατάστασης της νόσου του Peyronie Οδηγίες χρήσης
ES 11	Biodesign® Injerto de reparación de Peyronie de Instrucciones de uso
FR 12	Biodesign® Matériau de greffe pour réfection de Peyronie Mode d'emploi
HU 14	Biodesign® graft a Peyronie-kór kezelésére Használati utasítás
IT 16	Biodesign® Innesto per la riparazione della sindrome di Peyronie Istruzioni per l'uso
NL 17	Biodesign® Peyronies hersteltransplantaat Gebruiksaanwijzing
NO 19	Biodesign® Peyronies vevsgraft Bruksanvisning
PT 21	Biodesign® Enxerto para a correção da doença de Peyronie Instruções de utilização
SV 22	Biodesign® Peyronies reparationstransplantat Bruksanvisning

FP0046-01H

BIODESIGN® PEYRONIE'S REPAIR GRAFT

INTENDED USE

The Cook® Biodesign® Peyronie's Repair Graft is intended for implantation to reinforce soft tissue where weakness exists in the urological anatomy, including but not limited to urethroplasty, repair of tunica albuginea defects, and reinforcement in the repair of Peyronie's disease. The graft is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

PEYRONIE'S REPAIR GRAFT This symbol means the following: Peyronie's Repair Graft

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used for patients with known sensitivity to porcine materials.

PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- Graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard graft if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Ensure that the graft is rehydrated prior to suturing.
- Graft performance has not been evaluated with suture spacing greater than 2 mm.
- Ensure that all layers of the graft are secured when suturing.

POTENTIAL COMPLICATIONS

The following complications are possible when surgical graft materials are used in the penis. If any of these conditions occur and cannot be managed medically, careful graft removal should be considered:

- Infection
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

NOTE: Graft failure resulting in a lack of soft tissue reinforcement may lead to reoccurrence of the Peyronies plaque. If this occurs, additional surgical intervention to remove the plaque may be necessary.

NOTE: *De novo* post-operative erectile dysfunction is a known complication associated with the use of graft materials in the penis. If post-operative erectile dysfunction occurs, graft removal may exacerbate erectile changes. Therefore, graft removal is not recommended. Alternative medical and surgical treatments should be considered.

STORAGE

This graft should be stored in a clean, dry location at room temperature.

STERILIZATION

This graft has been sterilized with ethylene oxide.

USE OF ANTIMICROBIALS

Because the graft is at times used in surgical fields where sterility cannot be assured, the use of antimicrobials is common practice and may prevent infectious complications.¹ In these cases both antibiotic prophylaxis of the patient and antimicrobial soaking of the graft have been used. Typical flora can be expected to include a variety of aerobic and facultative anaerobic organisms, including, but not limited to, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Pseudomonas aeruginosa*, and *Escherichia coli*. Therefore the following points should be considered:

- Antimicrobials, if used topically or systemically, should provide coverage against a wide spectrum of aerobic and anaerobic organisms.
- Antibacterial prophylaxis, if chosen, should be started prior to surgery and continued post-operatively.¹

The presence of certain antimicrobials may inhibit revascularization and/or infiltration of cells into the graft. For example, gentamicin is known to hinder neovascularization, epithelialization, and keratinocyte growth,²⁻⁴ while povidone iodine,⁵ bacitracin,^{3,5} polymyxin B,⁶ and vancomycin⁷ have all been reported to slow or inhibit wound healing. However, no studies have been conducted to evaluate the combination of antimicrobials with respect to graft placement.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede the institutional protocols or professional clinical judgment concerning patient care.

NOTE: Handle graft sheets using aseptic technique, minimizing contact with latex gloves.

REQUIRED MATERIALS

- A sterile dish (kidney dish or other bowl)
 - Sterile forceps
 - Rehydration fluid: at least 50 mL of room temperature sterile saline or sterile lactated Ringer's solution for each graft sheet.
1. Using aseptic technique, remove the graft inner pouch from its outer bag, and place the inner pouch in the sterile field.
 2. Open the inner pouch carefully, and aseptically remove the graft sheet with the sterile forceps.
 3. Place the graft sheet into the sterile dish in the sterile field. (Multiple graft sheets may be rehydrated simultaneously in the same dish.)
 4. Add enough rehydration fluid to the dish to immerse the graft sheet.
 5. Allow graft sheet to rehydrate for at least ten minutes.
 6. Prepare the graft site using standard surgical techniques.
 7. Using aseptic technique, trim the graft sheet to fit the site. User feedback suggests that the graft should be sized approximately 30% larger than the site.

NOTE: An alternative method is to cut the graft sheet to size prior to rehydration. If this method is selected, be sure to rehydrate the graft sheet prior to suturing it into place. See step 5.

8. Using aseptic technique, transfer the graft sheet to the graft site and suture into place using resorbable suture, while avoiding excess tension.

NOTE: Surgical experience indicates that suturing graft sheets with close tissue approximation produces better outcomes.

9. Complete the standard surgical procedure.
10. Discard any unused portions of the graft sheet following standard technique for disposal of medical waste.

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West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.