COOK®
MEDICAL

CE 0123

**EN
3**

**Biodesign®
Fistula Plug**

Instructions for Use

**DA
8**

**Biodesign®
fistelprop**

Brugsanvisning

**DE
13**

**Biodesign®
Fistelverschluss**

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**EL
18**

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
BIODESIGN® FISTULA PLUG

INTENDED USE

The Biodesign® Fistula Plug is for implantation to reinforce soft tissue for repair of recto-vaginal or anorectal fistulas. The device is supplied sterile and is intended for one time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

- This device is derived from a porcine source and should not be used for patients sensitive to porcine materials.
- **Not for vascular use.**

PRECAUTIONS

- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- Device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard device if mishandling has caused possible damage or contamination or if the device is past its expiration date.
- Do not implant the device in a grossly infected or abscessed fistula tract.
- **In fistula cases involving evidence of acute inflammation, purulence, or excessive drainage, a draining seton should be used to allow the tract to mature and stabilize for six to eight weeks before placing the plug.**
- Ensure that the device is rehydrated prior to placement, cutting, or suturing.
- Placement of the plug in tracts less than 1 cm in length can result in incomplete incorporation and/or expulsion of the device.

GENERAL

- Users should be familiar with surgical technique for recto-vaginal or anorectal fistula repair.
- Users should exercise good surgical practice for the management of clean-contaminated, contaminated or infected fields.
- The potential for infection of the graft material following implantation may be reduced by the use of prophylactic antibiotics and cleaning of the fistula tract. **(See Use of Antimicrobials)**
- The Fistula Plug should always be introduced through the rectal (primary) fistula opening.
- The Fistula Plug should be drawn into the fistula tract completely, until the button is set against the rectal wall.
- The external/vaginal (secondary) fistula opening should remain open to allow drainage to occur.
- **IMPORTANT:** It is recommended that users counsel patients on abstaining from strenuous physical activity for at least 2 weeks following anorectal fistula repair, or 6 weeks following recto-vaginal fistula repair. Refer to the section on Post-Operative Care.

POTENTIAL COMPLICATIONS

Complications that can occur with the Fistula Plug include, but are not limited to:

- Inflammation
- Induration
- Migration
- Erosion
- Extrusion
- Seroma formation
- Infection
- Abscess
- Fistula recurrence
- Delayed or failed incorporation of the device

If any of the following conditions occur and cannot be resolved, device removal should be considered:

- Infection
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

STORAGE

This device should be stored in a clean, dry location at room temperature.

STERILIZATION

This device has been sterilized with ethylene oxide and should not be resterilized.

INSTRUCTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

NOTE: Handle device using aseptic technique. Minimize contact with latex gloves.

REQUIRED MATERIALS

- A sterile dish (kidney dish or other bowl)
- Rehydration fluid: at least 200ml of room temperature sterile saline or sterile lactated Ringer's solution
- Suitable resorbable suture, such as: 2-0 or 0 polydioxanone (PDO or PDS) or coated polyglycolic acid suture (coated PGA)
- A short 5Fr catheter
- Hydrogen peroxide or saline for flushing
- 10cc Syringe for flushing
- Cook Fistula Brush for tract debridement

PREPARATION

1. Select the proper Fistula Plug size according to PROCEDURAL STEP 6.
2. Remove the packaging containing the Fistula Plug from the box.
3. Remove the inner pouch containing the Fistula Plug from the outer package using aseptic technique. Place the inner pouch into the sterile field.
4. Using sterile gloved hands carefully open the inner pouch and remove the tray containing the Fistula Plug. Place the tray into the sterile field.
5. Open the tray and place the Fistula Plug into a sterile dish using aseptic technique. Add enough rehydration fluid to the sterile dish to submerge the Fistula Plug. Allow the Fistula Plug to rehydrate in the sterile dish until the desired handling characteristics are achieved, but do not allow the plug to rehydrate any longer than 2 minutes.
6. Carefully remove the hydrated Fistula Plug from the sterile dish using forceps.
7. Attach a suitable resorbable suture (approximately 30cm in length) around the narrow end of the Fistula Plug to use as a tether for pulling the plug into the fistula tract.

NOTE: The recommended practice for preoperative bowel preparation in elective colorectal surgery includes mechanical bowel cleansing through the use of enemas and cathartic agents, and administration of prophylactic oral or intravenous antimicrobial agents. Insufficient cleansing or inadequate antibacterial prophylaxis can predispose the patient to infections.¹⁻³ (See Use of Antimicrobials)

PROCEDURAL

1. Perform under local, regional, or general anesthesia.
2. If a seton is not already in place, identify the rectal (primary) fistula opening by inserting a sterile probe into the secondary opening and navigating it through the fistula tract. Alternatively, inject an appropriate sterile fluid (hydrogen peroxide or saline) into the secondary opening of the fistula tract to identify the site of emergence at the primary opening.
NOTE: Failure to locate the rectal opening may lead to persistence of the fistula. If the rectal opening cannot be reliably identified, an alternative method of treatment should be considered.
3. If a seton is in place, cut the seton and attach the Cook Fistula Brush to the seton at the rectal (primary) opening.
4. Pull the Cook Fistula Brush into the fistula tract using a suitable suture and gently clean/debride the fistula tract to remove non-vascularized tissue being careful not to enlarge the tract.
NOTE: Bleeding should be noted on the bristles of the Cook Fistula Brush and at both the primary and secondary openings after adequate debridement.
5. Flush the tract as thoroughly as possible with hydrogen peroxide or sterile saline using the 10cc syringe and 5Fr catheter.
6. Use a sterile probe or other appropriate measuring instrument to approximate the diameter of the fistula tract at the primary opening. Select the proper size of Fistula Plug according to the table below.

Primary Opening Diameter (mm)	Fistula Plug size (cm)
>1 to 2	0.2 (2mm)
>2 to 4	0.4 (4mm)
>4 to 7	0.7 (7mm)

7. Prepare appropriately sized plug according to instructions in the preparation section.

8. Attach the tether connected to the Fistula Plug to the end of the Cook Fistula Brush located at the primary opening.


IMPORTANT: The button of the Fistula Plug should closely approximate the rectal wall or be fully implanted below the mucosa in order to minimize contact of the Fistula Plug material with bowel contents.


IMPORTANT: The rectal opening is the high-pressure zone of the fistula, as well as the site of ingress of fecal debris. The button end of the plug must therefore be securely snug against the rectal opening or fully implanted below the mucosa to prevent ingress of fecal debris. In addition, the higher pressures within the rectum and anal canal assist in maintaining the plug in the fistula tract by simple mechanical force.

Superficial Button Placement

9. Draw the tether attached to the Fistula Plug into the rectal (primary) opening and through the fistula tract until the Fistula Plug is drawn completely into the fistula tract and the Biodesign button is against the rectal wall.


Recto-vaginal fistula repair


Anorectal fistula repair

10. When the Fistula Plug is properly positioned, grip the tether used to draw the Fistula Plug into the fistula tract and apply gentle tension to it. This will hold the Fistula Plug in place during fixation of the device on the rectal side. Refer to the illustration below.


11. Using a 2-0 PDO, PDS, coated PGA or comparable suture, secure the button to adjacent tissue by stitching through the center portion of the button and directly through the rectal wall and mucosa. Obtain an adequate bite of sphincter muscle with the suture and secure the button to the tissue by tying the two suture ends together. This should be repeated four times and the sutures should be equally spaced around the button.

Superficial Button Placement


Submucosal Button Placement

12. When implanting the Biodesign button under the rectal mucosa, mobilize a mucosal flap or circumferentially dissect a pocket of sufficient size to accommodate the Biodesign button.
13. Draw the tether attached to the plug into the rectal (primary) opening and through the fistula tract until the plug is drawn completely into the fistula tract and the Biodesign button is underneath the mobilized mucosal flap or dissected pocket.
14. Suture the plug according to PROCEDURAL STEP 11.

Submucosal Button Placement

15. Reapproximate the mobilized mucosal flap or dissected pocket and suture the tissue securely over the Biodesign button after the button has been secured to the adjacent tissue using a 2-0 PDO, PDS, coated PGA or comparable suture.

Submucosal Button Placement



16. If the secondary opening is external (skin level), a sterile dressing may be placed to prevent soiling.

IMPORTANT: Do not close the secondary opening when suturing the plug to tissue. This is to allow continued drainage of the fistula tract. If necessary, consider widening the secondary opening to facilitate drainage. Complete obstruction of the fistula opening may result in accumulation of fluid, infection, or abscess.

NOTE: Continued drainage of the fistula tract can be expected to last several days up to 12 weeks after plug placement.


OPTIONAL FIXATION FOR RECTO-VAGINAL FISTULAS ONLY

17. Using a 2-0 coated PGA or comparable suture, place a minimum of two stitches through the middle of the Fistula Plug at the vaginal fistula opening. Obtain an adequate bite of both the plug and tissue with the suture, allowing 2mm between each stitch.

IMPORTANT: Do not close the secondary opening completely when suturing the plug to tissue. This is to allow continued drainage of the fistula tract. Complete obstruction of the fistula opening may result in accumulation of fluid, infection, or abscess.

18. When the Fistula Plug is properly secured, trim away and discard any remaining portion of the plug that is not implanted within the fistula tract. Leave at least 2mm of material between the end of the trimmed Fistula Plug and suture to provide the appropriate suture retention strength.

19. Discard any unused portions according to institutional guidelines for medical waste.

POST-OPERATIVE CARE

To provide the best environment for tissue integration into the Fistula Plug, patient activity should be minimized. Provide patients with a list of post-procedural care recommendations. The following patient guidelines should be considered:

1. Patients should avoid any strenuous physical activity beyond a gentle walk for at least 2 weeks following anorectal fistula repair, or 6 weeks following recto-vaginal fistula repair.
2. Patients should avoid any heavy lifting over 10lbs (5kg) for at least two weeks following anorectal fistula repair, or at least 6 weeks following recto-vaginal fistula repair.
3. Patients should use a stool softener for at least 2 weeks after surgery.
4. Patients should observe a liquid diet for the first 2 days after surgery.
5. Patients should begin a high fiber diet 2 days after surgery.
6. Patients should abstain from sexual intercourse and other forms of vaginal/rectal insertion (e.g. tampons) for at least 2 weeks following anorectal fistula repair, or 6 weeks following recto-vaginal fistula repair.

USE OF ANTIMICROBIALS WITH THE FISTULA PLUG

Because the Fistula Plug is used in surgical fields where sterility cannot be assured, the use of antimicrobials is common practice and may prevent infectious complications.¹⁻³ Both mechanical bowel cleansing and antibiotic prophylaxis of the patient have been used successfully, and the U.S. Centers for Disease Control recommends bowel cleansing and administration of prophylactic enteric and intravenous antimicrobial agents before elective colorectal operations.³ (See reference for specific dose, timing, and drug choices.)

Typical gastrointestinal flora can be expected to include a variety of aerobic and anaerobic organisms. Therefore the following points should be considered:

- Antimicrobials, if used topically or systemically, should provide coverage against a wide spectrum of aerobic and anaerobic organisms.⁴
- Mechanically prepare the colon by appropriate use of enemas and cathartic agents.
- A dose of prophylactic antimicrobial agent should be given intravenously and timed such that a bactericidal concentration of the drug is established in serum and tissues when the device is implanted.
- Maintain therapeutic levels of the agent in serum and tissues throughout the operation.³

The presence of certain antimicrobials may inhibit revascularization and/or infiltration of cells into the repairing tissue.⁵⁻⁷ For example, gentamicin is known to hinder neovascularization, epithelialization, and keratinocyte growth,⁶ while povidone iodine,⁸ bacitracin,^{5,8} polymyxin B,⁹ and vancomycin¹⁰ have all been reported to slow or inhibit wound healing. Careful consideration is required before using any antimicrobial or antiseptic (topical or systemic) that has not been proven compatible with surgical implantation and wound healing. However, no studies have been conducted to evaluate the combination of antimicrobials with Fistula Plug placement.

REFERENCES

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis* 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci* 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. Guideline for Prevention of Surgical Site Infection, 1999. Centers for Disease Control and Prevention (CDC).

4. Aldridge KE, et al. Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of Bacteroides fragilis group, Prevotella, Fusobacterium, Porphyromonas, and Peptostreptococcus species. *Antimicrob Agents Chemother* 2001;45(4):1238-43.
5. Petroutsos G, et al. Antibiotics and corneal epithelial wound healing. *Arch Ophthalmol* 1983;101(11):1775-8.
6. Bang K, et al. Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice. *Scand J Plast Reconstr Surg Hand Surg* 1998;32(2):129-33.
7. Nelson JD, et al. Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics. *Curr Eye Res* 1990;9(3):277-85.
8. Kjolseth D, et al. Comparison of the effects of commonly used wound agents on epithelialization and neovascularization. *J Am Coll Surg* 1994;179(3):305-12.
9. Nakamura M, et al. Effects of antimicrobials on corneal epithelial migration. *Curr Eye Res* 1993;12(8):733-40.
10. Petroutsos G, Guimaraes R, and Pouliquen Y. The effect of concentrated antibiotics on the rabbit's corneal epithelium. *Int Ophthalmol* 1984;7(2):65-9.

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LITTERATUR

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis* 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci* 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. Guideline for Prevention of Surgical Site Infection, 1999. Centers for Disease Control and Prevention

(CDC).

4. Aldridge KE, et al. Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of Bacteroides fragilis group, Prevotella, Fusobacterium, Porphyromonas, and Peptostreptococcus species. *Antimicrob Agents Chemother* 2001;45(4):1238-43.
5. Petroustos G, et al. Antibiotics and corneal epithelial wound healing. *Arch Ophthalmol* 1983;101(11):1775-8.
6. Bang K, et al. Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice. *Scand J Plast Reconstr Surg Hand Surg* 1998;32(2):129-33.
7. Nelson JD, et al. Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics. *Curr Eye Res* 1990;9(3):277-85.
8. Kjolseth D, et al. Comparison of the effects of commonly used wound agents on epithelialization and neovascularization. *J Am Coll Surg* 1994;179(3):305-12.
9. Nakamura M, et al. Effects of antimicrobials on corneal epithelial migration. *Curr Eye Res* 1993;12(8):733-40.
10. Petroustos G, Guimaraes R, and Pouliquen Y. The effect of concentrated antibiotics on the rabbit's corneal epithelium. *Int Ophthalmol* 1984;7(2):65-9.

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QUELLEN

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis* 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci* 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. Guideline for Prevention of Surgical Site Infection, 1999. Centers for Disease Control and Prevention (CDC).
4. Aldridge KE, et al. Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of *Bacteroides fragilis* group, *Prevotella*, *Fusobacterium*, *Porphyromonas*, and *Peptostreptococcus* species. *Antimicrob Agents Chemother* 2001;45(4):1238-43.
5. Petroutsos G, et al. Antibiotics and corneal epithelial wound healing. *Arch Ophthalmol* 1983;101(11):1775-8.
6. Bang K, et al. Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice. *Scand J Plast Reconstr Surg Hand Surg* 1998;32(2):129-33.
7. Nelson JD, et al. Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics. *Curr Eye Res* 1990;9(3):277-85.
8. Kjolseth D, et al. Comparison of the effects of commonly used wound agents on epithelialization and neovascularization. *J Am Coll Surg* 1994;179(3):305-12.
9. Nakamura M, et al. Effects of antimicrobials on corneal epithelial migration. *Curr Eye Res* 1993;12(8):733-40.
10. Petroutsos G, Guimaraes R, and Pouliquen Y. The effect of concentrated antibiotics on the rabbit's corneal epithelium. *Int Ophthalmol* 1984;7(2):65-9

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1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis* 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci* 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. Guideline for Prevention of Surgical Site Infection, 1999. Centers for Disease Control and Prevention (CDC).
4. Aldridge KE, et al. Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of *Bacteroides fragilis* group, *Prevotella*, *Fusobacterium*, *Porphyromonas*, and *Peptostreptococcus* species. *Antimicrob Agents Chemother* 2001;45(4):1238-43.
5. Petroutsos G, et al. Antibiotics and corneal epithelial wound healing. *Arch Ophthalmol* 1983;101(11):1775-8.
6. Bang K, et al. Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice. *Scand J Plast Reconstr Surg Hand Surg* 1998;32(2):129-33.
7. Nelson JD, et al. Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics. *Curr Eye Res* 1990;9(3):277-85.
8. Kjolseth D, et al. Comparison of the effects of commonly used wound agents on epithelialization and neovascularization. *J Am Coll Surg* 1994;179(3):305-12.
9. Nakamura M, et al. Effects of antimicrobials on corneal epithelial migration. *Curr Eye Res* 1993;12(8):733-40.
10. Petroutsos G, Guimaraes R, and Pouliquen Y. The effect of concentrated antibiotics on the rabbit's corneal epithelium. *Int Ophthalmol* 1984;7(2):65-9.

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REFERENCIAS

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis* 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci* 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. Guideline for Prevention of Surgical Site Infection, 1999. Centers for Disease Control and Prevention (CDC).
4. Aldridge KE, et al. Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of *Bacteroides fragilis* group, *Prevotella*, *Fusobacterium*, *Porphyromonas*, and *Peptostreptococcus* species. *Antimicrob Agents Chemother* 2001;45(4):1238-43.
5. Petroutsos G, et al. Antibiotics and corneal epithelial wound healing. *Arch Ophthalmol* 1983;101(11):1775-8.
6. Bang K, et al. Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice. *Scand J Plast Reconstr Surg Hand Surg* 1998;32(2):129-33.

7. Nelson JD, et al. Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics. *Curr Eye Res* 1990;9(3):277-85.
8. Kjolseth D, et al. Comparison of the effects of commonly used wound agents on epithelialization and neovascularization. *J Am Coll Surg* 1994;179(3):305-12.
9. Nakamura M, et al. Effects of antimicrobials on corneal epithelial migration. *Curr Eye Res* 1993;12(8):733-40.
10. Petroustos G, Guimaraes R, and Pouliquen Y. The effect of concentrated antibiotics on the rabbit's corneal epithelium. *Int Ophthalmol* 1984;7(2):65-9.

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BIBLIOGRAFIA

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis* 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci* 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. Guideline for Prevention of Surgical Site Infection, 1999. Centers for Disease Control and Prevention (CDC).
4. Aldridge KE, et al. Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of *Bacteroides fragilis* group, *Prevotella*, *Fusobacterium*, *Porphyromonas*, and *Peptostreptococcus* species. *Antimicrob Agents Chemother* 2001;45(4):1238-43.
5. Petroustos G, et al. Antibiotics and corneal epithelial wound healing. *Arch Ophthalmol* 1983;101(11):1775-8.
6. Bang K, et al. Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice. *Scand J Plast Reconstr Surg Hand Surg* 1998;32(2):129-33.
7. Nelson JD, et al. Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics. *Curr Eye Res* 1990;9(3):277-85.
8. Kjolseth D, et al. Comparison of the effects of commonly used wound agents on epithelialization and neovascularization. *J Am Coll Surg* 1994;179(3):305-12.
9. Nakamura M, et al. Effects of antimicrobials on corneal epithelial migration. *Curr Eye Res* 1993;12(8):733-40.

10. Petroutsos G, Guimaraes R, and Pouliquen Y. The effect of concentrated antibiotics on the rabbit's corneal epithelium. *Int Ophthalmol* 1984;7(2):65-9.

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

REFERENSER

1. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis* 1997;Apr24(4):609-19.
2. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci* 1997;Dec44(4):75-80.
3. Mangram A, Horan TC, Pearson ML, et al. Guideline for Prevention of Surgical Site Infection, 1999. Centers for Disease Control and Prevention (CDC).
4. Aldridge KE, et al. Multicenter survey of the changing in vitro antimicrobial susceptibilities of clinical isolates of *Bacteroides fragilis* group, *Prevotella*, *Fusobacterium*, *Porphyromonas*, and *Peptostreptococcus* species. *Antimicrob Agents Chemother* 2001;45(4):1238-43.
5. Petroustos G, et al. Antibiotics and corneal epithelial wound healing. *Arch Ophthalmol* 1983;101(11):1775-8.
6. Bang K, et al. Gentacoll hampers epithelialisation and neovascularisation in excisional wounds in hairless mice. *Scand J Plast Reconstr Surg Hand Surg* 1998;32(2):129-33.
7. Nelson JD, et al. Corneal epithelial wound healing: a tissue culture assay on the effect of antibiotics. *Curr Eye Res* 1990;9(3):277-85.
8. Kjolseth D, et al. Comparison of the effects of commonly used wound agents on epithelialization and neovascularization. *J Am Coll Surg* 1994;179(3):305-12.
9. Nakamura M, et al. Effects of antimicrobials on corneal epithelial migration. *Curr Eye Res* 1993;12(8):733-40.
10. Petroustos G, Guimaraes R, and Pouliquen Y. The effect of concentrated antibiotics on the rabbit's corneal epithelium. *Int Ophthalmol* 1984;7(2):65-9.

Keep dry

Manufacturer

Use-by date

Attention, see instructions for use

Do not re-use

Temperature limit

Sterilized using ethylene oxide

MANUFACTURER

COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.
Phone: 812 339-2235
Toll Free: 800 457-4500
Toll Free Fax: (800) 554-8335

COOK (CANADA) INC.
111 Sandiford Drive
Stouffville, Ontario
L4A 7X5 CANADA
Phone: 905 640-7110
Toll Free: 800 668-0300

EC REPRESENTATIVE

COOK IRELAND
O'Halloran Road
National Technological Park
Limerick, IRELAND
Phone: +353 61 334440

WILLIAM A. COOK
AUSTRALIA PTY. LTD.
Brisbane Technology Park
95 Brandl Street
Eight Mile Plains
Brisbane, QLD 4113 Australia
Phone: +61 7 38 41 11 88