**COOK**

BIOTECH

EN  
3

**Biodesign®  
Fistula Plug**

Instructions for Use

ES  
7

**Tapón para fístulas  
Biodesign®**

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
PT-BR  
12

**Plugue para fístula  
Biodesign®**

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FP0073-01H

## ENGLISH

### BIODESIGN® FISTULA PLUG

#### INTENDED USE

The Biodesign® Fistula Plug is for implantation to reinforce soft tissue for repair of recto-vaginal or anorectal fistulas.

The plug is supplied sterile and is intended for one time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**MR** This symbol means the following: Magnetic Resonance Safe

This product is intended for use by medical professionals trained in the repair of recto-vaginal and anorectal fistulas.

#### CONTRAINDICATIONS

- This device is derived from a porcine source and should not be used for patients sensitive to porcine materials.
- **Not for vascular use.**

#### PRECAUTIONS

- The plug is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions.
- The plug is sterile if the package is dry, unopened, and undamaged. Do not use if the package seal is broken.
- Discard the plug if mishandling has caused possible damage or contamination or if the plug is past its expiration date.
- Do not implant the plug in a grossly infected or abscessed fistula tract.
- **A seton should be used until there is no evidence of acute inflammation, purulence, or excessive drainage. Allow the tract to stabilize and mature for at least six to eight weeks before placing the plug.**
- Ensure that the plug is hydrated prior to placement, cutting, or suturing.
- Placement of the plug in fistula tracts less than 1 cm in length can result in incomplete incorporation and/or expulsion.

#### GENERAL

- Exercise good surgical practice for the management of clean-contaminated, contaminated, or infected fields.
  - Risk of infection following implantation may be reduced by debriding the fistula tract prior to plug insertion.
  - Always introduce the Fistula Plug through the rectal (primary) fistula opening.
  - Draw the Fistula Plug into the fistula tract completely, until the button is set against the rectal wall.
  - To allow drainage to occur, do not completely occlude the external/vaginal (secondary) fistula opening.
- IMPORTANT:** Users should instruct patients to abstain from strenuous physical activity for at least two weeks following anorectal fistula repair, or six weeks following recto-vaginal fistula repair. Refer to the section on Post-Operative Care.
- Expect some drainage for two to four weeks. Drainage can continue for up to 12 weeks after the procedure as the plug incorporates and the fistula tract closes.

#### POTENTIAL COMPLICATIONS

Complications that can occur with the Fistula Plug include, but are not limited to:

- |   |                                    |
|---|------------------------------------|
| • Abscess                                       | • Incontinence to stool or flatus* |
| • Allergic reaction                             | • Induration                       |
| • Delayed or failed incorporation of the device | • Infection                        |
| • Erosion                                       | • Inflammation                     |
| • Extrusion                                     | • Migration                        |
| • Fistula recurrence                            | • Seroma formation                 |

\*Note: Reported incontinence rates are very low and are typically temporary due to the fact that sphincter muscles are not divided. Decreased/minimal risk of incontinence is generally recognized as a benefit of using fistula plugs. Fistula plugs have lower reported incontinence rates when compared to other fistula repair procedures like advancement flaps.<sup>1</sup>

#### STORAGE

Store the plug in a clean, dry location at room temperature.

#### STERILIZATION

The plug has been sterilized with ethylene oxide. Do not resterilize.

## INSTRUCTIONS FOR USE

### REQUIRED MATERIALS

- A sterile basin
- Hydration fluid: room temperature sterile saline or sterile lactated Ringer's solution
- Resorbable suture, such as 2-0 or 0 polydioxanone (PDO or PDS) or coated polyglycolic acid suture (coated PGA)
- Silk suture
- 5 Fr catheter
- Hydrogen peroxide or saline
- 10 cc syringe
- Brush for tract debridement (such as the Cook Fistula Brush)

**NOTE:** Handle the plug using aseptic technique. Minimize contact with latex gloves.

### PREPARATION

1. Use a sterile probe or other appropriate measuring instrument to approximate the diameter of the fistula tract at the internal (rectal) opening. Select the proper plug size according to the table below.

Diameter of Internal Fistula Opening	Plug Size	Plug Part Number
> 1 mm to 2 mm	2 mm	C-FP-0.2
> 2 mm to 4 mm	4 mm	C-FP-0.4
> 4 mm to 7 mm	7 mm	C-FP-0.7

2. Using aseptic technique, open the outer pouch and pass the inner pouch containing the plug onto the sterile field.
3. Open the inner pouch and remove the tray containing the plug. Open the tray and place the plug into a sterile basin.
4. Add enough hydration fluid to the basin to fully submerge the plug. Allow the plug to hydrate until the desired handling characteristics are achieved. Do not hydrate the plug for more than one minute.
5. Remove the hydrated plug from the basin using forceps.
6. Attach a silk suture (approximately 30 cm in length) around the tail (narrow end) of the plug to use as a tether to pull the plug into the fistula tract.
7. Prepare the patient using standard surgical techniques appropriate for fistula repair.

**NOTE:** Bowel preparation and/or rectal enemas may be performed at the surgeon's discretion. A single dose of systemic antibiotics is recommended before the procedure.<sup>2,3</sup>

### PROCEDURE

1. Perform under local, regional, or general anesthesia.
2. If a seton is not in place, identify the internal (rectal) fistula opening.
 

**NOTE:** Failure to locate the internal opening may lead to persistence of the fistula. If the internal opening cannot be reliably identified, an alternative treatment method should be considered.
3. Insert a fistula probe, seton, or suitable instrument through the fistula tract, entering through the external opening and exiting via the internal opening.
4. Attach the brush to the instrument at the internal opening.
5. If a seton is in place, cut the seton and attach the brush to the seton at the internal (rectal) opening.
6. Pull the brush into the tract and gently debride to remove non-vascularized tissue. Do not enlarge the tract.
 

**NOTE:** Blood should be noted on the bristles of the brush and at both the internal and external fistula openings after adequate debridement.
7. Flush the tract thoroughly with hydrogen peroxide or sterile saline using the 10 cc syringe and 5 Fr catheter.
8. Attach the tether connected to the Fistula Plug to the end of the brush at the internal opening.


### Superficial Button Placement

- Use the tether to pull the plug into the internal opening and through the fistula tract. Ensure the plug is drawn completely into the tract. Do not use the entire plug unless the tract is the full length of the plug.


**IMPORTANT:** The internal opening is the high-pressure zone of the fistula. The button end of the plug must be securely fastened to the internal opening or fully implanted below the mucosa to prevent ingress of fecal debris. The higher pressures within the rectum and anal canal assist in maintaining the plug in the fistula tract.


**Recto-vaginal fistula repair**

**Anorectal fistula repair**

- Once the plug is properly positioned, gently pull the tether to hold the plug in place during fixation on the internal side.


- Using a suitable resorbable suture, secure the button to the adjacent tissue by suturing through the center portion of the button and directly through the rectal wall and mucosa. Obtain an adequate bite of sphincter muscle and secure the button to the tissue by tying the two suture ends together. Equally space four interrupted sutures around the button to evenly distribute force.

**Superficial Button Placement**

### Submucosal Button Placement


- Create a mucosal flap of sufficient size to accommodate the 16 mm diameter button.
- Use the tether to pull the plug into the internal opening and through the fistula tract. Ensure the plug is drawn completely into the tract and that the button is positioned within the mucosal flap.
- Suture the plug according to PROCEDURAL STEP 11.

**Submucosal Button Placement**

- Reapproximate the mucosal flap and securely suture the tissue over the button.

**Submucosal Button Placement**

16. Trim the excess external plug flush with the skin, and place a sterile dressing over the external opening.




**NOTE:** Do not close the secondary opening. Complete obstruction of the fistula opening may result in accumulation of fluid, infection, or abscess.

**NOTE:** For anorectal fistulas, fixation of the plug at the external opening is not recommended. Fixation at the vaginal opening can be considered when treating recto-vaginal fistulas.

**NOTE:** Continued drainage of the fistula tract can be expected to last several days up to 12 weeks after plug placement.

**OPTIONAL FIXATION FOR RECTO-VAGINAL FISTULAS ONLY**


17. Using a suitable resorbable suture, place a minimum of two stitches through the middle of the Fistula Plug at the vaginal fistula opening. Obtain an adequate bite of both the plug and tissue with the suture, allowing 2 mm between each stitch.

**Cross-sectional View of Vaginal Wall**

**IMPORTANT:** Do not close the secondary opening completely when suturing the plug to tissue. Complete obstruction of the fistula opening may result in accumulation of fluid, infection, or abscess.

18. Once the plug is properly secured, trim and discard any remaining portion of the plug. Leave at least 2 mm of material between the end of the trimmed Fistula Plug and suture to prevent suture pull-through.

19. Discard any unused portions according to institutional guidelines for disposal of medical waste.

**POST-OPERATIVE CARE**

To provide the best environment for tissue integration into the Fistula Plug, post-operative patient activity should be minimized. Provide patients with a list of post-operative care recommendations. The following patient guidelines should be considered:

1. Patients should avoid any strenuous physical activity beyond a gentle walk for at least two weeks following anorectal fistula repair, or six weeks following recto-vaginal fistula repair.
2. Patients should avoid any heavy lifting over 10 lbs (5 kg) for at least two weeks following anorectal fistula repair, or at least six weeks following recto-vaginal fistula repair.

3. Patients should use a stool softener for at least two weeks after their procedure.
4. Patients should observe a liquid diet for 48 hours after their procedure.
5. Patients should begin a high fiber diet 48 hours after their procedure.
6. Patients should abstain from sexual intercourse and other forms of vaginal/rectal insertion (e.g. tampons) for at least two weeks following anorectal fistula repair, or six weeks following recto-vaginal fistula repair.
7. Patients should take over-the-counter pain medication as needed.
8. Patients should avoid the use of topical steroids.

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


  
**LOT****Batch code****REF****Catalogue number****Consult instructions for use****Contains a biological material  
of animal origin****Country of manufacture****Date of manufacture****Do not resterilize****Do not re-use****Do not use if package is  
damaged and consult  
instructions for use****Keep dry****Manufacturer****MD****Medical device****Double sterile barrier system****STERILE****Sterile****STERILE EO****Sterilized using ethylene oxide****UDI****Unique device identifier****Use-by date****MANUFACTURER**

COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.