Biodesign®

SINONASAL REPAIR GRAFT

FP0089-01G

COOK®
MEDICAL

LOT Batch code

REF Catalogue number

i Consult instructions for use

BIO Contains biological material of animal origin

US Country of manufacture

STERILIZE Do not re-sterilize

2 Do not re-use

Do not use if package is damaged and consult instructions for use

 Keep dry

 Manufacturer

MD Medical device

 Double sterile barrier system

STERILE EO Sterilized using ethylene oxide

UDI Unique device identifier

 Use-by date

BIODESIGN® SINONASAL REPAIR GRAFT

INTENDED USE

The Biodesign® Sinonasal Repair Graft is intended to separate tissue or structures compromised by surgical trauma and act as an adjunct to aid in the natural healing process. The device is indicated for use where an open wound dressing material is required in the nasal and/or sinus cavities following nasal and/or sinus surgery where separation of tissues or structures is desired during nasal septal perforation repair and donor site repair.

The device is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

MR This symbol means the following: Magnetic Resonance Safe

SINONASAL REPAIR GRAFT This symbol means the following: Sinonasal Repair Graft

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This graft is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.

COMPOSITION

The Biodesign Sinonasal Repair Graft is an extracellular membrane derived from the Small Intestinal Submucosa (SIS) of pigs from qualified animal production facilities. SIS is obtained from the intestine using a process that retains the natural composition of matrix molecules such as collagen (Types I, III, IV, VI), glycosaminoglycans (hyaluronic acid, chondroitin sulfate A and B, heparin, and heparan sulfate), proteoglycans, growth factors (FGF-2, TGF-β), and fibronectin.^{1,2,3}

PROPERTIES

The Biodesign Sinonasal Repair Graft conforms to the defect once hydrated. After hydration at the site of implant, the material remains in place and is not disrupted by the closure of the surgical site. If desired, the graft can be affixed with a suture.

SAFETY

A viral inactivation study using the Biodesign Sinonasal Repair Graft material demonstrated at least a 6.1-log reduction of the PrV virus, a 5.5-log reduction in Reo-3, 4.3-log reduction in A-MuLV and 3.6 log in PPV. This represented a reduction in viral burden to the level of detection in the test. Data were published in the peer reviewed article cited in the references.⁴

PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, re-sterilize, and/or reuse may lead to device failure and/or transmission of disease.
- Do not re-sterilize.** Discard all open and unused portions of the graft.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- The graft should not be applied until excessive exudate, bleeding, acute swelling, and infection is controlled.
- Ensure that all layers of the graft are secured during fixation.

POTENTIAL COMPLICATIONS

The complications listed below are possible with the use of surgical graft materials in nasal/sinus surgery. If any of these conditions occur, the graft should be removed if possible. However, if graft removal is not possible due to graft incorporation, then appropriate care should be given to treat the complication.

- Allergic reaction
- Bleeding
- Blistering
- Excessive redness
- Hematoma
- Infection/Toxic Shock Syndrome
- Inflammation
- Nasal obstruction
- Pain
- Persistence of perforation or recurrence

STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

STERILIZATION

The graft has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

NOTE: Handle the graft using aseptic technique, minimizing contact with latex gloves.

Preparation

- Prepare the wound bed using standard methods to ensure it is free of exudate and devitalized tissue. Debride the wound as necessary to ensure the wound edges contain

viable tissue.

- Ensure that any bleeding is controlled before applying the graft.
- Irrigate the wound thoroughly with sterile saline.

Procedure

- Choose the graft size based on the size of the defect, ensuring that the graft covers the entire defect and extends slightly beyond its margins.
- Remove the packaging containing the graft from the box.
- Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
- Open the inner pouch and remove the graft.
- If necessary, trim the graft.
- Place the graft in a dry state and smooth into place to ensure complete contact with the defect.
- If multiple grafts are used, slightly overlap the edges of the grafts.
- Thoroughly hydrate the graft in place using sterile saline.
- Securely affix the graft using sutures, nasal packing, splinting, or other appropriate method.

Reassessment

NOTE: Successful integration of the sinonasal graft may result in the formation of a carmel-colored or off-white gel that looks like slough. **DO NOT DISRUPT THIS GEL DURING REASSESSMENT.** This gel is part of the wound healing process and will eventually resolve without intervention.

REFERENCES

- Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med.* 2007;18:537-543.
- Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL. Glycosaminoglycan content of small intestinal submucosa: A bioscaffold for tissue replacement. *Tissue Eng.* 1996;2:209-217.
- Data on File, Cook Biotech, Inc.
- Hodde J, Hiles M. Virus safety of a porcine-derived medical device: Evaluation of a viral inactivation method. *Biotechnol Bioeng* 2002;79:211-216.

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RÉFÉRENCES

1. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med.* 2007;18:537-543.
2. Hodde JP, Badylak SF, Brightman AO, Voytik-Harbin SL. Glycosaminoglycan content of small intestinal submucosa: A bioscaffold for tissue replacement. *Tissue Eng.* 1996;2:209-217.
3. Data on File, Cook Biotech, Inc.
4. Hodde J, Hiles M. Virus safety of a porcine-derived medical device: Evaluation of a viral inactivation method. *Biotechnol Bioeng* 2002;79:211-216.

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© COOK BIOTECH 2024