FP0095-01F

	Batch code		Keep dry
	Catalogue number		Manufacturer
	Consult instructions for use		Medical device
	Contains biological material of animal origin		Double sterile barrier system
	Country of manufacture		Sterilized using ethylene oxide
	Do not re-sterilize		Unique device identifier
	Do not re-use		Use-by date
	Do not use if package is damaged and consult instructions for use		

MANUFACTURER
 COOK BIOTECH
 INCORPORATED
 1425 Innovation Place
 West Lafayette, IN 47906 U.S.A.

ENGLISH

BIODESIGN® DURAPLASTY GRAFT

DESCRIPTION

The Biodesign Duraplasty Graft is a nonporous sheet of extracellular collagen matrix derived from porcine small intestinal submucosa. The graft is non-pyrogenic and has sufficiently low endotoxin levels to make it suitable for applications where it will contact cerebrospinal fluid during repair of dura mater.

INDICATIONS FOR USE / INTENDED USE

The Biodesign Duraplasty Graft is intended for use as a dura substitute for the repair of dura mater. This graft is supplied sterile in peel-open packages and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

MR This symbol means the following: Magnetic Resonance Safe

DURAPLASTY GRAFT This symbol means the following: Duraplasty Graft

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

The Biodesign Duraplasty Graft is not designed, sold or intended for use except as described in the indications for use and is contraindicated:

- For use in patients with a known sensitivity to porcine materials;
- For repair of spinal neural tube defects; and
- For anterior spinal surgery with dural resection (e.g. transoral surgery)

Additionally:

- Use with caution in infected regions.
- Covering defects involving mastoid cells is not recommended.
- It is not recommended for large defects at the skull base following surgery; however the Biodesign Duraplasty Graft may be used to augment other forms of specific repair (e.g. pedicled flaps or vascularized pedicled flaps).

PRECAUTIONS

• The Biodesign Duraplasty Graft is designed to augment skull base repair where layering techniques such as bony buttressing, pedicled flaps or vascularized pedicled flaps and packing are currently used. The graft should not replace standard layering techniques or be implanted as a stand-alone repair.

- Peer reviewed literature¹ has reported the use of this material as a dural substitute at the skull base in defects up to 4.5 cm².
- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- Do not use if the product package is damaged or opened.
- Discard the graft if damage or contamination is observed, or if the graft is past its expiration date.
- Prior to touching the graft, rinse surgical gloves thoroughly to remove the glove powder.
- Ensure that graft sterility is maintained during preparation and implantation at repair site.
- Suturing is not required but if the graft is to be sutured, tensionless suturing technique is recommended
- Discard all open and unused portions of the graft.
- The graft should be cut to size ensuring an overlap to cover the existing dura.
- Ensure that the graft is hydrated and all layers of the graft are secured if suturing the graft into place.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials in neurosurgical procedures.

- Adhesion
- Allergic reaction
- Calcification
- CSF leak
- Delayed hemorrhage
- Infection
- Inflammation

STORAGE

Store in a clean place at room temperature. Do not place in freezer. Avoid excessive heat. Keep dry.

STERILIZATION

The graft has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

NOTE: The Biodesign Duraplasty Graft is an onlay graft and does not require suturing into place.

Required Materials

- Sterile basin
- Smooth forceps

- Hydration fluid: sterile, room temperature saline or lactated Ringer's solution

NOTE: Handle the graft using aseptic technique, minimizing contact with latex gloves.

Procedure

1. Choose the graft size based on the size of the dura mater defect, ensuring that the graft covers the entire defect and extends slightly beyond the margins.
2. Remove the packaging containing the graft from the box.
3. Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
4. Open the inner pouch and remove the graft.
5. If necessary, trim the graft to fit the patient's defect, providing a small allowance for overlap of the margins.
6. Place the graft into the sterile basin on the sterile field. Multiple grafts may be hydrated simultaneously in the same basin.
7. Add enough hydration fluid to the basin to completely cover the graft(s).
8. Hydrate the graft in the fluid until the desired handling characteristics are achieved. A hydration time of at least two minutes is recommended.
9. Transfer the graft to the defect site and affix into place. Suturing is not required, but tensionless, atraumatic stay sutures may be used if desired.
10. Complete the surgical procedure.
11. Discard unused portions of the graft according to institutional guidelines for disposal of medical waste.

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