**COOK**  
MEDICAL

**CE**  
2797

**EN**  
3

**Biodesign®  
Duraplasty Graft**

Instructions for Use

**CS**  
4

**Biodesign®  
Graft pro duraplastiku**

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**DA**  
6

**Biodesign®  
Duraplastik-implantat**

Brugsanvisning

**DE**  
8

**Biodesign®  
Duraplastik-Implantat**

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**EL**  
9

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
# BIODESIGN® DURAPLASTY GRAFT

### DESCRIPTION

Cook® Biodesign Duraplasty Graft is a nonporous sheet of extracellular collagen matrix derived from porcine small intestinal submucosa. The graft is non-pyrogenic and has sufficiently low endotoxin levels to make it suitable for applications where it will contact cerebrospinal fluid during repair of dura mater. Biodesign® Duraplasty Graft is comprised of an extracellular matrix (ECM) and is fully remodeled during the healing process. Animal studies have demonstrated that complete resorption of the graft, when used as indicated, is evident by 60 days.<sup>1,2</sup>

### INDICATIONS FOR USE / INTENDED USE

The Biodesign® Duraplasty Graft is intended for use as a dura substitute for the repair of dura mater. This graft is supplied sterile in peel-open packages and is intended for one-time use.

Biodesign® Duraplasty Graft is designed and labeled specifically for the above indication. Use of this product for indications other than specified above is not recommended and may result in increased complications.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**MR**

This symbol means the following: Magnetic Resonance Safe

**DURAPLASTY GRAFT** This symbol means the following: Duraplasty Graft

This product is intended for use by trained medical professionals.

### CONTRAINDICATIONS

The Biodesign Duraplasty Graft is not designed, sold or intended for use except as described in the indications for use and is contraindicated:

- For use in patients with a known history of hypersensitivity to porcine derived materials;
  - For repair of spinal neural tube defects; and
  - For anterior spinal surgery with dural resection (e.g. transoral surgery)
- Additionally:
- Use with caution in infected regions.
  - It is not recommended to cover defects involving mastoid cells.
  - It is not recommended for large defects at the skull base following surgery; however the Biodesign Duraplasty Graft may be used to augment other forms of specific repair (e.g. pedicled flaps or vascularized pedicled flaps).

### PRECAUTIONS

- The Biodesign Duraplasty Graft is designed to augment skull base repair where layering techniques such as bony buttressing, pedicled flaps or vascularized pedicled flaps and packing are currently used. The graft should not replace standard layering techniques or be implanted as a stand-alone repair.
- Peer reviewed literature<sup>3</sup> has reported the use of this material as a dural substitute at the skull base in defects up to 4.5 cm<sup>2</sup>.
- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- Do not use if the product package is damaged or opened.
- Discard graft if damage or contamination is observed, or if the graft is past its expiration date.
- Prior to touching the graft, wash surgical gloves thoroughly to remove the glove powder.
- Ensure that graft sterility is maintained during preparation and implantation at repair site.
- Suturing is not required but if the graft is to be sutured, tensionless suturing technique must be used.
- Discard all open and unused portions of the graft sheets.
- The graft should be cut to size ensuring an overlap to cover the existing dura.
- Ensure that graft is rehydrated and all layers of the graft are secured if fixation with suture is employed.

### POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials in neurosurgical procedures.

- Infection, adhesion, CSF leak, delayed hemorrhage and calcification
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation.)
- Allergic reaction

## STORAGE

This graft should be stored in a clean, dry location at controlled room temperature. Do not freeze this graft.

## STERILIZATION

This graft has been sterilized with ethylene oxide.

## INSTRUCTIONS FOR USE

Biodesign Duraplasty Graft is an onlay graft and does not require sutures.

### Required Materials

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Rehydration fluid: at least 100 mL of sterile, room temperature, saline or sterile, room temperature lactated Ringer's solution for each graft sheet

**NOTE: Always handle the graft using aseptic technique, minimizing contact with latex gloves.**

1. Using aseptic technique, remove the graft inner pouch from its outer pouch, and place the inner pouch in the sterile field.
2. Open the inner pouch carefully, and aseptically remove the graft sheet with the sterile forceps.
3. Place the graft sheet into the sterile dish in the sterile field. (Multiple graft sheets may be rehydrated simultaneously in the same dish.)
4. Add to the dish at least 100 mL of the rehydration fluid for each graft sheet.

5. Allow graft sheets to rehydrate for at least two minutes.
6. Using aseptic technique, trim the graft sheet to fit the site, providing a small allowance for overlap.

**NOTE:** An alternative method is to cut the graft sheet to size prior to rehydration. If this method is selected, be sure to rehydrate the graft sheet prior to fixing it into place. See step 5.

7. Using aseptic technique, transfer the graft sheet to the graft site and fix into place. Suturing is not required, but tensionless, atraumatic stay sutures may be used if desired.

**NOTE:** Surgical experience indicates that fixing graft sheets with close tissue approximation produces better outcomes.

8. Complete the surgical procedure.
9. Closed suction wound drainage is recommended for 1-3 days postoperatively.
10. Discard any unused portions according to institutional guidelines for medical waste.

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© COOK BIOTECH 2024
West Lafayette, IN 47906 U.S.A.