EN 3	<b>Biodesign® Fistula Plug</b> Instructions for Use
CS 6	<b>Píštělová zátka Biodesign®</b> Návod k použití
DA 10	<b>Biodesign® fistelprop</b> Brugsanvisning
DE 14	<b>Biodesign® Fistelverschluss</b> Gebrauchsanweisung
EL 18	<b>Βύσμα συριγγίου Biodesign®</b> Οδηγίες χρήσης
ES 22	<b>Tapón para fístulas Biodesign®</b> Instrucciones de uso
FR 26	<b>Obturateur de fistule Biodesign®</b> Mode d'emploi
HR 30	<b>Biodesign® čep za fistulu</b> Upute za uporabu
IT 33	<b>Tampone per fistola Biodesign®</b> Istruzioni per l'uso
MK 37	<b>Biodesign® Фистулен чеп</b> Упатство за употреба
NL 41	<b>Biodesign® fistelplug</b> Gebruiksaanwijzing
NO 45	<b>Biodesign® fistelplugg</b> Bruksanvisning
PT-BR 49	<b>Plugue para fístula Biodesign®</b> Instruções de uso
SL 52	<b>Čep za fistulo Biodesign®</b> Navodila za uporabo
SR 56	<b>Čep za fistulu Biodesign®</b> Uputstvo za upotrebu

**SV  
60**

**Biodesign®  
fistelplugg**

Bruksanvisning

**TR  
63**

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Kullanma Talimatı

## BIODESIGN® FISTULA PLUG

### DEVICE DESCRIPTION

The Biodesign® Fistula Plug is a tapered small intestinal submucosa (SIS) cylinder with an SIS button. The SIS button is fastened to the wide end of the cylinder with resorbable polyglycolic acid (PGA) suture. The plug is used to reinforce soft tissue for the repair of anorectal fistulas. Once implanted, the natural composition of SIS allows for the patient's healing mechanisms to deposit cells and collagen during cellular and extracellular matrix metabolism. Publications by Ellis et al.<sup>1</sup> and Jayne et al.<sup>2</sup> report no evidence of residual fistula tract or fluid observed in 75% (6/8) and 49% (54/110) of patients, respectively, under MRI visualization 12-months after plug implantation. This is consistent with histological evidence in which patients' tissues completely incorporate SIS material over 7-17 months, resulting in new, remodeled tissue.<sup>3-5</sup> The plug can be cut to size to accommodate the patient's anatomy and is provided sterile for single use only.

### INTENDED USE

The Biodesign Fistula Plug is for implantation to reinforce soft tissue for repair of anorectal fistulas. The plug is supplied sterile and is intended for one-time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**MR** This symbol means the following: Magnetic Resonance Safe

This product is intended for use by medical professionals trained in the repair of anorectal fistulas.

### CONTRAINDICATIONS

- The plug is derived from a porcine source and should not be used for patients sensitive to porcine materials.
- **Not for vascular use.**

### PRECAUTIONS

- The plug is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to failure and/or transmission of disease.
- Do not resterilize. Discard all open and unused portions.
- The plug is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the plug if mishandling has caused possible damage or contamination or if the plug is past its expiration date.
- Do not implant the plug in a grossly infected or abscessed fistula tract.
- **A seton should be used until there is no evidence of acute inflammation, purulence, or excessive drainage. Allow the tract to stabilize and mature for at least six to eight weeks before placing the plug.**
- Ensure that the plug is hydrated prior to placement, cutting, or suturing.
- Placement of the plug in fistula tracts less than 1 cm in length can result in incomplete incorporation and/or expulsion.
- No studies have been conducted to evaluate performance of the plug in pediatric populations.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the plug.

### GENERAL

- Users should be familiar with surgical technique for anorectal fistula repair.
- Users should exercise good surgical practice for the management of clean-contaminated, contaminated or infected fields.
- Risk of infection following implantation may be reduced by debriding the fistula tract prior to plug insertion.
- **IMPORTANT:** Users should instruct patients to abstain from strenuous physical activity for at least six weeks after plug placement.
- **Expect some drainage for two to four weeks. Drainage can continue for up to 12 weeks after the procedure as the plug incorporates and the fistula tract closes.**

### POTENTIAL COMPLICATIONS

Complications that can occur with the use of the plug include:

- |   |                |
|---|----------------|
| • Abscess                                     | • Infection    |
| • Allergic reaction                           | • Inflammation |
| • Bleeding                                    | • Migration    |
| • Delayed or failed incorporation of the plug | • Pain         |
| • Erosion                                     | • Redness      |
| • Extrusion                                   | • Sepsis       |

- Fistula persistence or recurrence
- Induration
- Seroma Formation
- Swelling

## STORAGE

Store the plug in a clean, dry location at room temperature.

## STERILIZATION

The plug has been sterilized with ethylene oxide. Do not resterilize.

## INSTRUCTIONS FOR USE

### Required Materials

- Sterile basin
- Hydration fluid: room temperature sterile saline or sterile lactated Ringer's solution
- Resorbable suture, such as 2-0 or 0 polydioxanone (PDO or PDS) or coated polyglycolic acid suture (coated PGA)
- Silk suture
- 5 Fr catheter
- Hydrogen peroxide or saline
- 10 cc syringe
- Brush for debridement (such as the Cook® Fistula Brush)

**NOTE:** Handle the plug using aseptic technique, minimizing contact with latex gloves.

## PREPARATION

1. Use a sterile probe or other appropriate measuring instrument to approximate the diameter of the fistula tract at the internal (rectal) opening. Select the proper plug size according to the table below.


Diameter of Internal Fistula Opening	Plug Size	Plug Part Number
> 1 mm to 2 mm	2 mm	C-FP-0.2-2
> 2 mm to 4 mm	4 mm	C-FP-0.4-2
> 4 mm to 7 mm	7 mm	C-FP-0.7-2

2. Using aseptic technique, open the outer pouch and pass the inner pouch containing the plug onto the sterile field.
3. Open the inner pouch and remove the tray containing the plug. Open the tray and place the plug into a sterile basin.
4. Add enough hydration fluid to the basin to fully submerge the plug. Allow the plug to hydrate until the desired handling characteristics are achieved. Do not hydrate the plug for more than 1 minute.
5. Remove the hydrated plug from the basin using forceps.
6. Attach a silk suture (approximately 30 cm in length) around the tail (narrow end) of the plug to use as a tether to pull the plug into the fistula tract.
7. Prepare the patient using standard surgical techniques appropriate for anal fistula repair.

**NOTE:** The recommended practice for preoperative bowel preparation in elective colorectal surgery includes mechanical bowel cleansing using enemas and cathartic agents, and administration of prophylactic oral or intravenous antimicrobial agents. Insufficient cleansing or inadequate antibacterial prophylaxis can predispose the patient to infections.<sup>6-8</sup>


## PROCEDURE

1. Perform under local, regional, or general anesthesia.
2. **If a seton is not in place**, identify the internal (rectal) fistula opening. **NOTE:** Failure to locate the internal opening may lead to persistence of the fistula. If the internal opening cannot be reliably identified, an alternative treatment method should be considered.
3. Insert a fistula probe, seton, or suitable instrument through the fistula tract, entering through the external opening and exiting via the internal opening.
4. Attach the brush to the instrument at the internal opening.
5. **If a seton is in place**, cut the seton and attach the brush to the seton at the internal (rectal) opening.
6. Pull the brush into the tract and gently debride to remove non-vascularized tissue. Do not enlarge the tract. **NOTE:** Blood should be noted on the bristles of the brush and at both the internal and external fistula openings after adequate debridement.
7. Flush the tract thoroughly with hydrogen peroxide or sterile saline using the 10 cc syringe and 5 Fr catheter.
8. Attach the tether connected to the tail of the plug to the end of the brush at the internal opening.


**Superficial Button Placement**

9. Use the tether to pull the plug into the internal opening and through the fistula tract. Ensure the plug is drawn completely into the tract. Do not use the entire plug unless the tract is the full length of the plug.  
**IMPORTANT:** The internal opening is the high-pressure zone of the fistula. The button end of the plug must be securely fastened to the internal opening or fully implanted below the mucosa to prevent ingress of fecal debris. The higher pressures within the rectum and anal canal assist in maintaining the plug in the fistula tract.
10. Once the plug is properly positioned, gently pull the tether to hold the plug in place during fixation on the internal side.
11. Using a 2-0 PDO, PDS, coated PGA or comparable suture, secure the button to the adjacent tissue by suturing through the center portion of the button and directly through the rectal wall and mucosa. Obtain an adequate bite of sphincter muscle and secure the button to the tissue by tying the two suture ends together. Equally space four interrupted sutures around the button to evenly distribute force.


**Submucosal Button Placement**

12. Create a mucosal flap of sufficient size to accommodate the 16 mm diameter button.
13. Use the tether to pull the plug into the internal opening and through the fistula tract. Ensure the plug is drawn completely into the tract and that the button is positioned within the mucosal flap.
14. Suture the plug according to PROCEDURE STEP 11.

15. Reapproximate the mucosal flap and securely suture the tissue over the button.

16. Trim the excess external plug flush with the skin. Do not suture the plug to the external opening. The external opening may be enlarged to facilitate drainage.  
**NOTE:** Complete obstruction of the external opening may result in accumulation of fluid, infection, or abscess.

17. Place a sterile dressing over the external opening.
18. Discard any unused portions of the plug according to institutional guidelines for disposal of medical waste.

## POST-OPERATIVE CARE

To provide the best environment for tissue integration into the plug, post-operative patient activity should be minimized. Provide patients with a list of post-operative care recommendations. The following patient guidelines should be considered:

1. Patients should avoid any strenuous physical activity beyond a gentle walk for at least six weeks.
2. Patients should avoid any heavy lifting over 5 kg (10 lbs) for at least six weeks.
3. Patients should refrain from inserting anything into the vagina or rectum for at least six weeks.
4. Patients should shower standing up to keep the area clean.
5. Patients should observe a liquid diet for 48 hours after surgery.
6. Patients should begin a high fiber diet after 48 hours.
7. Patients should use a stool softener for at least two weeks after surgery.
8. Patients should take over-the-counter pain medication as needed.
9. Patients should avoid the use of topical steroids.

## REFERENCES

1. Ellis CN, Rostas JW, Greiner FG. Long-term outcomes with the use of bioprosthetic plugs for the management of complex anal fistulas. *Dis Colon Rectum*. 2010;53(5):798-802.
2. Jayne DG, Scholefield J, Tolan D, et al. Anal fistula plug versus surgeon's preference for surgery for trans-sphincteric anal fistula: the FIAT RCT. *Health Technol Assess*. 2019;23(21):1-76.
3. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol*. 2004;172(1):215-218.
4. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestinesubmucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
5. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinesubmucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.
8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg*. 2017;152(8):784-791.

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- Jayne DG, Scholefield J, Tolan D, et al. Anal fistula plug versus surgeon's preference for surgery for trans-sphincteric anal fistula: the FIAT RCT. *Health Technol Assess.* 2019;23(21):1-76.
- Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
- Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestinesubmucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
- Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinesubmucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
- Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
- Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery: comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
- Berrios-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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## QUELLEN

1. Ellis CN, Rostas JW, Greiner FG. Long-term outcomes with the use of bioprosthetic plugs for the management of complex anal fistulas. *Dis Colon Rectum*. 2010;53(5):798-802.
2. Jayne DG, Scholefield J, Tolan D, et al. Anal fistula plug versus surgeon's preference for surgery for trans-sphincteric anal fistula: the FIAT RCT. *Health Technol Assess*. 2019;23(21):1-76.
3. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol*. 2004;172(1):215-218.
4. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestinesubmucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
5. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinesubmucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.

7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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5. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinals submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery: comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.
8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg*. 2017;152(8):784-791.

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### REFERENCIAS

1. Ellis CN, Rostas JW, Greiner FG. Long-term outcomes with the use of bioprosthetic plugs for the management of complex anal fistulas. *Dis Colon Rectum*. 2010;53(5):798-802.
2. Jayne DG, Scholefield J, Tolan D, et al. Anal fistula plug versus surgeon's preference for surgery for trans-sphincteric anal fistula: the FIAT RCT. *Health Technol Assess*. 2019;23(21):1-76.
3. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol*. 2004;172(1):215-218.
4. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestinesubmucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
5. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinesubmucosa as a prosthetic device for laparoscopic repair of hernias incontaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.

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### REFERANSER

1. Ellis CN, Rostas JW, Greiner FG. Long-term outcomes with the use of bioprosthetic plugs for the management of complex anal fistulas. *Dis Colon Rectum*. 2010;53(5):798-802.
2. Jayne DG, Scholefield J, Tolan D, et al. Anal fistula plug versus surgeon's preference for surgery for trans-sphincteric anal fistula: the FIAT RCT. *Health Technol Assess*. 2019;23(21):1-76.
3. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol*. 2004;172(1):215-218.
4. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestinesubmucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
5. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinesubmucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery: comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.
8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg*. 2017;152(8):784-791.

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#### REFERENCE

1. Ellis CN, Rostas JW, Greiner FG. Long-term outcomes with the use of bioprosthetic plugs for the management of complex anal fistulas. *Dis Colon Rectum*. 2010;53(5):798-802.
2. Jayne DG, Scholefield J, Tolan D, et al. Anal fistula plug versus surgeon's preference for surgery for trans-sphincteric anal fistula: the FIAT RCT. *Health Technol Assess*. 2019;23(21):1-76.
3. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol*. 2004;172(1):215-218.
4. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestinesubmucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
5. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinesubmucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.

8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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## REFERANSLAR

1. Ellis CN, Rostas JW, Greiner FG. Long-term outcomes with the use of bioprosthetic plugs for the management of complex anal fistulas. *Dis Colon Rectum*. 2010;53(5):798-802.
2. Jayne DG, Scholefield J, Tolan D, et al. Anal fistula plug versus surgeon's preference for surgery for trans-sphincteric anal fistula: the FIAT RCT. *Health Technol Assess*. 2019;23(21):1-76.
3. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol*. 2004;172(1):215-218.
4. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestinesubmucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
5. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinesubmucosa as a prosthetic device for laparoscopic repair of hernias incontaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.
8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg*. 2017;152(8):784-791.

Keep dry

Manufacturer

Use-by date

Caution

Do not re-use

Date of Manufacture

Magnetic resonance safe

Sterilized using ethylene oxide

**MANUFACTURER**

COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.

**EC REPRESENTATIVE**

Cook Biotech Europe ApS  
Sandet 6, DK-4632  
Bjaeverskov, DENMARK  
Tel: +45 56868686