EN 3	<b>Biodesign® Rectopexy Graft</b> Instructions for Use
CS 6	<b>Biodesign® Graft pro rektopexi</b> Návod k použití
DA 9	<b>Biodesign® Rektopeksiimplantat</b> Brugsanvisning
DE 12	<b>Biodesign® Rektopexie-Implantat</b> Gebrauchsanweisung
EL 15	<b>Biodesign® Μόσχευμα πρωκτοπηξίας</b> Οδηγίες χρήσης
ES 19	<b>Biodesign® Injerto para rectopexia</b> Instrucciones de uso
FR 22	<b>Biodesign® Implant biologique pour rectopexie</b> Mode d'emploi
HR 26	<b>Biodesign® transplantat za rektopeksiju</b> Upute za uporabu
IT 29	<b>Biodesign® Innesto per rettopenesi</b> Istruzioni per l'uso
MK 32	<b>Biodesign® Графт за ректопексија</b> Упатство за употреба
NL 36	<b>Biodesign® Rectopexieprothese</b> Gebruiksaanwijzing
NO 39	<b>Biodesign® Rektopeksi-implantat</b> Bruksanvisning
PT-BR 42	<b>Biodesign® Enxerto de retopexia</b> Instruções de uso
SL 46	<b>Biodesign® Presadek za rektopeksijo</b> Navodila za uporabo
SR 49	<b>Biodesign® Mrežica za rektopeksiju</b> Uputstvo za upotrebu
SV 52	<b>Biodesign® Rektopexigraft</b> Bruksanvisning
TR 55	<b>Biodesign® Rektopeksi Grefti</b> Kullanma Talimatı

## BIODESIGN® RECTOPEXY GRAFT

### DEVICE DESCRIPTION

The Biodesign® Rectopexy Graft is a dried multi-layered small intestinal submucosa (SIS) sheet. The graft is used to reinforce soft tissue for the repair of rectal prolapse and rectal intussusception. Once implanted, the natural composition of SIS allows for the patient's healing mechanisms to deposit cells and collagen during cellular and extracellular matrix metabolism. Over the course of approximately 7 to 17 months, the patient's tissues completely incorporate the SIS until it can no longer be detected, resulting in new, remodeled tissue.<sup>1,3</sup> The graft can be cut to size to accommodate the patient's anatomy and is provided sterile for single use only. A study using the Biodesign Rectopexy Graft for rectopexy procedures revealed high symptom improvement and patient satisfaction scores.<sup>4</sup>

### INTENDED USE

The Biodesign Rectopexy Graft is intended to support/reinforce soft tissue in surgical procedures for open and laparoscopic repair of rectal prolapse/rectal intussusception. The graft is supplied sterile and is intended for one time use.

**Rx ONLY** This symbol means the following:

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**RECTOPEXY GRAFT** This symbol means the following: Rectopexy Graft

**MR** This symbol means the following: Magnetic Resonance Safe

This product is intended for use by medical professionals trained in colorectal surgery.

### CONTRAINDICATIONS

The graft is derived from a porcine source and should not be used for patients sensitive to porcine material.

### PRECAUTIONS

- The graft is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to graft failure and/or transmission of disease.
- **Do not resterilize.** Discard all open and unused portions of the graft.
- The graft is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard the graft if mishandling has caused possible damage or contamination, or if the graft is past its expiration date.
- Ensure that the graft is hydrated prior to suturing, stapling, or loading the graft laparoscopically.
- Place the graft in maximum possible contact with healthy, well-vascularized tissue to encourage cell ingrowth and tissue remodeling.
- Suturing, stapling, or tacking more than one graft together may decrease graft performance.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the graft.
- Extended hydration or excessive handling could lead to partial delamination of superficial layers of the graft.
- Care should be taken to avoid damage to the graft when loading laparoscopically. Loading through a 5 mm or larger port is recommended.
- Avoid implantation of the graft in a grossly infected or abscessed surgical field.
- No studies have been conducted to evaluate the performance of the graft in pediatric populations or in patients with severe intra-abdominal adhesions.

### GENERAL

- Users should be familiar with surgical technique for rectal prolapse repair, specifically the ventral rectopexy procedure.
- Users should exercise good surgical practice for the management of clean-contaminated, contaminated or infected fields.

### POTENTIAL COMPLICATIONS

Complications that can occur with the use of any surgical graft may include:

- Abscess
- Adhesions
- Allergic reaction
- Bowel obstruction
- Bowel or vaginal erosion
- Bowel or vaginal perforation
- Constipation
- De novo stress urinary incontinence
- Fever
- Impaction
- Infection
- Inflammation
- Migration
- Nerve damage
- Osteomyelitis
- Pain
- Recurrent prolapse
- Periostitis
- Seroma formation
- Spondylodiscitis

- Incisional herniation
- Induration
- Urinary retention

If any of these conditions occur and cannot be resolved, graft removal should be considered.

## STORAGE

The graft should be stored in a clean, dry location at room temperature.

## STERILIZATION

The graft has been sterilized with ethylene oxide and should not be resterilized.

## INSTRUCTIONS FOR USE

### Required Materials

- Sterile basin
- Sterile smooth forceps
- Hydration fluid: room temperature, sterile saline or sterile lactated Ringer's solution
- Suitable resorbable suture, such as: 2-0 or 0 polydioxanone (PDS) or coated polyglycolic acid suture (coated PGA), tacks, and/or screws
- Measuring instrument

**NOTE:** Handle the graft using aseptic technique, minimizing contact with latex gloves.


## PREPARATION

1. Remove the packaging containing the graft from the envelope.
2. Using aseptic technique, open the outer pouch and pass the inner pouch containing the graft onto the sterile field.
3. Open the pouch and shape the graft to accommodate the patient's anatomy; typically this is 4 cm wide by 18-20 cm long.<sup>5</sup>


**NOTE:** The recommended practice for preoperative bowel preparation in elective colorectal surgery includes mechanical bowel cleansing through the use of enemas and cathartic agents, and administration of prophylactic oral or intravenous antimicrobial agents. Insufficient cleansing or inadequate antibacterial prophylaxis may predispose the patient to infections.<sup>6-8</sup>

## PROCEDURAL


1. Perform under regional or general anesthesia.
2. Position the patient in the modified lithotomy or Lloyd Davies position with both arms secured on either side.
3. If performing the procedure laparoscopically, place the ports according to the surgeon's preference.
4. Retract uterus, sigmoid colon, and vagina.

5. Begin dissection at the sacral promontory. Take care to avoid damage to the left hypogastric nerve and left iliac vein.

6. Continue dissection distally to the right side of the mesorectum to the pouch of Douglas. Take care to avoid damage to the right iliac vessels and ureter.

7. Incise the pouch of Douglas to develop a plane from the rectovaginal septum to the pelvic floor.
8. In males, limit the ventral dissection in the rectovesical pouch to the apex of the prostate and avoid lateral dissection around the seminal vesicles.
9. Hydrate the graft in sterile saline or sterile lactated Ringer's solution until the desired handling characteristics are achieved. A hydration time of greater than 1 minute is not required.
10. Deliver the distal end of the graft to the ventral rectum. Secure the graft to the anterolateral rectum with interrupted deep muscle sutures on each side of the graft. Absorbable sutures, such as 2-0 PDS, are recommended.

**NOTE:** Avoid full thickness suturing of the rectum.

11. Place additional sutures on the left side of the intraperitoneal rectum.
12. Secure the graft to the sacral promontory under minimal tension with sutures, tacks, or screws.
13. Close the peritoneal tissue over the graft. Ensure maximum contact between the graft and tissue to produce better surgical outcomes.

14. Following confirmation of hemostasis, release the sigmoid colon and pelvic structures from retraction, and close the port sites.
15. Discard any unused portion of the graft according to institutional guidelines for medical waste.

### POST-OPERATIVE CARE

To provide the best environment for tissue integration into the graft, patient activity should be minimized. Provide patients with a list of post-procedural care recommendations. The following patient guidelines should be considered:

1. Patients should avoid any strenuous physical activity beyond a gentle walk for at least 2 weeks following rectal prolapse repair.
2. Patients should avoid any heavy lifting over 5 kg (10 lbs) for at least 4 weeks following rectal prolapse repair.
3. Patients should use a stool softener for at least 4 weeks after surgery.
4. The patient may resume sexual activity when comfortable.

### REFERENCES

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
3. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
4. Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis.* 2018;33(4):449-457.
5. Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis.* 2014;16(2):82-88.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.

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### LITERATURA

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
3. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.

- Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis.* 2018;33(4):449-457.
- Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis.* 2014;16(2):82-88.
- Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
- Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
- Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.

- Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
- Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
- Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis*. 2018;33(4):449-457.
- Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis*. 2014;16(2):82-88.
- Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
- Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery—comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.
- Berrios-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg*. 2017;152(8):784-791.

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### REFERENCIAS

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol*. 2004;172(1):215-218.
2. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology*. 2003;62(5):805-809.
3. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
4. Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis*. 2018;33(4):449-457.
5. Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis*. 2014;16(2):82-88.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery—comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.
8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg*. 2017;152(8):784-791.

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## RÉFÉRENCES

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
3. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
4. Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis.* 2018;33(4):449-457.
5. Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis.* 2014;16(2):82-88.

- Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
- Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
- Berrios-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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## LITERATURA

- Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
- Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
- Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
- Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis.* 2018;33(4):449-457.
- Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis.* 2014;16(2):82-88.
- Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
- Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
- Berrios-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

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### **REFERÊNCIAS**

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.

- Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia*. 2004;8(3):186-189.
- Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis*. 2018;33(4):449-457.
- Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis*. 2014;16(2):82-88.
- Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis*. 1997;24(4):609-619.
- Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci*. 1997;44(4):75-80.
- Berrios-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg*. 2017;152(8):784-791.

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### REFERANSLAR

1. Wiedemann A, Otto M. Small intestinal submucosa for pubourethral sling suspension for the treatment of stress incontinence: First histopathological results in humans. *J Urol.* 2004;172(1):215-218.
2. Rutner AB, Levine SR, Schmaelzle JF. Processed porcine small intestine submucosa as a graft material for pubovaginal slings: durability and results. *Urology.* 2003;62(5):805-809.
3. Franklin ME, JR., Gonzalez JJ, Jr., Glass JL. Use of porcine small intestinal submucosa as a prosthetic device for laparoscopic repair of hernias in contaminated fields: 2-year follow-up. *Hernia.* 2004;8(3):186-189.
4. Brunner M, Roth H, Gunther K, Grutzmann R, Matzel KE. Ventral rectopexy with biological mesh: short-term functional results. *Int J Colorectal Dis.* 2018;33(4):449-457.
5. Mercer-Jones MA, D'Hoore A, Dixon AR, et al. Consensus on ventral rectopexy: Report from a panel of international experts. *Colorectal Dis.* 2014;16(2):82-88.
6. Nichols RL, Smith JW, Garcia RY, et al. Current practices of preoperative bowel preparation among North American colorectal surgeons. *Clin Infect Dis.* 1997;24(4):609-619.
7. Yabata E, Okabe S, Endo M. A prospective, randomized clinical trial of preoperative bowel preparation for elective colorectal surgery-comparison among oral, systemic, and intraoperative luminal antibacterial preparations. *J Med Dent Sci.* 1997;44(4):75-80.
8. Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for Disease Control and Prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surg.* 2017;152(8):784-791.

**Keep dry**

**Manufacturer**

**Use-by date**

**Caution**

**Do not re-use**

**Magnetic resonance safe**

**Date of manufacture**

**Sterilized using ethylene oxide**

---

**MANUFACTURER**

COOK BIOTECH  
INCORPORATED  
1425 Innovation Place  
West Lafayette, IN 47906 U.S.A.

**EC REPRESENTATIVE**

Cook Biotech Europe ApS  
Sandet 6, DK-4632  
Bjaeverskov, DENMARK  
Tel: +45 56868686