- EN**
3 **Amplatz Ultra Stiff Wire Guides**
Instructions for Use
- BG**
4 **УлтраТВърди телени водачи Amplatz**
Инструкции за употреба
- CS**
6 **Ultra tuhé vodící dráty Amplatz**
Návod k použití
- DA**
7 **Amplatz ultrastive kateterledere**
Brugsanvisning
- DE**
9 **Amplatz ultrasteife Führungsdrähte**
Gebrauchsanweisung
- EL**
10 **Εξαιρετικά δύσκαμπτοι συρμάτινοι οδηγοί Amplatz**
Οδηγίες χρήσης
- ES**
12 **Guías ultrarrígidas Amplatz**
Instrucciones de uso
- ET**
14 **Ülijäigad juhtetraadid Amplatz**
Kasutusjuhend
- FR**
15 **Guides d'Amplatz ultra rigides**
Mode d'emploi
- HR**
17 **Ultrakrute žice vodilice Amplatz**
Upute za uporabu
- HU**
18 **Amplatz ultramerev vezetődrótok**
Használati utasítás
- IT**
20 **Guide ultra rigide Amplatz**
Istruzioni per l'uso
- LT**
22 **„Amplatz“ ultrastandieji vieliniai kreipikliai**
Naudojimo nurodymai
- LV**
23 **Amplatz īpaši stingras vadītājstīgas**
Lietošanas instrukcija
- NL**
25 **Amplatz ultrastugge voerdraden**
Gebruiksaanwijzing
- NO**
26 **Amplatz ultrastive ledevaiere**
Bruksanvisning
- PL**
28 **Ultrasztywne przewodniki Amplatz**
Instrukcja użycia

PT 29	Fios guia ultrarrígidos Amplatz Instruções de utilização
RO 31	Fire de ghidaj Amplatz Ultra Rigid Instrucțiuni de utilizare
RU 33	Проводники повышенной жесткости Amplatz Инструкция по применению
SK 34	Ultra tuhé vodiace drôty Amplatz Návod na použitie
SL 36	Ultra toga žična vodila Amplatz Navodila za uporabo
SR 37	Izuzetno krute žice vodiči Amplatz Uputstvo za upotrebu
SV 39	Amplatz ultrastyva ledare Bruksanvisning
TR 40	Amplatz Ultra Sert Kılavuz Teller Kullanma Talimatı
UK 42	Ультразорсткі провідники Amplatz Інструкція з застосування

AMPLATZ ULTRA STIFF WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Amplatz Ultra Stiff Wire Guides straight (THSF/-AUS1/AUS2) and curved (THSCF/-AUS1/AUS2) are PTFE-coated stainless steel wire guides with a length of 60-260 cm and a 3 cm (AUS1) or 7 cm (AUS2) flexible tip. The J-curve radius for the curved wire guides is 3 mm. Amplatz Ultra Stiff Wire Guide consists of a core mandril tapered on the distal end with a spring coil anchored over the entire mandril allowing utilization of an ultra-stiff mandril while maintaining tip flexibility.

The Amplatz Ultra Stiff Wire Guide is used both to assist in anatomical access for other devices (not included) and to support the delivery of medical devices. The wire guide is introduced into the target vessel; other devices, such as a sheath, catheter, stent, or endovascular graft can then be passed over the wire guide to be positioned or manipulated within the vascular system.

INTENDED USE

The Amplatz Ultra Stiff Wire Guides are intended to facilitate catheterization and/or placement of devices during vascular diagnostic procedures and vascular interventional procedures.

The Amplatz Ultra Stiff Wire Guides are intended for use in the major vessels, the aorta and vena cava, including their access vessels and adjacent vessels.

CONTRAINDICATIONS

Known allergic reactions to one of the materials in the device (e.g. nickel) should always be considered a contraindication.

WARNINGS

- Possible allergic reactions should always be considered.
- Wire guide stiffness and diameter should be considered.
- Not intended for vessels in the central nervous system or the coronary vessels.
- Insert wire guide from the flexible end. Inserting from the rigid end may cause damage to tissue or device.
- If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take appropriate steps to relieve the resistance to avoid

the risk of vessel perforation. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

- This single use device is not designed for re-use. If the device is re-used it could pose risks of cross contamination with microbiological agents. Attempts to reprocess, (resterilize) and/or re-use may lead to device failure and/or transmission of disease.

PRECAUTIONS

- Manipulation of wire guides requires fluoroscopic visualization.
- End-hole size and length of the catheter must be taken into consideration to ensure proper fit between wire guide and catheter.
- Altering tip configuration of the wire guide may result in damage.
- Withdrawal or manipulation of wire guide through a needle tip may result in damage or breakage.
- The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- If increased stiffness is needed, use e.g. Lunderquist Extra Stiff Wire Guides.

POTENTIAL ADVERSE EVENTS

Adverse events that may occur and/or require intervention include:

- Allergic reaction
- Arrhythmia or cardiac arrest
- Cardiac perforation
- Death
- Embolism
- Foreign body/wire fracture
- Hemorrhage or hematoma
- Infection or sepsis
- Organ injury
- Pseudoaneurysms
- Stroke or other neurological event
- Valvular/retrograde regurgitation and hypotension/ sudden drop in blood pressure
- Vascular access site complications
- Vascular thrombus/thromboembolism
- Vessel occlusion
- Vessel spasm
- Vessel trauma, dissection, perforation, rupture or injury
- Wire entrapment/entanglement
- Wound complications

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the wire guide from its outer packaging and place in the sterile field.
2. Rinse or wipe the wire guide with heparinized saline or sterile water prior to use.
3. Carefully insert the distal tip of the wire guide.
4. Under fluoroscopy and while maintaining position of the peripheral vascular access catheter, advance the wire guide to the targeted site.

NOTE: Under fluoroscopy, observe all wire guide movement in the vessel. The wire guide should be advanced only when visualizing the tip of the wire guide fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

5. With the wire guide secured in place, advance the percutaneous catheter or other therapeutic device to the target site.

NOTE: Always maintain a continuous flush while exchanging wire guides to prevent air entry into the system.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred. Do not use the device if damage has occurred or if the sterilization barrier has been damaged or broken; instead, return the product to Cook.

DISPOSAL OF DEVICE

This device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2022