EN
2**Frova Intubating Introducer**

Instructions for Use

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DA
6**Frova intubationsindfører**

Brugsanvisning

DE
8**Frova Intubationshilfe**

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

ENGLISH

FROVA INTUBATING INTRODUCER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Frova Intubating Introducer is a 14 French, 70 cm long, radiopaque catheter introducer with centimetre markings and a blunt, curved tip that can be passed blindly into the trachea. The catheter has a through lumen design with two distal sideports to ensure adequate airflow. A stiffening stylet and two Rapi-Fit® Adapters for connection to a ventilatory device are supplied as indicated below in **Table 1**.

Table 1

ORDER NUMBER	COMPONENTS
C-CAE-14.0-70-FII	Catheter introducer, stiffening stylet and two Rapi-Fit Adapters (15 mm connector and Luer lock connector)
C-CAE-14.0-70-FI	Catheter introducer and two Rapi-Fit Adapters (as above)
C-CAE-14.0-70-FIC	Catheter introducer only, supplied as a box of 10 catheters
C-CAE-14.0-70-FIC-SPOPS	Catheter introducer only, supplied in single units in a slip peel-pouch

Refer to label for further information.

INTENDED USE

The Frova Intubating Introducer is intended to facilitate endotracheal intubation in patients where the visualisation of the glottis is inadequate. The 14 French catheter introducer has been designed for placement of a single lumen endotracheal tube whose inner diameter is 6 mm or larger.

CONTRAINDICATIONS

Do not use if the epiglottis cannot be visualised when performing laryngoscopy, i.e., Grade IV Cormack & Lehane laryngoscopy classification.

WARNINGS

- Do not use the Frova Intubating Introducer with double lumen endotracheal or endobronchial tubes.
- To avoid barotrauma and/or pneumothorax, examine the patient’s anatomy to help determine the optimal placement for the Frova Intubating Introducer. Ensure the catheter is in the physician-preferred location relative to the carina by referencing its centimeter markings.
- Care must be taken not to provoke injury to the epiglottis and glottis, perforation of the sinus pyriformis, trachea or bronchus.
- The use of removable Rapi-Fit Adapters permits the utilisation of a high- and low-pressure oxygen source, if necessary, during the procedure. Use of the Rapi-Fit Adapter for oxygenation may be associated with a risk of barotrauma.
- Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume.
- If a high-pressure oxygen source is used for insufflation (e.g., jet ventilator), begin at a lower pressure and work up gradually. Rising chest wall, pulse oximetry and oral air flow should be carefully monitored. **Table 2** provides information on catheter oxygenation.
- Ensure that the Rapi-Fit Adapter is securely connected to the Frova catheter prior to oxygen delivery. Failure to properly secure the adapter to the catheter introducer may result in hypoxia, hypoxaemia and serious adverse events.

- Lubricate the catheter introducer and endotracheal tube before use.
- Ensure proper sizing of the endotracheal tube to be used in combination with the Frova Intubating Introducer.
- Care must be taken when introducing/removing the catheter introducer from the endotracheal tube; contact with sharp edges on the internal surface of the endotracheal tube may cause small fragments to be shaved off the catheter introducer during introduction/removal.
- Possible allergic reactions (e.g., to butyl rubber) should be considered.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in airway management.
- Sizing of the Rapi-Fit Adapter and the Frova Intubating Introducer must match in order to provide oxygenation through the catheter introducer.
- Not for intravascular use.

POTENTIAL ADVERSE EVENTS

- Allergic reaction (e.g., to butyl rubber)
- Barotrauma
- Hypoxia or hypoxaemia
- Injury to the epiglottis and glottis, perforation of the sinus pyriformis, trachea or bronchus.
- Pneumothorax

INSTRUCTIONS FOR USE

1. Lubricate an appropriately sized single lumen endotracheal tube and the distal tip of the catheter introducer.
2. Preload the endotracheal tube onto the proximal part of the catheter introducer.
3. Using laryngoscopy, introduce the tip of the catheter introducer beyond the epiglottis and advance it in a straight line towards the glottis.
4. Advance the catheter introducer into the trachea, approximately 2-3 cm. If resistance is encountered, do not force the catheter introducer; instead, it should be gently rotated and advanced.

5. Remove the stiffening stylet (not supplied with C-CAE-/-FI; C-CAE-/-FIC; C-CAE-/-FIC-SPOPS) and continue to advance the catheter introducer into the trachea for not more than 10 cm, the distance being determined by the size of the patient.

NOTE: During introduction you may feel the tracheal rings.

6. Confirm catheter position using standard methods (e.g., capnography, breath sounds, X-ray).
7. Whilst maintaining position of the catheter introducer, advance the preloaded endotracheal tube into the trachea to the appropriate distance.
8. Whilst maintaining position of the endotracheal tube, remove the catheter introducer and the laryngoscope.
9. Confirm correct position of the endotracheal tube in accordance with the endotracheal tube device labeling.

NOTE: If necessary, oxygen may be given during the procedure, by utilising the Rapi-Fit Adapters (not supplied with C-CAE-/-FIC; C-CAE-/-FIC-SPOPS).

USE OF RAPI-FIT ADAPTERS FOR OXYGENATION

Rapi-Fit Adapters should only be used when oxygen requirements are high and intubation is unsuccessful. Use of an oxygen source should only be considered if the patient has sufficient egression of the insufflated gas volume. If an oxygen source is used for insufflation, begin at a lower pressure and work up gradually. Observe the chest for outward and inward movements to confirm oxygen insufflation and egression. Pulse oximetry and oral air flow should be carefully monitored as well. In cases of upper airway obstruction, gas discharge from the patient's lungs may require more time. See **Table 2** for information on catheter oxygenation.

1. To attach the Rapi-Fit Adapter, position the adapter on the catheter introducer, then push the white collar forward and lock into position. (**Fig. 1**)
2. To remove the adapter, pull the white collar back to release, and then remove from the catheter introducer. (**Fig. 2**)

CATHETER OXYGENATION

The Frova Intubating Introducer is designed for both positive airway pressure ventilation (15 mm CPAP adapter) and jet ventilation (Luer lock adapter). In the table below, the I:E ratio, delivered minute volume and measured average maximum airway pressure are given for jet ventilation in adult patients with healthy lung tissue.

Table 2

Catheter RPN	I:E Ratio	Delivered Minute Volume ¹ (L/min)	Measured Average Maximum Airway Pressure ¹ (cm H ₂ O)
C-CAE-14.0-70-FII	1:4	8.6	12.1 ²

¹ The presented numbers are mean values and the following test conditions were used in an active model:
Input pressure set to: 50 psi, Lung Compliance: 100 mL/cm H₂O, Resistance: 3 cm H₂O/(L/s), tested with ASL 5000, Ingmar Medical, Ltd.

² Measured Average Maximum Airway Pressure did not exceed 14 cm H₂O in the testing performed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019