
EN	Microcoils™-18 Coil Pusher	Instructions for Use	3
CS	Pružinový tlačný mandrén Microcoils™-18	Návod k použití	3
DA	Microcoils™-18 coil-skubber	Brugsanvisning	4
DE	Microcoils™-18 Coil Pusher	Gebrauchsanweisung	5
EL	Διάταξη ώθησης σπειραμάτων Microcoils™-18	Οδηγίες χρήσης	6
ES	Empujador de espirales Microcoils™-18	Instrucciones de uso	7
ET	Microcoils™-18 spiraali tõukur	Kasutusjuhised	8
FR	Pousse-spirale Microcoils™-18	Mode d'emploi	8
HR	Alat za potiskivanje zavojnice Microcoils™-18	Upute za uporabu	9
HU	Microcoils™-18 tekerestoló	Használati utasítás	10
IT	Spingitore per spirale Microcoils™-18	Istruzioni per l'uso	11
LT	Ritės stūmiklis „Microcoils™-18“	Naudojimo nurodymai	12
LV	Microcoils™-18 spirāles bīdītājs	Lietošanas instrukcija	12
NL	Microcoils™-18 coilpusher	Gebruiksaanwijzing	13
NO	Microcoils™-18 spiralskyver	Bruksanvisning	14
PL	Popychacz spirali Microcoils™-18	Instrukcja użycia	15
PT	Empurrador de coil Microcoils™-18	Instruções de utilização	16
RO	Împingător de bobină Microcoils™-18	Instrucțiuni de utilizare	17
RU	Толкатель спирали Microcoils™-18	Инструкция по применению	17
SK	Nástroj na zatlačenie cievky Microcoils™-18	Návod na použitie	18
SR	Microcoils™-18 potiskivač spirale	Uputstvo za upotrebu	19
SV	Microcoils™-18 spiralspusher	Bruksanvisning	20
TR	Microcoils™-18 Sarmal İtici	Kullanma Talimatı	21

MICROCOILS™-18 COIL PUSHER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Illustrations

The Microcoils-18 Coil Pusher is a .018 inch diameter PTFE coated stainless steel mandril shaft with a 49 cm very soft distal spring coil tip, a 1 cm radiopaque marker is placed 1 cm from the distal tip. (Fig 1)

The product also contains a wire guide inserter, and a pin vise handle for manipulation.

Refer to label for further information.

INTENDED USE

The Microcoils-18 Coil Pusher is intended for placement of .018 inch microcoils through a microcatheter with an ID of .018 to .025 inches.

CONTRAINDICATIONS

None known.

WARNINGS

- Possible allergic reactions should be considered.
- Manipulation of products requires fluoroscopic control.
- To determine correct catheter position, perform an angiogram prior to embolization.
- Never advance the coil pusher outside the tip of the microcatheter after the embolization coil has been placed. This might displace the coil and damage the vessel.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Prior to introduction of the coil, flush the microcatheter with saline.

INSTRUCTIONS FOR USE

1. For coil loading instructions, please see IFU supplied with the coil.
2. Introduce the coil pusher through the wire guide inserter into the catheter.
3. Slowly advance the coil pusher under fluoroscopic control until the coil has reached the target lesion at the distal tip of the catheter. The coil must be advanced completely outside the microcatheter.

NOTE: Never advance the coil pusher outside the microcatheter tip.

4. To facilitate coil pusher control the radiopaque marker (1 cm from the tip) should remain inside the microcatheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. The coil pusher may be used to deliver multiple coils in one procedure. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

DISPOSAL OF DEVICE

This device is potentially contaminated with infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2021