EN

18

CS

22

DA

27

DE

32

EL

37

ES

42

FI

47

FR

52

HR

57

HU

62

IT

67

NL

72

Günther Tulip® Vena Cava Filter Retrieval Set for Jugular Vein Approach

Instructions for Use

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Bruksanvisning

TABLE OF CONTENTS

1. DEVICE DESCRIPTION	18
1.1 Construction of Device	18
1.2 Performance Characteristics	18
1.3 Specifications	18
1.4 Device Compatibility	18
1.5 Patient Population	18
1.6 Intended User	18
1.7 Contact with Body Tissue	18
1.8 Operating Principle	18
2. INTENDED USE	19
3. INDICATIONS FOR USE	19
4. CLINICAL BENEFITS	19
5. CONTRAINDICATIONS	19
6. WARNINGS	19
6.1 General Warnings	19
6.2 Sterile and Single Use	19
6.3 Malfunctions/Changes in Performance	19
7. PRECAUTIONS	19
7.1 General Precautions	19
7.2 Malfunctions/Changes in Performance	19
7.3 CMR 1A/1B and/or Endocrine Disrupting Substances	20
7.4 Allergic Reaction	20
7.5 Training Requirements for User	20
7.6 Requirements for Facilities	20
8. POTENTIAL ADVERSE EVENTS	20
9. CLINICAL STUDIES	20
10. STERILIZATION	20
11. HOW SUPPLIED	21
12. INSPECTION OF DEVICE	21
13. INSTRUCTIONS FOR USE	21
13.1 Patient Preparation	21
13.2 Step-by-Step Guide	21
13.2.1 Preparation	21
13.2.2 Filter Retrieval	21
14. REQUIRED EQUIPMENT AND ACCESSORIES	22
14.1 Materials Required	22
15. DISPOSAL OF DEVICE	22
16. REFERENCES	22

17. PATIENT COUNSELING	22
17.1 Symptoms	22
17.2 Follow-Up Information for the Patient	22
17.3 Summary of Safety and Clinical Performance (SSCP)	22
18. SERIOUS INCIDENT REPORTING	22

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







10

GÜNTHER TULIP® VENA CAVA FILTER RETRIEVAL SET FOR JUGULAR VEIN APPROACH

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

STERILE – DO NOT RESTERILIZE – SINGLE USE ONLY.

1. DEVICE DESCRIPTION

1.1 Construction of Device

The Günther Tulip® Vena Cava Filter Retrieval Set for Jugular Vein Approach (hereinafter referred to as the GTRS) consists of a retrieval loop system with a braided platinum retrieval loop, an 11 French coaxial retrieval sheath system, an 18 gage entry needle, a 0.035 inch wire guide, a 12 French pre-dilator, and a three-way stopcock. (Fig. 1) The coaxial retrieval sheath system consists of the outer sheath, inner sheath, and introducer dilator. The tip of the outer sheath has a radiopaque band.

1.2 Performance Characteristics

The GTRS is intended for percutaneous introduction

into the inferior vena cava (IVC) via a jugular vein access. The GTRS is designed for retrieval of implanted Günther Tulip, Cook Celect®, and Celect Platinum® Vena Cava Filters in patients who no longer require a filter. The loop snare is used to catch the filter by first expanding around the ball-tipped hook of the filter and then retracting until it engages the hook. The pin vise and Y-fitting are used to manipulate the retrieval loop during the procedure. The side-arm fitting is used to collapse the filter into the coaxial retrieval sheath system and to prevent blood loss.

The radiopaque band on the outer sheath of the coaxial retrieval sheath system is used to identify the precise location of the distal tip of the sheath for positioning accuracy. Likewise, the catheter for the angled wire loop retriever, the retrieval loop, and the inner sheath are radiopaque to facilitate fluoroscopic visualization of the device during the procedure. The inner sheath is used for injection of contrast media. The stopcock is used for injection of contrast media. If necessary, the pre-dilator may be used to widen the puncture site prior to insertion of the coaxial sheath system to minimize vessel trauma.

1.3 Specifications

The specifications for the GTRS are shown in Table 1.

Table 1 – Sizes of Parts Included in the Retrieval Set

Entry needle (Fig. 1a)	18 gage
Wire guide (Fig. 1b)	0.035 inch
Pre-dilator (Fig. 1c)	12 Fr, 20 cm long
Coaxial retrieval sheath system (Fig. 1d)	11 Fr

1.4 Device Compatibility

The GTRS is compatible with the Günther Tulip, Cook Celect, and Celect Platinum Vena Cava Filters, and a 0.035 inch wire guide.

1.5 Patient Population

The target population for the GTRS are patients with an implanted Günther Tulip, Cook Celect, or Celect Platinum Vena Cava Filter, who have been assessed as no longer requiring a filter for the prevention of pulmonary embolism.

1.6 Intended User

The product is intended for use by physicians trained and experienced in diagnostic and interventional endovascular techniques.

1.7 Contact with Body Tissue

The GTRS is an external communicating device, which is in direct contact with circulating blood for a limited contact duration (≤24 hours).

1.8 Operating Principle

The GTRS is indicated for retrieval of implanted Günther Tulip, Cook Celect, or Celect Platinum Vena Cava Filters in patients who no longer require a filter. The GTRS is used to capture the filter by engaging the retrieval loop on the filter. The filter is secured to the wire loop and then collapsed into the coaxial retrieval sheath system which is then retrieved.

Fluoroscopy should be used prior to insertion of the GTRS to evaluate positioning of the filter and be used during and after retrieval of the filter.

Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed. Local anesthesia, sedation, and analgesia will be used as required according to institution standard of care.

2. INTENDED USE

The Günther Tulip Vena Cava Filter Retrieval Set has been designed for retrieval of implanted Günther Tulip, Cook Celect, and Celect Platinum Vena Cava Filters. Retrieval of the filter can be performed only by jugular approach.

3. INDICATIONS FOR USE

The Günther Tulip Vena Cava Filter Retrieval Set is indicated for the retrieval of an implanted Günther Tulip, Cook Celect, or Celect Platinum Vena Cava Filter in patients assessed as no longer requiring a filter for prevention of pulmonary embolism.

4. CLINICAL BENEFITS

Filter retrieval is a patient-specific, clinically complex decision that should be based on the individual patient's risk/benefit profile. The clinical benefit of the GTRS is to facilitate endovascular retrieval of indwelling Günther Tulip, Cook Celect, and Celect Platinum Vena Cava Filters to remove the risk of filter related complications when the filter has been assessed as no longer needed for the prevention of PE.

5. CONTRAINDICATIONS

- Retrieval of the filter with significant amounts of trapped thrombus (greater than 25% of the volume of the cone).
- Retrieval of the filter for patients with an ongoing high risk for pulmonary embolism (PE).

6. WARNINGS

6.1 General Warnings

- An inferior vena caval imaging evaluation for residual captured thrombus should be performed prior to an attempted retrieval.
- Excessive force should not be exerted to retrieve the filter.
- The Cook Celect and the Celect Platinum Vena Cava Filters should be completely collapsed in the inner sheath before subsequently advancing the outer sheath.
- Never attempt to re-deploy a retrieved filter.

6.2 Sterile and Single Use

- Do not use the device if the sterile packaging is damaged or unintentionally opened before use. If a contaminated device is introduced to the bloodstream, it may cause an infection and lead to severe harm.

- This single-use device is not designed for reuse. If the device is reused it could pose risks of cross contamination with microbiological agents. Attempts to reprocess (resterilize) and/or to reuse may lead to device failure and/or transmission of disease.

6.3 Malfunctions/Changes in Performance

- For the Günther Tulip Vena Cava Filter, do not advance the inner sheath of the coaxial retrieval sheath over the anchors of the filter as doing so may cause particles to scratch off the sheath.

7. PRECAUTIONS

7.1 General Precautions

- Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- Manipulation of the GTRS requires imaging control.
- Follow the instructions thoroughly to ensure successful retrieval, and to avoid any harm to the patient or damage to the device.
- For filter retrieval, the right jugular vein is usually preferred due to its straighter route to the vena cava.
- Filter tilt, vena cava wall penetration/perforation, filter fracture, filter migration, and/or filter/fragment embolization have been reported. These observations may be symptomatic or asymptomatic and may contribute to the complexity of the attempted filter retrieval procedure.
- Failure to store the device correctly may result in material degradation and/or damage to the device.

7.2 Malfunctions/Changes in Performance

- Device modification or alteration is not recommended, as the device's safety and effectiveness has not been established following any modifications.
- Do not pull on the filter beyond what is required to keep tension on the loop. Doing so may cause damage to the caval wall.
- Do not retract the retrieval loop system before the tip of the coaxial retrieval sheath is at the filter anchors. Doing so may cause damage to the caval wall.
- Before the final retrieval step, ensure that the entire filter is covered by the outer sheath to prevent the filter anchors damaging the caval wall.

7.3 CMR 1A/1B and/or Endocrine Disrupting Substances

- This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and carcinogenic (Category 1B). However, the device contains cobalt as part of a stainless steel alloy at a very low concentration (up to 0.4% w/w), it has limited direct exposure (≤ 24 hours) and does not release cobalt at levels which cause an increased risk of cancer or adverse reproductive effect to the patient.

- This symbol on the label indicates that the device contains lead (Pb) at a level above 0.1% w/w, which is a substance that is toxic to reproduction (Category 1A). However, the device contains lead as part of a brass collet at a concentration up to 2.5% w/w, it has limited indirect exposure (≤ 24 hours) and does not release lead at levels which cause an increased risk of adverse reproductive effects to the patient.

7.4 Allergic Reaction

- Possible allergic reactions to cobalt, chromium, nickel, and platinum should be considered.
- Ensure that the patient does not have impaired tolerance to general, regional, or local anesthesia to avoid adverse reactions associated with the anesthetic procedure.
- Ensure that the patient is not allergic/sensitive to contrast media since the use of contrast media during the procedure and/or during postoperative imaging may cause an allergic reaction and/or other contrast-induced harms.

7.5 Training Requirements for User

The GTRS is intended for use by physicians trained and experienced in diagnostic and interventional endovascular techniques.

7.6 Requirements for Facilities

The device is intended for use at hospital operating rooms. Users must wear standard operating room personal protective equipment (PPE; gloves, face mask, sterile gown, etc.) and follow local guidelines for sterile procedures.

8. POTENTIAL ADVERSE EVENTS

- Access site related complications
- Air embolism
- Arrhythmia
- Arteriovenous fistula
- Cardiac damage
- Death
- Device failure
- Embolization
- Hematoma
- Hemorrhage
- Infection at vascular access site
- Pain
- Pneumothorax
- Pulmonary embolism
- Radiation exposure
- Thrombosis
- Vena cava stenosis
- Vessel injury

9. CLINICAL STUDIES

Physician practice guidelines and published guidance from regulatory agencies recommend that patients with indwelling filters undergo routine follow-up. The risks/benefits of filter retrieval should be considered for each patient during follow-up. Refer to **Section 16, REFERENCES** for citations that include recommendations related to filter follow-up and retrieval.

Once protection from PE is no longer necessary, filter retrieval should be considered. Filter retrieval should be attempted when feasible and clinically indicated. Filter retrieval is a patient-specific, clinically complex decision; the decision to remove a filter should be based on each patient's individual risk/benefit profile (e.g., a patient's continued need for protection from PE compared to their experience with and/or ongoing risk of experiencing filter-related complications). For all retrievable IVC filters, retrieval becomes more challenging with time, and this is commonly due to encapsulation of the filter legs or hook (in a tilted filter) by tissue ingrowth.

Please refer to the **CLINICAL STUDIES** section of the Instructions for Use for Günther Tulip, Cook Celect, or Celect Platinum Vena Cava Filters for further data on filter retrieval; clinical outcomes are described for the Günther Tulip, Cook Celect, or Celect Platinum Vena Cava Filters.

10. STERILIZATION

All devices are sterilized using ethylene oxide (EO) gas.

11. HOW SUPPLIED

Supplied in peel-open packages. Keep the device dry and away from sunlight. Do not use after the expiration date printed on the label.

12. INSPECTION OF DEVICE

Inspect the device thoroughly including all levels of the packaging to verify that there is no damage prior to use. Inspect and confirm that the sterile barrier has not been compromised in any way. Inspect and confirm that the device corresponds to the label and the IFU.

Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

13. INSTRUCTIONS FOR USE

13.1 Patient Preparation

The physician should follow institutional protocols and local guidelines.

13.2 Step-by-Step Guide

13.2.1 Preparation

1. Flush the pre-dilator, coaxial retrieval sheath system, and retrieval loop system.
2. On the retrieval loop system, ensure that the retrieval loop is fully expanded and tighten the pin vise. (Fig. 2)
3. Hold the Y-fitting and pull back the plastic pin vise to retract the retrieval loop into the catheter. Tighten the screw of the Y-fitting to keep the retrieval loop inside the catheter. (Fig. 3)

13.2.2 Filter Retrieval

4. Access the chosen jugular vein using the Seldinger technique.
5. Place a supportive 0.035 inch wire guide in the IVC.
6. If necessary, dilate the puncture site with the 12 French pre-dilator over the wire guide.
7. Remove the pre-dilator and advance a flush catheter (not included in the set).
8. Remove the wire guide.
9. Perform diagnostic imaging to evaluate the filter position and check for thrombus.
10. After confirming the plan for filter retrieval on diagnostic imaging, insert the wire guide.
11. Exchange the flush catheter for the coaxial retrieval sheath system, advancing it over the wire guide to a point approximately 3-5 cm above the hook of the filter.
12. Remove the introducer dilator and the wire guide from the coaxial retrieval sheath system. Verify the position of the coaxial retrieval sheath by diagnostic imaging. The three-way stopcock should be attached to the inner sheath hub during injection of contrast media. (Fig. 4)

13. Introduce the retrieval loop system through the coaxial retrieval sheath, then advance and connect the side-arm fitting of the loop system to the sheath. Tighten the side-arm fitting adaptor around the catheter to prevent blood loss. (Fig. 5)

14. Loosen the screw of the Y-fitting. Hold the Y-fitting and advance the pin vise until the retrieval loop has fully expanded inside the vena cava and surrounds the filter. (Fig. 6)

15. Pull back the retrieval loop until it engages the hook of the filter. (Fig. 7)

CAUTION: Do not pull on the filter beyond what is required to keep tension on the retrieval loop. Doing so may cause damage to the caval wall.

16. Hold the retrieval loop steady with the pin vise, then push the Y-fitting with the catheter forward, still ensuring that the filter hook is caught in the retrieval loop. Lock the Y-fitting to secure the snare around the filter hook. (Fig. 8)

NOTE: If the retrieval loop loses its shape during the retrieval attempt, it can be removed and gently reshaped. After reshaping, clean the retrieval loop and proceed from step 13.

17. While holding the retrieval loop and Y-fitting steady, advance the side-arm fitting to collapse the filter into the coaxial retrieval sheath and disengage the anchors of the filter from the caval wall. (Fig. 9)

CAUTION: Do not retract the retrieval loop system at this step. Doing so may cause damage to the caval wall.

18. When the tip of the coaxial retrieval sheath is at the filter anchors, loosen the hub of the outer sheath and advance the outer sheath forward to cover the entire filter. (Fig. 10)

WARNING: Specifically for the Günther Tulip Vena Cava Filter do not advance the inner sheath over the anchors of the filter; doing so may cause particles to scratch off the sheath.

CAUTION: Ensure that the entire filter is covered by the outer sheath to prevent the filter anchors from damaging the caval wall.

19. Retrieve the entire assembly leaving only the outer sheath in place.
20. Perform diagnostic imaging to assess the inferior vena cava for potential injury following filter retrieval.

NOTE: Hospital standard of care should be followed for removing the complete assembly and providing hemostasis to prevent bleeding at the vascular access site.

14. REQUIRED EQUIPMENT AND ACCESSORIES

14.1 Materials Required

- Flush catheter
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Flushing solution
- Sterile gauze pads

15. DISPOSAL OF DEVICE

After the procedure, this device and its accessories may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

16. REFERENCES

These instructions for use are based on experience from physicians and/or their published literature, IVC filter guidelines, ISO 25539-3, and regulatory safety communications regarding IVC filters. Refer to your local Cook Medical sales representative for information on available literature.

The published clinical literature includes descriptions of alternative techniques for filter retrieval (e.g., techniques that vary from that described in this instruction for use); use of these techniques varies according to physician experience, patient anatomy, and filter position. **The safety or effectiveness of these alternative retrieval techniques has not been established.**

17. PATIENT COUNSELING

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

Furthermore, the patient should be informed about the residual risks when treated with this device and the potential adverse events related to the use of this device.

17.1 Symptoms

The physician must advise the patient to contact their healthcare provider if he/she experiences changes in their medical condition or in case of concern.

17.2 Follow-Up Information for the Patient

Follow-up will be individual, so the physician must advise the patient on the recommended follow-up.

17.3 Summary of Safety and Clinical Performance (SSCP)

You can use the basic unique device identifier (BUDI)

for this device (0827002-WCE-MDF0000638-H7) to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: ec.europa.eu/tools/eudamed.

18. SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, and the competent authority and/or regulatory authority of the country where the device was used.

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© COOK 2025
Australian Sponsor**William A. Cook Australia Pty Ltd**