EN
2

Lunderquist® Extra Stiff Wire Guides

Instructions for Use

CS
3

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Bruksanvisning

LUNDERQUIST® EXTRA STIFF WIRE GUIDES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Lunderquist® Extra Stiff Wire Guide (TS(C)MG-/-LES) is a TFE-coated stainless steel wire guide. Lengths ≥ 260 cm have 4 cm of tip flexibility including an inner gold coil for enhanced visibility. Lengths < 260 cm the distal tip have 4 or 7 cm of tip flexibility, and the TS(C)MG-/-LES1 3, 5, 7 cm flexible tip. The Lunderquist Extra Stiff DC Wire Guide (TSCMG-/-E)-LESDC is a TFE-coated stainless steel wire guide with a double curved tip design. The double curved distal tip has 4 cm of tip flexibility including an inner gold coil for enhanced visibility. LES-DC has a primary/secondary curve radius on 75/15 mm, and E-LESDC has a primary/secondary curve radius on 55/15 mm.

The wire guide is used both to assist anatomical access for other devices (not included), and to support the delivery of medical devices. The wire guide is introduced into the target vessel; other devices, such as a sheath, catheter, stent, or endovascular graft can then be passed over the wire guide to be positioned or manipulated within the vascular system.

INTENDED USE

The device is intended for complex diagnostic and interventional procedures where increased body, flexibility, and low surface friction of the wire guide is needed.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known.

WARNINGS

- Possible allergic reactions should always be considered.
- Wire guide stiffness should be considered.
- Insert wire guide from the flexible end. (Inserting from the rigid end may cause damage to tissue or device).
- If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take appropriate steps to relieve the resistance to avoid the risk of vessel perforation. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

PRECAUTIONS

- Manipulation of wire guides always requires fluoroscopic visualization.
- End hole size and length of the catheter must be taken into consideration to ensure proper fit between wire guide and catheter.
- Altering tip configuration of the wire guide may result in damage.
- Withdrawal or manipulation of wire guide through a needle tip may result in damage or breakage.
- The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

POTENTIAL ADVERSE EVENTS

Using this product may cause the following adverse events:

- Allergic reaction
- Arrhythmia or cardiac arrest
- Hemorrhage or hematoma
- Infection
- Organ injury
- Product failures (broken or damaged wire guide)
- Pseudoaneurysms
- Vascular thrombus
- Vessel spasm
- Vessel trauma, dissection, perforation, rupture or injury

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. Attach a syringe with heparinized saline solution or sterile water to the Luer Lock fitting of the wire guideholder.
3. Inject enough solution to fill the holder and bathe the wire in solution.
4. Carefully insert the distal tip of the wire guide in the peripheral vascular access catheter.
5. Under fluoroscopy and while maintaining position of the peripheral vascular access catheter advance the wire to the targeted site.

NOTE: Under fluoroscopy, observe all wire guide movement in the vessel. The wire guide should be advanced only when visualizing the tip of the wire guide fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance to avoid the risk of perforation. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

6. With the wire guide secured in place, advance the therapeutic device to the target site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light.

Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2013