- EN**
3 **Lunderquist® Extra Stiff Wire Guides**
Instructions for Use
- BG**
4 **Екстра твърди телени водачи Lunderquist®**
Инструкции за употреба
- CS**
6 **Lunderquist® extra tuhé vodící dráty**
Návod k použití
- DA**
8 **Lunderquist® ekstrastive kateterledere**
Brugsanvisning
- DE**
10 **Lunderquist® extrasteife Führungsdrähte**
Gebrauchsanweisung
- EL**
12 **Εξαιρετικά δύσκαμπτοι συρμάτινοι οδηγοί Lunderquist®**
Οδηγίες χρήσης
- ES**
14 **Guías extrarrígidas Lunderquist®**
Instrucciones de uso
- ET**
15 **Eriti jäigad juhtetraadid Lunderquist®**
Kasutusjuhend
- FR**
17 **Guides Lunderquist® extra rigides**
Mode d'emploi
- HR**
19 **Ekstra krute žice vodilice Lunderquist®**
Upute za uporabu
- HU**
21 **Lunderquist® extramerev vezetődrótok**
Használati utasítás
- IT**
22 **Guide extra rigide Lunderquist®**
Istruzioni per l'uso
- LT**
24 **„Lunderquist®“ itin standūs vieliniai kreipikliai**
Naudojimo nurodymai
- LV**
26 **Lunderquist® īpaši stingras vadītājstīgas**
Lietošanas instrukcija
- NL**
28 **Lunderquist® extra stugge voerdraden**
Gebruiksaanwijzing
- NO**
30 **Lunderquist® ekstra stive ledevaiere**
Bruksanvisning
- PL**
31 **Ekstrasztwyne przewodniki Lunderquist®**
Instrukcja użycia
- PT**
33 **Fios guia extrarrígidos Lunderquist®**
Instruções de utilização
- RO**
35 **Fire de ghidaj Lunderquist® Extra Rigid**
Instrucţiuni de utilizare
- RU**
37 **Сверхжесткие проводники Lunderquist®**
Инструкция по применению
- SK**
39 **Extra tuhé vodiace drôty Lunderquist®**
Návod na použitie
- SL**
40 **Zelo toga žična vodila Lunderquist®**
Navodila za uporabo
- SR**
42 **Posebno krute žice vodiči Lunderquist®**
Uputstvo za upotrebu
- SV**
44 **Lunderquist® extra styva ledare**
Bruksanvisning

TR
46

Lunderquist® Ekstra Sert Kılavuz Teller

Kullanma Talimatı

UK
47

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LUNDERQUIST® EXTRA STIFF WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

This group of wire guides contain the following products:

Straight (TSMG-/LES1) and Curved (TSCMG-/LES1) Lunderquist Extra Stiff Wire Guide. PTFE-coated stainless steel wire guides with a length of 90 cm and a 3 cm, 5 cm or 7 cm flexible tip. The J-curve radius for the curved wire guides is 3 mm.

Straight (TSMG-/LES) and Curved (TSCMG-/LES) Lunderquist Extra Stiff Wire Guide. PTFE-coated stainless steel wire guides with lengths ranging from 90 cm to 300 cm and either a 4 cm or 7 cm flexible tip. For the 260 cm and 300 cm lengths, the flexible tip includes an inner gold coil for enhanced visibility. The J-curve radius for the curved wire guides is either 3 mm or 7 mm. Double curved (TSCMG-/LESDC) and extended double curved (TSCMG-/E-LESDC) Lunderquist Extra Stiff DC Wire Guide. PTFE-coated stainless steel wire guides with a 4 cm flexible tip that includes an inner gold coil for enhanced visibility. TSCMG-/LESDC has a primary/secondary curve radius on 75/15 mm, and TSCMG-/E-LESDC has an extended primary/secondary curve radius on 55/15 mm.

The Lunderquist Extra Stiff Wire Guide is used both to assist in anatomical access for other devices (not included) and to support the delivery of medical devices. The wire guide is introduced into the target vessel; other devices, such as a sheath, catheter, stent, or endovascular graft can then be passed over the wire guide to be positioned or manipulated within the vascular system.

The Lunderquist Extra Stiff Wire Guides are available in different lengths, shapes, etc. to accommodate different anatomical applications.

INTENDED USE

The Lunderquist Extra Stiff Wire Guides are intended to facilitate catheterization and/or placement of devices during vascular diagnostic procedures and vascular interventional procedures.

The Lunderquist Extra Stiff Wire Guides are intended for use in the major vessels, the aorta and vena cava, including their access vessels and adjacent vessels.

CONTRAINDICATIONS

Known allergic reactions to one of the materials in the device (e.g. nickel) should always be considered a contraindication.

WARNINGS

- Possible allergic reactions should always be considered.
- Wire guide stiffness and diameter should be considered.
- Not intended for vessels in the central nervous system or the coronary vessels.
- Insert wire guide from the flexible end. Inserting from the rigid end may cause damage to tissue or device.
- If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take appropriate steps to relieve the resistance to avoid the risk of vessel perforation. Advancement and withdrawal of the wire guide should be performed slowly and carefully.
- This single use device is not designed for re-use. If the device is re-used it could pose risks of cross contamination with microbiological agents. Attempts to reprocess, (resterilize) and/or re-use may lead to device failure and/or transmission of disease.

PRECAUTIONS

- Manipulation of wire guides always requires fluoroscopic visualization.
- End-hole size and length of the catheter must be taken into consideration to ensure proper fit between the wire guide and catheter.
- Altering the tip configuration of the wire guide may result in damage.
- Withdrawal or manipulation of the wire guide through a needle tip may result in damage or breakage.
- The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques.
- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- If a less stiff wire guide is needed, use e.g. Amplatz Ultra Stiff Wire Guides.

POTENTIAL ADVERSE EVENTS

Adverse events that may occur and/or require intervention include:

- Allergic reaction
- Arrhythmia or cardiac arrest
- Cardiac perforation
- Death
- Embolism
- Foreign body/wire fracture

- Hemorrhage or hematoma
- Infection or sepsis
- Organ injury
- Pseudoaneurysm
- Stroke or other neurological event
- Valvular/retrograde regurgitation and hypotension/sudden drop in blood pressure
- Vascular access site complications
- Vascular thrombus/thromboembolism
- Vessel occlusion
- Vessel spasm
- Vessel trauma, dissection, perforation, rupture or injury
- Wire entrapment/entanglement
- Wound complications

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the wire from its outer packaging and place it in the sterile field.
2. Attach a syringe with heparinized saline solution or sterile water to the Luer lock fitting of the wire guide holder.
3. Inject enough solution to fill the holder and bathe the wire in solution.
4. Carefully insert the distal tip of the wire guide in the peripheral vascular access catheter.
5. Under fluoroscopy, maintain the position of the peripheral vascular access catheter and advance the wire guide to the targeted site.

NOTE: Under fluoroscopy, observe all wire guide movement in the vessel. Do not advance the wire guide if you cannot visualize its tip fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance, to avoid the risk of perforation. Be sure to advance and withdraw the wire guide slowly and carefully.

6. With the wire guide secured in place, advance the therapeutic device to the target site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred. Do not use the device if damage has occurred or if the sterilization barrier has been damaged or broken; instead, return the product to Cook.

DISPOSAL OF DEVICE

This device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2022
William A. Cook Australia Pty Ltd