- EN 13** **Zenith® Dissection Endovascular Stent with the Z-Trak® Plus Introduction System**
Instructions for Use
- CS 17** **Endovaskulární stent Zenith® Dissection se zaváděcím systémem Z-Trak® Plus**
Návod k použití
- DA 21** **Zenith® Dissection endovaskulær stent med Z-Trak® Plus indføringssystemet**
Brugsanvisning
- DE 25** **Zenith® Dissection endovaskulärer Stent mit Z-Trak® Plus Einführsystem**
Gebrauchsanweisung
- EL 30** **Ενδαγγειακή ενδοπρόσθεση Zenith® Dissection με το σύστημα εισαγωγής Z-Trak® Plus**
Οδηγίες χρήσης
- ES 35** **Stent endovascular Zenith® Dissection con el sistema de introducción Z-Trak® Plus**
Instrucciones de uso
- FR 40** **Endoprothèse vasculaire Zenith® Dissection à système d'introduction Z-Trak® Plus**
Mode d'emploi
- HU 45** **Z-Trak® Plus felvezetőrendszerrel ellátott Zenith® Dissection endovaszkuláris sztent**
Használati utasítás
- IT 50** **Stent endovascolare Zenith® Dissection con sistema di introduzione Z-Trak® Plus**
Istruzioni per l'uso
- NL 55** **Zenith® Dissection endovasculaire stent met het Z-Trak® Plus introductiesysteem**
Gebruiksaanwijzing
- NO 60** **Zenith® Dissection endovaskulær stent med Z-Trak® Plus innføringssystem**
Bruksanvisning
- PL 64** **Stent wewnętrznożylniowy Zenith® Dissection z systemem wprowadzającym Z-Trak® Plus**
Instrukcja użycia
- PT 69** **Stent endovascular Zenith® Dissection com o sistema de introdução Z-Trak® Plus**
Instruções de utilização
- SV 74** **Zenith® Dissection endovaskulärt stent med Z-Trak® Plus införingssystem**
Bruksanvisning

TABLE OF CONTENTS

1 DEVICE DESCRIPTION	13
1.1 The Zenith Dissection Endovascular Stent	13
Table 1 Zenith Dissection Endovascular Stent.....	13
1.2 Zenith Dissection Endovascular Stent with the Z-Trak Plus Introduction System	13
2 INTENDED USE	13
Introduction System Sizes.....	13
3 CONTRAINDICATIONS	13
4 WARNINGS AND PRECAUTIONS	13
4.1 General	13
4.2 Patient Selection, Treatment and Follow-Up	13
4.3 Implant Procedure	13
4.4 MRI Information	14
5 POTENTIAL ADVERSE EVENTS	14
6 POTENTIAL RISKS AND BENEFITS	14
7 PATIENT SELECTION AND TREATMENT	14
7.1 Individualization of Treatment	14
8 PATIENT COUNSELING INFORMATION	14
9 HOW SUPPLIED	15
10 CLINICAL USE INFORMATION	15
10.1 Physician Training	15
Patient selection	15
A multidisciplinary team that has combined procedural experience with.....	15
10.2 Inspection Prior to Use	15
10.3 Materials Recommended (Not included with the Zenith Dissection Endovascular Stent)	15
10.4 Device Length Selection Guidelines	15
10.5 Device Diameter Selection Guidelines	15
11 DIRECTIONS FOR USE	15
Pre-Implant Determinants	15
Device Preparation	15
11.1 The Zenith Dissection Endovascular Stent System	15
11.1.1 Preparation/Flush	15
11.1.2 Placement	15
11.2 Final Angiogram	15
12 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP	15
12.1 General	15
Table 2 Recommended Imaging Schedule for Endograft Patients.....	16
12.2 Contrast and Non-Contrast CT Recommendations	16
Table 3 Acceptable Imaging Protocols.....	16
12.3 Thoracic and Abdominal Device Radiographs	16
12.4 MRI Information	16
12.5 Additional Surveillance and Treatment	16
13 REFERENCES	16

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
7b

8

ZENITH® DISSECTION ENDOVASCULAR STENT WITH THE Z-TRAK® PLUS INTRODUCTION SYSTEM

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.
- CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- CAUTION: All contents of the inner pouch (including the introduction system and endovascular stent) are supplied sterile, for single use only.

Table 1 Zenith Dissection Endovascular Stent

Cook Reorder Number	Stent Diameter	Introducer Sheath Size ID Fr/OD mm	Stent Length (at nominal diameter)	Stent Length Maximum (at 20/28 mm diameter)	Introducer Sheath Length
ZDES-36-80	36 mm	16 Fr/6.0 mm	80 mm	91 mm at 20 mm	100 cm
ZDES-36-120	36 mm	16 Fr/6.0 mm	120 mm	136 mm at 20 mm	100 cm
ZDES-36-180	36 mm	16 Fr/6.0 mm	180 mm	201 mm at 20 mm	100 cm
ZDES-46-80	46 mm	16 Fr/6.0 mm	80 mm	93 mm at 28 mm	100 cm
ZDES-46-120	46 mm	16 Fr/6.0 mm	120 mm	137 mm at 28 mm	100 cm
ZDES-46-185	46 mm	16 Fr/6.0 mm	185 mm	208 mm at 28 mm	100 cm

NOTE: The Zenith Dissection Endovascular Stent is typically used together with Zenith® TX2® Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component, whose purpose is to seal entry tear(s). For information regarding the use and deployment of these components, please refer to the Zenith TX2 Dissection Endovascular Graft Instructions for Use.

1.2 Zenith Dissection Endovascular Stent with the Z-Trak Plus Introduction System

The Zenith Dissection Endovascular Stent is shipped preloaded in a 16 French introducer sheath that is 100 cm long and has a radiopaque band at the tip. The introduction system has a sequential deployment method with built-in features to provide continuous control of the stent throughout the deployment procedure. The introduction system enables precise positioning of the device prior to and throughout deployment to its final release. To keep the gray positioner from rotating in relation to the sheath, an anti-torque device is mounted on the introduction system at the level of the Captor® Hemostatic Valve. For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. All introduction systems feature Flexor introducer sheaths, which resist kinking and are hydrophilically coated. The system uses a single trigger-wire release mechanism to secure the endovascular stent onto the introduction system until released by the physician. The introduction system is compatible with a .035 inch wire guide (Fig. 2).

2 INTENDED USE

The Zenith Dissection Endovascular Stent with the the Z-Trak Plus Introduction System is indicated for the endovascular treatment of patients with symptomatic aortic dissection distal to the left subclavian artery (can be used to provide support to delaminated segments of the aorta typical in dissection) having vascular morphology suitable for endovascular repair (Fig. 3), including:

- Adequate iliac/femoral access compatible with the required introduction systems; when used together with the Zenith TX2 Dissection Endovascular Graft with Pro-Form.

Introduction System Sizes

TX2 Straight Component 22-34 mm	OD = 7.7 mm
TX2 Straight Component 36-42 mm	OD = 8.5 mm
TX2 Tapered Component 32, 34 mm	OD = 7.7 mm
TX2 Tapered Component 36-42 mm	OD = 8.5 mm
Dissection Stent	OD = 6.0 mm

- Radius of curvature greater than 35 mm along the entire length of aorta intended to be treated by the Zenith Dissection Endovascular Stent,
- Localized angulation less than 45 degrees, and
- Diameter at intended implant site (measured outer wall to outer wall) of no greater than 38 mm (true lumen) and no less than 20 mm (total aortic diameter).

The Zenith Dissection Endovascular Stent is intended to be used as a distal component to a Zenith TX2 Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component, to provide support to delaminated segments of aortic dissection distal to the left subclavian artery.

Standard techniques for placement of the descending thoracic aortic access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Dissection Endovascular Stent with the Z-Trak Plus Introduction System.

3 CONTRAINDICATIONS

The Zenith Dissection Endovascular Stent with the Z-Trak Introduction System is contraindicated in:

- Patients with known sensitivities or allergies to nitinol, polyester or gold.
- Patients with a systematic infection who may be at increased risk of endovascular stent infection.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.
- The Zenith Dissection Endovascular Stent should be used only by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in Section 10.1, Physician Training.
- The long-term performance of this device has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their device. Patients with specific clinical findings (e.g., changes in the structure or position of the device) should receive enhanced follow-up. Specific follow-up guidelines are described in Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP.
- After device placement, patients should be regularly monitored for flow in the false lumen of the dissection, or changes in the structure or position of the device. At a minimum, annual imaging is required, including: 1) thoracic and abdominal device radiographs to examine device integrity (e.g., stent fracture); and 2) contrast and non-contrast CT to examine flow in the false lumen of

1.1 The Zenith Dissection Endovascular Stent

The Zenith Dissection Endovascular Stent is a one-piece cylindrical device constructed from self-expanding nitinol stent segments sewn together with polyester suture (Fig. 1). The stents are uncovered to avoid blockage of dominant spinal cord intercostal arteries and allow for deployment of the stents across branch vessel origins as needed to treat the dissection. The Zenith Dissection Endovascular Stent is available in multiple lengths (4, 6 or 8 stent segments) and in two diameters (36 mm and 46 mm). Given the nature of the uncovered stent design, overall device length will vary in vivo with vessel diameter (Table 1).

the dissection, patency, tortuosity, device position and progressive disease. If renal complications or other factors preclude the use of image contrast media, use of other imaging modalities (e.g., TEE, duplex ultrasound, IVUS) should be considered.

- The Zenith Dissection Endovascular Stent is not recommended in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation studies as described in Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP.
- Following initial endovascular repair, additional endovascular interventions or conversion to standard open surgical repair should be considered for patients experiencing continued flow in the false lumen of the dissection which may lead to rupture. Further intervention should be considered for patients exhibiting compromise of organ vessel flow.
- Always have a qualified surgical team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- Interventions such as defibrillation, cardioversion, or CPR, although not specifically evaluated in studies, may have the potential to disrupt device position, and should be followed by imaging to confirm continued device function.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

4.2 Patient Selection, Treatment and Follow-Up

- Access vessel diameter (measured inner wall to inner wall) and morphology (tortuosity, occlusive disease, and/or calcification) should be compatible with vascular access techniques and a 20 or 22 French introducer system as the device is typically used together with Zenith TX2 Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude femoral introduction of the endovascular graft and/or stent and may increase the risk of embolization.
- Key anatomic elements that may affect successful treatment of dissection include severe angulation (radius of curvature < 35 mm and localized angulation > 45 degrees).
- The Zenith Dissection Endovascular Stent is not recommended for patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging.
- The Zenith Dissection Endovascular Stent is not recommended for patients whose weight or size would compromise or prevent the necessary imaging requirements.
- Patients with connective tissue disorders have not been evaluated.

4.3 Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital- and physician-preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be used.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during introduction system insertion.
- Do not bend or kink the introduction system. Doing so may cause damage to the introduction system and the Zenith Dissection Endovascular Stent.
- Always use fluoroscopy for guidance, delivery, and observation of the Zenith Dissection Endovascular Stent within the vasculature.
- The use of the Zenith Dissection Endovascular Stent with the Z-Trak Plus Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure.
- As the sheath is withdrawn, do not advance the introduction system. Doing so can cause the stent to become inverted.
- Overlapping of bare stent(s) or overlap with the Zenith TX2 Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component is left to the discretion of the implanting physician. Factors affecting whether or not to overlap, such as locations of reentries or expanded false lumen, should be judged by individual patient anatomy. When overlapping the bare stent within the stent graft component, no more than one-half of a partially overlapped bare stent body should be un-overlapped, so as to prevent flaring of the bare stent.
- If the distal end of the stent will be deployed in a funnel-shaped or angulated section of the aorta, or if the distal end of the stent appears conical in shape upon deployment, it is recommended to extend the treated segment distally with an additional stent, or choose a longer device so it ends in a straight part of the aorta. Similarly, if the distal end of the stent will be deployed at the level of the diaphragm, or in a segment adjacent to the origin of the Celiac Trunk, Superior Mesenteric Artery and /or Renal Arteries, it is also recommended to extend the treated segment distally with an additional stent or choose a longer device.
- Use of a molding balloon inside a section of aorta treated with the Zenith Dissection Endovascular Stent is not recommended.

- To avoid twisting the device, never rotate the introduction system during the procedure. Allow the device to conform naturally to the curves and tortuosity of the aorta.
- Avoid twisting or rotating the gray positioner against the introducer sheath assembly. Doing so may cause the loaded stent to become entangled and to deploy in a twisted state, or not to release from the introduction system.
- If necessary to increase hemostasis, an uninflated molding balloon or introduction system dilator may be placed within the Captor Valve.
- As the sheath is withdrawn, anatomy and stent position may change. Constantly monitor device position and perform angiography to check position as necessary.
- Do not continue advancing the wire guide or any portion of the introduction system if resistance is felt. Stop and assess the cause of resistance; vessel, catheter, or stent damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis, or calcified or tortuous vessels.
- Use caution during manipulation of catheters, wires and sheaths within a dissection. Significant disturbances may cause further injury or dislodge fragments of thrombus, which can cause distal or cerebral embolization.
- Avoid damaging the device or disturbing device positioning after placement in the event reinstrumentation (secondary intervention) of the stent is necessary. Exercise caution when manipulating a wire guide through an in situ Zenith Dissection Endovascular Stent; the wire guide may become entangled with the stent.
- Care should be taken not to advance the sheath while the device is still within it. Advancing the sheath at this stage may cause malpositioning of the device.
- Do not attempt to re-sheath the stent after partial or complete deployment.

4.4 MRI Information

Nonclinical testing has demonstrated that the Zenith TX2 Dissection Endovascular Graft overlapped with the nitinol Zenith Dissection Endovascular Stent (ZDES) is **MR Conditional** according to ASTM F2503. A patient with these devices may be safely scanned after placement under the following conditions.

- Static magnetic fields of 1.5 or 3.0 Tesla
- Maximum spatial magnetic gradient of 720 Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode) for 15 minutes of scanning or less (i.e., per scanning sequence)

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MR-Related Heating

1.5 Tesla Temperature Rise

In nonclinical testing, the Zenith TX2 Dissection Endovascular Graft overlapped with the Zenith Dissection Endovascular Stent (ZDES) produced a maximum temperature rise of 1.2°C (scaled to an SAR of 2.0 W/kg) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Magnetom, Software Numaris/4).

3.0 Tesla Temperature Rise

In nonclinical testing, the Zenith TX2 Dissection Endovascular Graft overlapped with the Zenith Dissection Endovascular Stent (ZDES) produced a maximum temperature rise of 1.4°C (scaled to an SAR of 2.0 W/kg) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (General Electric Excite, HDx, Software 14X.M5).

Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 80 mm of the position of the Zenith TX2 Dissection Endovascular Graft overlapped with the Zenith Dissection Endovascular Stent (ZDES) as found during nonclinical testing using T1-weighted, spin echo and gradient echo pulse sequence in a MR 3.0 Tesla system (General Electric Excite, HDx, Software 14X.M5). Therefore, it may be necessary to optimize MR imaging parameters for the presence of these devices.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail:	MedicAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382
Phone:	888-633-4298 (toll free) 209-668-3333 from outside the US
Fax:	209-669-2450
Web:	www.medicalert.org

5 POTENTIAL ADVERSE EVENTS

Adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Aorto-bronchial fistula
- Aorto-esophageal fistula
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma, or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, tamponade, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoprostheses: improper component placement; incomplete component deployment; component migration and/or separation; suture break; occlusion; infection; stent fracture; stent corrosion; dilatation; erosion and puncture
- Entry-Flow
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g.,

ischemia, erosion, fistula, urinary incontinence, hematuria, infection)

- Hepatic failure
- Impotence
- Infection of the dissection, device or access site, including abscess formation, transient fever and pain
- Local or systemic neurologic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, spinal cord shock, paraparesis, paralysis)
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula, lymphocele)
- Occlusion of device or native vessel
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)

6 POTENTIAL RISKS AND BENEFITS

The device is an implantable endoprosthesis intended to provide support to delaminated segments of aortic dissection and, with coverage of primary entry tear, reduce the risk of rupture and treat disturbances of end organ perfusion.

The hazards associated can be categorized as device-related (e.g., lack of sterility, toxicity, biodegradation of the device), deployment-related (e.g., failure to traverse the iliac arteries, misdeployment), performance-related (e.g., migration, stent fracture), and disease-related (e.g., persistent false lumen flow and/or end organ malperfusion). The consequent risks to the patient depend on the incidence and effects of each hazard, which have been explored in a number of experimental and clinical insertions. These risks of endovascular repair must be weighed against the risks associated with the current alternative forms of descending thoracic aortic dissection management.

Implantation of the Zenith Dissection Endovascular Stent is likely a less invasive procedure than open surgical repair. Therefore, potential clinical benefits to patients treated with the Zenith Dissection Endovascular Stent may include a suitable thoracic aorta dissection repair with less risk and fewer complications than those treated with open surgical repair. Zenith Dissection Endovascular Stent patients may benefit from a reduced risk of serious treatment-related complications, shorter anesthesia times, shorter procedure times, reduced procedural blood loss and reduced need for blood products. Additional potential risks and benefits are unknown.

7 PATIENT SELECTION AND TREATMENT

(See Section 4.2, Patient Selection, Treatment and Follow-Up)

7.1 Individualization of Treatment

Cook recommends that the Zenith Dissection Endovascular Stent component lengths described in **Table 1** be selected to correspond to the length of dissection to be treated. All device lengths necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. The risks and benefits previously described in **Section 6, POTENTIAL RISKS AND BENEFITS** should be carefully considered for each patient before use of the Zenith Dissection Endovascular Stent. Additional considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- The risk of rupture or organ failure compared to the risk of treatment with the Zenith Dissection Endovascular Stent
- Ability to tolerate general, regional, or local anesthesia
- Ilio-femoral access vessel size and morphology (thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of 20 or 22 French as the device is typically used together with Zenith TX2 Dissection Endovascular Graft with Pro-Form - Straight Component or Tapered Component, including:
 - Adequate iliac/femoral access compatible with the required introduction systems,
 - Radius of curvature greater than 35 mm along the entire length of aorta intended to be treated by the Zenith Dissection Endovascular Stent,
 - Localized angulation less than 45 degrees,
 - Diameter at intended implant site (measured outer wall to outer wall) of no greater than 38 mm and no less than 20 mm,
 - The ends of the device should not terminate in a curvature < 35 mm and localized angulation > 45 degrees.

If the distal end of the stent will be deployed in a funnel-shaped or angulated section of the aorta , or if the distal end of the stent appears conical in shape upon deployment, it is recommended to extend the treated segment distally with an additional stent, or choose a longer device so it ends in a straight part of the aorta. Similarly, if the distal end of the stent will be deployed at the level of the diaphragm, or in a segment adjacent to the origin of the Celiac Trunk, Superior Mesenteric Artery and/or Renal Arteries, it is also recommended to extend the treated segment distally with an additional stent or choose a longer device. The final treatment decision is at the discretion of the physician and patient.

8 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the potential risks and benefits when discussing this endovascular device and procedure, including:

- Risks and differences between endovascular repair and open surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- Potential advantages of medical therapy
- The possibility that subsequent interventional or open surgical repair of the dissection may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment to and compliance with postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- The long-term performance of this device has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their device. Patients with specific clinical findings (e.g., changes in the structure or position of the device) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment

of thoracic aortic dissections. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements is required and should be considered a life-long commitment to the patient's health and well-being.

- The patient should be told that successful dissection repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise every patient that it is important to seek prompt medical attention if he/she experiences signs of rupture. Rupture may be asymptomatic, but usually presents as pain, numbness, weakness in the legs, any back or chest pain, persistent cough, dizziness, fainting, rapid heartbeat, or sudden weakness.

The physician should complete the Patient Card and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

9 HOW SUPPLIED

The Zenith Dissection Endovascular Stent is supplied sterilized using ethylene oxide, preloaded in a 16 French introducer sheath, and packaged in a sterile, peel-open pouch. This device is intended for one-time use. Do not use the product if package is opened or damaged. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

- Do not use after the expiration date printed on the label.
- Store in a dark, cool, dry place.

The Zenith Dissection Endovascular Stent is available in the following lengths and diameters.

NOTE: The length of the device varies with vessel diameter (Table 1).

NOTE: The loaded stent is compressed lengthwise. Movement applied to the positioner as the stent is unsheathed may allow the deployed stent to lengthen.

10 CLINICAL USE INFORMATION

10.1 Physician Training

Familiarity with the Zenith® family of products is advantageous.

CAUTION: Always have a qualified surgical team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Dissection Endovascular Stent should only be used by physicians and teams trained in vascular interventional techniques (endovascular and surgical) and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Dissection Endovascular Stent are outlined below:

Patient selection:

- Knowledge of the natural history of thoracic aortic dissection and co-morbidities associated with thoracic aortic dissection repair.
- Knowledge of radiographic image interpretation, patient selection, device selection, planning and sizing.

A multidisciplinary team that has combined procedural experience with:

- Femoral and brachial cutdown, arteriotomy, and repair or conduit technique
- Percutaneous access and closure techniques
- Nonselective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

10.2 Inspection Prior to Use

- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

10.3 Materials Recommended (Not included with the Zenith Dissection Endovascular Stent)

- Zenith TX2 Dissection Endovascular Graft with Pro-Form
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Power injector
- Syringe
- Heparinized saline solution
- Sterile gauze pads
- .035 inch (0.89 mm) extra stiff wire guide, 260/300 cm; for example:
 - Cook Lunderquist™ Double Curved Extra Stiff Wire Guides (LESDC)
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide; for example:
 - Cook .035 inch Wire Guides
 - Cook .035 inch Benton Wire Guide
 - Cook Nimble® Wire Guides
- Introducer sets; for example:
 - Cook Check-Flo® Introducer Sets
- Sizing catheter; for example:
 - Cook Auros® Centimeter Sizing Catheters
- Angiographic radiopaque marker catheters; for example:
 - Cook Beacon® Tip Angiographic Catheters
 - Cook Beacon® Tip Royal Flush Catheters
- Entry needles; for example:
 - Cook Single Wall Entry Needles

10.4 Device Length Selection Guidelines

The choice of length should be determined from the pre-implant examinations, taking into consideration the fact that device length varies with vessel diameter, the components may be overlapped, and degree of tortuosity.

10.5 Device Diameter Selection Guidelines

The Zenith Dissection Endovascular Stent is intended for use as a distal component in combination with the Zenith TX2 Dissection Endovascular Graft with Pro-Form. Therefore, the diameter of the Zenith Dissection Endovascular Stent should be selected with consideration to the distal diameter of the Zenith TX2 Dissection Endovascular Graft with Pro-Form. The 36 mm diameter stent is intended for use in conjunction with distal graft diameters ranging from 22 to 34 mm. The 46 mm diameter component is intended for use in conjunction with distal graft diameter ranging from 36 to 42 mm. Additional considerations may affect the choice of stent diameter.

11 DIRECTIONS FOR USE

The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the introduction system(s)
2. Angulation of aorta, aneurysm (if present) and iliac arteries
3. Quality of the implantation sites (e.g., angulation)
4. Diameters of the implantation sites

Device Preparation

For instruction on the Zenith TX2 Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component deployment, consult the Zenith TX2 Dissection Endovascular Graft with Pro-Form Instructions for Use (IFU). **Leave the sheath and wire guide from the Zenith TX2 Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component in place, as the introducer for the Zenith Dissection Endovascular Stent is introduced through it coaxially.** The ID of the Zenith TX2 Dissection Endovascular Graft Straight Component or Tapered Component introducer sheath will accommodate introduction of the Zenith Dissection Endovascular Stent Introducer Sheath.

11.1 The Zenith Dissection Endovascular Stent System

11.1.1 Preparation/Flush

1. Remove yellow-hubbed shipping stylet and remove cannula protector tube (at the handle). Remove Peel-Away® sheath from back of valve assembly (Fig. 4).
2. Elevate distal tip of system and flush through the hemostatic valve until fluid emerges from the sideport near the tip of the introduction sheath (Fig. 5). Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube. **NOTE:** Graft flushing solution of heparinized saline is often used.
3. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the distal dilator tip (Fig. 6).
4. Soak sterile gauze pads in saline solution and use to wipe the Flexor Introducer Sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

11.1.2 Placement

1. Perform angiography at the appropriate level. If using radiopaque markers, adjust position as necessary and repeat angiography.
2. Ensure system has been flushed with heparinized saline (appropriate flush solution), and all air has been removed.
3. Give systemic heparin. Flush all catheters and wet all wire guides with a heparin solution. This should be repeated following each exchange.
4. Remove pigtail flush catheter and leave the sheath and wire guide in place.
5. Introduce the Zenith Dissection Endovascular Stent introduction system over the wire guide through the Zenith TX2 Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component and advance until the desired device position is reached. Make sure that the valve assembly of the Zenith Dissection Endovascular Stent sheath docks with the previously placed sheath.

During coaxial introduction of the Zenith Dissection Endovascular Stent Introducer Sheath inside of the Zenith TX2 Dissection Endovascular Graft with Pro-Form Straight Component or Tapered Component sheath, take care not to inadvertently advance the outer sheath. Dislodgement of the in situ Graft Component can occur.

CAUTION: To avoid twisting the device, never rotate the introduction system during the procedure. Allow the device to conform naturally to the curves and tortuosity of the aorta.

NOTE: The dilator tip will soften at body temperature.

6. Verify wire guide position in the aortic arch. Ensure correct stent position.
7. Ensure that the Captor Hemostatic Valve on the introduction sheath is turned to the open position (Fig. 7a).
8. Just before withdrawing the sheath to deploy the stent, unlock the black cap on the anti-torque device by rotating it counter-clockwise. The anti-torque device is now released from the gray dilator and attached only to the Captor Hemostatic Valve (Fig. 7b).
9. Stabilize the gray positioner (introduction system shaft) and begin withdrawing the sheath until the stent is fully expanded and the valve assembly docks with the control handle.
10. Loosen the safety lock from the green trigger-wire release mechanism. Withdraw the trigger-wire in a continuous movement until the proximal end of the device opens. Do not rotate the green trigger-wire knob (Fig. 8). The distal end is still attached. Continue to withdraw the trigger-wire in a continuous movement until the distal end opens. Withdraw the trigger-wire completely.

As the distal end of the stent is still attached to the introduction system do not move the gray positioner until both ends of the stent are fully released.

NOTE: Check to make sure that the trigger-wire is removed prior to withdrawal of the introduction system.

NOTE: When using the sheath as a conduit through which other devices will be inserted, stabilize the sheath and remove the inner introduction system entirely, leaving sheath and wire guide in position. Remove the anti-torque device from the Captor Hemostatic Valve by twisting and removing it (Fig. 8). Close the Captor Hemostatic Valve by turning it clockwise until it stops. Before any secondary procedure, open the Captor Hemostatic Valve by turning it counter-clockwise until it stops.

11. Remove the introduction system, leaving the wire guide in the graft.

11.2 Final Angiogram

1. Position angiographic catheter just above the level of the endovascular graft. Perform angiography to verify correct positioning. Verify patency of vessels inside the stented area.
2. Repair vessels and close in standard surgical fashion.

12 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP

12.1 General

The long-term performance of this device has not yet been established.

All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and performance of their device. Patients with specific clinical findings (e.g., changes in the structure or position of the device) should receive additional follow-up. Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of thoracic aortic dissections.

Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in Table 2. This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.

Annual imaging follow-up should include thoracic and abdominal device radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, thoracic and abdominal device radiographs and non-contrast CT may be used.

- The combination of contrast and non-contrast CT imaging provides information on flow within the false lumen of the dissection, patency, tortuosity, progressive disease, stent-to-vessel apposition, and other morphological changes.

- The thoracic and abdominal device radiographs provide information on device integrity (e.g., stent separation, stent fracture).

Table 2 lists the minimum requirements for imaging follow-up for patients with the Zenith Dissection Endovascular Stent. Patients requiring enhanced follow-up should have interim evaluations.

Table 2 Recommended Imaging Schedule for Endograft Patients

	Angiogram	CT (Contrast and non-contrast)	Thoracic and Abdominal Device Radiographs
Pre-procedure		X ¹	
Procedural	X		
Pre-discharge (within 7 days)		X ²	X
1 month		X ²	X
6 month		X ²	X
12 month (annually thereafter)		X ²	X

¹ Imaging should be performed within 6 months before the procedure.

² If flow persists within the false lumen of the dissection, prompt intervention and additional follow-up post-intervention is recommended.

See Section 12.5, Additional Surveillance and Treatment.

12.2 Contrast and Non-Contrast CT Recommendations

- Film sets should include all sequential images at lowest possible slice thickness (≤ 3 mm). Do NOT perform large slice thickness (> 3 mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14" x 17" sheets if film is used.
- Both non-contrast and contrast runs are required, with matching or corresponding table positions.

- Pre-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.

Non-contrast and contrast enhanced baseline and follow-up imaging are important for optimal patient surveillance. It is important to follow acceptable imaging protocols during the CT exam. **Table 3** lists examples of acceptable imaging protocols.

Table 3 Acceptable Imaging Protocols

	Non-Contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral capable of > 40 seconds	Spiral capable of > 40 seconds
Injection volume	n/a	150 cc
Injection rate	n/a	> 2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E. or equivalent
Coverage - start	Neck	Proximal to the left subclavian artery
Coverage - finish	Diaphragm	Profunda femoris origin
Collimation	< 3 mm	< 3 mm
Reconstruction	2.5 mm throughout - soft algorithm	2.5 mm throughout - soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

12.3 Thoracic and Abdominal Device Radiographs

The following views are required:

- Two films: supine-frontal (AP) and cross-table lateral.
- Record the table-to-film distance and use the same distance at each subsequent examination.
- Middle photo cell should be used to fully penetrate mediastinum and allow visualization of device.

Ensure entire device is captured on each single image format lengthwise. If devices extend into abdomen, additional imaging is required.

If there is any concern about the device integrity (e.g., stent breaks), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length) using 2-4X magnification visual aid.

12.4 MRI Information

Nonclinical testing has demonstrated that the Zenith TX2 Dissection Endovascular Graft overlapped with the nitinol Zenith Dissection Endovascular Stent (ZDES) is **MR Conditional** according to ASTM F2503. A patient with these devices may be safely scanned after placement under the following conditions.

- Static magnetic fields of 1.5 or 3.0 Tesla
- Maximum spatial magnetic gradient of 720 Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode) for 15 minutes of scanning or less (i.e., per scanning sequence)

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Temperature Rise

In nonclinical testing, the Zenith TX2 Dissection Endovascular Graft overlapped with the Zenith Dissection Endovascular Stent (ZDES) produced a maximum temperature rise of 1.2°C (scaled to an SAR of 2.0 W/kg) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Magnetom, Software Numaris/4).

3.0 Tesla Temperature Rise

In nonclinical testing, the Zenith TX2 Dissection Endovascular Graft overlapped with the Zenith Dissection Endovascular Stent (ZDES) produced a maximum temperature rise of 1.4°C (scaled to an SAR of 2.0 W/kg) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (General Electric Excite, HDx, Software 14X.M5).

Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 80 mm of the position of the Zenith TX2 Dissection Endovascular Graft overlapped with the Zenith Dissection Endovascular Stent (ZDES) as found during nonclinical testing using T1-weighted, spin echo and gradient echo slice sequence in a MR 3.0 Tesla system (General Electric Excite, HDx, Software 14X.M5). Therefore, it may be necessary to optimize MR imaging parameters for the presence of these devices.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail:	MedicAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382
Phone:	888-633-4298 (toll free) 209-668-3333 from outside the US
Fax:	209-669-2450
Web:	www.medicalert.org

12.5 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Growth or extension of the false lumen
- Flow in false lumen of the dissection
- Obstruction/compromise of flow to end organs
- Migration
- Inadequate stent-to-vessel apposition

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy, and the patient's personal choices. Patients should be counseled that subsequent reinterventions, including catheter-based and open surgical conversion, are possible following endograft placement.

13 REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019