<b>EN</b> 26	<b>Zenith® Fenestrated</b> AAA Endovascular Graft with the H&L-B One-Shot™ Introduction System Instructions for Use
<b>CS</b> 35	Endovaskulární graft <b>Zenith® Fenestrated</b> AAA se zaváděcím systémem H&L-B One-Shot™ Návod k použití
<b>DA</b> 44	<b>Zenith® Fenestrated</b> AAA endovaskulær protese med H&L-B One-Shot™ indføringssystem Brugsanvisning
<b>DE</b> 53	<b>Zenith® Fenestrated</b> endovaskuläre AAA-Prothese mit dem Einführsystem H&L-B One-Shot™ Gebrauchsweisung
<b>EL</b> 64	Ενδαγγειακό μόσχευμα AAA <b>Zenith® Fenestrated</b> με το σύστημα εισαγωγής H&L-B One-Shot™ Οδηγίες χρήσης
<b>ES</b> 75	Endoprótesis vascular para AAA <b>Zenith® Fenestrated</b> con el sistema de introducción H&L-B One-Shot™ Instrucciones de uso
<b>FI</b> 86	Endovaskulaarinen <b>Zenith® Fenestrated</b> AAA -grafti ja H&L-B One-Shot™ -sisäänvientijärjestelmä Käyttöohjeet
<b>FR</b> 95	Endoprothèse vasculaire <b>Zenith® Fenestrated</b> AAA à système d'introduction H&L-B One-Shot™ Mode d'emploi
<b>IT</b> 106	Endoprotesi addominale <b>Zenith® Fenestrated</b> AAA con sistema di introduzione H&L-B One-Shot™ Istruzioni per l'uso
<b>NL</b> 117	<b>Zenith® Fenestrated</b> AAA endovasculaire prothese met H&L-B One-Shot™ introductiesysteem Gebruiksaanwijzing
<b>PL</b> 128	Stent-graft wewnętrznożylniowy <b>Zenith® Fenestrated</b> AAA z systemem wprowadzającym H&L-B One-Shot™ Instrukcja użycia
<b>PT</b> 139	Prótese endovascular <b>Zenith® Fenestrated</b> AAA com o sistema de introdução H&L-B One-Shot™ Instruções de utilização
<b>SK</b> 150	Endovaskulárny štep <b>Zenith® Fenestrated</b> AAA so zavádzacím systémom H&L-B One-Shot™ Návod na použitie
<b>SV</b> 159	<b>Zenith® Fenestrated</b> AAA endovaskulärt graft med H&L-B One-Shot™ införingssystem Bruksanvisning

## TABLE OF CONTENTS

Illustrations .....	11-25
<b>1 DEVICE DESCRIPTION .....</b>	<b>26</b>
1.1 Proximal Body Graft .....	26
1.2 Proximal Body Graft Delivery System .....	26
1.3 Distal Bifurcated Body Graft .....	26
1.4 Distal Bifurcated Body Graft Delivery System .....	26
1.5 Iliac Leg .....	26
1.6 Ancillary Components .....	26
<b>2 INTENDED USE .....</b>	<b>26</b>
<b>3 CONTRAINDICATIONS .....</b>	<b>26</b>
<b>4 WARNINGS AND PRECAUTIONS .....</b>	<b>26</b>
4.1 General Use Information .....	26
4.2 Patient Selection, Treatment and Follow-Up .....	26
4.3 Implant Procedure .....	27
4.4 Molding Balloon Use .....	27
4.5 MRI Safety and Compatibility .....	27
<b>5 POTENTIAL ADVERSE EVENTS .....</b>	<b>27</b>
<b>6 PATIENT SELECTION AND TREATMENT .....</b>	<b>27</b>
6.1 Individualisation of Treatment .....	27
<b>7 PATIENT COUNSELING INFORMATION .....</b>	<b>28</b>
<b>8 HOW SUPPLIED .....</b>	<b>28</b>
8.1 Proximal Body Graft .....	28
8.2 Distal Bifurcated Body Graft .....	29
<b>9 CLINICAL USE INFORMATION .....</b>	<b>29</b>
9.1 Physician Training .....	29
9.2 Inspection Prior to Use .....	29
9.3 Materials Required .....	29
9.4 Materials Recommended .....	29
9.5 Device Diameter Sizing Guidelines .....	29
9.5.1 Proximal Body Graft Diameter Sizing .....	30
9.5.2 Distal Ipsilateral Graft Diameter Sizing .....	30
<b>10 INSTRUCTIONS FOR USE .....</b>	<b>30</b>
General Use Information .....	30
Pre-Implant Determinants .....	30
Patient Preparation .....	30
10.1 Fenestrated System .....	30
10.1.1 Bifurcated Main Body Preparation/Flush .....	30
10.1.2 Distal Bifurcated Body Graft Preparation/Flush .....	30
10.1.3 Vascular Access and Angiography .....	30
10.1.4 Proximal Body Placement .....	30
10.1.5 Docking of Top Cap .....	31
10.1.6 Fenestration Stent Placement and Deployment .....	31
10.1.7 Distal Bifurcated Body Placement .....	31
10.1.8 Contralateral Iliac Wire Guide Placement .....	31
10.1.9 Distal Bifurcated Body Deployment .....	31
10.1.10 Iliac Leg (Contralateral) Placement .....	31
10.1.11 Distal Bifurcated Body Deployment .....	32
10.1.12 Molding Balloon Insertion .....	32
Final Angiogram .....	32
<b>11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP .....</b>	<b>32</b>
11.1 General .....	32
11.2 Contrast and Non-Contrast CT Recommendations .....	32
11.3 Abdominal Radiographs .....	33
11.4 Ultrasound .....	33
11.5 MRI Safety and Compatibility .....	33
11.6 Additional Surveillance and Treatment .....	33
<b>12 TRIGGER-WIRE RELEASE TROUBLESHOOTING .....</b>	<b>33</b>
12.1 Alternate Proximal Body Deployment .....	33
12.2 Docking of Top Cap .....	33
<b>13 SUPRARENAL STENT DEPLOYMENT TROUBLESHOOTING .....</b>	<b>33</b>
13.1 Proximal Body Placement with distal attachment .....	33
13.2 Proximal Body Placement without distal attachment .....	34

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
Figure 69

## ENGLISH

### ZENITH® FENESTRATED AAA ENDOVASCULAR GRAFT WITH THE H&L-B ONE-SHOT™ INTRODUCTION SYSTEM

Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

#### 1 DEVICE DESCRIPTION

The Zenith Fenestrated AAA Endovascular Graft is a modular system consisting of three components, a proximal body graft, a distal bifurcated body graft and one iliac leg. (Figure 1) The graft modules are constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel Cook-Z® stents with braided polyester and monofilament polypropylene suture.

The modules are fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Additionally, the Cook-Z stents provide the necessary attachment and seal of the graft to the vessel wall.

Ancillary devices such as main body extensions, iliac leg extensions, converters, and iliac plugs may also be required.

Each individual device has its own separate delivery system.

Each component comes in a range of lengths and diameters which allows the physician to tailor the device to individual patient anatomies and select the best proximal and distal fixation sites.

##### 1.1 Proximal Body Graft

The bare suprarenal stent at the proximal end of the proximal body graft contains barbs that are placed at 3 mm increments for additional fixation of the device. This graft contains up to three precisely located holes (fenestration(s)), and cut-outs from the proximal margin (scallop(s)) of the graft material. (Figure 2)

The purpose of these scallops and fenestrations is to allow the proximal margin of the device to sit higher than standard AAA devices and allow uninterrupted blood flow to branch vessels of the aorta such as the renal and superior mesenteric arteries.

To facilitate fluoroscopic visualisation of the stent graft, gold radiopaque markers are positioned as follows: one on the lateral aspect of the most distal stent and four in a circumferential orientation within 1 mm of the most superior aspect of the graft material.

To facilitate orientation of the graft during deployment gold radiopaque markers are positioned on the anterior and posterior surfaces of the graft, below the level of the fenestrations. Additional gold markers are positioned anteriorly in the shape of a tick or check marker.

##### 1.2 Proximal Body Graft Delivery System

The Zenith Fenestrated AAA Endovascular Proximal Body Graft is shipped preloaded onto the H&L-B One-Shot Introduction System. (Figure 3) It has sequential deployment method with built-in features to provide continuous control of the graft throughout the deployment procedure.

The graft is reduced in diameter by an independent wire tied to diameter reducing ties, which allows the graft to be manipulated within the aorta to allow accurate positioning of the graft, which enables the fenestration(s) to line up with the desired arteries.

The bare suprarenal stent is constrained within a top cap and held by a trigger-wire. The distal end of the graft is also attached to the delivery system and held by an independent wire. The H&L-B One-Shot Introduction System enables precise positioning and allows readjustment of the final graft position before deployment of the bare barbed suprarenal stent.

The delivery system uses a 6.7 mm I.D. (20 French) or 7.3 mm I.D. (22 French) H&L-B One-Shot Introduction System. All Systems are compatible with a .035 inch (0.89 mm) wire guide.

For added hemostasis, the Captor® Hemostatic Valve can be opened or closed for the introduction and/or removal of ancillary devices into and out of the sheath. The proximal body graft delivery system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

##### 1.3 Distal Bifurcated Body Graft

The Zenith Fenestrated AAA Endovascular Distal Bifurcated Body Graft has one long ipsilateral iliac limb and one short contralateral limb. To facilitate fluoroscopic visualisation of the stent graft, there is a radiopaque marker at the graft bifurcation and a contralateral limb radiopaque marker at the distal end of the contralateral limb. (Figure 4)

##### 1.4 Distal Bifurcated Body Graft Delivery System

The Zenith Fenestrated AAA Endovascular Distal Bifurcated Body Graft is shipped preloaded onto the H&L-B One-Shot Introduction System. (Figure 5)

The bifurcated graft has one long ipsilateral limb and one short contralateral limb. There is a radiopaque marker at the graft bifurcation and a radiopaque marker at the distal end of the contralateral limb. It has a sequential deployment method with built-in features to provide continuous control of the graft throughout the deployment procedure. Both the proximal and distal segments of the graft are attached to the delivery system and held by independent wires.

The H&L-B One Shot Introduction System enables precise positioning and allows readjustment of the graft position before deployment of the graft.

The delivery system uses a 6.7 mm I.D. (20 French) H&L-B One-Shot Introduction System. All systems are compatible with a .035 inch (0.89 mm) wire guide.

For added hemostasis, the Captor Hemostatic Valve can be opened or closed for the introduction and/or removal of ancillary devices into and out of the sheath. The distal bifurcated body graft delivery system features a Flexor introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

##### 1.5 Iliac Leg

The Zenith Fenestrated AAA Endovascular Graft includes a Zenith iliac leg graft component. Examples include the Zenith Spiral-Z AAA Iliac Leg (ZSLE). Refer to the respective Instructions for Use.

All iliac leg delivery systems use a 4.7 mm I.D. (14 Fr) or 5.3 mm I.D. (16 French) Introduction System. All systems are compatible with a .035 inch (0.89 mm) wire guide system and are designed for ease of use with minimal preparation.

#### 1.6 Ancillary Components

Additional ancillary components (main body extensions, iliac leg extensions, converters, and iliac plugs) are available. (Figure 6)

Zenith ancillary components are constructed from the same polyester fabric, self-expanding stainless steel Cook-Z stents, and polypropylene suture used in constructing the principal graft modules. Refer to Instructions for Use enclosed in device packaging.

#### 2 INTENDED USE

The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is indicated for the treatment of infrarenal abdominal aortic aneurysms (AAAs) in high-risk patients who are not suitable for conventional open surgical repair or who are not suitable for repair with a standard Zenith Endovascular Graft, when the aneurysmal disease extends up to the level of the renal arteries.

The patient must have aneurysmal morphology suitable for endovascular repair, including:

- Adequate femoral/iliac access compatible with the required delivery systems;
- A non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm with:
  - A length of at least 4 mm, i.e. at least 4 mm circumferential wall contact around fenestrations and/or limits of scallop;
  - A diameter measured outer wall to outer wall of no greater than 31 mm and no less than 19 mm;
  - An angle less than 45 degrees relative to the long axis of the aneurysm; and
  - An angle less than 45 degrees relative to the axis of the suprarenal aorta.
- An ipsilateral iliac artery distal fixation site greater than 10 mm in length and 9–21 mm in diameter (measured outer wall to outer wall); and
- A contralateral iliac artery distal fixation site greater than 10 mm in length and 7–21 mm in diameter (measured outer wall to outer wall).

#### 3 CONTRAINDICATIONS

The Zenith Fenestrated AAA Endovascular Graft is contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, nitinol (nickel), polyester, solder (tin, silver), polypropylene, urethane or gold
- Patients with a systemic infection may be at increased risk of endovascular graft infection.
- Patients who cannot tolerate contrast agents necessary for intra-operative and post-operative follow-up imaging.
- Patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.

#### 4 WARNINGS AND PRECAUTIONS

##### 4.1 General Use Information

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Fenestrated grafts are made to a customised design to a specification requested by the responsible Physician, and are tailored to a specific patient's anatomy.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System must only be used by physicians and teams trained in vascular interventional techniques and in the use of this device.
- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Preprocedure imaging reconstruction thickness > 3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenosis from CT
- **The long-term performance and safety of endovascular grafts has not yet been established.** As a result, life-long, regular follow-up must be undertaken in all patients to assess the ongoing performance of the Zenith Fenestrated AAA Endovascular Graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in Section 11.
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft.
  - At a minimum, annual imaging is required, including:
    - 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease.
  - If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is not recommended in patients unable to undergo, or who will not be compliant with the necessary preoperative and post-operative imaging and implantation studies as described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up.**
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak.
- An increase in aneurysm size and/or persistent endoleak may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the haemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the grey positioner. After the grey positioner has been removed, if blood loss is excessive, consider placing an uninflated moulding balloon or an introduction system dilator within the valve, restricting flow.

##### 4.2 Patient Selection, Treatment and Follow-Up

- Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease, and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 5.5 mm O.D. (14 French) to 8.6 mm O.D. (22 French) vascular introducer sheath.
- Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolisation.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (> 45 degrees for infrarenal neck to axis of AAA or > 45 degrees for suprarenal neck relative

- to the immediate infrarenal neck); short proximal aortic neck (< 4 mm); greater than 10% increase in diameter over 15 mm of proximal aortic neck length; and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface.
- Irregular calcification and/or plaque may compromise the fixation and sealing of the implantation sites. Necks exhibiting these key anatomic elements may be more conducive to graft migration.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.

#### 4.3 Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with 4x4 gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Maintain wire guide position during delivery system insertion.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Fenestrated AAA Endovascular Graft.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System requires administration of intravascular contrast.

Patients with pre-existing renal insufficiency may have an increased risk of renal failure post-operatively. Care should be taken to limit the amount of contrast media used during the procedure.

- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
- Inaccurate placement and/or incomplete sealing of the Zenith Fenestrated AAA Endovascular Graft within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.
- Inadequate fixation of the Zenith Fenestrated AAA Endovascular Graft may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation bars.

- Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system.
- Stop and assess the cause of resistance. Vessel or catheter damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Unless medically indicated, do not deploy the Zenith Fenestrated AAA Endovascular Graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (**exception is the inferior mesenteric artery only if the superior mesenteric artery and the anastomosis of Riolan are patent and functioning correctly**) with the endoprosthesis. Vessel occlusion may occur. Clinical studies have identified that with use of the Zenith Fenestrated AAA Endovascular Graft there is a risk of renal impairment/loss (which could lead to a need for dialysis).
- Clinical studies have shown that this risk is reduced when the renal arteries are stented as part of the procedure to implant a Zenith Fenestrated AAA Endovascular Graft.
- Stenting of the renal arteries is therefore **strongly** recommended with use of the Zenith Fenestrated AAA Endovascular Graft to reduce the risk of renal impairment/loss.
- Take care during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolisation.
- Before deployment of the suprarenal stent, verify that the position of the access wire guide extends just distal to the aortic arch.
- Verify that the predetermined contralateral iliac leg is selected for insertion on the contralateral side of the patient before implantation.
- Care should be taken not to damage the graft or disturb graft positioning after graft placement in the event reinstrumentation of the graft is necessary.

#### 4.4 Molding Balloon Use

- Prior to molding in the vicinity of any fenestration stent(s) confirm that the aortic section of the stent has been flared.
- Confirm complete deflation of balloon prior to repositioning.
- Do not inflate balloon in iliac vessel outside of graft.

#### 4.5 MRI Safety and Compatibility

Non-clinical testing conducted on the standard Zenith AAA Endovascular Graft has demonstrated that the graft is MR Conditional.

The Zenith Fenestrated AAA Graft contains the same metal as the standard Zenith AAA device (stainless steel), but also contains a small amount of nitinol. It is not expected that this small amount of Nitinol would alter the MRI Conditional rating that was observed for the standard Zenith AAA device.

The standard Zenith AAA Endovascular Graft can be scanned safely under the following conditions:

##### 1.5 Tesla Systems:

- Static magnetic field strength of 1.5 Tesla
- Spatial gradient field of 450 Gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the standard Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.4°C at a maximum whole body averaged specific absorption rate (SAR) of 2.8 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla static magnetic field strength Magnetom, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole body averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

##### 3.0 Tesla Systems:

- Static magnetic field strength of 3.0 Tesla
- Spatial gradient field of 720 Gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the standard Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9°C at a maximum whole body averaged specific absorption rate (SAR) of 3.0 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, G3.0-052B Software, MR Scanner. The maximum whole body averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

MR image quality of the standard Zenith AAA Endovascular Graft may be compromised if the area of interest is in the exact same area or within approximately 20 cm of the device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo, in a 3.0 Tesla static magnetic field strength, Excite, GE Electric Healthcare, with G3.0-052B software, MR system with body radiofrequency coil. Therefore it may be necessary to optimize MR imaging parameters for this metallic implant.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

**NOTE:** For the Zenith Fenestrated AAA Graft, the clinical benefit of an MRI scan should be balanced with the potential risk of the procedure.

#### 5 POTENTIAL ADVERSE EVENTS

Use of the Zenith Fenestrated AAA Endovascular Graft poses the following potential risks in addition to those associated with use of a standard Zenith Endovascular Graft:

- Organ impairment/loss due to side-branch vessel occlusion (in particular, renal and/or gastrointestinal impairment/loss).

Clinical studies have identified that with use of the Zenith Fenestrated AAA Endovascular Graft there is a risk of renal impairment/loss (which could lead to a need for dialysis). Clinical studies have shown that this risk is reduced when the renal arteries are stented as part of the procedure to implant a Zenith Fenestrated AAA Endovascular Graft.

Stenting of the renal arteries is therefore **strongly** recommended with use of the Zenith Fenestrated AAA Endovascular Graft to reduce the risk of renal impairment/loss.

Other adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g. buttock, lower limb)
- Death
- Edema
- Embolisation (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Fever and localised inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)

#### 6 PATIENT SELECTION AND TREATMENT

(See Warnings and Precautions)

##### 6.1 Individualisation of Treatment

Each patient must be evaluated on an individual basis with careful consideration given to both the potential benefits and specific risks associated with the procedure.

Considerations regarding the use of the Zenith Fenestrated AAA Endovascular Graft (see Warnings and Precautions) include:

- Risk of aneurysm rupture
- Morbidity and mortality associated with conventional surgical repair
- Comorbidities
- Size of aneurysm
- History of renal failure
- Life expectancy
- Anaesthetic risk
- Age of patient
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 14 French (5.5 mm O.D.) to 22 French (8.6 mm O.D.) vascular introducer sheath.
- Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
  - with a length of at least 4 mm,
  - with a diameter measured outer wall to outer wall of no greater than 31 mm and no less than 19 mm,
  - with an angle less than 45 degrees relative to the long axis of the aneurysm, and
  - with an angle less than 45 degrees relative to the axis of the suprarenal aorta.

- Ipsilateral iliac artery distal fixation site greater than 10 mm in length and 9-21 mm in diameter (measured outer wall to outer wall).
- Contralateral iliac artery distal fixation site greater than 10 mm in length and 7-21 mm in diameter (measured outer wall to outer wall).
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.

#### 7 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure including:

- Risks and differences between endovascular repair and surgical repair.
- Potential advantages of traditional open surgical repair.
- Potential advantages of endovascular repair.

The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance and safety of endovascular grafts has not yet been established.** As a result, life-long, regular follow-up must be undertaken in all patients to assess the ongoing performance of the Zenith Fenestrated AAA Endovascular Graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up.**

#### 8 HOW SUPPLIED

The Zenith Fenestrated AAA Endovascular Graft is supplied sterile (100% ethylene oxide) and pre-loaded in peel-open packages.

This device is intended for single patient use only. Do not reuse, reprocess or resterilize any part of this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilisation barrier has been damaged or broken. If damage has occurred, do not use the product and return to your Cook representative or your nearest Cook office.

Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient. Do not use after the expiration date printed on the label. Store in a cool dry place.

The proximal body and distal bifurcated body grafts are loaded into either a 20 French or 22 French Flexor introducer sheath. The sheath's surface is treated with a hydrophilic coating that, when activated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a 4x4 gauze pad soaked in saline solution.

The Zenith Fenestrated AAA Endovascular Graft and delivery systems are available in the lengths and diameters detailed on **pages 30 to 31.**

Reorder Number	Proximal Body Diameter	8.1 Proximal Body Graft		Introduction Sheath Length
		Proximal Body Length	Introduction Sheath French Size (I.D. / O.D.)	
ZFEN-P-1-24-76	24 mm	76 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-26-76	26 mm	76 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-28-76	28 mm	76 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-30-76	30 mm	76 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-32-76	32 mm	76 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-24-91	24 mm	91 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-26-91	26 mm	91 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-28-91	28 mm	91 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-30-91	30 mm	91 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-32-91	32 mm	91 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-24-106	24 mm	106 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-26-106	26 mm	106 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-28-106	28 mm	106 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-30-106	30 mm	106 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-32-106	32 mm	106 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-24-121	24 mm	121 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-26-121	26 mm	121 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-28-121	28 mm	121 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-30-121	30 mm	121 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-32-121	32 mm	121 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-34-84	34 mm	84 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-36-84	36 mm	84 mm	22 Fr 7.3 mm / 8.5 mm	50 cm
ZFEN-P-1-34-99	34 mm	99 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-36-99	36 mm	99 mm	22 Fr 7.3 mm / 8.5 mm	50 cm
ZFEN-P-1-34-114	34 mm	114 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-36-114	36 mm	114 mm	22 Fr 7.3 mm / 8.5 mm	50 cm
ZFEN-P-1-34-129	34 mm	129 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-1-36-129	36 mm	129 mm	22 Fr 7.3 mm / 8.5 mm	50 cm
ZFEN-P-2-24-94	24 mm	94 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-26-94	26 mm	94 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-28-94	28 mm	94 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-30-94	30 mm	94 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-32-94	32 mm	94 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-24-109	24 mm	109 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-26-109	26 mm	109 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-28-109	28 mm	109 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-30-109	30 mm	109 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-32-109	32 mm	109 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-24-124	24 mm	124 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-26-124	26 mm	124 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-28-124	28 mm	124 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-30-124	30 mm	124 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-32-124	32 mm	124 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-34-107	34 mm	107 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-36-107	36 mm	107 mm	22 Fr 7.3 mm / 8.5 mm	50 cm
ZFEN-P-2-34-122	34 mm	122 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-36-122	36 mm	122 mm	22 Fr 7.3 mm / 8.5 mm	50 cm
ZFEN-P-2-34-137	34 mm	137 mm	20 Fr 6.7 mm / 7.7 mm	50 cm
ZFEN-P-2-36-137	36 mm	137 mm	22 Fr 7.3 mm / 8.5 mm	50 cm

Reorder Number	Ipsilateral Leg Diameter	8.2 Distal Bifurcated Body Graft		Introduction Sheath French Size	Introduction Sheath (I.D. / O.D.)	Introduction Sheath Length
		Distal Body Length	Distal Leg Length			
ZFEN-D-12-28-76	12 mm	76 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-28-76	16 mm	76 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-28-76	20 mm	76 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-28-94	12 mm	94 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-28-94	16 mm	94 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-28-94	20 mm	94 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-28-109	12 mm	109 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-28-109	16 mm	109 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-28-109	20 mm	109 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-28-124	12 mm	124 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-28-124	16 mm	124 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-28-124	20 mm	124 mm	28 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-45-76	12 mm	76 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-45-76	16 mm	76 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-45-76	20 mm	76 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-45-76	24 mm	76 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-45-94	12 mm	94 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-45-94	16 mm	94 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-45-94	20 mm	94 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-45-94	24 mm	94 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-45-109	12 mm	109 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-45-109	16 mm	109 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-45-109	20 mm	109 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-45-109	24 mm	109 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-45-124	12 mm	124 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-45-124	16 mm	124 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-45-124	20 mm	124 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-45-124	24 mm	124 mm	45 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-62-76	12 mm	76 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-62-76	16 mm	76 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-62-76	20 mm	76 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-62-76	24 mm	76 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-62-94	12 mm	94 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-62-94	16 mm	94 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-62-94	20 mm	94 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-62-94	24 mm	94 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-62-109	12 mm	109 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-62-109	16 mm	109 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-62-109	20 mm	109 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-62-109	24 mm	109 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-12-62-124	12 mm	124 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-16-62-124	16 mm	124 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-20-62-124	20 mm	124 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
ZFEN-D-24-62-124	24 mm	124 mm	62 mm	20 Fr	6.7 mm / 7.7 mm	40 cm

## 9 CLINICAL USE INFORMATION

### 9.1 Physician Training

**CAUTION: Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.**

**CAUTION: The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System MUST only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System are outlined below:**

#### Patient Selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair.
- Knowledge of radiographic image interpretation, device selection and sizing.

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolisation
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimise radiation exposure
- Expertise in necessary patient follow-up modalities

### 9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilisation barrier has been damaged or broken. If damage has occurred, do not use the product and return to the nearest Cook office or Representative. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

### 9.3 Materials Required

(Not included in 3-piece modular system)

- Zenith AAA Endovascular Graft Ancillary Kit
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinised saline solution

### 9.4 Materials Recommended

(Not included in 3-piece modular system)

The following products are recommended:

- 0.035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
  - Cook Amplatz Ultra Stiff Wire Guides (AUS)
  - Cook Lunderquist Extra Stiff Wire Guides (LES)
- 0.035 inch (0.89 mm) standard wire guide; for example:
  - Cook 0.035 inch wire guides
  - Cook Nimble™ Wire Guides
- Molding Balloons
- Introducer sets; for example:
  - Cook Check-Flo® Introducer Sets
  - Cook Extra Large Check-Flo Introducer Sets
  - Cook Flexor Balkin Up & Over® Contralateral Introducers
  - Cook Flexor Guiding Sheaths
- Sizing catheter; for example:
  - Cook Aurores® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
  - Cook Beacon® Tip Angiographic Catheters
  - Cook Royal Flush® Catheters with Beacon Tip
- Entry needles; for example:
  - Cook single wall entry needles

### 9.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

### 9.5.1 Proximal Body Graft Diameter Sizing

Intended Aortic Vessel Diameter	Main Body Diameter	Introduction Sheath	
		French Size	(I.D. / O.D.)
19 mm	24 mm	20 Fr	6.7 mm / 7.7 mm
20 mm	24 mm	20 Fr	6.7 mm / 7.7 mm
21 mm	24 mm	20 Fr	6.7 mm / 7.7 mm
21 mm	26 mm	20 Fr	6.7 mm / 7.7 mm
22 mm	26 mm	20 Fr	6.7 mm / 7.7 mm
23 mm	28 mm	20 Fr	6.7 mm / 7.7 mm
24 mm	28 mm	20 Fr	6.7 mm / 7.7 mm
24 mm	30 mm	20 Fr	6.7 mm / 7.7 mm
25 mm	30 mm	20 Fr	6.7 mm / 7.7 mm
26 mm	30 mm	20 Fr	6.7 mm / 7.7 mm
26 mm	32 mm	20 Fr	6.7 mm / 7.7 mm
27 mm	32 mm	20 Fr	6.7 mm / 7.7 mm
28 mm	32 mm	20 Fr	6.7 mm / 7.7 mm
28 mm	34 mm	20 Fr	6.7 mm / 7.7 mm
29 mm	34 mm	20 Fr	6.7 mm / 7.7 mm
29 mm	36 mm	22 Fr	7.3 mm / 8.5 mm
30 mm	36 mm	22 Fr	7.3 mm / 8.5 mm
31 mm	36 mm	22 Fr	7.3 mm / 8.5 mm

### 9.5.2 Distal Ipsilateral Graft Diameter Sizing

Intended Iliac Vessel Diameter	Ipsilateral Fenestrated Leg	Introduction Sheath	
		French Size	(I.D. / O.D.)
9 mm	12 mm	20 Fr	6.7 mm / 7.7 mm
10 mm	12 mm	20 Fr	6.7 mm / 7.7 mm
11 mm	12 mm	20 Fr	6.7 mm / 7.7 mm
12 mm	16 mm	20 Fr	6.7 mm / 7.7 mm
13 mm	16 mm	20 Fr	6.7 mm / 7.7 mm
14 mm	16 mm	20 Fr	6.7 mm / 7.7 mm
15 mm	16 mm	20 Fr	6.7 mm / 7.7 mm
15 mm	20 mm	20 Fr	6.7 mm / 7.7 mm
16 mm	20 mm	20 Fr	6.7 mm / 7.7 mm
17 mm	20 mm	20 Fr	6.7 mm / 7.7 mm
18 mm	20 mm	20 Fr	6.7 mm / 7.7 mm
18 mm	24 mm	20 Fr	6.7 mm / 7.7 mm
19 mm	24 mm	20 Fr	6.7 mm / 7.7 mm
20 mm	24 mm	20 Fr	6.7 mm / 7.7 mm
21 mm	24 mm	20 Fr	6.7 mm / 7.7 mm

## 10 INSTRUCTIONS FOR USE

### General Use Information

Prior to use of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System, review this Suggested Instructions for Use booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary.

These instructions are intended to help guide the physician and do not take the place of physician judgment.

### Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the main body system, (i.e., define respective contralateral and ipsilateral iliac arteries).
2. Angulation of aortic neck, aneurysm and iliacs.
3. Quality of the aortic neck.
4. Diameters of infrarenal aortic neck and distal iliac vessels.
5. Distance from renal arteries to the aortic bifurcation.
6. Distance from the renal arteries to the hypogastric (internal iliac) arteries/attachment site(s).
7. Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
8. Consider the degree of vascular calcification.

**NOTE:** Each respective vessel diameter and length (aorta, ipsilateral iliac and contralateral iliac) provides the necessary criteria for choosing the appropriate endovascular graft.

### Patient Preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualisation from the aortic arch to the femoral bifurcations.
3. Expose both common femoral arteries using standard surgical technique.
4. Establish adequate proximal and distal vascular control of both femoral vessels.

### 10.1 Fenestrated System

#### 10.1.1 Bifurcated Main Body Preparation/Flush

1. Remove black-hubbed shipping stilet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away® sheath from back of the hemostatic valve. (Figure 7) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove in the tip. (Figure 8) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

**NOTE:** Graft flushing solution of Heparinised saline is always used.

2. Attach syringe with normal Heparinised saline to the hub on the inner cannula. Flush until fluid exits the distal tip. (Figure 9)

**NOTE:** When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak 4x4 gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator tip liberally.

#### 10.1.2 Distal Bifurcated Body Graft Preparation/Flush

1. Remove black-hubbed shipping stilet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. (Figure 10) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove in the tip. (Figure 11) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

**NOTE:** Graft flushing solution of Heparinised saline is always used.

2. Attach syringe with normal Heparinised saline to the hub on the inner cannula. Flush until fluid exits the distal tip. (Figure 9)

**NOTE:** When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak 4x4 gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator tip liberally.

### 10.1.3 Vascular Access and Angiography

1. Puncture the selected common femoral arteries using standard technique with an 18 or 19UT gage arterial needle. Upon vessel entry, insert:
  - Wire guides – standard .035 inch (0.89 mm) diameter, 145 cm long, J tip or Bentson Wire Guide
  - Appropriate size sheaths (e.g., 6 French (2.0 mm I.D.) or 8 French (2.7 mm I.D.))
  - Flush catheter (often radiopaque sizing catheters – e.g., Centimeter Sizing Catheter or straight flush catheter)
2. Perform angiography to identify level(s) of renals, aortic bifurcation and iliac bifurcations.

**NOTE:** If fluoroscope angulation is used with an angulated neck it may be necessary to perform angiograms using various projections.

**NOTE:** A previous planning exercise will have determined which side will be used to introduce the proximal and distal bodies.

### 10.1.4 Proximal Body Placement

**CAUTION:** Verify that the predetermined access site is chosen for the introduction and placement of the proximal body.

1. Ensure the delivery system has been flushed with Heparinised saline and that all air is removed from the system.
2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

**NOTE:** Monitor the patient's coagulation status throughout the procedure.

3. On ipsilateral side, replace J wire with stiff wire guide (AUS or LES) .035 inch (0.89 mm), 260 cm long and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.

**NOTE:** A straight angiographic catheter may be inserted up the contralateral side to aid in placement of graft.

4. Before insertion, position proximal body delivery system on patient's abdomen under fluoroscopy to assist with orientation and positioning. Rotate to a position where the anterior markers are situated in the most anterior (12:00 o'clock) position. The sidearm of the hemostatic valve may serve as an external reference to the fenestration(s) and/or scallop(s), anterior and posterior markers and body side markers.

**CAUTION: Maintain wire guide position during delivery system insertions.**

**CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).**

5. Advance the delivery system until the radiopaque markers indicating the fenestration(s) and/or scallop(s) are at the level of the appropriate arteries. Check that the distal end of the graft is in a satisfactory position above the aortic bifurcation and that the anterior and posterior markers indicate that the graft is in satisfactory orientation. (Figure 15a)
6. The 'tick' marker can be used to assist with orientation of the graft during deployment and prevent 180-degree mis-orientation of the graft.
  - A marker formation of '✓' indicates an anterior position of the anterior markers. (Figure 15a, Illustration A)
  - A marker formation of '∨' indicates a posterior position of the anterior markers. (Figure 15a, Illustration B)

**NOTE:** Angiography should be performed as needed throughout deployment, to confirm correct placement of the graft.

7. Verify position of the wire guide in the thoracic aorta. Ensure that fenestration(s) and/or scallops are at the level of the appropriate arteries and the anterior markers are in the most anterior (12:00 o'clock) position.

**NOTE:** The vertical anterior markers, and the horizontal posterior markers should form a cross, on the fluoroscopic image, when correctly oriented. (Figure 15b)

**NOTE:** The fenestration/scallop markers should be in close apposition to the appropriate side branch vessels.

**NOTE:** Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (Figure 16)

Clear identification of fenestration position(s) may not be possible until the graft has been fully unsheathed.

8. Stabilise the grey positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first two (2) covered stents by withdrawing the sheath while monitoring device location.
9. Perform angiography, and adjust graft placement as necessary. Continue to withdraw the sheath making positional adjustments as necessary.

**NOTE:** Techniques to ensure that the fenestration(s) and/or scallop(s) will accurately align with their respective vessels will vary, and will depend upon vessel anatomy, graft design, and physician preferences.

10. Proceed with deployment until the graft has been fully unsheathed. (Figure 17)

11. When a satisfactory graft position has been achieved, withdraw the angiographic catheter and wire guide, then exchange to selective wire guide/selective catheter to below the level of the proximal body. Cannulate the partially deployed proximal main body.

**NOTE:** If a small fenestration is being utilised, care should be taken to properly align the fenestration with the respective vessel.

12. Utilising contralateral access sheath and wire guide, cannulate and advance a guiding catheter into each small fenestration and its respective vessel. (Figure 18)

**NOTE:** Non-compliant angioplasty balloons may be used as an alternative to guiding catheters.

**NOTE:** Cannulation of the scallop and its respective vessel may also be achieved using similar techniques.

**NOTE:** It is not recommended to use balloons or guiding catheters to guide final placement of large fenestrations as stent struts across fenestration may interfere.

**CAUTION:** Before release of the diameter reducing ties, verify that the position of the ipsilateral access wire extends just distal to the aortic arch.

**CAUTION:** During proximal trigger-wire removal, top cap advancement, and subsequent suprarenal stent deployment, verify that the position of the main body wire guide extends just distal to the aortic arch and that support of the system is maximized.

13. Verify proper position of proximal body. Remove the safety lock from the gold trigger-wire release mechanism. Withdraw and remove the trigger-wire to release diameter reducing ties by sliding the gold trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 19)

**NOTE:** At this point, the proximal main body graft should be fully expanded with the proximal bare stent contained within the top cap.

14. Remove the safety lock from the black trigger-wire release mechanism. Under fluoroscopy, withdraw and remove the trigger-wire to unlock the suprarenal stent from the top cap by sliding the black trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 20)

**NOTE:** If resistance is felt or system bowing is noticed, the trigger-wire is under tension. Excessive force may cause the graft position to be altered. If excessive resistance or delivery system movement is noted, stop and assess the situation.

If unable to remove the black trigger-wire release mechanism from the top cap, perform the following steps under fluoroscopy:

- a. Remove tension on the trigger-wire by loosening the pin vise and slightly pulling the inner cannula to move the top cap down over the suprarenal stent. Avoid compressing the Zenith Fenestrated proximal body.
- b. Retighten the pin vise.
- c. Remove the black trigger-wire release mechanism.
- d. Continue with (15) in Section 10.1.4 Proximal Body Placement.

**NOTE:** If still unable to remove the black trigger-wire release mechanism from the top cap, see Section 12 Trigger-Wire Release Troubleshooting.

15. Loosen the pin vise. (Figure 21) Control the position of the graft by stabilizing the grey positioner of the introducer.

**CAUTION:** Before deployment of suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement and, if required, re-position the access wire guide into the aortic arch to accommodate.

16. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed. (Figures 22a and 22b) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise to avoid contact with the deployed suprarenal stent.

**WARNING:** The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

17. Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire to detach the distal end of the endovascular graft from the delivery system by sliding the white trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Figure 23)

**NOTE:** Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

#### 10.1.5 Docking of Top Cap

1. Loosen the pin vise. (Figure 24)
2. Secure sheath and inner cannula to avoid any movement of these components.
3. Advance the grey positioner over the inner cannula until it docks with the top cap. (Figures 25a, 25b and 25c)

**NOTE:** If resistance occurs, slightly rotate grey positioner and continue to gently advance.

4. Retighten the pin vise and withdraw the entire top cap and grey positioner through the graft and through the sheath by pulling on the inner cannula. (Figure 26) Leave the sheath and wire guide in place.
5. Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops.

#### 10.1.6 Fenestration Stent Placement and Deployment

##### General Use Information

In the event that small fenestrations are being utilized, stents may be placed to secure positive alignment.

Standard techniques for placement of arterial stents should be employed during use of stents.

1. Return to the guide catheter and wire guide which cannulate the small fenestration and respective vessel.
2. Introduce appropriately sized balloon expandable stent and advance to the ostium of the fenestration/vessel. Advance into the vessel, leaving approximately 5 mm of stent in the aorta. (Figure 27)

**NOTE:** Fluoroscopic views tangential to the fenestration will optimise visualisation of the stent position relative to the stent graft.

3. Expand stent.
4. Remove the balloon and replace with an oversized angioplasty balloon. Advance the balloon until the proximal tip is positioned at the ostium.
5. Inflate the balloon to flare the intra-aortic segment of the stent. (Figure 28)

**CAUTION:** This technique requires high quality imaging. Mobile image intensifiers provide less than adequate imaging quality.

6. Remove the angioplasty balloon.

**NOTE:** In the event that there is more than one fenestration, repeat the preceding steps for each additional small fenestration.

7. Withdraw renal access sheaths, catheters and wire guides in the contralateral side to a level just above the aortic bifurcation.

#### 10.1.7 Distal Bifurcated Body Placement

1. Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.
2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.
3. Before insertion, position distal bifurcated body delivery system on patient's abdomen under fluoroscopy to determine the orientation of the contralateral limb. The side arm of the hemostatic valve may serve as an external reference to the contralateral limb radiopaque marker.

**NOTE:** Distal bifurcated body delivery system will not pass through the sheath used to deliver the proximal body.

**NOTE:** The proximal body delivery sheath must be removed prior to insertion of the distal bifurcated body delivery system.

4. Insert Distal Bifurcated Body delivery system over the wire, into the femoral artery with attention to sidearm reference.

**CAUTION:** Maintain wire guide position during delivery system insertion.

**CAUTION:** To avoid any twist in the endovascular graft, during any rotation of the delivery, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

5. Advance delivery system until the contralateral limb is positioned above and anterior to the origin of the contralateral iliac. (Figure 29) If the contralateral limb radiopaque marker is not properly aligned, rotate the entire system until it is correctly positioned half way between a lateral and an anterior position on the contralateral side.
6. Repeat angiogram to verify:
  - The degree of overlap with proximal body (no less than 2 stents)
  - The position of the contralateral limb
  - The position of the ipsilateral iliac cuff with respect to the common iliac bifurcation.

Reposition distal bifurcated body as required.

**CAUTION:** When introducing distal bifurcated body, observe proximal body closely to avoid any disruption to its position.

**NOTE:** Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (Figure 16)

7. Stabilise the grey positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first two (2) covered stents by withdrawing the sheath while monitoring device location. Proceed with deployment until contralateral limb is fully deployed. (Figure 30)

**NOTE:** Tick marker on the contralateral limb of the distal bifurcated body is used to determine anterior/posterior orientation of the contralateral limb. It is not intended to line up with the anterior tick mark on the proximal body.

#### 10.1.8 Contralateral Iliac Wire Guide Placement

1. Advance the contralateral catheter and wire guide into the common iliac artery to a level below the short contralateral limb and then manipulate the wire guide into the contralateral limb and into the Distal Bifurcated Body. (Figure 31) AP and oblique fluoroscopic views can aid in verification of device cannulation.
2. Advance the angiographic catheter into the body of the graft. Perform angiography to confirm correct position inside the Distal Bifurcated Body. Advance the catheter to where the proximal end of the Distal Bifurcated Body is attached to the introducer.

#### 10.1.9 Distal Bifurcated Body Deployment

1. Perform angiography to confirm proper position of the iliac leg with respect to the internal iliac (hypogastric) artery. Adjust position if necessary.
2. Withdraw sheath until the iliac leg is fully deployed.
3. Remove the safety lock from the black trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Figure 32) Stop withdrawing sheath.

#### 10.1.10 Iliac Leg (Contralateral) Placement

1. Position the image intensifier to show both the contralateral internal iliac artery and contralateral common iliac artery.
2. Prior to the introduction of the contralateral limb delivery system, inject contrast through the contralateral femoral sheath to locate the contralateral internal iliac artery.
3. Introduce the contralateral iliac leg delivery system into the artery. Advance slowly until the iliac leg graft overlaps at least one full iliac leg stent (i.e., proximal stent of iliac leg graft) inside the contralateral limb of the main body. (Figure 33) If there is any tendency for the

distal bifurcated body graft to move during this maneuver, hold it in position by stabilizing the grey positioner on the distal bifurcated body component (on the ipsilateral side).

**NOTE:** If difficulty is encountered advancing the iliac leg delivery system, exchange to a more supportive wire guide. In tortuous vessels the anatomy may alter significantly with the introduction of the rigid wires and sheath systems.

4. Confirm position of distal end of the iliac leg graft. Reposition the iliac leg graft if necessary to ensure both internal iliac patency and a minimum overlap of one full iliac leg stent (i.e., proximal stent of iliac leg graft, maximum overlap of 1.5 stents) within the main body endovascular graft.
5. To deploy, hold the iliac leg graft in position with the grey positioner while withdrawing the sheath. (Figures 34a and 34b)

Ensure one stent overlap is maintained.

6. Stop withdrawing the sheath as soon as the distal end of the iliac leg graft is released.
7. Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise, retract inner cannula to dock tapered dilator to grey positioner. Tighten pin vise. Maintain sheath position while withdrawing grey positioner with secured inner cannula. (Figure 35)
8. Re-check the position of the wire guide.

#### 10.1.11 Distal Bifurcated Body Deployment

1. Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the white trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Figure 36)
2. Under fluoroscopy and after verification of iliac leg graft position, loosen pin vise, retract inner cannula to dock tapered dilator to grey positioner. Tighten pin vise. Maintain sheath position while withdrawing grey positioner with secured inner cannula.
3. Re-check the position of the wire guides. Leave sheath and wire guide in place.
4. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning in a clockwise direction until hemostasis is achieved. (Figure 37)

#### 10.1.12 Molding Balloon Insertion

1. Prepare Molding balloon as follows:
  - Flush wire lumen with Heparinised saline.
  - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning it counter-clockwise.
3. Advance the Molding balloon over the wire guide and through the Hemostatic Valve of the distal bifurcated body introduction system to level of renal arteries. Maintain proper sheath positioning.

**NOTE:** Captor Hemostatic Valve may be utilized to assist with hemostasis by turning in a clockwise rotation to the "close" position.

**NOTE:** Captor Hemostatic Valve should always be in the "open" position when repositioning of molding balloon.

4. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the suprarenal stent and the infrarenal neck, starting proximally and working in the distal direction. (Figure 38)

**CAUTION: Prior to molding in the vicinity of any Fenestration stent(s) confirm that the aortic section of the stent has been flared.**

**CAUTION: Confirm complete deflation of balloon prior to repositioning.**

5. Withdraw the Molding balloon to the ipsilateral limb distal fixation site and expand.

**CAUTION: Do not inflate balloon in iliac vessel outside of graft.**

6. Deflate and remove molding balloon. Transfer the molding balloon onto the contralateral wire guide and into the contralateral iliac leg introduction system. Advance molding balloon to the contralateral limb overlap and expand.

**CAUTION: Confirm complete deflation of balloon prior to repositioning.**

7. Withdraw the molding balloon to the contralateral iliac leg/vessel distal fixation and expand. (Figure 38)

**CAUTION: Do not inflate balloon in iliac vessel outside of graft.**

8. Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
9. Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

#### Final Angiogram

1. Position angiographic catheter just above the level of the renal arteries. Perform angiography to verify that the renal arteries are patent and that there are no endoleaks. Verify patency of internal iliac arteries.
2. Confirm there are no endoleaks or kinks and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.

**NOTE:** If endoleaks or other problems are observed, refer to Section 1.6, Ancillary Components.

3. Repair vessels and close in standard surgical fashion.

### 11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

#### 11.1 General

**The long-term performance and safety of endovascular grafts has not yet been established.** As a result, life-long, regular follow-up must be undertaken in all patients to assess the ongoing performance of the Zenith Fenestrated AAA Endovascular Graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up.

Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.

Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in Table 11.1. This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.

Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs, non-contrast CT and duplex ultrasound may be used.

- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (separation between components, stent fracture and barb separation).
- Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound. Ultrasound may be a less reliable and sensitive diagnostic method compared to CT. Table 11.1 lists the minimum requirements for imaging follow-up for patients with the Zenith Fenestrated AAA Endovascular Graft. Patients requiring enhanced follow-up should have interim evaluations.

Table 11.1 Recommended Imaging Schedule for Endograft Patients

	Angiogram	CT (Contrast and non-contrast)	Abdominal Radiographs
Pre-procedure	X <sup>1</sup>	X <sup>1</sup>	
Procedural	X		
Pre-discharge (within 7 days)		X <sup>2,3,4</sup>	X
1 month		X <sup>2,3,4</sup>	X
3 month		X <sup>2,4,5</sup>	
6 month		X <sup>2,4</sup>	X
12 month (annually thereafter)		X <sup>2,4</sup>	X

<sup>1</sup>Imaging should be performed within 6 months before the procedure.

<sup>2</sup>Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan. With ultrasound, non-contrast CT is still recommended.

<sup>3</sup>Either pre-discharge or 1 month CT recommended.

<sup>4</sup>If Type I or III endoleak, prompt intervention and additional follow-up post-intervention recommended. See Section 11.6, Additional Surveillance and Treatment.

<sup>5</sup>Recommended if endoleak reported at pre-discharge or 1 month.

#### 11.2 Contrast and Non-Contrast CT Recommendations

- Film sets should include all sequential images at lowest possible slice thickness (≤3 mm). DO NOT perform large slice thickness (>3 mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14 inch x 17 inch sheets if film is used.
- Both non-contrast and contrast runs are required, with matching or corresponding table positions.
- Non-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.

Non-contrast and contrast enhanced baseline and follow-up imaging are important for optimal patient surveillance. It is important to follow acceptable imaging protocols during the CT exam. Table 11.2 lists examples of acceptable imaging protocols.

Table 11.2 Acceptable Imaging Protocols

	Non-Contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral capable of >40 seconds	Spiral capable of >40 seconds
Injection volume	n/a	150 cc
Injection rate	n/a	>2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E. or equivalent
Coverage - start	Diaphragm	1 cm superior to celiac axis
Coverage - finish	Proximal femur	Profunda femoris origin
Collimation	<3 mm	<3 mm
Reconstruction	2.5 mm throughout – soft algorithm	2.5 mm throughout – soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

### 11.3 Abdominal Radiographs

The following views are required:

- Four films: supine-frontal (AP), cross-table lateral, 30 degree LPO and 30 degree RPO views centered on umbilicus.
- Record the table-to-film distance and use the same distance at each subsequent examination.

Ensure entire device is captured on each single image format lengthwise.

If there is any concern about the device integrity (e.g., kinking, stent breaks, barb separation, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length including components) using 2-4X magnification visual aid.

### 11.4 Ultrasound

Ultrasound imaging may be performed in place of contrast CT when patient factors preclude the use of image contrast media. Ultrasound may be paired with non-contrast CT. A complete aortic duplex is to be videotaped for maximum aneurysm diameter, endoleaks, stent patency and stenosis. Included on the videotape should be the following information as outlined below:

- Transverse and longitudinal imaging should be obtained from the level of the proximal aorta demonstrating mesenteric and renal arteries to the iliac bifurcations to determine if endoleaks are present utilizing color flow and color power angiography (if accessible).
- Spectral analysis confirmation should be performed for any suspected endoleaks.
- Transverse and longitudinal imaging of the maximum aneurysm should be obtained.

### 11.5 MRI Safety and Compatibility

Non-clinical testing conducted on the standard Zenith AAA Endovascular Graft has demonstrated that the graft is MR Conditional.

The Zenith Fenestrated AAA Graft contains the same metal as the standard Zenith AAA device (stainless steel), but also contains a small amount of Nitinol. It is not expected that this small amount of nitinol would alter the MRI Conditional rating that was observed for the standard Zenith AAA device.

The standard Zenith AAA Endovascular Graft can be scanned safely under the following conditions:

#### 1.5 Tesla Systems:

- Static magnetic field strength of 1.5 Tesla
- Spatial gradient field of 450 Gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the standard Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.4°C at a maximum whole body averaged specific absorption rate (SAR) of 2.8 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla static magnetic field strength Magnetom, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole body averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

#### 3.0 Tesla Systems:

- Static magnetic field strength of 3.0 Tesla
- Spatial gradient field of 720 Gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the standard Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9°C at a maximum whole body averaged specific absorption rate (SAR) of 3.0 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, G3.0-052B Software, MR Scanner. The maximum whole body averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

MR image quality of the standard Zenith AAA Endovascular Graft may be compromised if the area of interest is in the exact same area or within approximately 20 cm of the device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo, in a 3.0 Tesla static magnetic field strength, Excite, GE Electric Healthcare, with G3.0-052B software, MR system with body radiofrequency coil. Therefore it may be necessary to optimize MR imaging parameters for this metallic implant.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

**NOTE:** For the Zenith Fenestrated AAA Graft, the clinical benefit of an MRI scan should be balanced with the potential risk of the procedure.

### 11.6 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, 5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's comorbidities, life expectancy and the patient's personal choices.

Patients should be counseled that subsequent reinterventions including catheter based and open surgical conversion are possible following endograft placement.

### 12 TRIGGER-WIRE RELEASE TROUBLESHOOTING

**CAUTION:** The following steps should be performed only if unable to remove the proximal trigger-wire as described in Section 10.1.4, Proximal Body Placement (14).

**NOTE:** Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

#### 12.1 Alternate Proximal Body Deployment

1. Cut the exposed suprarenal stent trigger-wire between the white and black release mechanisms (Figure 39) and remove the black release mechanism from the handle.
2. Remove the safety lock from the white (distal) trigger-wire release mechanism.
3. Withdraw the white trigger-wire release mechanism and then remove via its slot over the device inner cannula.

**NOTE:** This will detach the distal end of the graft from the grey positioner.

4. Using locking forceps, clamp and secure the cut end of the suprarenal stent trigger-wire. (Figure 40)

5. Loosen the pin vise and, while maintaining inner cannula and trigger-wire position, advance the grey positioner and sheath into the graft until the tip of the grey positioner is approximately 2 cm from the gold markers on the proximal edge of the proximal body (Figure 41). The advanced grey positioner provides added support to the inner cannula.

**NOTE:** Take care when advancing the grey positioner as there will be sheaths and wire guides positioned in target vessels. Ensure that the tip of the grey positioner is not advanced into the top cap.

**NOTE:** Maintain gentle tension on the suprarenal stent trigger-wire to remove any slack in the wire as the grey positioner and sheath are being advanced.

6. Lock the pin vise. Confirm that the suprarenal stent trigger-wire is secured by the forceps.
7. Stabilize the grey positioner and slowly advance the sheath until the sheath tip is 2 mm from the gold markers. (Figure 42)

**NOTE:** Take care when advancing the sheath as there will be other sheaths and wire guides positioned in target vessels. Take care not to advance the graft itself during sheath advancement.

8. Stabilize the sheath and slightly retract the grey positioner with inner cannula to move the top cap down over the suprarenal stent. (Figure 43)

**NOTE:** Avoid compressing the body of the graft.

9. Ensure the fenestrations are positioned correctly. Test the resistance on the trigger wire and make minor adjustments as necessary to reduce trigger wire resistance (Section 12.1, 8).
10. Remove the suprarenal stent trigger-wire.
11. Withdraw the sheath until the tapered tip of the grey positioner is exposed.
12. If there are multiple fenestrations, withdraw all but one of the fenestration guiding catheters. It is the physician's preference as to which guiding catheter to withdraw. It is recommended that the decision be made based on the ease of cannulation of the fenestrations and their respective vessels. (Figure 44)

**NOTE:** Leave the wire guides in place when removing the guiding catheter(s).

12. (a) If there is a single fenestration, then the proximal body needs to be cannulated with a suitable wire guide from the same side as the in-situ guide catheter. This is to allow placement of a molding balloon within the proximal body.
13. Advance a molding balloon along the now available wire guide into the proximal body and position it just superior to the distal-most end of the graft.

**NOTE:** When using the molding Balloon, ensure the appropriate sheath is used – either through the in-situ 20 Fr contralateral sheath or, if direct punctures have been used, through a 14 Fr introducer sheath. This will ensure the safe retrieval of the molding balloon.

14. Inflate the balloon to the full diameter of the graft. (Figure 45)
15. Loosen the pin vise (Figure 46). Control the position of the graft by stabilizing the grey positioner and balloon catheter.

**CAUTION:** Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

16. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed (Figures 47 and 48). Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (Figure 49) to avoid contact with the deployed suprarenal stent.

**WARNING:** The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation bars. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

17. If there are multiple fenestrations (Section 12.1, 12), deflate the balloon then withdraw it leaving the wire guide in place
17. (a) If there is a single fenestration (Section 12.1, 12a), then the molding balloon and wire guide can be safely removed.

**NOTE:** Care should be taken during removal to not disturb the guide sheath and wire guide, which remain in the target vessel

18. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (Figure 50)

**NOTE:** Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

### 12.2 Docking of Top Cap

1. Loosen the pin vise. (Figure 51)
2. Secure sheath and inner cannula to avoid any movement of these components.
3. Advance the grey positioner over the inner cannula until it docks with the top cap. (Figures 52, 53 and 54)

**NOTE:** If resistance occurs, slightly rotate grey positioner and continue to gently advance.

4. Retighten the pin vise and withdraw the entire top cap and grey positioner through the graft and through the sheath by pulling on the inner cannula. (Figure 55) Leave the sheath and wire guide in place.
5. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning it in a clockwise direction until it stops.
6. Return to Section 10.1.6, Fenestration Stent Placement and Deployment.

### 13 SUPRARENAL STENT DEPLOYMENT TROUBLESHOOTING

**CAUTION:** The following steps should be performed only if unable to deploy the suprarenal stent by advancing the top cap as described in Section 10.1.4, Proximal Body Placement (15).

**NOTE:** Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

#### 13.1 Proximal Body Placement with distal attachment

If the suprarenal stent cannot be fully deployed by advancing the top cap inner cannula, perform the following steps under fluoroscopy.

1. Tighten the pin-vise. If there are multiple fenestrations, withdraw all but one of the fenestration guiding catheters. It is the physician's preference as to which guiding catheter to withdraw. It is recommended that the decision be made based on the ease of cannulation of the fenestrations and their respective vessels. (Figure 56)

- (a) If there is a single fenestration, then the proximal body needs to be cannulated with a suitable wire guide from the same side as the in-situ guide catheter. This is to allow placement of a molding balloon within the proximal body.

**NOTE:** Leave the wire guides in place when removing the guiding catheter(s).

2. Advance a molding balloon along the now available wire guide into the proximal body and position it just superior to the distal-most end of the graft.

**NOTE:** When using the molding Balloon, ensure the appropriate sheath is used – either through the in-situ 20 Fr Contralateral sheath or, if direct punctures have been used, through a 14 Fr introducer sheath. This will ensure the safe retrieval of the molding balloon.

3. To add support to the inner cannula, inflate the balloon to the full diameter of the graft. (Figure 57)
4. Loosen the pin vise. (Figure 58)
5. Control the position of the graft by stabilising the grey positioner and balloon catheter.

**CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.**

6. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed. (Figures 59 and 60) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (Figure 61) to avoid contact with the deployed suprarenal stent.

**NOTE:** Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

**If the suprarenal stent is fully deployed:**

7. a) If there are multiple fenestrations (Section 13.1, 1), deflate the balloon then withdraw it leaving the wire guide in place. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (Figure 62)
- b) If there is a single fenestration (Section 13.1, 1a), then the molding balloon and wire guide can be safely removed.

**NOTE:** Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

**WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.**

8. Return to Section 10.1.5 Docking of Top Cap.

### 13.2 Proximal Body Placement without distal attachment

If still unable to fully deploy the suprarenal stent, perform the following steps:

1. Tighten the pin vise and deflate the balloon, while maintaining balloon position.
2. Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire to detach the distal end of the endovascular graft from the delivery system by sliding the trigger-wire release mechanism off the handle and remove via its slot over the device inner cannula. (Figure 63)
3. Loosen the pin vise (Figure 64) and, while maintaining inner cannula position, advance the grey positioner and sheath into the graft until the tip of the grey positioner is approximately 2 cm from the gold markers on the proximal edge of the proximal body (Figure 65). The advanced grey positioner provides added support to the inner cannula.

**NOTE:** Take care when advancing the grey positioner as there will be sheaths and wire guides positioned in target vessels. Ensure that the tip of the grey positioner is not advanced into the top cap.

4. Lock the pin vise.
5. Verify position of the gold markers and ensure the fenestrations are positioned correctly.
6. To add support to the inner cannula, inflate the balloon to the full diameter of the graft. (Figure 66)
7. Loosen the pin vise (Figure 67). Control the position of the graft by stabilising the grey positioner and balloon catheter.

**CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.**

8. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed (Figures 68 and 69). Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (Figure 70) to avoid contact with the deployed suprarenal stent.

**NOTE:** Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

9. a) If there are multiple fenestrations (Section 13.1, 1), deflate the balloon then withdraw it leaving the wire guide in place. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (Figure 71)
- b) If there is a single fenestration (Section 13.1, 1a), then the molding balloon and wire guide can be safely removed.

**NOTE:** Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

**WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.**

10. Return to Section 10.1.5, Docking of Top Cap.

**NOTE:** Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

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© COOK 2020
WILLIAM A. COOK AUSTRALIA PTY. LTD.
COOK MEDICAL EUROPE LTD.