


Zenith® Fenestrated AAA Endovascular Graft

Instructions for Use

TABLE OF CONTENTS

Illustrations	4-14
1 DEVICE DESCRIPTION	15
1.1 Proximal Body Graft	
1.2 Proximal Body Graft Delivery System	
1.3 Distal Bifurcated Body Graft	
1.4 Distal Bifurcated Body Graft Delivery System	
1.6 Ancillary Components and Delivery System	
1.7 Adjunctive Zenith Alignment Stent and Delivery System	15
2 INTENDED USE	
3 CONTRAINDICATIONS	
4 WARNINGS AND PRECAUTIONS	
4.1 General Use Information	
4.2 Patient Selection, Treatment and Follow-Up	
4.4 Molding Balloon Use.	
4.5 MRI Safety and Compatibility	
5 ADVERSE EVENTS	
6 SUMMARY OF CLINICAL STUDIES	
7 PATIENT SELECTION AND TREATMENT	
7.1 Individualization of Treatment	
8 PATIENT COUNSELING INFORMATION	25
9 HOW SUPPLIED	25
10 CLINICAL USE INFORMATION	25
10.1 Physician Training	25
10.2 Inspection Prior to Use	
10.3 Materials Required	
10.4 Materials Recommended	
10.6 Device Length Sizing Guidelines	
10.7 Graft Fenestration/Scallop Guidelines	
11 INSTRUCTIONS FOR USE	
11.1 General Use Information	
11.2 Pre-Implant Determinants	
11.3 Patient Preparation	
11.4 Fenestrated System	26
12 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP	28
12.1 General	28
12.2 Angiography	28
12.3 Computed Tomography (CT)	
12.4 Device X-Ray	28
12.5 Renal Duplex Ultrasound	
12.6 MRI Safety and Compatibility	
13 TRIGGER-WIRE RELEASE TROUBLESHOOTING	
13.1 Alternate Proximal Body Deployment	
14 SUPRARENAL STENT DEPLOYMENT TROUBLESHOOTING 14.1 Proximal Body Placement with distal attachment	
14.1 Proximal Body Placement with distal attachment	29 20



- Iliac Leg
 Distal Bifurcated Body
 Proximal Body






Figure 2

- Gold Radiopaque Marker
 Vertical Gold Markers
 Suprarenal Stent
 Scallop
 Gold Radiopaque Markers (4)
 Fenestration(s)
 Proximal Body




- 1. Hub 2. Pin Vise 3. Safety Locks 4. Peel-Away® Sheath 5. Stopcock 6. Connecting Tube

- 7. Sheath 8. Flushing Groove 9. Dilator Tip 10. Hemostatic Valve 11. Captor® Hemostatic Valve 12. Grey Positioner
- 13. White Trigger-Wire Release Mechanism 14. Black Trigger-Wire Release Mechanism 15. Gold Trigger-Wire Release Mechanism 16. Top Cap Inner Cannula 17. Proximal Body Graft




- Gold Radiopaque Markers
 Distal Bifurcated Body
 Gold Radiopaque Marker



- 1. Hub 2. Pin Vise 3. Safety Locks 4. Peel-Away® Sheath 5. Stopcock 6. Connecting Tube

- 7. Sheath 8. Flushing Groove 9. Dilator Tip 10. Hemostatic Valve 11. Captor® Hemostatic Valve 12. Grey Positioner
- 13. White Trigger-Wire Release Mechanism 14. Black Trigger-Wire Release Mechanism 15. Top Cap Inner Cannula 16. Distal Bifurcated Body Graft







Figure 7







Figure 8











Figure 10

Figure 11












Figure 22




Figure 23




Figure 24




Figure 25




Figure 26




Figure 27




Figure 28




Figure 29




Figure 30




Figure 31




Figure 32

1. Balloon Expansion/Graft Sealing Sites




Figure 33




Figure 34




Figure 35




Figure 36




Figure 37




Figure 38




Figure 39




Figure 41




Figure 40

Figure 42



Figure 43




Figure 44




Figure 45




Figure 46



Figure 47




Figure 48




Figure 49




Figure 50




Figure 51



Figure 52




Figure 53




Figure 54




Figure 55




Figure 56




Figure 57




Figure 58




Figure 59




Figure 60



Figure 61




Figure 62




Figure 63



Figure 64




Figure 65

ZENITH® FENESTRATED AAA ENDOVASCULAR GRAFT WITH THE H&L-B ONE-SHOT™ INTRODUCTION SYSTEM

Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

CAUTION: U.S Federal law restricts this device to sale by or on the order of a physician.

1 DEVICE DESCRIPTION

The Zenith Fenestrated AAA Endovascular Graft is a modular system consisting of three components, a proximal body graft, a distal bifurcated body graft and one iliac leg. (Figure 1) The graft modules are constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel Cook-Z® stents with braided polyester and monofilament polypropylene suture. The modules are fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Additionally, the Cook-Z stents provide the necessary attachment and seal of the graft to the vessel wall. Ancillary devices such as main body extensions, iliac leg extensions, converters, and iliac plugs may also be required. Each individual device has its own separate delivery system. Each component comes in a range of lengths and diameters which allows the physician to tailor the device to individual patient anatomies and select the best proximal and distal fixation sites.

1.1 Proximal Body Graft

The bare suprarenal stent at the proximal end of the proximal body graft contains barbs that are placed at 3 mm increments for additional fixation of the device. This graft contains up to three precisely located holes (fenestration(s)), and cut-outs from the proximal margin (scallop(s)) of the graft material. (Figure 2) The fenestrations are either small (fit entirely between struts of the seal stent) or large (cross struts of the seal stent). The purpose of these scallops and fenestrations is to allow the proximal margin of the device to sit higher than standard AAA devices and allow uninterrupted blood flow to branch vessels of the aorta such as the renal and superior mesenteric arteries. It is recommended that all vessels accommodated by a small fenestration be stented in order to secure positive alignment of the graft fenestration with the vessel origin. Stenting is optional for vessels accommodated by a scallop and not recommended for vessels accommodated by a large fenestration. To facilitate fluoroscopic visualization of the stent graft, gold radiopaque markers are positioned as follows; one on the lateral aspect of the most distal stent and four in a circumferential orientation within 1 mm of the most superior aspect of the graft material. The proximal body graft also has vertically-aligned gold markers on the anterior side (at the 12:00 o'clock position) that should form a cross (+) with the horizontallyaligned gold markers on the posterior side (180 degrees opposite the vertical markers) when the device is properly oriented.

1.2 Proximal Body Graft Delivery System

The Zenith Fenestrated AAA Endovascular Proximal Body Graft is shipped preloaded onto the H&L-B One-Shot Introduction System. (Figure 3) It has a sequential deployment method with built-in features to provide continuous control of the graft throughout the deployment procedure. The graft is reduced in diameter by an independent wire tied to diameter reducing ties, which allows the graft to be manipulated within the aorta to allow accurate positioning of the graft, which enables the fenestration(s) to line up with the desired arteries. The bare suprarenal stent is constrained within a top cap and held by a trigger-wire. The distal end of the graft is also attached to the delivery system and held by an independent wire. The H&L-B One-Shot Introduction System enables precise positioning and allows readjustment of the final graft position before deployment of the bare barbed suprarenal stent. The delivery system uses a 6.7 mm I.D. (20 French) or 7.3 mm I.D. (22 French) H&L-B One-Shot Introduction System. All Systems are compatible with a .035 inch wire guide. For added hemostasis, the Captor™ Hemostatic Valve can be loosened or tightened for the introduction and/or removal of accessory/ancillary devices into and out of the sheath. The proximal body graft delivery system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.3 Distal Bifurcated Body Graft

The Zenith Fenestrated AAA Endovascular Distal Bifurcated Body Graft has one long ipsilateral iliac limb and one short contralateral limb. To facilitate fluoroscopic visualization of the stent graft, there is a radiopaque marker at the graft bifurcation, at the distal end of the contralateral limb, and at the proximal end (contralateral side) of the graft. (Figure 4)

1.4 Distal Bifurcated Body Graft Delivery System

The Zenith Fenestrated AAA Endovascular Distal Bifurcated Body Graft is shipped preloaded onto the H&L-B One-Shot Introduction System. (Figure 5) It has a sequential deployment method with built-in features to provide continuous control of the graft throughout the deployment procedure. Both the proximal and distal segments of the graft are attached to the delivery system and held by independent wires. The H&L-B One Shot Introduction System enables precise positioning and allows readjustment of the graft position before deployment of the graft. The delivery system uses a 6.7 mm I.D. (20 French) H&L-B One-Shot Introduction System. All systems are compatible with a .035 inch wire guide. For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of accessory/ancillary devices into and out of the sheath. The distal body graft delivery system features a Flexor introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.5 Iliac Leg Graft and Delivery System

The Zenith Fenestrated AAA Endovascular Graft utilizes the same iliac leg graft as is available for the standard Zenith Flex AAA Endovascular Graft. Zenith iliac leg grafts are constructed from polyester fabric, self-expanding stainless steel and nitinol Z-stents, and polypropylene suture. Refer to the iliac leg graft instructions for Use enclosed in device packaging for more information.

1.6 Ancillary Components and Delivery System

The Zenith Fenestrated AAA Endovascular Graft utilizes the same ancillary components (main body extensions, iliac leg extensions,

converters, and iliac plugs) as are available for the standard Zenith Flex AAA Endovascular Graft.

Zenith ancillary components are constructed from the same polyester fabric, self-expanding stainless steel Z-stents, and polypropylene suture. Refer to the ancillary component Instructions for Use enclosed in device packaging for more information.

1.7 Adjunctive Zenith Alignment Stent and Delivery System

It is recommended that all vessels accommodated by a small fenestration be stented in order to secure positive alignment of the graft fenestration with the vessel origin (stenting optional for scallops and not recommended for large fenestrations). The Zenith Alignment Stent is available for this purpose. The Zenith Alignment Stent is a balloon-expandable stent that can be deployed through scallops or fenestrations in a Zenith Fenestrated AAA Endovascular Graft into branch vessels of the aorta. Refer to the Zenith Alignment Stent Instructions for Use for more information.

2 INTENDED USE

The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is indicated for the endovascular treatment of patients with abdominal aortic or aorto-iliac aneurysms having morphology suitable for endovascular repair, including:

- Adequate iliac/femoral access compatible with the required introduction systems,
- Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
 - with a length that is at least 4 mm and unsuitable for a nonfenestrated graft,
 - with a diameter measured outer wall to outer wall of no greater than 31 mm and no less than 19 mm,
 - with an angle less than 45 degrees relative to the long axis of the aneurysm, and
 - with an angle less than 45 degrees relative to the axis of the suprarenal aorta.
- Ipsilateral iliac artery distal fixation site greater than 30 mm in length and 9-21 mm in diameter (measured outer wall to outer wall).
- Contralateral Iliac artery distal fixation site greater than 30 mm in length and 7-21 mm in diameter (measured outer wall to outer wall).

3 CONTRAINDICATIONS

The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is contraindicated in the following:

- Patients with known sensitivities or allergies to stainless steel, polyester, nitinol, solder (tin, silver), polypropylene or gold
- Patients with systemic or local infection that may increase the risk of endovascular graft infection.

4 WARNINGS AND PRECAUTIONS

4.1 General Use Information

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Fenestrated grafts are made to a customized design to a specification requested by the responsible Physician, and are tailored to a specific patient's anatomy.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device, which requires precise planning/sizing as well as accurate longitudinal positioning and rotational orientation during placement.
- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Preprocedure imaging reconstruction thickness > 3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenosis from CT
- Implantation of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System requires high quality imaging. Some types of mobile image intensifiers may not provide adequate imaging quality.
- The long-term performance of fenestrated endovascular grafts, including the stents placed in fenestrations/scallops, has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, changes in the structure or position of the endovascular graft, or stenosis/occlusion of vessels accommodated by fenestrations) should receive enhanced follow-up. Specific follow-up guidelines are described in Section 12.
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth, patency of vessels accommodated by a fenestration/scallop, or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is recommended, including: 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is not recommended in patients unable to undergo, or who will not be compliant with the necessary preoperative and post-operative imaging and implantation studies as described in

Section 12, Imaging Guidelines and Post-Operative Follow-Up.

- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb/ fenestration and/or leaks may be required to undergo secondary interventions or surgical procedures.

- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

4.2 Patient Selection, Treatment and Follow-Up

- Inappropriate patient selection may result in poor performance of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System.
- Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease, and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 14 French to 22 French vascular introducer sheath. Iliac conduits may be used to ensure the safe insertion of the introduction system. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization/trauma.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (> 45 degrees for infrarenal neck to axis of AAA or > 45 degrees for suprarenal neck relative to the immediate infrarenal neck); short proximal aortic neck (<4 mm); greater than 10% increase in diameter over 15 mm of proximal aortic neck length; and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. Irregular calcification and/or plaque may compromise the fixation and sealing of the implantation sites. Necks exhibiting these key anatomic elements may be more conducive to graft migration.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is not recommended in patients who cannot tolerate contrast agents necessary for intra-operative and postoperative follow-up imaging.
- The use of this device requires administration of radiographic agents.
 Patients with pre-existing renal insufficiency may have an increased risk of post-operative renal failure.
- The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System is not recommended in patients of excessive weight and/or size that would limit, compromise, or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- Patients with recurrent aortic aneurysmal disease or with disease above the renal arteries may be prone to further aortic dilation in the renal/ visceral segment, which could compromise device integrity/fixation.
- The Zenith Fenestrated AAA Endovascular Graft has not been evaluated in the following patient populations:
- Less than 18 years of age
- Females who are pregnant or breast-feeding
- Leaking/ruptured or symptomatic aneurysms
- Patients with connective tissue disorders
- Patients with previous stent placement in vessels to be accommodated by fenestrations

4.3 Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with 4X4 gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Maintain wire guide position during delivery system insertion.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Fenestrated AAA Endovascular Graft.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure post-operatively. Care should be taken to limit the amount of contrast media used during the procedure.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
- Inaccurate placement and/or incomplete sealing of the Zenith Fenestrated AAA Endovascular Graft within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications. It is recommended that all vessels accommodated by a small fenestration be stented in order to secure positive alignment of the graft fenestration with the vessel origin.
- Inadequate fixation of the Zenith Fenestrated AAA Endovascular Graft may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance. Vessel or catheter damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Unless medically indicated, do not deploy the Zenith Fenestrated AAA Endovascular Graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (exception is the inferior mesenteric artery) with the endoprosthesis.
- Take care during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization.
- Care should be taken not to damage the graft or disturb graft positioning after graft placement in the event reinstrumentation of the graft is necessary.

4.4 Molding Balloon Use

- Prior to molding in the vicinity of any fenestration stent(s) confirm that the aortic section of the stent has been flared.
- Confirm complete deflation of balloon prior to repositioning
- Do not inflate balloon in the vessel outside of graft, as doing so could result in damage to the vessel (e.g., rupture).

4.5 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the Zenith Fenestrated AAA Endovascular Graft is MR Conditional. A patient with this endovascular graft in place for at least 6 months can be scanned safely under the following conditions:

- · Static magnetic field of 3.0 Tesla or 1.5 Tesla
- Maximum spatial magnetic gradient of 720 Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning or less (i.e., per scanning sequence)
- · Normal operating mode.

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Systems:

In non-clinical testing, the Zenith AAA Endovascular Graft (similar construction as the Zenith Fenestrated AAA Endovascular Graft) produced a temperature rise of less than or equal to 1.4 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.8 W.kg, for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

3.0 Tesla Systems:

In non-clinical testing, the Zenith AAA Endovascular Graft (similar construction as the Zenith Fenestrated AAA Endovascular Graft) produced a temperature rise of less than or equal to 1.9 °C at a maximum whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, G3.0-052B Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

Image Artifact

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in nonclinical testing using the sequence: Fast spin echo, in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B Software, MR system with body radiofrequency coil. For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest. Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners.

Mail: MedicAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382 Phone: 888-633-4298 (toll free)

209-668-3333 from outside the US

Fax: 209-669-2450 Web: www.medicalert.org

5 ADVERSE EVENTS

Potential adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- · Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g. buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction

- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Fever and localized inflammation
- Fistula (e.g., aortoenteric, arteriovenous)
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., confusion, stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- · Occlusion of device or native vessel
- Organ impairment/loss due to side-branch vessel occlusion (in particular, renal and/or gastrointestinal impairment/loss)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery stenosis or occlusion, contrast toxicity, infarct, insufficiency, failure)
- · Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula, dissection
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- · Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection, hematoma, seroma, cellulitis).

6 SUMMARY OF CLINICAL STUDIES

The Zenith Fenestrated AAA Endovascular Graft US clinical study is a nonrandomized, multi-center study that was conducted to help evaluate the safety and effectiveness of the Zenith Fenestrated AAA Endovascular Graft in the treatment of abdominal aortic aneurysms in patients with short infrarenal neck lengths (>4 mm and <15 mm). A total of 42 patients were enrolled among 7 investigational sites between January 6, 2005 and August 18, 2010. Each patient was treated with an individually tailored Fenestrated Graft. The study was initially approved for 30 patients and use of available balloon-expandable stents in combination with the Fenestrated Graft. Following completion of the initial 30 patient enrollment, the study hypothesis and requirements for approval were agreed upon. The study was then expanded to include the Zenith Alignment Stent and enrollment of 12 additional patients, thus providing 42 total patients, which was the pre-specified study sample size. The adjunctive Zenith Alignment Stent was used in 11 total patients.

The primary safety and effectiveness endpoint was based on treatment success, which was defined as technical success (i.e., successful access of the aneurysm site and deployment of the Zenith Fenestrated AAA Endovascular Graft in the intended location, with all vessels targeted by fenestrations patent at the completion of the procedure) plus freedom from the following at 6 months: Type I and III endoleak, aneurysn growth >0.5 cm, any AAA-related serious adverse event (death, rupture, conversion), and any AAA-related major complication (Q-wave MI; congestive heart failure; cardiac ischemia requiring intervention; renal failure requiring permanent dialysis; bowel obstruction, ischemia, or $fistula; stroke\ with\ permanent\ deficit;\ paralysis).\ The\ study\ results\ for$ the primary endpoint, treatment success, were explored in patients treated with the Zenith Fenestrated AAA Endovascular Graft and in matched patients treated with the standard Zenith AAA Endovascular Graft. Propensity score methods with a pre-specified one-to-one global matching algorithm were used to select patients from the Zenith AAA Endovascular Graft multicenter study. The covariates for matching were pre-specified and included the following relevant demographic, comorbid, and anatomic characteristics, all of which were comparable between the two groups: age, gender, height, weight, arrhythmia, cancer, chronic heart failure, chronic obstructive pulmonary disease cerebrovascular disease, diabetes, hypertension, previous diagnosis of systemic infection, previous myocardial infarction, peripheral vascular disease, previous surgeries at access site, thromboembolic event, maximum aneurysm diameter, minimum aneurysm diameter, and neck

Additional measures assessed in the cohort of patients treated with the Fenestrated Graft included mortality, pre-specified morbid events, change in aneurysm size, endoleak, migration, device integrity, and secondary interventions. The patients were to be seen for clinical and imaging (CT and X-ray) follow-up at pre-discharge, 1 month, 6 months, 12 months, and yearly thereafter through 5 years.

Table 6.1 reports the patient availability for follow-up. Of 42 patients enrolled in the clinical study, 95% (40) were evaluable for the primary endpoint analysis. The 42 patient cohort, combines patient data from the feasibility study (n = 30; implanted device between January 2005 – January 2006) with pivotal study data (n = 12).

Table 6.1 Follow-up availability

					lable	J.I FUIIUW	-up avana	Dility					
		Percent	t of Data Ava	ailable¹	Adequate	Imaging to	Assess the Pa	arameter ²		Events Oc	curring Before	Next Interva	I
Follow-up Visit	Patients Eligible for Follow-up ³	Clinical	X-ray	СТ	Size Increase	Endoleak	Migration	Fracture	Death	Conversion	Lost to Follow- up (LTF) or Withdrawal		Not Due for Next Visit
Pre-discharge	42 (0)	100.0% (42/42)	95.2% (40/42)	95.2% (40/42)	95.2% (40/42)	92.9% (39/42)	95.2% (40/42)	100.0% (42/42)	0	0	0	0	0
30-day	42 (0)	97.6% (41/42)	88.1% (37/42)	97.6% (41/42)	97.6% (41/42)	92.9% (39/42)	95.2% (40/42)	97.6% (41/42)	1	0	1	0	0
6-month	40 (0)	97.5% (39/40)	95.0% (38/40)	95.0% (38/40)	95.0% (38/40)	87.5% (35/40)	95.0% (38/40)	97.5% (39/40)	0	0	1	0	4
12-month	35 (3)	91.4% (32/35)	82.9% (29/35)	88.6% (31/35)	82.9% (29/35)	80.0% (28/35)	82.9% (29/35)	82.9% (29/35)	1	0	0	0	7
24-month	27 (0)	96.3% (26/27)	85.2% (23/27)	96.3% (26/27)	96.3% (26/27)	81.5% (22/27)	96.3% (26/27)	96.3% (26/27)	1	0	0	6	0
3-year	20 (0)	90.0% (18/20)	70.0% (14/20)	90.0% (18/20)	75.0% (15/20)	70.0% (14/20)	80.0% (16/20)	75.0% (15/20)	0	0	1	0	0
4-year	19 (0)	94.7% (18/19)	73.7% (14/19)	89.5% (17/19)	84.2% (16/19)	63.2% (12/19)	84.2% (16/19)	89.5% (17/19)	0	0	2	0	0
5-year	17 (0)	100.0% (17/17)	70.6% (12/17)	82.4% (14/17)	64.7% (11/17)	52.9% (9/17)	64.7% (11/17)	64.7% (11/17)	0	0	0	0	0

¹Site submitted data.

Table 6.2 summarizes the demographics and patient characteristics of patients implanted with the Zenith Fenestrated AAA Endovascular Graft.

	Table 6.2 Demographics and patient characteristics
Demographic	Result ¹
Age (years)	75.3 ± 7.4 (58 - 86), 42
Gender	
Male	78.6% (33/42)
Female	21.4% (9/42)
Ethnicity	
White	92.9% (39/42)
Hispanic or Latino	2.4% (1/42)
Black or African American	0.0% (0/42)
American Indian or Alaska Native	2.4% (1/42)
Asian	2.4% (1/42)
Native Hawaiian or other Pacific Islander	0.0% (0/42)
Other	0.0% (0/42)
Height (in)	67.5 ± 4.4 (51 - 74), 41
Weight (lbs)	190.6 ± 46.9 (110 - 342), 42
Body mass index	28.7 ± 4.8 (19.5 - 40.8), 41

Mean values +/- the standard deviation, with the range of values shown in parentheses, followed by the number of patients evaluated

²Based on core lab analysis – does not include imaging exams received by the core lab for analysis, but that have not yet been analyzed.

³Number in parenthesis indicates the number of patients without submitted data who are still eligible for follow-up.

Initial cohort of 30 patients consented only for 2-year follow-up and therefore were asked to reconsent for 3-5 year follow-up.

Table 6.3 Pre-existing comorbid medical conditions

Medical History	Percent Patients (number/total number)
Cardiovascular	
Previous myocardial infarction	23.8% (10/42)
Previous diagnosis of symptomatic congestive heart failure	9.5% (4/42)
Previous diagnosis of coronary artery disease	52.4% (22/42)
Previous diagnosis of cardiac arrhythmia	40.5% (17/42)
Vascular	
Thromboembolic event	11.9% (5/42)
Peripheral vascular disease	23.8% (10/42)
Family history of aneurysmal disease	14.3% (6/42)1
Hypertension	92.9% (39/42)
Pulmonary	
Chronic obstructive pulmonary disease	33.3% (14/42)
Renal	
Diagnosis of renal failure requiring dialysis	0.0% (0/42)
Renal insufficiency	9.5% (4/42)
GFR ≤ 60 ml/min/1.73 m ²	21.4% (9/42)
Endocrine	
Diabetes	26.2% (11/42)
Infectious disease	
Previous diagnosis of sepsis	7.1% (3/42)
Gastrointestinal	
Gastrointestinal disease	40.5% (17/42)
Hepatobiliary	
Previous diagnosis of liver disease	2.4% (1/42)
Neoplasms	
Previous diagnosis of cancer	35.7% (15/42)
Neurologic	
Previous diagnosis of cerebrovascular disease	16.7% (7/42)
Previous endarterectomy	2.4% (1/42)
Substance use	
Excessive alcohol use	0.0% (0/42)
Tobacco use: currently smokes	28.6% (12/42)
quit smoking	57.1% (24/42)
never smoked	14.3% (6/42)
Access site	
Previous surgery at the intended access site	11.9% (5/42)2
Previous radiation at the intended access site	0.0% (0/42)

 $^{^{1}}$ In 11.9% (5/42) of patients family history of an eurysmal disease was reported as unknown.

Table 6.4 lists the anatomical characteristics of the subject population for this study, as assessed by the core lab.

Table 6.4 Presenting anatomical dimensions, as assessed by core lab

Measure	Mean ± S.D. (range), N=42
Aortic diameters (mm)	
Diameter at celiac artery	28.2 ± 3.2 (21.2 - 35.9)
Diameter at SMA	28 ± 3.5 (22.3 - 39.8)
Diameter at lowest patent renal artery	25.7 ± 3.2 (19.2 - 33.2)
Diameter at midpoint of renal arteries	25.5 ± 5.1 (0.0 - 32.2)
Maximum aneurysm diameter – long axis	61.1 ± 10.9 (45.2 - 94.2)
Maximum aneurysm diameter – short axis	56.8 ± 10.3 (43.4 - 90.4)
Proximal neck length (mm)	9.7 ± 3.5 (2.4 - 19.1)
Angles (°)	
Angle between immediate suprarenal neck and immediate infrarenal neck	15.9 ± 9.6 (2 - 40)
Angle between the proximal neck and the longitudinal axis of the aneurysm	34 ± 14.2 (7 - 57)
Diameter of renal artery ostia (mm)	
Right renal artery	$6.5 \pm 1.2 (4.6 - 8.9)$
Left renal artery	$6.8 \pm 1.3 (4.0 - 9.4)$

Table 6.5 reports the type of stent-graft components that were deployed during the index procedure. All but one patient received the standard 3-piece system (proximal graft, distal graft, contralateral leg) – one patient received only a proximal graft, which landed in a previous open surgical graft and thus did not require a distal graft or contralateral leg.

Table 6.5 Stent-graft components deployed

Type	Percent Patients (number/total number)
Proximal graft	100% (42/42)
Distal graft	97.6% (41/42)*
Contralateral leg	97.6% (41/42)*
Ancillary components	
Main body extension	0.0% (0/42)
Additional iliac leg	7.1% (3/42)
Ipsilateral iliac leg extension	2.4% (1/42)
Contralateral iliac leg extension	7.1% (3/42)
Occluder	0.0% (0/42)
Converter	0.0% (0/42)

^{*} One patient that had undergone prior open surgical AAA repair received only the proximal fenestrated component

Table 6.6 reports the sizes (diameters and lengths) of the proximal grafts used during the initial implant procedure. The full range of available graft diameters and lengths was utilized.

Table 6.6 Proximal graft sizes used

Diameter (mm)			Length (mm)			Total
	97	107	109	122	124	
24	0	0	0	0	2	2
26	0	0	4	0	0	4
28	1	0	7	0	4	12
30	1	0	6	0	9	16
32	0	0	1	0	2	3
34	0	2	0	2	0	4
36	0	0	0	1	0	1
Total	2	2	18	3	17	42

Table 6.7 reports the sizes (diameters and lengths) of the distal grafts used during the initial implant procedure. The full range of available graft diameters and lengths was utilized.

Table 6.7 Distal graft sizes used

Diameter (mm)				Length (mm)				Total
	119	121	136	138	151	153	168	
12	0	2	2	1	0	3	0	8
16	0	2	3	4	2	2	4	17
20	1	0	1	3	1	5	3	14
24	0	0	1	0	1	0	0	2
Total	1	4	7	8	4	10	7	41

² In 2.4% (1/42) of patients previous surgery at intended access site was reported as unknown.

 Table 6.8 reports the sizes (diameters and lengths) of the contralateral leg grafts used during the initial implant procedure.

Table 6.8 Contralateral leg sizes used

Diameter (mm)	Length (mm)						Total
	54	56	71	73	88	90	
12	3	1	2	0	0	0	6
14	1	0	10	0	1	0	12
16	3	0	3	1	0	0	7
18	6	0	3	0	1	0	10
20	0	0	3	0	0	1	4
22	0	0	1	0	0	0	1
24	0	0	1	0	0	0	1
Total	13	1	23	1	2	1	41

The location of the most proximal graft margin relative to the renal arteries, SMA, and celiac artery is provided in **Table 6.9**. The proximal margin of the graft was above the renal arteries in all patients.

Table 6.9 Graft Location

Location of proximal graft margin relative to specified vessel		Percent Patients (number/total number)		
Renal arteries	Above	100.0% (42/42)		
Renaranteries	Below	0.0% (0/42)		
SMA	Above	66.7% (28/42)		
	Below	33.3% (14/42)		
C-li	Above	0.0% (0/42)		
Celiac	Polow	100 00/- (42 /42)		

The specific graft fenestration/scallop configurations that were utilized to accommodate the vessels intended to remain patent are provided in **Table 6.10**. The most commonly used configuration was 2 fenestrations and 1 scallop.

Table 6.10 Fenestrated configurations used

Configurations	% (n/N)
1 scallop	9.5% (4/42)
1 small fenestration and 1 scallop	11.9% (5/42)
2 small fenestrations	4.8% (2/42)
2 small fenestrations and 1 scallop	69.0% (29/42)
2 small fenestrations and 1 large fenestration	4.8% (2/42)

Table 6.11 provides the total number of each stent type used during the initial implant procedure. Eleven (11) patients received a Zenith Alignment Stent.

Table 6.11 Type and number of fenestration stents used

Stent type/description	(n)
Zenith Alignment Stent (uncovered, balloon-expandable, 316 L stainless steel)	22
Uncovered, balloon-expandable, 316L stainless steel biliary stent	28
Uncovered, balloon-expandable, 316L stainless steel biliary/iliac stent	20
Uncovered, balloon-expandable, 316L stainless steel biliary/renal stent	8
Covered, balloon-expandable, 316L stainless steel tracheobronchial stent	2

Table 6.12 indicates which vessels were targeted by either a fenestration or scallop and were either stented or unstented. All vessels accommodated by a small fenestration were stented.

Table 6.12 Fenestration and vessel stenting

Vessel	Small fenestration		Large fe	nestration	Sca	Total	
	Stented	Unstented	Stented	Unstented	Stented	Unstented	
Celiac	0	0	0	0	0	0	0
SMA	0	0	0	2	0	29	31
Right renal	35	0	0	0	3	0	38
Left renal	36	0	0	0	4	2	42
Accessory	0	0	0	0	0	0	0
Total	71	0	0	2	7	31	111

Primary Endpoint

Table 6.13 reports the 6-month treatment success for the Zenith Fenestrated AAA Endovascular Graft, as compared to the matched patients treated with the standard Zenith AAA Endovascular Graft. Of 42 patients enrolled in the clinical study, 40 were evaluable for the primary endpoint analysis (two patients were lost to follow-up). The 6-month treatment success rate was 97.5% in the Fenestrated endovascular treatment group compared to 95% in the matched Zenith AAA cohort.

Table 6.13 Results for 6-month treatment success

Measure	Zenith Fenestrated	Zenith AAA
Treatment success	97.5% (39/40) ¹	95.0% (38/40) ²

¹Failure due to bowel ischemia

Safety Data

Table 6.14 reports the technical success results for the Zenith Fenestrated AAA Endovascular Graft, which was defined as successful access of the aneurysm site and deployment of the Zenith Fenestrated AAA Endovascular Graft in the intended location, with all vessels targeted by fenestrations patent at the completion of the procedure. Technical success was 100%.

Table 6.14 Technical success

Measure	Percent (n/N)	
Technical success	100.0% (42/42)	

 $^{{}^2\}textit{Failure due to congestive heart failure in one and congestive heart failure as well as cardiac is chemia requiring intervention in another.}$

Table 6.15 reports the Kaplan-Meier survival estimates for freedom from major adverse events (MAE) within 30 days (death, Q-wave MI, bowel ischemia, paralysis, stroke, reintubation, renal failure requiring dialysis). There were no MAEs in the Zenith AAA group within 30 days.

Table 6.15 Results from Kaplan-Meier analysis for freedom from 30-day MAE

Event	Parameter	Zenith Fenestrated	
Any MAE	Number at risk	41	
•	Cumulative events	1	
	Cumulative censored	0	
	Kaplan-Meier estimate	0.976	
	Standard error	0.024	
Death	Number at risk	42	
	Cumulative events	0	
	Cumulative censored	0	
	Kaplan-Meier estimate	1.000	
	Standard error	0.000	
Q-wave MI	Number at risk	42	
	Cumulative events	0	
	Cumulative censored	0	
	Kaplan-Meier estimate	1.000	
	Standard error	0.000	
Bowel ischemia	Number at risk	41	
	Cumulative events	1	
	Cumulative censored	0	
	Kaplan-Meier estimate	0.976	
	Standard error	0.024	
Paralysis	Number at risk	42	
•	Cumulative events	0	
	Cumulative censored	0	
	Kaplan-Meier estimate	1.000	
	Standard error	0.000	
Stroke	Number at risk	42	
	Cumulative events	0	
	Cumulative censored	0	
	Kaplan-Meier estimate	1.000	
	Standard error	0.000	
Re-intubation	Number at risk	42	
	Cumulative events	0	
	Cumulative censored	0	
	Kaplan-Meier estimate	1.000	
	Standard error	0.000	
Renal failure requiring dialysis	Number at risk	42	
3 - 7 - 7	Cumulative events	0	
	Cumulative censored	0	
	Kaplan-Meier estimate	1.000	
	Standard error	0.000	

Table 6.16 provides the Kaplan-Meier estimates for freedom from serious adverse events (death [all-cause and AAA-related], rupture, and conversion), as compared to Zenith AAA. AAA-related death was defined as any death occurring within 30 days of the initial implant procedure (or secondary intervention) or any death determined by the independent clinical events to be related. The cause of death was unknown in one patient from the Zenith Fenestrated group, which the CEC was therefore unable to adjudicate – this was the only patient death counted as AAA-related in the Zenith Fenestrated group. No aneurysm ruptures or conversions to open repair where reported in the Zenith Fenestrated group.

 $Table\,6.16\,Results\,from\,Kaplan-Meier\,analysis\,for\,serious\,adverse\,events\,(death,\,rupture,\,conversion)$

Event	Group	Parameter	30 Days	365 Days	730 Days	1095 Days	1460 Days	1825 Days
All-cause mortality	Zenith Fenestrated	Number at risk	42	36	27	20	19	14
		Cumulative events	0	1	2	4	4	4
		Cumulative censored	0	5	13	18	19	24
		Kaplan-Meier estimate	1.000	0.976	0.943	0.861	0.861	0.861
		Standard error	0.000	0.024	0.040	0.066	0.066	0.066
	Zenith AAA	Number at risk	33	32	30	30	28	17
		Cumulative events	0	0	2	2	2	2
		Cumulative censored	0	1	1	1	3	14
		Kaplan-Meier estimate	1.000	1.000	0.938	0.938	0.938	0.938
		Standard error	0.000	0.000	0.043	0.043	0.043	0.043
AAA-related mortality	Zenith Fenestrated	Number at risk	42	36	27	20	19	14
•		Cumulative events	0	0	0	1*	1	1
		Cumulative censored	0	6	15	21	22	27
		Kaplan-Meier estimate	1.000	1.000	1.000	0.955	0.955	0.955
		Standard error	0.000	0.000	0.000	0.044	0.044	0.044
	Zenith AAA	Number at risk	33	32	30	30	28	17
		Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
Rupture	Zenith Fenestrated	Number at risk	42	37	28	21	21	18
•		Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	5	14	21	21	24
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
	Zenith AAA	Number at risk	33	32	30	30	28	17
		Cumulative events	0	1	1	1	1	1
		Cumulative censored	0	0	2	2	4	15
		Kaplan-Meier estimate	1.000	0.970	0.970	0.970	0.970	0.970
		Standard error	0.000	0.030	0.030	0.030	0.030	0.030
Conversion	Zenith Fenestrated	Number at risk	42	37	28	21	21	18
		Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	5	14	21	21	24
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
	Zenith AAA	Number at risk	33	32	30	30	28	17
		Cumulative events	0	1	1	1	1	1
		Cumulative censored	0	0	2	2	4	15
		Kaplan-Meier estimate	1.000	0.970	0.970	0.970	0.970	0.970
		Standard error	0.000	0.030	0.030	0.030	0.030	0.030

^{*} 1 case of death that the CEC was unable to adjudicate, which was conservatively counted as AAA-related for the purpose of analysis.

Table 6.17 reports the Kaplan-Meier survival estimates for freedom from any pre-specified cardiovascular, pulmonary, renal, GI, wound, neurologic, and vascular event reported by the investigative sites, as compared to Zenith AAA. The procedure-related incidence (i.e., within 30 days) of cardiovascular, pulmonary, gastrointestinal, wound, and neurologic events appeared comparable between Zenith Fenestrated and Zenith AAA, and the occurrence of events in these categories beyond 30 days was not surprising given the pre-existing comorbid conditions of the patient populations. The percent of patients experiencing renal events or vascular events within 30 days trended higher for Zenith Fenestrated compared to Zenith AAA patients, the details of which are provided in **Tables 6.18** and **6.19**, respectively.

Table 6.17 Kaplan-Meier estimates (freedom from morbidity	, by category)
-------------------------------------	------------------------	----------------

Event	Group	Parameter	30 Days	365 Days	730 Days	1095 Days	1460 Days	1825 Days
Cardiovasculara	Zenith Fenestrated	Number at risk	40	30	21	15	15	10
		Cumulative events	2	6	7	8	8	10
		Cumulative censored	0	8	14	19	19	22
		Kaplan-Meier estimate	0.952	0.854	0.817	0.766	0.766	0.656
		Standard error	0.033	0.055	0.064	0.078	0.078	0.098
	Zenith AAA	Number at risk	29	27	24	24	23	14
		Cumulative events	4	5	6	6	6	6
		Cumulative censored	0	1	3	3	4	13
		Kaplan-Meier estimate	0.879	0.848	0.817	0.817	0.817	0.817
		Standard error	0.057	0.062	0.068	0.068	0.068	0.068
Pulmonary ^b	Zenith Fenestrated	Number at risk	41	32	24	18	17	12
		Cumulative events	1	4	4	5	5	6
		Cumulative censored	0 0.976	6 0.903	14 0.903	19 0.855	20	24
		Kaplan-Meier estimate Standard error	0.976	0.903	0.903	0.064	0.855 0.064	0.798 0.081
	Zenith AAA	Number at risk	33	32	29	29	28	17
	Zelliti AAA	Cumulative events	0	0	1	1	1	1
		Cumulative censored	0	1	3	3	4	15
		Kaplan-Meier estimate	1.000	1.000	0.969	0.969	0.969	0.969
		Standard error	0.000	0.000	0.031	0.031	0.031	0.031
Renal ^c	Zenith Fenestrated	Number at risk	37	30	21	15	14	10
.c.iui	Zeman enestrateu	Cumulative events	5	6	8	9	10	10
		Cumulative censored	0	6	13	18	18	22
		Kaplan-Meier estimate	0.881	0.856	0.791	0.742	0.692	0.692
		Standard error	0.050	0.054	0.067	0.079	0.088	0.088
	Zenith AAA	Number at risk	33	31	29	29	27	16
		Cumulative events	0	1	1	1	1	1
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	0.970	0.970	0.970	0.970	0.970
		Standard error	0.000	0.030	0.030	0.030	0.030	0.030
GI ^d	Zenith Fenestrated	Number at risk	40	33	25	20	19	14
		Cumulative events	2	2	2	2	2	2
		Cumulative censored	0	7	15	20	21	26
		Kaplan-Meier estimate	0.952	0.952	0.952	0.952	0.952	0.952
		Standard error	0.033	0.033	0.033	0.033	0.033	0.033
	Zenith AAA	Number at risk	33	32	30	30	28	17
		Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
14/ 10	7 21 5	Standard error	0.000	0.000	0.000	0.000	0.000	0.000
Wound ^e	Zenith Fenestrated	Number at risk	41	32	24	19	18	13
		Cumulative events	1 0	3 7	3 15	3 20	3 21	3 26
		Cumulative censored Kaplan-Meier estimate	0.976	0.927	0.927	0.927	0.927	0.927
		Standard error	0.024	0.927	0.927	0.927	0.927	0.927
	Zenith AAA	Number at risk	32	31	29	29	26	15
	Zellitil AAA	Cumulative events	1	1	1	1	2	2
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	0.970	0.970	0.970	0.970	0.936	0.936
		Standard error	0.030	0.030	0.030	0.030	0.044	0.044
Neurologic ^f	Zenith Fenestrated	Number at risk	42	35	26	19	18	13
		Cumulative events	0	0	0	1	1	1
		Cumulative censored	0	7	16	22	23	28
		Kaplan-Meier estimate	1.000	1.000	1.000	0.950	0.950	0.950
		Standard error	0.000	0.000	0.000	0.049	0.049	0.049
	Zenith AAA	Number at risk	33	32	28	28	26	15
		Cumulative events	0	0	2	2	2	2
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	1.000	0.934	0.934	0.934	0.934
		Standard error	0.000	0.000	0.045	0.045	0.045	0.045
Vascularg	Zenith Fenestrated	Number at risk	34	28	20	17	16	10
		Cumulative events	8	8	8	8	8	9
		Cumulative censored	0	6	14	17	18	23
			0.810	0.810	0.810	0.810	0.810	0.759
		Kaplan-Meier estimate						
		Standard error	0.061	0.061	0.061	0.061	0.061	0.075
	Zenith AAA	Standard error Number at risk	0.061 30	0.061 28	26	26	23	13
	Zenith AAA	Standard error Number at risk Cumulative events	0.061 30 3	0.061 28 4	26 4	26 4	23 5	13 5
	Zenith AAA	Standard error Number at risk Cumulative events Cumulative censored	0.061 30 3 0	0.061 28 4 1	26 4 3	26 4 3	23 5 5	13 5 15
	Zenith AAA	Standard error Number at risk Cumulative events	0.061 30 3	0.061 28 4	26 4	26 4	23 5	13 5

^a Cardiovascular pre-specified events: cardiac ischemia requiring intervention, inotropic support, Q-wave MI, non-Q-wave MI, congestive heart failure (CHF), arrhythmia requiring intervention or new treatment, and medically intractable hypertension.

Given the involvement of the renal arteries in the repair with a fenestrated graft, there is an expected higher risk for renal adverse events as compared to use of a standard, non-fenestrated endograft. Renal morbidity was therefore closely monitored during the study by evaluating several pre-specified events (renal infarct, renal insufficiency, renal failure requiring dialysis, renal artery occlusion). **Table 6.18** reports the Kaplan-Meier estimates for freedom from the individual pre-specified renal morbid events, as compared to Zenith AAA. Also included in **Table 6.18** is the Kaplan-Meier estimate for freedom from stenosis/compression events that required reintervention, as also counted (as device/renal stenosis) in **Table 6.24** (Reasons for Secondary Intervention).

There were five incidental findings of renal infarct on imaging (without an associated clinical event). Each occurred in a patient from the Zenith Fenestrated group that had some degree of calcification/thrombus in the sealing zone (one also with a history of infarct and coverage of an accessory renal artery at the time of the procedure).

Renal insufficiency was observed in Zenith Fenestrated (n=3) as well as Zenith AAA (n=1) patients. One of the Zenith Fenestrated patients with renal insufficiency was also the only patient in either group to require dialysis, which the CEC judged to be unrelated to AAA repair due to underlying renal dysfunction. Renal insufficiency in one of the other patients from the Zenith Fenestrated group was also judged unrelated to AAA repair by the CEC due to underlying renal dysfunction.

There were two reports of renal occlusion in the Zenith Fenestrated group, neither of which was associated with graft migration. One required reintervention and occurred in a patient with suboptimal placement of the renal stent in the middle/upper portion of the fenestration. There were seven patients with stenosis/compression events requiring secondary intervention (one associated with migration), four of which had a peak systolic velocity <280 cm/s prior to reintervention.

b Pulmonary pre-specified events: pneumonia requiring antibiotics, supplemental oxygen at discharge, ventilation (>24 hours and >72 hours), and re-intubation.

^c Renal pre-specified events: renal failure requiring dialysis, renal insufficiency, renal infarct, and occlusion of fenestrated renal vessel.

^d GI pre-specified events: bowel obstruction, bowel ischemia/mesenteric ischemia, paralytic ileus >4 days, and aorto-enteric fistula.

^{*} Wound pre-specified events: incisional hernia, wound infection requiring antibiotic treatment, wound complication requiring return to the operating room (OR), seroma requiring treatment, lymph fistula, and wound breakdown requiring debridement.

f Neurologic pre-specified events: transient ischemic attack (TIA)/reversible ischemic neurological deficit (RIND), stroke, spinal cord ischemia/paralysis.

⁹ Vascular pre-specified events: embolization resulting in tissue loss or requiring intervention; limb thrombosis, aneurysm leak/rupture, pseudoaneurysm, increase in aneurysm size by >0.5 cm, vascular injury, and post-procedure transfusion.

Table 6.18 Kaplan-Meier estimates for freedom from pre-specified renal events (regardless of whether determined by the clinical events committee to be related or unrelated to AAA repair)

		committee to be relat						
Event	Group	Parameter	30 Days	365 Days	730 Days	1095 Days	1460 Days	1825 Days
Renal infarct*	Zenith Fenestrated	Number at risk	37	31	24	19	18	14
		Cumulative events	5 ^{a,b,c,d,e}	5	5	5	5	5
		Cumulative censored	0	6	13	18	19	23
		Kaplan-Meier estimate	0.881	0.88	0.881	0.881	0.881	0.881
		Standard error	0.050	0.050	0.050	0.050	0.050	0.050
	Zenith AAA	Number at risk	33	32	30	30	28	17
		Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
Renal insufficiency**	Zenith Fenestrated	Number at risk	42	35	25	18	16	11
on two or more		Cumulative events	0	0	1 ^f	2 ⁹	3 ^h	3
follow-up tests)		Cumulative censored	0	7	16	22	23	28
·		Kaplan-Meier estimate	1.000	1.000	0.963	0.912	0.862	0.862
		Standard error	0.000	0.000	0.036	0.060	0.075	0.075
	Zenith AAA	Number at risk	33	31	29	29	27	16
		Cumulative events	0	1	1	1	1	1
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	0.970	0.970	0.970	0.970	0.970
		Standard error	0.000	0.030	0.030	0.030	0.030	0.030
Dialysis***	Zenith Fenestrated	Number at risk	42	35	26	20	18	13
3141/313		Cumulative events	0	0	0	0	19	1
		Cumulative censored	0	7	16	22	23	28
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	0.947	0.947
		Standard error	0.000	0.000	0.000	0.000	0.051	0.051
	Zenith AAA	Number at risk	33	32	30	30	28	17
	Zellitil AAA	Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
Renal occlusion	Zenith Fenestrated	Number at risk	42	34	24	18	18	13
relial occiusion	Zenitn renestrateu	Cumulative events	0	1 ⁱ	2 ⁴	2	2	2
		Cumulative events Cumulative censored	0	7	16	22	22	27
		Kaplan-Meier estimate	1.000	0.975	0.945	0.945	0.945	0.945
		Standard error	0.000	0.975	0.943	0.943	0.945	0.943
	Zenith AAA	Number at risk	33	32	30	30	28	17
	Zenith AAA	Cumulative events	0	0	0	0	28 0	0
			0	1	3	3	5	16
		Cumulative censored	-		-	-	-	
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
Stenosis/compression	Zenith Fenestrated	Number at risk	41	33	24	18	16	11
requiring		Cumulative events	1 ^k	3 ^{l,m}	4 ⁿ	5°	6 ^b	7º
reintervention		Cumulative censored	0	6	14	19	20	24
		Kaplan-Meier estimate	0.976	0.927	0.897	0.850	0.800	0.747
		Standard error	0.024	0.040	0.049	0.065	0.078	0.089

^{*} As reported by sites, regardless of whether confirmed by core lab.

Given the expected longer procedure times for Zenith Fenestrated compared to Zenith AAA (252.2 \pm 75.5 minutes for Zenith Fenestrated vs. 160.6 \pm 60.6 minutes for Zenith AAA) and correspondingly greater amount of procedural blood loss (537.4 \pm 498.5 cc for Zenith Fenestrated vs. 281.2 \pm 192.4 cc for Zenith AAA), the need for post-procedure transfusion was also greater, as shown in **Table 6.19**, which provides the Kaplan-Meier estimates for freedom from individual pre-specified vascular events occurring in either the Zenith Fenestrated or Zenith AAA groups.

^{**} Creatinine rise >2 mg/dl and >30% from baseline

^{***} Although dialysis in patients with a normal pre-operative renal function was pre-specified, the analysis was performed with consideration to dialysis in any patient.

e' (0111011) Incidental finding on imaging without an associated clinical event at the time of reporting; patient noted to have (mild) thrombus and calcification in the seal zone on pre-procedure imaging.

b (0511009) Incidental finding on imaging without an associated clinical event at the time of reporting; patient noted to have (mild) thrombus and calcification in the seal zone on pre-procedure imaging; patient also with bilateral renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) treated with bilateral angioplasty and stenting. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

⁽⁰⁵¹¹⁰¹⁰⁾ Incidental finding on imaging without an associated clinical event at the time of reporting; patient noted to have (moderate) thrombus and (mild) calcification in the seal zone on pre-procedure imaging; patient also with hydronephrosis on POD# 237.

^a (1111002) Incidental finding on imaging without an associated clinical event at the time of reporting; patient noted to have (moderate) calcification in the seal zone as well as renal infarct on pre-procedure imaging, and also underwent intentional coverage of an accessory renal artery at the time of aneurysm repair.

e (1111007) Incidental finding on imaging without an associated clinical event at the time of reporting; patient noted to have (mild) calcification in the seal zone on pre-procedure imaging; patient also with renal insufficiency (creatinine rise >2 mg/dl and >30% from baseline) on a single follow-up (at POD# 424).

f (021101) Patient with a decrease in GFR >30% at the 24-month follow-up, but not on subsequent follow-up at 36 months. The patient underwent secondary intervention to treat a Type II endoleak (on POD# 239) and hospitalization for congestive heart failure treated with Lasix (on POD# 314), but there were no reports of renal artery stenosis or occlusion at any time point.

g (0111006) Patient also with renal calculi noted on POD# 214; all events (renal calculi, renal insufficiency, and dialysis) were determined un-related to AAA repair by the CEC.

h (0421001) Patient also with stenosis of an unstented renal artery proximal to the graft margin, which underwent stenting on POD# 1221; all events (renal insufficiency and stenosis) was determined up related to AAA renair by the CEC

insufficiency and stenosis) were determined un-related to AAA repair by the CEC.

(0211008) No evidence of graft migration, but with compression of the fenestration stent (uncovered, balloon-expandable 316L stainless steel biliary/renal stent), due likely to suboptimal deployment of the renal stent into the middle/upper portion of the fenestration; patient underwent secondary intervention (ilio-renal

¹ (0611003) No evidence of graft migration or fenestration stent compression (uncovered, balloon-expandable 316L stainless steel biliary/iliac stent), suggesting occlusion likely resulted from the development and progression of thrombus or intimal hyperplasia within the stented vessel; patient did not undergo secondary intervention; patient also with site-reported atrophy of kidney (on POD# 177).

^{* (0211011)} Angiography revealed that the right renal artery was severely stenosed. Attempted cannulation was unsuccessful, as the fenestration stent (Zenith* Alignment Stent) was not flared at the time of the initial implant procedure. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

⁽⁰¹¹¹⁰⁰⁸⁾ Right renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) was treated by angioplasty and additional stent placement. No evidence of gardt migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

⁽⁰¹¹¹⁰¹⁴⁾ Bilateral renal artery stenoses (uncovered, balloon-expandable 316L stainless steel biliary stent) were treated by angioplasty and additional stent placement. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

[&]quot; (0211007) Right renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary/renal stent) was treated by angioplasty and additional stent placement. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

^{° (0511006)} Right renal artery stent compression and subsequent stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) treated by angioplasty and stent placement. Compression of fenestration stent associated with graft migration (approximately 12 mm by 24 months) due likely to longitudinal progression of disease with further a

p (0511003) Left renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) from slight compression of fenestration stent (with no measurable graft movement > 5 mm) treated by angioplasty and stent placement.

Table 6.19 Kaplan-Meier estimates for freedom from pre-specified vascular events occurring in either Zenith Fenestrated or Zenith AAA

Event	Group	Parameter	30 Days	365 Days	730 Days	1095 Days	1460 Days	1825 Days
Embolization resulting	Zenith Fenestrated	Number at risk	41	34	25	19	18	13
in tissue loss or		Cumulative events	1	1	1	1	1	1
requiring intervention		Cumulative censored	0	7	16	22	23	28
		Kaplan-Meier estimate	0.976	0.976	0.976	0.976	0.976	0.976
		Standard error	0.024	0.024	0.024	0.024	0.024	0.024
	Zenith AAA	Number at risk	33	32	30	30	28	17
		Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
Limb thrombosis	Zenith Fenestrated	Number at risk	42	35	26	20	19	14
		Cumulative events	0	0	0	0	0	0
		Cumulative censored	0	7	16	22	23	28
		Kaplan-Meier estimate	1.000	1.000	1.000	1.000	1.000	1.000
		Standard error	0.000	0.000	0.000	0.000	0.000	0.000
	Zenith AAA	Number at risk	33	31	29	29	27	16
		Cumulative events	0	1	1	1	1	1
		Cumulative censored	0	1	3	3	5	16
		Kaplan-Meier estimate	1.000	0.970	0.970	0.970	0.970	0.970
		Standard error	0.000	0.030	0.030	0.030	0.030	0.030
Post-procedure	Zenith Fenestrated	Number at risk	34	28	20	17	16	10
transfusion		Cumulative events	8	8	8	8	8	9
		Cumulative censored	0	6	14	17	18	23
		Kaplan-Meier estimate	0.810	0.810	0.810	0.810	0.810	0.759
		Standard error	0.061	0.061	0.061	0.061	0.061	0.075
	Zenith AAA	Number at risk	30	29	27	27	24	14
		Cumulative events	3	3	3	3	4	4
		Cumulative censored	0	1	3	3	5	15
		Kaplan-Meier estimate	0.909	0.909	0.909	0.909	0.875	0.875
		Standard error	0.050	0.050	0.050	0.050	0.058	0.058

Effectiveness Data

Table 6.20 reports the percent of patients with an increase (> 5 mm), decrease (> 5 mm), or no change (≤ 5 mm) in aneurysm size at each follow-up time point, as compared to pre-discharge based on the results from core lab analysis. There were two cases of aneurysm expansion, both of which occurred in patients with a persistent Type II endoleak.

Table 6.20 Change in aneurysm size based on results from core lab analysis

Item	1-month	6-month	12-month	24-month	36-month	48-month	60-month
Increase (> 5mm)	0.0% (0/39)	0.0% (0/38)	0.0% (0/29)	0.0% (0/26)	6.7% (1/15)1	6.3% (1/16)2	0.0% (0/11)
Decrease (> 5mm)	2.6% (1/39)	50.0% (19/38)	69.0% (20/29)	69.2% (18/26)	73.3% (11/15)	75.0% (12/16)	72.7% (8/11)
No change (≤ 5mm)	97.4% (38/39)	50.0% (19/38)	31.0% (9/29)	30.8% (8/26)	20.0% (3/15)	18.8% (3/16)	27.3% (3/11)

¹ Patient 0511004 had a persistent Type II endoleak requiring secondary intervention at 1393 days post-procedure.

Table 6.21 reports endoleaks by type, as assessed by the core lab at each exam period. Except for two endoleaks of unknown type, all other reported endoleaks were Type II.

Table 6.21 Endoleak based on results from core lab analysis

Type	Pre-discharge	1-month	6-month	12-month	24-month	36-month	48-month	60- month
Any (new only)	32.5% (13/40)	2.4% (1/41)	5.3% (2/38)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)
Any (new and persistent)	32.5% (13/40)	22.0% (9/41)	23.7% (9/38)	27.6% (8/29)	15.4% (4/26)	12.5% (2/16)	12.5% (2/16)	0.0% (0/11)
Multiple	0.0% (0/40)	0.0% (0/41)	0.0% (0/38)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)
Proximal Type I	0.0% (0/40)	0.0% (0/41)	0.0% (0/38)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)
Distal Type I	0.0% (0/40)	0.0% (0/41)	0.0% (0/38)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)
Type II	30.0% (12/40)	22.0% (9/41)	21.1% (8/38)	27.6% (8/29)	15.4% (4/26)	12.5% (2/16)	12.5% (2/16)	0.0% (0/11)
Type III	0.0% (0/40)	0.0% (0/41)	0.0% (0/38)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)
Type IV	0.0% (0/40)	0.0% (0/41)	0.0% (0/38)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)
Unknown	2.5% (1/40)	0.0% (0/41)	2.6% (1/38)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)

Table 6.22 reports the percent of patients with CEC-confirmed radiographic migration (≥10 mm) or clinically significant migration (measurable movement of the stent-graft >5 mm and that developed a type I endoleak or renal stenosis/occlusion with demonstrable deformation of the mating renal stent by core lab) at each follow-up time point (date of first occurrence). There were two reports of migration, one of which required secondary intervention (due to associated renal stenosis). Neither case was associated with aneurysm growth or endoleak. Both cases of migration occurred in patients with evidence of disease progression at follow-up (without aneurysm pressurization).

Table 6.22 CEC-confirmed migration (date of first occurrence)

Item	1-montn	6-montn	12-montn	24-montn	36-montn	48-montn	60-montn
Radiographic migration	0.0% (0/40)	0.0% (0/38)	0.0% (0/30)	3.6%1 (1/28)	0.0% (0/16)	0.0% (0/16)	9.1%² (1/11)
Clinically significant migration	0.0% (0/40)	0.0% (0/38)	0.0% (0/30)	3.6%1 (1/28)	0.0% (0/16)	0.0% (0/16)	0.0% (0/11)

¹ Patient 0511006 with renal stenosis from associated stent compression (uncovered, balloon-expandable 316L stainless steel biliary stent) requiring secondary intervention. Longitudinal progression of disease with further aortic neck dilatation likely resulted in migration. There was no endoleak or increase in aneurysm size in this patient. The total amount of graft movement detected at the time of the clinically significant migration was approximately 12 mm (relative to the collision).

Device integrity observations are summarized in **Table 6.23**. Losses in device integrity included three patients with barb separation, one patient with possible fenestration stent fracture, and one patient with seal stent and fenestration stent fracture (who also had evidence of disease progression during follow-up in the absence of aneurysm pressurization). None of the integrity findings were associated with adverse clinical sequelae or the need for reintervention. Although not associated with a device integrity loss (i.e., fracture), other observations included 4 cases with fenestration stent deformation/compression (1 also with migration, 3 without migration), 2 of which underwent reintervention to treat stenosis.

² Patient 0211010 had a persistent Type II endoleak requiring secondary intervention at 239 days post-procedure, but the Type II endoleak was still evident on the 48-month exam.

² Patient 0511008 was without any associated renal stenosis requiring reintervention and additionally did not have any endoleak or increase in aneurysm size. Longitudinal progression of disease with further aortic neck dilatation likely resulted in migration. The total amount of graft movement was approximately 10 mm (relative to the celiac), which retrospectively occurred over 60 months. No interventions have been performed on this patient.

Table 6.23 Device integrity findings by core lab (time of first occurrence)

Finding	Pre-discharge	1-month	6- month	12-month	24-month	36-month	48-month	60- month
			Ste	nt-graft				
Barb separation	0.0% (0/42)	0.0% (0/41)	2.6% (1/39)1	3.4% (1/29)2	3.8% (1/26)3	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Stent fracture (single)	0.0% (0/42)	0.0% (0/41)	0.0% (0/39)	3.4% (1/29)4	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Stent fracture (multiple)	0.0% (0/42)	0.0% (0/41)	0.0% (0/39)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Component separation	0.0% (0/42)	0.0% (0/41)	0.0% (0/39)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Limb separation	0.0% (0/42)	0.0% (0/41)	0.0% (0/39)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Stent-to-graft separation	0.0% (0/42)	0.0% (0/41)	0.0% (0/39)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Other	0.0% (0/42)	0.0% (0/41)	0.0% (0/39)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Fenestration stent								
Fracture	0.0% (0/42)	0.0% (0/41)	2.6% (1/39)5	3.4% (1/29)4	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Separation	0.0% (0/42)	0.0% (0/41)	0.0% (0/39)	0.0% (0/29)	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)
Other	0.0% (0/42)	0.0% (0/41)	7.7% (3/39)6,7,8	3.4% (1/29)9	0.0% (0/26)	0.0% (0/16)	0.0% (0/17)	0.0% (0/11)

¹ Patient 0421003: Separation of a single fixation barb. No clinical sequelae related to the barb separation have been reported

Table 6.24 summarizes the site reported reasons for secondary intervention. Of the 11 patients who underwent a secondary intervention, 7 did so because of renal stenosis (1 associated with graft migration and stent deformation, 1 associated with stent deformation without migration). In 4 patients, the peak systolic velocity was <280 cm/s prior to reintervention. The other reported reasons for reintervention included renal occlusion in 1, Type II endoleak in 2, and suspected Type I endoleak in 1 (ruled out by angiogram).

Table 6.24 Reasons for secondary intervention (as reported by site)							
Finding	0-30 Days	31-365 Days	366-730 Days	731-1095 Days	1096-1460 Days	1461-1825 Days	
Aneurysm rupture	0	0	0	0	0	0	
Symptoms	0	0	0	0	0	0	
Device/renal stenosis	11	25,6	17	1 ⁸	19	111	
Device migration	0	0	0	0	0	0	
Device separation	0	0	0	0	0	0	
Occlusion	0	12	0	0	0	0	
Device kink	0	0	0	0	0	0	
Infection	0	0	0	0	0	0	
Endoleak							
Type I proximal	0	1 ³	0	0	0	0	
Type I distal	0	0	0	0	0	0	
Type IIA (vessel perfusion)	0	14	0	0	110	0	
Type IIB (vessel perfusion)	0	0	0	0	0	0	
Type III (graft overlap joint)	0	0	0	0	0	0	
Type IV (through graft body)	0	0	0	0	0	0	
unknown	0	0	0	0	0	0	
Other	0	0	0	18	110	0	

¹ Patient 0211011: Angiography revealed that the right renal artery was severely stenosed. Attempted cannulation was unsuccessful, as the fenestration stent (Zenith® Alignment Stent) was not flared at the time of the initial implant procedure. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

Summary

This study enrolled 42 patients treated with the Zenith Fenestrated AAA Endovascular Graft, a line extension of the Zenith AAA Endovascular Graft that is customized to the individual anatomy of patients having an infrarenal aortic neck length that is too short for the standard endograft. A variety of fenestration/scallop configurations were utilized, the most frequent of which was 2 small fenestrations and 1 scallop. A total of 111 vessels were targeted by either a fenestration or scallop, 78 of which received a fenestration stent (all stented vessels were main renal arteries accommodated by either a small fenestration or a scallop), including 21 vessels (11 patients) receiving 22 Zenith Alignment Stents. All devices deployed successfully in the intended location, and all graft components and vessels targeted by a fenestration were patent upon completion of deployment, yielding a technical success rate of 100%.

The primary safety and effectiveness data showed that the 6-month treatment success rate for Zenith Fenestrated (97.5%) was similar to that for matched patients treated with Zenith AAA (95%).

There were no ruptures or conversions following treatment with Zenith Fenestrated at any time point. Only one death was counted as AAA-related because the cause was unknown and the CEC was therefore unable to adjudicate it – all other deaths in the Zenith Fenestrated group (3) were determined unrelated to AAA-repair by the CEC.

Pre-specified renal adverse events included renal infarct, renal insufficiency, renal failure requiring dialysis, and occlusion of a fenestrated renal vessel. There were five patients with renal infarct (none were associated with a clinical event), each of which occurred in a patient with some degree of either thrombus or calcification in the seal zone (as well as a history of infarct and coverage of an accessory renal in one). Two of three patients with renal insufficiency in the Zenith Fenestrated group had renal dysfunction prior to treatment and were considered unrelated to AAA-repair by the CEC, one of which was also the only patient in the Zenith Fenestrated group requiring dialysis (also unrelated according to the CEC). Two patients developed occlusion of a fenestrated renal vessel (neither was associated with graft migration), one of which had evidence of fenestration stent compression (from suboptimal stent placement in the mid/upper portion of the fenestration) that required reintervention.

There were no reports of Type I or Type III endoleak, and the only reports of aneurysm growth (2) occurred in patients with a Type II endoleak. There were 2 reports of migration, both in patients with evidence of disease progression at follow-up (without aneurysm pressurization), one of which had associated fenestration stent compression requiring secondary intervention. One patient was noted to have fracture of a fenestration stent as well as the seal stent on the Fenestrated Graft, neither of which resulted in endoleak, a clinical renal event, or the need for secondary intervention. This patient also exhibited disease

² Patient 0111009: Separation of a single fixation barb. No clinical sequelae related to the barb separation have been reported.

³ Patient 0511008: Separation of two barbs. No clinical sequelae related to the barb separation have been reported, although radiographic migration (approximately 10 mm over 5 years) was observed and was due likely to longitudinal progression of disease with further aortic neck dilatation.

⁴ Patient 0411001: Fracture of sealing stent (at the distal edge of the scallop fenestration) and left renal fenestration stent (uncovered, balloon-expandable 316L stainless steel biliary/iliac stent), but in a patient with progressive aneurismal disease within and proximal to the treated segment, which likely resulted in uncharacteristic tension/loading of the stents. No subsequent renal events, endoleak, or secondary interventions reported in this patient.

⁵ Patient 0511010: Fracture of left renal fenestration stent (Zenith* Alignment Stent) not readily confirmed based on subsequent bench top CT imaging studies that showed the same appearance of fracture, but in an entirely intact stent.

⁶ Patient 1111011: Deformation of fenestration stent (Zenith® Alignment Stent) with no measurable graft movement > 5 mm and not requiring secondary intervention.

⁷ Patient 0511003: Slight compression of fenestration stent (uncovered, balloon-expandable 316L stainless steel biliary stent) with no measurable graft movement > 5 mm. Angioplasty and stent placement was performed 1539 days post-procedure to treat stenosis.

⁸ Patient 0511007: Slight compression of fenestration stent (uncovered, balloon-expandable 316L stainless steel biliary stent) with no measurable graft movement > 5 mm and not requiring secondary intervention.

⁹ Patient 0511006: Compression of fenestration stent (uncovered, balloon-expandable 316L stainless steel biliary stent) associated with graft migration (approximately 12 mm by 24 months) due likely to longitudinal progression of disease with further aortic neck dilatation. Angioplasty and stent placement were performed 883 days post-procedure to treat stenosis.

² Patient 0211008: Angiogram demonstrated occluded left renal artery with proximal compression of the left renal stent (uncovered, balloon-expandable 316L stainless steel biliary/renal stent), which was treated with iliorenal bypass. Compression without evidence of migration due likely to suboptimal deployment of the renal stent into the middle/upper portion of the fenestration

³ Patient 0411004: Selective left renal angiography was performed for suspected Type I endoleak. No type I endoleak was identified; however, Type II endoleak was identified but not treated.

Patient 0211010: Persistent Type II endoleak was treated by coil embolization.

⁵ Patient 0111008: Right renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) was treated by angioplasty and additional stent placement. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

⁶ Patient 0111014: Bilderal renal artery stenoses (uncovered, balloon-expandable 316L stainless steel biliary stent) were treated by angioplasty and additional stent placement. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented seament.

⁷ Patient 0211007: Right renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary/renal stent) was treated by angioplasty and additional stent placement. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

⁸ Patient 0511006: Right renal artery stent compression and subsequent stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) treated by angioplasty and stent placement. Compression of fenestration stent associated with graft migration (approximately 12 mm by 24 months) due likely to longitudinal progression of disease with further aortic neck dilatation.

^o Patient 0511009: Bilateral renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) was treated by bilateral angioplasty and stenting. No evidence of graft migration or stent deformation, suggesting the stenosis may have resulted from local thrombus formation or intimal hyperplasia within the stented segment.

¹⁰ Patient 0511004: Underwent diagnostic angiogram for suspect Type IIa and Type III endoleak, which were not detected at 1137 days post-procedure; additional intervention performed 1393 days post-procedure, involving laparotomy, suture ligation of IIMA, and exploration of aneurysm sac as treatment for Type II endoleak with aneurysm growth.

¹¹ Patient 0511003: Left renal artery stenosis (uncovered, balloon-expandable 316L stainless steel biliary stent) from slight compression of fenestration stent (with no measurable graft movement > 5 mm) treated by angioplasty and stent placement.

progression at follow-up in the absence of aneurysm pressurization. A possible second patient with fenestration stent fracture was identified without a subsequent clinical renal event or need for reintervention.

The majority of patients who underwent reintervention following treatment with the Zenith Fenestrated Graft (7 of 11) did so for renal stenosis. There was evidence of fenestration stent deformation in 2 of 7 patients that underwent reintervention for renal stenosis (1 from suboptimal stent placement in the mid/upper portion of the fenestration, and 1 from migration due to progression of disease at follow-up in the absence of aneurysm pressurization).

7 PATIENT SELECTION AND TREATMENT

(See Warnings and Precautions)

7.1 Individualization of Treatment

Each patient must be evaluated on an individual basis to determine the specific location of the graft fenestrations (refer to Planning and Sizing Sheet), with careful consideration also given to both the potential benefits and specific risks associated with the procedure.

Considerations regarding the use of the Zenith Fenestrated AAA Endovascular Graft (see Warnings) include:

- Risk of aneurysm rupture
- Morbidity and mortality associated with conventional surgical repair
- Comorbidities
- · Size of aneurysm
- · History of renal failure
- · Life expectancy
- · Anesthetic risk
- · Age of patient
- Iliofemoral access vessel size and morphology (minimal thrombus. calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 14 French to 22 French vascular introducer sheath.

NOTE: Iliac conduits may be used to ensure the safe insertion of the delivery system

- · Non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurvsm:
- with a length that is at least 4 mm and unsuitable for a nonfenestrated graft,
- with a diameter measured outer wall to outer wall of no greater than 31 mm and no less than 19 mm
- · with an angle less than 45 degrees relative to the long axis of the aneurysm, and
- with an angle less than 45 degrees relative to the axis of the suprarenal aorta.
- Ipsilateral iliac artery distal fixation site greater than 30 mm in length and 9-21 mm in diameter (measured outer wall to outer wall).
- Contralateral iliac artery distal fixation site greater than 30 mm in length and 7-21 mm in diameter (measured outer wall to outer wall).
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.

The final treatment decision is at the discretion of the physician and patient.

8 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this fenestrated endovascular device and procedure including:

- · Risks and differences between endovascular repair (fenestrated and non-fenestrated) and surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of fenestrated endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- $\bullet \ There \ have \ been \ limited \ numbers \ of \ patients \ treated \ with \ fenestrated$ endovascular grafts when compared to non-fenestrated endovascular
- · Long-term performance of fenestrated endovascular grafts and stents in the fenestrations/scallops has not yet been established
- All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft
- Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, changes in the structure or position of the endovascular graft, or stenosis/occlusion of vessels accommodated by fenestrations) should receive enhanced follow-up.

Specific follow-up guidelines are described in Section 12, Imaging Guidelines and Post-Operative Follow-up.

9 HOW SUPPLIED

The Zenith Fenestrated AAA Endovascular Graft is supplied sterile and pre-loaded in peel-open packages. The device is intended for single use only and the fenestration/scallop location is individually tailored for each patient. Do not re-sterilize the device. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to your Cook representative or your nearest Cook office. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient. Do not use after the expiration date printed on the label. Store in a cool dry place. The proximal body and distal bifurcated body grafts are loaded into either a 6.7 mm I.D. (20 French) or 7.3 mm I.D. (22 French) Flexor introducer sheath. The sheath's surface is treated with a hydrophilic coating that, when activate, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a 4X4 gauze pad soaked in saline solution

10 CLINICAL USE INFORMATION

10.1 Physician Training

CAUTION: Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device, which requires precise planning/sizing as well as accurate longitudinal positioning and rotational orientation during placement. The recommended skill/knowledge requirements for physicians using the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System are outlined below: **Patient Selection:**

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair
- $\bullet \ Knowledge \ of \ radio graphic \ image \ interpretation, device/fenestration$ selection and sizing.

A multi-disciplinary team that has combined procedural experience with:

- · Femoral cutdown, arteriotomy and repair
- · Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques, especially accessing visceral vessels (e.g., renal arteries)
- $\bullet \ Fluoroscopic \ and \ angiographic \ image \ interpretation$
- Embolization
- Angioplasty
- · Endovascular stent-graft placement
- · Renal/visceral stent placement
- Snare techniques
- · Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities.

10.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to your Cook representative or your nearest Cook office. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

10.3 Materials Required

(Not included in 3-piece modular system)

- · Zenith AAA Endovascular Graft Ancillary Kit
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- · Heparinized saline solution

10.4 Materials Recommended

(Not included in 3-piece modular system) The following products are recommended:

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
- Cook Amplatz Ultra Stiff Wire Guides (AUS)
- · Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide; for example:
- · Cook .035 inch wire guides
- Cook Nimble[™] Wire Guides
- · Cook Rosen Wire Guide • Molding Balloons (e.g., CODA)
- · Zenith Alignment Stents · Introducer sets; for example:
 - · Cook Check-Flo® Introducer Sets
 - Cook Extra Large Check-Flo Introducer Sets
 - Cook Flexor® Balkin Up & Over® Contralateral Introducers
- Sizing catheter; for example:
- Cook Aurous® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
 - Cook Beacon® Tip Angiographic Catheters
- · Cook Beacon Tip Royal Flush Catheters · Entry needles; for example
- · Cook single wall entry needles

10.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or excessive oversizing may result in incomplete sealing or comprised flow.

Table 10.5.1 Proximal Body Graft Diameter Sizing

Intended Aortic Vessel Diameter	Main Body Diameter	Introduction Sheath Fr	Introduction Sheath I.D./O.D.
19	24	20	6.7 mm/7.7 mm
20	24	20	6.7 mm/7.7 mm
21	24	20	6.7 mm/7.7 mm
21	26	20	6.7 mm/7.7 mm
22	26	20	6.7 mm/7.7 mm
23	28	20	6.7 mm/7.7 mm
24	28	20	6.7 mm/7.7 mm
24	30	20	6.7 mm/7.7 mm
25	30	20	6.7 mm/7.7 mm
26	30	20	6.7 mm/7.7 mm
26	32	20	6.7 mm/7.7 mm
27	32	20	6.7 mm/7.7 mm
28	32	20	6.7 mm/7.7 mm
28	34	20	6.7 mm/7.7 mm
29	34	20	6.7 mm/7.7 mm
29	36	22	7.3 mm/8.5 mm
30	36	22	7.3 mm/8.5 mm
31	36	22	7.3 mm/8.5 mm

Table 10.5.2 Distal Body (Ipsilateral Limb) Graft Diameter Sizing

Intended Iliac Vessel Diameter	Ipsilateral Limb Diameter	Introduction Sheath Fr	Introduction Sheath I.D./O.D.
9	12	20	6.7 mm/7.7 mm
10	12	20	6.7 mm/7.7 mm
11	12	20	6.7 mm/7.7 mm
12	16	20	6.7 mm/7.7 mm
13	16	20	6.7 mm/7.7 mm
14	16	20	6.7 mm/7.7 mm
15	16	20	6.7 mm/7.7 mm
15	20	20	6.7 mm/7.7 mm
16	20	20	6.7 mm/7.7 mm
17	20	20	6.7 mm/7.7 mm
18	20	20	6.7 mm/7.7 mm
18	24	20	6.7 mm/7.7 mm
19	24	20	6.7 mm/7.7 mm
20	24	20	6.7 mm/7.7 mm
21	24	20	6.7 mm/7.7 mm

10.6 Device Length Sizing Guidelines

The proximal body graft and distal body graft are available in multiple lengths. The chosen lengths should provide a minimum two-stent overlap should the graft components align completely along the greater curve of the aorta/aneurysm over time. Planning for a longer overlap length initially (e.g., 3-4 stents) is therefore preferable.

Table 10.6.1 Proximal Graft Lengths

Diameter	Body Length	
mm	mm	
24	76 / 91 / 94 / 106 / 109 / 121 / 124	
26	76 / 91 / 94 / 106 / 109 / 121 / 124	
28	76 / 91 / 94 / 106 / 109 / 121 / 124	
30	76 / 91 / 94 / 106 / 109 / 121 / 124	
32	76 / 91 / 94 / 106 / 109 / 121 / 124	_
34	84 / 99 / 114 / 129 / 107 / 122 / 137	_
36	84 / 99 / 114 / 129 / 107 / 122 / 137	

Table 10.6.2 Distal Graft Lengths

Ipsilateral Limb Diameter	Body Length	Ipsilateral Limb	
mm	mm	length mm	
12	76 / 91 / 106 / 121	28 / 45 / 62	
16	76 / 91 / 106 / 121	28 / 45 / 62	
20	76 / 91 / 106 / 121	28 / 45 / 62	
24	76 / 91 / 106 / 121	28 / 45 / 62	

10.7 Graft Fenestration/Scallop Guidelines

The proximal body graft may contain up to three precisely located holes (fenestration(s)), and cut-outs from the proximal margin (scallop(s)) of the graft material. The fenestration and/or scallop locations are individualized to the patient anatomy based on measurements from high resolution pre-operative CT imaging. Refer to the Planning and Sizing Sheet for information regarding how these locations are determined.

11 INSTRUCTIONS FOR USE

11.1 General Use Information

Prior to use of the Zenith Fenestrated AAA Endovascular Graft with the H&L-B One-Shot Introduction System, review this Suggested Instructions for Use booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary.

These instructions are intended to help guide the physician and do not

take the place of physician judgment. 11.2 Pre-Implant Determinants

Verify from pre-implant planning (refer to Planning and Sizing Sheet) that the correct device has been selected. Determinants include:

- Femoral artery selection for introduction of the main body system (i.e., define respective contralateral and ipsilateral iliac arteries).
- 2. Angulation of aortic neck, aneurysm and iliacs.
- 3. Quality of the aortic neck
- 4. Diameters of infrarenal aortic neck and distal iliac vessels.
- $5.\ Distance\ from\ renal\ arteries\ to\ the\ aortic\ bifurcation.$
- Distance from the renal arteries to the hypogastric (internal iliac) arteries/attachment site(s).
- Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
- 8. Consider the degree of vascular calcification.
- 9. Size and location/position of visceral vessel origins.

NOTE: Each respective vessel diameter and length (aorta, ipsilateral iliac and contralateral iliac) provides the necessary criteria for choosing the appropriate endovascular graft.

11.3 Patient Preparation

- 1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
- 2. Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
- ${\bf 3.}\ Expose both common femoral arteries using standard surgical technique.}$
- 4. Establish adequate proximal and distal vascular control of both femoral vessels.

11.4 Fenestrated System

11.4.1 Proximal Body Graft Preparation/Flush

1. Remove black-hubbed shipping stylet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away* sheath from back of the Captor hemostatic valve. (Figure 6) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove in the tip. (Figure 7) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

NOTE: Graft flushing solution of heparinized saline is often used.

Attach syringe with normal heparinized saline to the hub on the inner cannula. Flush until fluid exits the distal tip. (Figure 8)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

 Soak 4X4 gauze pad in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator tip liberally.

11.4.2 Distal Bifurcated Body Graft Preparation/Flush

1. Remove black-hubbed shipping stylet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the Captor hemostatic valve. (Figure 9) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove in the tip. (Figure 10) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

NOTE: Graft flushing solution of Heparinised saline is always used.

2. Attach syringe with normal Heparinized saline to the hub on the inner cannula. Flush until fluid exits the distal tip. (Figure 10)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air

Soak 4X4 gauze pad in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator tip liberally.

11.4.3 Iliac Leg (Contralateral) Preparation/Flush

Refer to the Instructions for Use enclosed in the device packaging for the iliac leg graft for instruction on preparation/flush.

11.4.4 Vascular Access and Angiography

- Puncture the selected common femoral arteries using standard technique with an 18 or 19UT gage arterial needle. Upon vessel entry, insert:
- Wire guides standard .035 inch diameter, 145 cm long, J tip or Bentson Wire Guide
- Appropriate size sheaths (e.g., 6.0 or 8.0 French)
- Flush catheter (often radiopaque sizing catheters e.g., Centimeter Sizing Catheter or straight flush catheter).
- 2. Perform angiography to identify level(s) of renals, aortic bifurcation and iliac bifurcations

NOTE: If fluoroscope angulation is used with an angulated neck it may be necessary to perform angiograms using various projections.

NOTE: A previous planning exercise will have determined which side will be used to introduce the proximal and distal bodies.

11.4.5 Proximal Body Placement

CAUTION: Verify that the predetermined access site is chosen for the introduction and placement of the proximal body.

- 1. Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.
- 2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

NOTE: Monitor the patient's coagulation status throughout the procedure.

 On ipsilateral side, replace J wire with stiff wire guide (LES) .035 inch, 260 cm long and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.

NOTE: A straight angiographic catheter should be inserted up the contralateral side to aid in placement of graft.

4. Before insertion, position proximal body delivery system on patient's abdomen under fluoroscopy to assist with orientation and positioning. Rotate to a position where the anterior markers are situated in the most anterior (12:00 o'clock) position. The sidearm of the hemostatic valve may serve as an external reference to the fenestration(s) and/or scallop(s), anterior and posterior markers and body side markers.

CAUTION: Maintain wire guide position during delivery system insertions

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

5. Advance the delivery system until the radiopaque markers indicating the fenestration(s) and/or scallop(s) are at the level of the appropriate arteries. Check that the distal end of the graft is in a satisfactory position above the aortic bifurcation and that the anterior and posterior markers indicate that the graft is in satisfactory orientation. (Figure 11)

NOTE: Angiography should be performed as needed throughout deployment, to confirm correct placement of the graft.

6. Verify position of the wire guide in the thoracic aorta. Ensure that fenestration(s) and/or scallops are at the level of the appropriate arteries and the anterior markers are in the most anterior (12:00 o'clock) position.

NOTE: The vertical anterior markers, and the horizontal posterior markers should form a cross, on the fluoroscopic image, when correctly oriented. (Figure 11)

NOTE: The fenestration/scallop markers should be in close apposition to the appropriate side branch vessels.

Clear identification of fenestration position(s) may not be possible until the graft has been partially unsheathed.

 $\label{NOTE:} \textbf{Ensure the Captor Hemostatic Valve on the introducer sheath is turned to the open position. (Figure 12)}$

- 7. Stabilize the grey positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first two (2) covered stents by withdrawing the sheath while monitoring device location.
- 8. Perform angiography, and adjust graft placement as necessary. Continue to withdraw the sheath making positional adjustments as necessary.

NOTE: Techniques to ensure that the fenestration(s) and/or scallop(s) will accurately align with their respective vessels will vary, and will depend upon vessel anatomy, graft design, and physician preferences.

NOTE: If a small fenestration is being utilized, care should be taken to properly align the fenestration with the respective vessel.

- Proceed with deployment until the graft has been fully unsheathed. (Figure 13)
- 10. When a satisfactory graft position has been achieved, exchange the contralateral angiographic catheter and wire guide with a selective wire guide/selective catheter positioned just below the level of the proximal body. Cannulate the partially deployed proximal main body. Advance the selective catheter over the selective wire guide into

- the renal artery. Exchange the selective wire with a Rosen wire or equivalent wire guide.
- Utilizing contralateral access sheath and wire guide, advance a guiding catheter into each small fenestration and its respective vessel. (Figure 14)

NOTE: Non-compliant angioplasty balloons may be used as an alternative to guiding catheters.

NOTE: Cannulation of the scallop and its respective vessel may also be achieved using similar techniques.

NOTE: It is not recommended to use balloons or guiding catheters to guide final placement of large fenestrations as stent struts across fenestration may interfere.

NOTE: To ensure renal stent placement in the lower portion of the fenestration, it may be necessary to slightly advance the graft after catheter/sheath access to the renal vessels, before removal of the diameter reducing ties.

CAUTION: Before release of the diameter reducing ties, verify that the position of the ipsilateral access wire extends just distal to the aprtic arch.

CAUTION: During proximal trigger-wire removal, top cap advancement, and subsequent suprarenal stent deployment, verify that the position of the main body wire guide extends just distal to the aortic arch and that support of the system is maximized.

12. Verify proper position of proximal body. Remove the safety lock from the first (distal) gold trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the gold trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 15)

NOTE: At this point, the proximal main body graft should be fully expanded with the proximal bare stent contained within the top cap.

13. Remove the safety lock from the top stent trigger-wire release mechanism. Under fluoroscopy, withdraw and remove the triggerwire to unlock the suprarenal stent from the top cap by sliding the black trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 16)

NOTE: If resistance is felt or system bowing is noticed, the trigger-wire is under tension. Excessive force may cause the graft position to be altered if excessive resistance or delivery system movement is noted, stop and assess the situation.

If unable to remove the black trigger-wire release mechanism from the top cap, perform the following steps under fluoroscopy: $\frac{1}{2} \int_{-\infty}^{\infty} \frac{1}{2} \int_{-\infty}^{\infty}$

- a. Remove tension on the trigger-wire by loosening the pin vise and slightly pulling the inner cannula to move the top cap down over the suprarenal stent. Avoid compressing the Zenith Fenestrated proximal body.
- b. Retighten the pin vise.
- c. Remove the black trigger-wire release mechanism

d. Continue with (14) in Section 11.4.5, Proximal Body Placement.

NOTE: If still unable to remove the black trigger-wire release mechanism from the top cap, see Section 13 Trigger-Wire Release Troubleshooting. **NOTE:** The distal stent is still secured by the trigger-wire.

 Loosen the pin vise. (Figure 17) Control the position of the graft by stabilizing the grey positioner of the introducer.

CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

15. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed. (Figures 18a and 18b) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise to avoid contact with the deployed suprarenal stent.

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

NOTE: If unable to deploy the suprarenal stent by advancing the top cap, see **Section 14, Suprarenal Stent Deployment Troubleshooting**.

16. Remove the safety lock from the second (proximal) white trigger-wire release mechanism. Withdraw and remove the trigger-wire to detach the distal end of the endovascular graft from the delivery system by sliding the white trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Figure 19)

NOTE: Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

11.4.6 Docking of Top Cap

- 1. Loosen the pin vise. (Figure 20)
- 2. Secure sheath and inner cannula to avoid any movement of these components.
- 3. Advance the grey positioner over the inner cannula until it docks with the top cap. (Figures 21a, 21b and 21c)

 $\begin{tabular}{ll} \textbf{NOTE:} If resistance occurs, slightly rotate grey positioner and continue to gently advance. \end{tabular}$

- 4. Retighten the pin vise and withdraw the entire top cap and grey positioner through the graft and through the sheath by pulling on the pink hub of the inner cannula. (Figure 22) Leave the sheath and wire guide in place.
- Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops.

11.4.7 Fenestration Stent Placement and Deployment General Use Information

In the event that small fenestrations are being utilized, stents may be placed to secure positive alignment. Standard techniques for placement of arterial stents should be employed during use of stents. The Zenith Alignment Stent is available for this purpose. Refer to the Zenith Alignment Stent Suggested Instructions for Use for details.

- Return to the guide catheter and wire guide which cannulate the small fenestration and respective vessel.
- Introduce appropriately sized balloon expandable stent and advance to the ostium of the fenestration/vessel. Advance stent partially into the vessel, leaving approximately 4-5 mm of stent in the aorta. (Figure 23)

NOTE: Fluoroscopic views tangential to the fenestration will optimize visualization of the stent position relative to the stent graft.

- 3. Expand stent
- Remove the balloon and replace with an oversized angioplasty balloon.
 Advance the balloon until the proximal tip is positioned at the ostium.
- Inflate the balloon to flare the intra-aortic segment of the stent (Figure 24)

CAUTION: This technique requires high quality imaging. Mobile image intensifiers may provide less than adequate imaging quality.

- 6. Remove the angioplasty balloon.
- 7. Withdraw renal access sheaths, catheters and wire guides in the contralateral side to a level just above the aortic bifurcation.

NOTE: In the event that there is more than one fenestration, repeat the preceding steps for each additional small fenestration.

11.4.8 Distal Bifurcated Body Placement

- Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.
- 2. Check flushing solutions. Flush after each catheter and/or wire guide exchange.
- 3. Before insertion, position distal bifurcated body delivery system on patient's abdomen under fluoroscopy to determine the orientation of the contralateral limb. The side arm of the hemostatic valve may serve as an external reference to the contralateral limb radiopaque marker.

NOTE: Distal bifurcated body delivery system will not pass through the sheath used to deliver the proximal body.

NOTE: The proximal body delivery sheath must be removed prior to insertion of the distal bifurcated body delivery system.

4. Insert Distal Bifurcated Body delivery system over the wire, into the femoral artery with attention to sidearm reference.

CAUTION: Maintain wire guide position during delivery system insertion.

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

- 5. Advance delivery system until the contralateral limb is positioned in suitable orientation above and anterior to the origin of the contralateral iliac. (Figure 25) If the contralateral limb radiopaque marker is not properly aligned, rotate the entire system until it is.
- 6. Repeat angiogram to verify:
- The degree of overlap with proximal body (no less than 2 stents)
- The position of the contralateral limb
- The position of the ipsilateral iliac limb with respect to the common iliac bifurcation.

Reposition distal bifurcated body as required.

CAUTION: When introducing distal bifurcated body, observe proximal body closely to avoid any disruption to its position.

NOTE: Ensure the Captor Hemostatic Valve on the introducer sheath is turned to the open position. **(Figure 26)**

7. Stabilize the grey positioner (the shaft of the delivery system) while withdrawing the sheath. Deploy the first two (2) covered stents by withdrawing the sheath while monitoring device location. Proceed with deployment until contralateral limb is fully deployed. (Figure 27)

11.4.9 Contralateral Iliac Wire Guide Placement

- Manipulate the wire guide from the contralateral side and into the contralateral limb and into the Distal Bifurcated Body. (Figure 28) AP and oblique fluoroscopic views can aid in verification of device cannulation.
- Advance the angiographic catheter into the body of the graft to the level of the overlap between the proximal and distal graft components. Perform angiography to confirm correct position inside the Distal Bifurcated Body. Advance the catheter to where the proximal end of the Distal Bifurcated Body is attached to the introducer.

11.4.10 Distal Bifurcated Body Deployment

- Perform angiography to confirm proper position of the ipsilateral iliac leg with respect to the internal iliac (hypogastric) artery. Adjust position if necessary.
- 2. Withdraw sheath until the iliac leg is fully deployed
- Remove the safety lock from the black trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Figure 29)
- Advance the contralateral catheter to above the level of the proximal graft, and pass a supportive wire guide (AUS or LES) through it to the level of the arch. Remove the catheter.

11.4.11 Iliac Leg (Contralateral) Placement

Refer to the Instructions for Use enclosed in the device packaging for the iliac leg graft. Once placement of the contralateral iliac leg is complete, continue with deployment of the distal bifurcated body. (Section 11.4.12)

11.4.12 Distal Bifurcated Body Deployment (Continued)

- Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the white trigger-wire release mechanism off the handle and then remove via its slot over the device inner cannula. (Figure 30)
- 2. Under fluoroscopy and after verification of iliac leg graft position, withdraw grey positioner with secured inner cannula.
- 3. Re-check the position of the wire guides. Leave sheath and wire guide in place.
- Close the Captor Hemostatic Valve on the introducer sheath by turning in a clockwise direction until hemostasis is achieved. (Figure 31)

11.4.13 Molding Balloon Insertion

- Prepare Molding balloon as follows:
- Flush wire lumen with heparinized saline.
- Remove all air from balloon.
- 2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning it counter-clockwise.
- Advance the Molding balloon over the wire guide and through the Captor Hemostatic Valve of the distal bifurcated body introduction system to level of renal arteries. Maintain proper sheath positioning

 $\begin{tabular}{ll} \textbf{NOTE:} Captor Hemostatic Valve may be utilized to assist with hemostasis by turning in a clockwise rotation to the "close" position. \end{tabular}$

NOTE: Captor Hemostatic Valve should always be in the "open" position when repositioning of molding balloon.

4. Expand the molding balloon with diluted contrast media (as directed by the manufacturer) in the area of the suprarenal stent and the infrarenal neck, starting proximally and working in the distal direction. (Figure 32)

CAUTION: Prior to molding in the vicinity of any Fenestration stent(s) confirm that the aortic section of the stent has been flared. CAUTION: Confirm complete deflation of balloon prior to repositioning.

5. Withdraw the Molding balloon to the ipsilateral limb distal fixation site and expand.

${\bf CAUTION: Do \ not \ inflate \ balloon \ in \ the \ vessel \ outside \ of \ graft, \ as \ doing \ so \ can \ result \ in \ vessel \ damage \ (e.g., \ rupture). }$

 Deflate and remove molding balloon. Transfer the molding balloon onto the contralateral wire guide and into the contralateral iliac leg introduction system. Advance molding balloon to the contralateral limb overlap and expand.

${\bf CAUTION: Confirm\ complete\ deflation\ of\ balloon\ prior\ to\ repositioning. }$

7. Withdraw the molding balloon to the contralateral iliac leg/vessel distal fixation and expand. (Figure 32)

CAUTION: Do not inflate balloon in the vessel outside of graft, as doing so can result in vessel damage (e.g., rupture).

- 8. Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
- Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

Final Angiogram

- Position angiographic catheter just above the level of the branch vessel(s) accommodated by a fenestration/scallop. Perform angiography to verify branch vessel (e.g., renal artery, superior mesenteric artery) patency and that there are no endoleaks. Verify patency of internal iliac arteries.
- Confirm there are no endoleaks or kinks and verify position of proximal gold radiopaque markers. Remove the sheaths, wires and catheters.

NOTE: If endoleaks or other problems are observed, refer to Section 1.6, Ancillary Components.

3. Repair vessels and close in standard surgical fashion.

12 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

12.1 General

All patients should be advised that endovascular treatment with this device requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) may require additional follow-up. Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be informed that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of abdominal aneurysms with this device. Physicians should evaluate patients on an individual basis and prescribe follow-up relative to the needs and circumstances of each individual patient.

Please refer to **Table 12.1** for recommended imaging and post-operative follow-up requirements. Precise planning and sizing of the Zenith Fenestrated AAA Endovascular Graft requires high resolution CT preprocedure in order to obtain accurate measurements for determining graft size (diameter and length) and fenestration/scallop location. Regular follow-up imaging is important to monitor the performance of the device, allowing for timely reintervention if necessary. The imaging recommended at follow-up (CT/X-ray) is the same as for a non-fenestrated device and is intended to similarly provide for an assessment of device integrity, endoleak, change in aneurysm size, and device position (migration, component overlap). Following placement of a fenestrated graft, it is also important to evaluate the patency of vessels accommodated by fenestrations, which can also be supported by high resolution CT imaging. Duplex ultrasound may also be a useful screening tool in assessing the patency of vessels accommodated by a fenestration, provided the results are interpreted using appropriate criteria.

Table 12.1 Recommended Imaging and Post-Operative Follow-up Schedule

	Pre-procedure	Procedure	30 Day	6 Month (optional)	12 Month ³
Clinical exam	X		Х	Χ	Х
СТ	X		Χ¹	X ¹	Χ¹
Device x-ray			Х	Χ	Х
Angiography	X ²	X			
Renal Duplex Ultrasound	X		Х	Х	Х
Blood Tests (Serum Creatinine, BUN)	X		Х	X	Х

¹ Duplex ultrasound along with a non-contrast CT may be used to assess the aneurysm for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan.

12.2 Angiography

Angiographic imaging is recommended during the procedure to evaluate anatomy and facilitate device placement. In addition, selective catheterization of the visceral vessels targeted by a fenestration is recommended. At the completion of the procedure, patency of the following arteries: aorta, celiac, superior mesenteric, right and left renals and/or accessory renals, and right and left internal external iliac arteries should be evaluated. Pre-procedural angiography may be required at the discretion of the implanting physician or film reviewer.

12.3 Computed Tomography (CT)

High resolution CT imaging is recommended pre-procedure in order to obtain accurate measurements for determining graft size (diameter and length) and fenestration/scallop location. CT images are also recommended at 30 days post-procedure, optionally at 6 months post-procedure and yearly thereafter to provide for an assessment of endoleak, change in aneurysm size, and device position (migration, component overlap). Following placement of a fenestrated graft, it is also important to evaluate the patency of vessels accommodated by fenestrations, which can also be supported by high resolution CT imaging.

Table 12.3.1 lists the general scan parameters for evaluation of patients with the Zenith Fenestrated AAA Endovascular Graft.

Table 12.3.1 Recommended CT imaging parameters

Scan Parameters	Recommendations for Optimal Imaging				
Acceptable machines	Spiral CT or high performance MDCT capable of >40 seconds				
Scan Parameters	Optimize the technique for body habitus and slice thickness				
Superior Extent	Above diaphragm				
Inferior Extent	Proximal femur				
Slice Thickness	≤1.0 mm				
Slice Spacing	At least 50% overlap				
Field of View (FOV)	Adjust for body habitus (include all anatomy / soft tissue)				
IV Contrast	100-200 cc, tailored to the needs or limitations of individual patient				
Reconstruction Algorithm	1 mm throughout				

12.4 Device X-Ray

X-rays are recommended at 30 days post-procedure, optionally at 6 months post-procedure and yearly thereafter as an adjunct to CT to assess device integrity. It is important to ensure that the entire device is captured on the images for device assessment. Recommended imaging parameters include:

- The patient be in the supine position
- 40 inch focal film distance (FFD)
- \bullet Obtain 4 views: AP, LAT, 30° RPO, and 30° LPO centered on umbilicus
- In order to properly penetrate and expose the abdomen utilize lumbar spine technique, center photo cell or manual technique.

12.5 Renal Duplex Ultrasound

Duplex ultrasound may be used as a screening tool for assessing the patency of vessels to be accommodated by a fenestration both pre- and post-procedure. Duplex ultrasound can be an important adjunct to noncontrast enhanced CT in patients with renal failure or who are otherwise contraindicated to receive iodinated contrast. The following information should be included in the evaluation if duplex ultrasound is performed:

- Transverse and longitudinal imaging from the level of the proximal abdominal aorta - demonstrating celiac, mesenteric and renal arteries - to the iliac bifurcations to verify if endoleaks are present and vessels are patent utilizing color flow and color power Doppler (if accessible)
- $\bullet \ Spectral \ analysis \ confirmation \ of \ any \ suspected \ endoleaks$
- Transverse and longitudinal imaging of the maximum aneurysm diameter.

12.6 MRI Safety and Compatibility

For MRI Safety and Compatibility refer to Section 4.5.

12.7 Supplemental Imaging

Additional radiological imaging may be necessary to further evaluate the endovascular graft in situ based on findings revealed by previous imaging assessments. The following recommendations may be considered.

- If there is evidence of poor or irregular position of the endovascular graft, severe angulation, kinking, or migration of the endovascular graft on X-rays, a spiral CT should be performed to assess aneurysm size and the presence or absence of an endoleak.
- If a new endoleak or increase in AAA size is observed by spiral CT, adjunctive studies such as 3-D reconstruction or angiographic assessment of the endovascular graft and native vasculature may be helpful in further evaluating any changes of the endovascular graft or aneurysm.
- Spiral CT without contrast or Renal Duplex Ultrasound may be considered in select patients who cannot tolerate contrast media or who have renal function impairment. For centers with appropriate expertise, gadolinium or CO₂ angiography may be considered in patients with renal function impairment requiring angiographic assessment.

13 TRIGGER-WIRE RELEASE TROUBLESHOOTING.

CAUTION: The following steps should be performed only if unable to remove the proximal trigger-wire as described in Section 11.4.5, Proximal Body Placement (13).

NOTE: Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

13.1 Alternate Proximal Body Deployment

- Cut the exposed suprarenal stent trigger-wire between the white and black release mechanisms (Figure 33) and remove the black release mechanism from the handle.
- 2. Remove the safety lock from the white (distal) trigger-wire release mechanism.
- 3. Withdraw the white trigger-wire release mechanism and then remove via its slot over the device inner cannula.

NOTE: This will detach the distal end of the graft from the grey positioner

4. Using locking forceps, clamp and secure the cut end of the suprarenal stent trigger-wire. **(Figure 34)**

 $^{{}^{2}\}operatorname{\textit{Pre-procedure angiography may be required at discretion of implanting physician or film reviewer.}$

³ Annually thereafter.

5. Loosen the pin vise and, while maintaining inner cannula and trigger-wire position, advance the grey positioner and sheath into the graft until the tip of the grey positioner is approximately 2 cm from the gold markers on the proximal edge of the proximal body (Figure 35). The advanced grey positioner provides added support to the inner cannula.

NOTE: Take care when advancing the grey positioner as there will be sheaths and wire guides positioned in target vessels. Ensure that the tip of the grey positioner is not advanced into the top cap.

NOTE: Maintain gentle tension on the suprarenal stent trigger-wire to remove any slack in the wire as the grey positioner and sheath are being advanced.

- 6. Lock the pin vise. Confirm that the suprarenal stent trigger-wire is secured by the forceps.
- Stabilize the grey positioner and slowly advance the sheath until the sheath tip is 2 mm from the gold markers. (Figure 36)

NOTE: Take care when advancing the sheath as there will be other sheaths and wire guides positioned in target vessels. Take care not to advance the graft itself during sheath advancement.

 Stabilize the sheath and slightly retract the grey positioner with inner cannula to move the top cap down over the suprarenal stent (Figure 37).

NOTE: Avoid compressing the body of the graft.

- Ensure the fenestrations are positioned correctly. Test the resistance on the trigger wire and make minor adjustments as necessary to reduce trigger wire resistance (Section 13.1, 8).
- 10. Remove the suprarenal stent trigger-wire.
- 11. Withdraw the sheath until the tapered tip of the grey positioner is exposed.
- 12. If there are multiple fenestrations, withdraw all but one of the fenestration guiding catheters. It is the physician's preference as to which guiding catheter to withdraw. It is recommended that the decision be made based on the ease of cannulation of the fenestrations and their respective vessels. (Figure 38)

NOTE: Leave the wire guides in place when removing the guiding catheter(s)

- 12. (a) If there is a single fenestration, then the proximal body needs to be cannulated with a suitable wire guide from the same side as the in-situ guide catheter. This is to allow placement of a molding balloon within the proximal body.
- Advance a molding balloon along the now available wire guide into the proximal body and position it just superior to the distal-most end of the graft.

NOTE: When using the molding Balloon, ensure the appropriate sheath is used – either through the in-situ 20 Fr contralateral sheath or, if direct punctures have been used, through a 14 Fr introducer sheath. This will ensure the safe retrieval of the molding balloon.

- 14. Inflate the balloon to the full diameter of the graft. (Figure 39)
- Loosen the pin vise (Figure 40). Control the position of the graft by stabilising the grey positioner and balloon catheter.

CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

16. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed (Figures 41 and 42). Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (Figure 43) to avoid contact with the deployed suprarenal stent.

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

- 17. If there are multiple fenestrations (Section 13.1, 12), deflate the balloon then withdraw it leaving the wire guide in place
- 17. (a) If there is a single fenestration (Section 13.1, 12a), then the molding balloon and wire guide can be safely removed.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide, which remain in the target vessel

 Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (Figure 44)

NOTE: Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.

13.2 Docking of Top Cap

- 1. Loosen the pin vise. (Figure 45)
- 2. Secure sheath and inner cannula to avoid any movement of these components.
- 3. Advance the grey positioner over the inner cannula until it docks with the top cap. (Figures 46, 47 and 48)

 $\begin{tabular}{ll} \textbf{NOTE:} If resistance occurs, slightly rotate grey positioner and continue to gently advance. \end{tabular}$

- Retighten the pin vise and withdraw the entire top cap and grey
 positioner through the graft and through the sheath by pulling on the
 inner cannula. (Figure 49) Leave the sheath and wire guide in place.
- 5. Close the Captor Hemostatic Valve on the Flexor introducer sheath by turning it in a clockwise direction until it stops.
- 6. Return to Section 11.4.7, Fenestration Stent Placement and Deployment.

14 SUPRARENAL STENT DEPLOYMENT TROUBLESHOOTING.

CAUTION: The following steps should be performed only if unable to deploy the suprarenal stent by advancing the top cap as described in Section 11.4.5, Proximal Body Placement (15).

 $\begin{tabular}{ll} \textbf{NOTE:} Technical assistance from a Cook product specialist may be obtained by contacting your local Cook representative. \end{tabular}$

14.1 Proximal Body Placement with distal attachment

If the suprarenal stent cannot be fully deployed by advancing the top cap inner cannula, perform the following steps under fluoroscopy.

Tighten the pin-vice. If there are multiple fenestrations, withdraw all but one of the fenestration guiding catheters. It is the physician's preference

- as to which guiding catheter to withdraw. It is recommended that the decision be made based on the ease of cannulation of the fenestrations and their respective vessels. (Figure 50)
- (a) If there is a single fenestration, then the proximal body needs to be cannulated with a suitable wire guide from the same side as the insitu guide catheter. This is to allow placement of a molding balloon within the proximal body.

NOTE: Leave the wire guides in place when removing the guiding catheter(s).

 Advance a molding balloon along the now available wire guide into the proximal body and position it just superior to the distal-most end of the graft.

NOTE: When using the molding Balloon, ensure the appropriate sheath is used – either through the in-situ 20 Fr Contralateral sheath or, if direct punctures have been used, through a 14 Fr introducer sheath. This will ensure the safe retrieval of the molding balloon.

- 3. To add support to the inner cannula, inflate the balloon to the full diameter of the graft. (Figure 51)
- 4. Loosen the pin vise. (Figure 52)
- 5. Control the position of the graft by stabilising the grey positioner and balloon catheter.

CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

6. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed. (Figures 53 and 54) Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (Figure 55) to avoid contact with the deployed suprarenal stent.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

If the suprarenal stent is fully deployed:

7. a) If there are multiple fenestrations (Section 14.1, 1), deflate the balloon then withdraw it leaving the wire guide in place. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (Figure 56) b) If there is a single fenestration (Section 14.1, 1a), then the molding balloon and wire guide can be safely removed.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

8. Return to Section 11.4.6, Docking of Top Cap.

If still unable to fully deploy the suprarenal stent, perform the following steps:

14.2 Proximal Body Placement without distal attachment

- 1. Tighten the pin vise and deflate the balloon, while maintaining balloon position.
- Remove the safety lock from the white trigger-wire release mechanism. Withdraw and remove the trigger-wire to detach the distal end of the endovascular graft from the delivery system by sliding the trigger-wire release mechanism off the handle and remove via its slot over the device inner cannula. (Figure 57)
- 3. Loosen the pin vise (Figure 58) and, while maintaining inner cannula position, advance the grey positioner and sheath into the graft until the tip of the grey positioner is approximately 2 cm from the gold markers on the proximal edge of the proximal body (Figure 59). The advanced grey positioner provides added support to the inner cannula.

NOTE: Take care when advancing the grey positioner as there will be sheaths and wire guides positioned in target vessels. Ensure that the tip of the grey positioner is not advanced into the top cap.

- 4. Lock the pin vise.
- 5. Verify position of the gold markers and ensure the fenestrations are positioned correctly.
- 6. To add support to the inner cannula, inflate the balloon to the full diameter of the graft. (**Figure 60**)
- 7. Loosen the pin vise (Figure 61). Control the position of the graft by stabilising the grey positioner and balloon catheter.

CAUTION: Before deployment of the suprarenal stent, verify that the position of the access wire extends just distal to the aortic arch. Ensure that the dilator tip will not extend beyond the end of the access wire guide during advancement, and if required re-position the access wire guide into the aortic arch to accommodate.

8. Deploy the suprarenal stent by advancing the top cap inner cannula 1 to 2 mm at a time while controlling the position of the proximal body until the top stent is fully deployed (Figures 62 and 63). Advance the top cap cannula an additional 1 to 2 cm and then retighten the pin vise (Figure 64) to avoid contact with the deployed suprarenal stent.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

- 9. a) If there are multiple fenestrations (Section 14.1, 1), deflate the balloon then withdraw it leaving the wire guide in place. Advance the access sheath and aligning stent, which was removed to facilitate the molding balloon, back over the wire guide, through the fenestration and into the respective vessel. (Figure 65)
- b) If there is a single fenestration (Section 14.1, 1a), then the molding balloon and wire quide can be safely removed.

NOTE: Care should be taken during removal to not disturb the guide sheath and wire guide(s), which remain in the target vessel(s).

WARNING: The Zenith Fenestrated AAA Endovascular Graft incorporates a suprarenal stent with fixation barbs. Exercise extreme caution when manipulating interventional devices in the region of the suprarenal stent.

10. Return to Section 11.4.6, Docking of Top Cap.

NOTE: Check to make sure that all trigger-wires are removed prior to withdrawal of the delivery system.



MR Conditional



Quantity per box



Keep dry



Keep away from sunlight



MANUFACTURER
WILLIAM A. COOK AUSTRALIA PTY. LTD.
95 Brandl Street, Brisbane Technology Park
Eight Mile Plains, Brisbane, QLD 4113 AUSTRALIA

www.cookmedical.com © COOK 2012

COOK INCORPORATED 750 Daniels Way Bloomington, IN 47404 U.S.A.