EN 16	Zenith® Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot™ Introduction System Instructions for Use
CS 21	Endovaskulární graft Zenith® Branch s iliakální bifurkací se zaváděcím systémem H&L-B One-Shot™ Návod k použití
DA 26	Zenith® Branch endovaskulær protese-iliaca bifurkation med H&L-B One-Shot™ indføringssystem Brugsanvisning
DE 31	Endovaskuläre iliakale Bifurkationsprothese Zenith® Branch mit dem Einführsystem H&L-B One-Shot™ Gebrauchsanweisung
EL 37	Ενδοαγγειακό μόσχευμα-διχασμός λαγονίου αρτηρίας Zenith® Branch με το σύστημα εισαγωγής H&L-B One-Shot™ Οδηγίες χρήσης
ES 43	Bifurcación endoprótesis vascular-iliaca Zenith® Branch con el sistema de introducción H&L-B One-Shot™ Instrucciones de uso
FI 48	H&L-B One-Shot™ -sisäänvientijärjestelmällä varustettu Zenith® Branch -endovaskulaarinen iliakaalinen bifurkaatiografti Käyttöohjeet
FR 53	Bifurcation iliaque de l'endoprothèse vasculaire Zenith® Branch à système d'introduction H&L-B One-Shot™ Mode d'emploi
IT 58	Biforcazione iliaca per protesi endovascolare Zenith® Branch con sistema di introduzione H&L-B One-Shot™ Istruzioni per l'uso
NL 63	Zenith® Branch endovasculair implantaat-iliacale bifurcatie met het H&L-B One-Shot™ introductiesysteem Gebruiksaanwijzing
PL 68	Rozwidlony biodrowy stent-graft wewnętrzznacyniowy Zenith® Branch z systemem wprowadzającym H&L-B One-Shot™ Instrukcja użycia
PT 73	Prótese endovascular Zenith® Branch com bifurcação iliaca e o sistema de introdução H&L-B One-Shot™ Instruções de utilização
SK 78	Endovaskulárny štep Zenith® Branch s iliackým rozdvojením so zavádzacím systémom H&L-B One-Shot™ Návod na použitie
SV 83	Zenith® Branch endovaskulärt iliabifurkationsgraft med H&L-B One-Shot™ införingssystem Bruksanvisning

ENGLISH

TABLE OF CONTENTS

Illustrations	11-15
1 DEVICE DESCRIPTION	16
1.1 Zenith Branch Endovascular Graft-Iliac Bifurcation	16
1.2 Zenith Branch Endovascular Graft-Iliac Bifurcation Delivery System	16
1.3 Additional Components	16
2 INDICATIONS FOR USE	16
3 CONTRAINDICATIONS	16
4 WARNINGS AND PRECAUTIONS	16
4.1 General	16
4.2 Patient Selection, Treatment and Follow-Up	16
4.3 Implant Procedure	16
Molding Balloon Use	16
4.4 MRI Safety and Compatibility	17
5 POTENTIAL ADVERSE EVENTS	17
6 PATIENT SELECTION AND TREATMENT	17
6.1 Individualization of Treatment	17
7 PATIENT COUNSELING INFORMATION	17
8 HOW SUPPLIED	17
9 CLINICAL USE INFORMATION	18
9.1 Physician Training	18
9.2 Inspection Prior to Use	18
9.3 Materials Required	18
9.4 Materials Recommended	18
9.5 Device Diameter Sizing Guidelines	18
10 DIRECTIONS FOR USE	18
General Use Information	18
General Deployment Information	18
Pre-Implant Determinants	18
Patient Preparation	18
10.1 Zenith Branch Endovascular Graft-Iliac Bifurcation Preparation/Flush	18
10.2 Vascular Access and Angiography	18
10.3 Zenith Branch Endovascular Graft-Iliac Bifurcation Placement	19
10.4 Through-and-Through Wire Guide Placement	19
10.5 Proximal Zenith Branch Endovascular Graft-Iliac Bifurcation Deployment	19
10.6 Up-and-Over Sheath Placement	19
10.7 Cannulation of Internal Iliac Artery	19
10.8 Balloon-Expandable Covered Peripheral Bridging Stent Placement	19
10.9 Distal Zenith Branch Endovascular Graft-Iliac Bifurcation Deployment	19
10.10 Balloon-Expandable Covered Peripheral Stent Deployment	19
10.11 Zenith AAA Graft and Iliac Leg Graft Deployment	19
10.12 Molding Balloon Insertion	20
10.13 Final Angiogram	20
11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP	20

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
Figure 20

ENGLISH

ZENITH® BRANCH ENDOVASCULAR GRAFT-ILIAC BIFURCATION WITH THE H&L-B ONE-SHOT™ INTRODUCTION SYSTEM

Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

For information on other Zenith components, please refer to the appropriate Instructions for Use.

1 DEVICE DESCRIPTION

1.1 Zenith Branch Endovascular Graft-Iliac Bifurcation

The Zenith Branch Endovascular Graft-Iliac Bifurcation is a bifurcated branch vessel graft with openings to connect the common iliac, sidebranch, and external iliac segments. (Figure 1)

The graft is constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel and nitinol Cook Z®-stents with braided polyester and monofilament polypropylene suture. The graft is fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Nitinol rings positioned at the proximal end of the graft and within the sidebranch help maintain lumen patency during access. Additionally, the Cook Z-stents provide the necessary seal of the graft to the vessel wall.

To facilitate fluoroscopic visualization of the stent graft, four radiopaque gold markers are positioned along the internal iliac side of the proximal portion of the graft to indicate the position of the sidebranch segment.

Additional components such as Zenith Iliac Leg grafts (e.g., TFLE, ZSLE, ZISL) and balloon-expandable covered peripheral stents are required. Each individual device has its own separate delivery system. Please refer to the relevant Instructions for Use.

1.2 Zenith Branch Endovascular Graft-Iliac Bifurcation Delivery System

The Zenith Branch Endovascular Graft-Iliac Bifurcation is shipped pre-loaded onto the H&L-B One-Shot Introduction System. (Figure 2) It has a sequential deployment method with built-in features to provide continuous control of the endovascular graft throughout the deployment procedure.

The proximal end of the graft is attached to the delivery system by two nitinol trigger wires. The distal end of the graft is also attached to the delivery system and held by an independent stainless steel trigger-wire. The H&L-B One-Shot Introduction System enables precise positioning and allows readjustment of the final graft position before full deployment. The delivery system uses a 20 French (6.7 mm I.D.) H&L-B One-Shot Introduction System that includes both a pre-loaded catheter, used to facilitate cannulation of the side branch, and wire guide, used to ensure that the lumen of the catheter is preserved during loading and shipping. All systems are compatible with a 0.035 inch (0.89 mm) wire guide. For added hemostasis, the Captor Hemostatic Valve can be opened or closed for the introduction and/or the removal of ancillary devices into and out of the sheath.

The Zenith Branch Endovascular Graft-Iliac Bifurcation delivery system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta. In addition, the pre-loaded catheter is curved in order to enhance the positioning capability of the wire guide during the snaring process.

1.3 Additional Components

Additional components (iliac leg graft, balloon-expandable covered peripheral stent) are also required.

Zenith Iliac Leg grafts are available (e.g., TFLE, ZSLE, ZISL) and are constructed from the same materials used in the construction of the Zenith Branch Endovascular Graft-Iliac Bifurcation. For all of the Zenith iliac leg grafts and delivery systems, refer to the respective Instructions for Use for deployment instructions. The common iliac section of the Zenith Branch Endovascular Graft-Iliac Bifurcation, positioned in the common iliac artery, will be connected to the short (contralateral) limb of a Zenith bifurcated graft (e.g., Zenith AAA Endovascular Graft) by a standard iliac leg of suitable length and with a distal diameter of 16 mm (e.g., TFLE-16-XX, ZSLE-16-XX-ZT, ZISL-16-XX).

A balloon-expandable covered peripheral stent will be deployed within the sidebranch segment and will extend into the internal iliac artery. Refer to the Instructions for Use of this device for deployment instructions.

2 INDICATIONS FOR USE

The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System is indicated for the endovascular treatment of patients with an aortoiliac or iliac aneurysm, an insufficient distal sealing site within the common iliac artery, and having morphology suitable for endovascular repair, including:

- Adequate iliac/femoral access compatible with a 20 French (7.7 mm O.D.) introduction system,
- Non-aneurysmal external iliac artery fixation segment distal to the aneurysm:
 - With a length of at least 20 mm,
 - With a diameter measured outer wall to outer wall of no greater than 11 mm and no less than 8 mm.
- Non-aneurysmal internal iliac artery segment distal to the aneurysm:
 - With a length of at least 10 mm (with 20-30 mm being preferred)
 - With a diameter acceptable for proper sealing.

3 CONTRAINDICATIONS

There are no known contraindications for these devices.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device.
- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thickness > 3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenoses from CT.
- The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up**.

- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft.
- At a minimum, annual imaging is required, including:
 - 1) abdominal radiographs to examine device integrity (separation between components or stent fracture) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease.
- If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.
- The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System is not recommended in patients unable to undergo, or who will not be compliant with the necessary preoperative and post-operative imaging and implantation studies as described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up**.
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the haemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the grey positioner. After the grey positioner has been removed, if blood loss is excessive, consider placing an uninflated moulding balloon or an introduction system dilator within the valve, restricting flow.

4.2 Patient Selection, Treatment and Follow-Up

- Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease, and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 20 French (7.7 mm O.D.) vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous, or thrombus lined may preclude placement of the endovascular graft and/or may increase the risk of embolization.
- Key anatomic elements that may affect successful exclusion of the aneurysm include tortuosity of any or all of the vessels involved, undersized or oversized iliac arteries, circumferential thrombus, aneurysm of the internal and/or external iliac artery and/or calcification of the arterial implantation sites. Irregular calcification and/or plaque may compromise the fixation and sealing of the implantation sites and the ability to advance the introducer systems.
- The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System is not recommended in patients who cannot tolerate contrast agents necessary for intra-operative and post-operative follow-up imaging.
- The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System is not recommended in patients exceeding weight and/or size limits, which compromise or prevent the necessary imaging requirements.
- The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System is not recommended in patients with known sensitivities or allergies to stainless steel, polyester, nitinol (nickel, titanium), solder (tin, silver), polypropylene or gold.
- Patients with a systemic infection may be at increased risk of endovascular graft infection.
- Inability to maintain patency of at least one internal iliac artery may increase the risk of pelvic/bowel ischemia.

4.3 Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulation should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with 4x4 gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Maintain wire guide position during delivery system insertion.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Branch Endovascular Graft-Iliac Bifurcation.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure post-operatively. Care should be taken to limit the amount of contrast media used during the procedure.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
- Inaccurate placement and/or incomplete sealing of the Zenith Branch Endovascular Graft-Iliac Bifurcation within the vessel or with additional components may result in increased risk of endoleak, migration or inadvertent occlusion of internal iliac artery.
- Inadequate fixation of the Zenith Branch Endovascular Graft-Iliac Bifurcation may result in increased risk of migration. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance. Vessel or catheter damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Take care during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus which can cause distal embolization.
- The balloon-expandable covered peripheral stent must be of adequate diameter to seal with both the sidebranch segment and the internal iliac artery, and of sufficient length to overlap 10-14 mm with the sidebranch segment and 10 mm (with 20-30 mm being preferred) with the distal fixation site in the internal iliac artery.
- Care should be taken not to damage any previously placed graft or disturb its position in the event reinstrumentation of the graft is necessary.

Molding Balloon Use

- Confirm complete deflation of balloon prior to repositioning.
- Do not inflate the balloon in vessel outside of graft.
- For added hemostasis, the Captor Hemostatic Valve can be opened or closed to accommodate the insertion and subsequent withdrawal of a molding balloon.

4.4 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the Zenith Branch Endovascular Graft-Iliac Bifurcation used with a balloon expandable covered stent is MR Conditional according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 720 gauss/cm (7.2 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Zenith Branch Endovascular Graft-Iliac Bifurcation used with a balloon expandable covered stent is not expected to result in a temperature rise of more than 2.4°C after 15 minutes of continuous scanning.

The image artifact extends approximately 122 mm from the Zenith Branch Endovascular Graft-Iliac Bifurcation used with a balloon expandable covered stent as found during non-clinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. The image artifact obscures the device lumen.

NOTE: For the Zenith Branch Endovascular Graft-Iliac Bifurcation, the clinical benefit of an MRI scan should be balanced with the potential risk of the procedure.

5 POTENTIAL ADVERSE EVENTS

Adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g. buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprostheses: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow and corrosion
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)

6 PATIENT SELECTION AND TREATMENT

(See section 4, Warnings and Precautions)

6.1 Individualization of Treatment

Cook recommends that the Zenith Branch Endovascular Graft-Iliac Bifurcation component diameters be selected as described in Table 8.1.

The Zenith Branch Endovascular Graft-Iliac Bifurcation and delivery systems are available in the following lengths and diameters:

Table 8.1 Branch Grafts-Iliac Bifurcation							
Reorder Number ¹	Proximal Diameter	Distal Diameter	Sidebranch Diameter ²	Introduction Sheath French Size	Introduction Sheath (I.D. / O.D.)	Iliac Segment Length Common	Total Graft Length
ZBIS-10-45-41	12 mm	10 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	45 mm	86 mm
ZBIS-10-61-41	12 mm	10 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	61 mm	102 mm
ZBIS-10-45-58	12 mm	10 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	45 mm	103 mm
ZBIS-10-61-58	12 mm	10 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	61 mm	119 mm
ZBIS-12-45-41	12 mm	12 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	45 mm	86 mm
ZBIS-12-61-41	12 mm	12 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	61 mm	102 mm
ZBIS-12-45-58	12 mm	12 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	45 mm	103 mm
ZBIS-12-61-58	12 mm	12 mm	8 mm	20 Fr	(6.7 mm / 7.7 mm)	61 mm	119 mm

¹ZBIS-XX-YY-ZZ is the Iliac Bifurcation graft where XX is distal diameter, YY is the Common Iliac Segment Length (length from the proximal graft edge to the tip of the sidebranch), and ZZ is the External Iliac Segment Length (length from the tip of the sidebranch to distal edge of the graft).

²Sidebranch contains a 6 mm diameter nitinol z-stent, but should be expanded to 8 mm during balloon-expandable covered peripheral stent deployment.

The Zenith iliac legs and delivery systems for use with the Zenith Branch Endovascular Graft-Iliac Bifurcation are available in the following lengths and diameters:

Table 8.2 Iliac Legs			
Reorder Number ¹	Iliac Leg Diameter	Introduction Sheath French Size	Iliac Leg Working Length ²
TFLE-16-YY	16 mm	14 Fr	(4.7 mm / 5.4 mm)
ZSLE-16-YY-ZT	16 mm	14 Fr	(4.7 mm / 5.4 mm)
ZISL-16-YY	16 mm	12 Fr	(4.0 mm / 4.7 mm)

¹Example: TFLE-XX-YY is the iliac leg where XX is the diameter and YY is the iliac leg working length.

²Overall leg length = working length + 22 mm docking stent (for TFLE and ZSLE) and 29.1 mm docking stents (for ZISL).

The length of the Zenith Branch Endovascular Graft-Iliac Bifurcation is generally chosen to extend from the proximal portion of the common iliac artery to the external iliac artery. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. The potential risks previously described in **Section 5, Potential Adverse Events** should be carefully considered for each patient before use of the Zenith Branch Endovascular Graft-Iliac Bifurcation. Additional considerations for patient selection include but are not limited to:

- Patient's age and life expectancy.
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity).
- Patient's suitability for open surgical repair.
- Patient's anatomical suitability for endovascular repair.
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith Branch Endovascular Graft-Iliac Bifurcation.
- Ability to tolerate general, regional or local anesthesia.
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery with vascular access techniques and accessories of the delivery profile of a 12 French (4.7 mm O.D.) to 20 French (7.7 mm O.D.) vascular introducer sheath.
- Non-aneurysmal external iliac artery fixation segment distal to the aneurysm:
 - With a length of at least 20 mm,
 - With a diameter measured outer wall to outer wall of no greater than 11 mm and no less than 8 mm.
- Non-aneurysmal internal iliac artery segment distal to the aneurysm:
 - With a length of at least 10 mm (with 20-30 mm being preferred)
 - With a diameter acceptable for proper sealing.
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.

The final treatment decision is at the discretion of the physician and patient.

7 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure including:

- Risks and differences between endovascular repair and surgical repair.
- Potential advantages of traditional open surgical repair.
- Potential advantages of endovascular repair.
- Potential advantages of alternative endovascular procedures.
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms, or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, Imaging Guidelines and Post-Operative Follow-up**.

8 HOW SUPPLIED

The Zenith Branch Endovascular Graft-Iliac Bifurcation is supplied sterile (100% ethylene oxide) and pre-loaded in peel-open packages.

This device is intended for single patient use only. Do not reuse, reprocess or resterilize any part of this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to Cook.

Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient. Do not use after the expiration date printed on the label. Store in a cool dry place.

The Zenith Branch Endovascular Graft-Iliac Bifurcation is loaded into a 20 Fr Flexorintroducer sheath. The sheath's surface is treated with a hydrophilic coating that, when activated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with 4x4 gauze pads soaked in saline solution.

9 CLINICAL USE INFORMATION

9.1 Physician Training

CAUTION: Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System should only be used by physicians and teams trained in vascular interventional techniques, general Zenith training, and training in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System are outlined below:

Patient Selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair.
- Knowledge of radiographic image interpretation, device selection and sizing.

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to COOK. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

9.3 Materials Required

(Not included in iliac bifurcated system)

- Zenith Iliac Leg graft (e.g. TFLE-16-(37,54,71,88), ZSLE-16-(39,56,74,90)-ZT, ZISL-16-(42,59,77,93))
- Zenith bifurcated graft (e.g. Zenith AAA Endovascular Graft)
- Zenith AAA Ancillary Kit
- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media

9.5.1 Zenith Branch Endovascular Graft-Iliac Bifurcation Diameter Sizing Guide* (External Iliac Segment)

Intended External Iliac Vessel Diameter ^{1,2}	Branch External Iliac Leg Diameter ³	Introduction Sheath French Size (I.D. / O.D.)
8 mm	10 mm	20 Fr (6.7 mm / 7.7 mm)
9 mm	10 mm	20 Fr (6.7 mm / 7.7 mm)
10 mm	12 mm	20 Fr (6.7 mm / 7.7 mm)
11 mm	12 mm	20 Fr (6.7 mm / 7.7 mm)

¹Maximum diameter along the distal fixation site.

²Round measured iliac diameter to nearest mm.

³Additional considerations may affect choice of diameter.

*All dimensions are nominal.

The diameter of the iliac leg graft (e.g. TFLE, ZSLE, ZISL) used to connect a Zenith AAA graft to the Zenith Branch Endovascular Graft-Iliac Bifurcation will be 16 mm.

The balloon-expandable covered peripheral stent must expand to 8 mm in the sidebranch of the Zenith Branch Endovascular Graft-Iliac Bifurcation, but should be sized to the internal iliac artery according to the manufacturer's instructions.

10 DIRECTIONS FOR USE

Prior to use of the Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System, review this Suggested Instructions for Use booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

General Use Information

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System. The Zenith Branch Endovascular Graft-Iliac Bifurcation with the H&L-B One-Shot Introduction System is compatible with 0.035 inch (0.89 mm) diameter wire guides.

General Deployment Information

The Zenith Branch Endovascular Graft-Iliac Bifurcation is designed to be placed via the common femoral artery on the side of the common iliac artery aneurysm to be treated. A Zenith AAA graft will be introduced on the side opposite the Zenith Branch Endovascular Graft-Iliac Bifurcation.

The suggested order of placement of the multiple components will be:

1. Deployment of the Zenith Branch Endovascular Graft-Iliac Bifurcation, including placement of the balloon-expandable covered peripheral stent (through the iliac artery) into the internal iliac artery.
2. Deployment of the Zenith AAA graft (introduced through the opposite iliac artery).
3. Deployment of a Zenith Iliac Leg graft (e.g. TFLE, ZSLE, ZISL) via the Zenith Branch Endovascular Graft-Iliac Bifurcation access site into the contralateral limb (i.e. short) of the abdominal graft, connecting the lumen of the abdominal graft to the Zenith Branch Endovascular Graft-Iliac Bifurcation.
4. Deployment of a Zenith Iliac Leg graft (e.g. TFLE, ZSLE, ZISL) into the ipsilateral limb (i.e. long) of the Zenith AAA graft (if necessary).

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the delivery system, (i.e., define relevant iliac artery).
2. Angulation of iliac arteries.
3. Diameters of relevant distal iliac vessel(s).
4. Distance from the aortic bifurcation to the hypogastric (internal iliac) artery and attachment site(s).
5. Consider the degree of vascular calcification.

- Syringe
- Heparinized saline solution
- Balloon-expandable covered peripheral stent
- Balloons of appropriate sizes

9.4 Materials Recommended

(Not included in iliac bifurcated system)

The following products are recommended.

Wire Guides

- 0.035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
 - Cook Amplatz Ultra Stiff Wire Guides (AUS2)
- 0.035 inch (0.89 mm) standard wire guide; for example:
 - Cook 0.035 inch (0.89 mm) wire guides
- 0.018 inch (0.46 mm) hydrophilic wire guide
- 0.035 inch (0.89 mm) hydrophilic wire guide, 260 cm; for example:
 - Cook Nimble™ Wire Guides
- 0.035 inch (0.89 mm) exchange wire guide, 260 cm; for example:
 - Cook Roadrunner® PC Wire Guide
 - Cook Amplatz Extra Stiff Wire Guide (AES)
 - Cook Rosen Wire Guide

Molding balloons

- Cook CODA® Balloon Catheter (32 mm)

Introducer sets (7Fr [2.3 mm I.D.], 8Fr [2.7 mm I.D.], 10Fr [3.3 mm I.D.] and 12Fr [4.0 mm I.D.])

- Cook Check-Flo® Introducer Sets
- Cook Flexor Check-Flo Introducers
- Cook Flexor Check-Flo Introducers – Raabe Modification
- Cook Extra Large Check-Flo Introducer Sets
- Cook Flexor Balkin Up & Over® Contralateral Introducers
- Cook Flexor Guiding Sheaths

Catheters

- Sizing catheter; for example:
 - Cook Aorous® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
 - Cook Beacon® Tip Angiographic Catheters
 - Cook Royal Flush® Catheters with Beacon Tip

Entry needles

- Cook single wall entry needles

Snare

- Dotter Intravascular Retriever Set

9.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Patient Preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
3. Expose common femoral arteries as necessary using standard surgical technique.
4. Establish adequate proximal and distal vascular control of femoral vessels, as necessary.

10.1 Zenith Branch Endovascular Graft-Iliac Bifurcation Preparation/Flush

1. Remove black-hubbed shipping stylet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. **(Figure 3)** Push the introduction sheath forward on the dilator tip until the curved catheter is no longer exposed. Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the sideport near the tip of the introduction sheath. **(Figure 4)** Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.

NOTE: Graft flushing solution of heparinized saline is always used.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. **(Figure 5)**

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the distal dilator tip. **(Figure 6)**

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

3. Remove the pre-loaded wire guide from the delivery system.
4. Attach syringe with heparinized saline to the pre-loaded catheter Luer fitting. Flush until fluid exits the distal catheter tip.
5. Insert a suitable hydrophilic nitinol core wire guide, 0.018 inch (0.46 mm) or 0.035 inch (0.89 mm) compatible and 260 cm in length, into the pre-loaded catheter.
6. Soak 4x4 gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator tip liberally.

10.2 Vascular Access and Angiography

1. Puncture the selected common femoral artery(s) using standard technique with an 18 or 19UT gauge arterial needle. Upon vessel entry, insert:
 - Wire guides – standard 0.035 inch (0.89 mm) diameter, 145 cm long, J tip or Bentson Wire Guide
 - Appropriate size sheaths (e.g., 5 French (1.7 mm I.D.) and 8 French (2.7 mm I.D.))
 - Flush catheter (often radiopaque sizing catheters – e.g., centimeter sizing catheter or straight flush catheter)

NOTE: Use an 8 French (2.7 mm I.D.) sheath advanced into the lower abdominal aorta on the contralateral side.

2. Perform angiography to determine the position of the aortic bifurcation and internal iliac artery on the side of Zenith Branch Endovascular Graft-Iliac Bifurcation implantation.

10.3 Zenith Branch Endovascular Graft-Iliac Bifurcation Placement

1. Ensure the delivery system has been flushed with heparinized saline and that all air is removed from the system.
2. Give systemic heparin and check flushing solutions. Flush after each catheter and/or wire guide exchange.

NOTE: Monitor the patient's coagulation status throughout the procedure.

3. On side of Zenith Branch Endovascular Graft-Iliac Bifurcation implantation, replace J wire with a stiff wire guide (AUS or LES) 0.035 inch (0.89 mm), 260 cm long and advance through catheter and up into the upper thoracic aorta. Remove the flush catheter and sheath. Maintain wire guide position.

NOTE: The location and orientation of the internal iliac artery origin as well as the aortic bifurcation should be identified from pre-operative or intraoperative imaging prior to device insertion.

4. Before insertion, check orientation of the sidebranch segment using fluoroscopy to ensure its position is appropriately aligned to the anatomy.

NOTE: The device will be completely covered by the outer sheath at this stage, but this image with the sheath withdrawn demonstrates the four gold radiopaque markers which are aligned with the most lateral aspect of the side branch, and the groove that accommodates the in dwelling catheter. (Figure 7)

NOTE: The second marker from the distal end (2D) is the "Aiming" marker for the most proximal location of the proximal end of the covered stent that extends the side branch into the internal iliac artery. The second marker from the proximal end (2P) is the "Aiming" marker for the most distal location of the distal end of the leg extension that will connect the iliac side branch to the bifurcated stent graft that will be deployed above the aortic bifurcation.

NOTE: Orientation can be confirmed by visual inspection since the groove in the tapered dilator tip (where the pre-loaded wire guide emerges from the system) is aligned with the sidebranch.

5. Insert Zenith Branch Endovascular Graft-Iliac Bifurcation delivery system over the wire guide, into the iliac artery until the tip of the in dwelling catheter, sitting in the groove in the dilator tip beneath the outer sheath, is just above the aortic bifurcation and in rotational alignment with the opposite common iliac artery origin. (Figure 8)

CAUTION: Maintain wire guide position during delivery system insertion.

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

NOTE: Slight advancement of the hydrophilic wire guide through the tip of the catheter may aid in visualization of the catheter tip position.

10.4 Through-and-Through Wire Guide Placement

1. Advance the hydrophilic wire guide and catheter into the aorta. Position the wire guide for snaring by advancing it through the tip of the curved catheter as required.
2. Advance a suitable snare through the access sheath in the opposite iliac artery, snare the tip of the hydrophilic wire guide, and pull it through to form a through-and-through wire guide. (Figure 9)

CAUTION: During this manoeuvre, the hydrophilic through-and-through wire guide should be advanced so that some slack is maintained at the aortic bifurcation.

10.5 Proximal Zenith Branch Endovascular Graft-Iliac Bifurcation Deployment

1. Advance the 8 French (2.7 mm I.D.) sheath in the opposite limb until it is above the level of the aortic bifurcation, and perform angiography through it to determine the location of the internal iliac artery on the side of the device implantation. (Figure 10)
2. Ensure the graft system is oriented such that the distal end of the sidebranch (as indicated by the most distal radiopaque marker) is positioned 10 mm above the origin of the internal iliac artery with rotational orientation such as to allow easy access into the internal iliac artery. If the radiopaque markers are not in the appropriate position, rotate the entire system until the markers are in line with the internal iliac artery orifice and/or advance or withdraw the system as necessary. (Figure 11)

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

3. When the position of the device is satisfactory, stabilize the grey positioner (the shaft of the delivery system) while withdrawing the outer introduction sheath. Deploy the Zenith Branch Endovascular Graft-Iliac Bifurcation until the distal end of the side branch is exposed. (Figure 12)
4. Stop withdrawing sheath.

WARNING: Do not withdraw the sheath any further, or the external iliac component will expand and prevent further positional changes.

10.6 Up-and-Over Sheath Placement

1. Withdraw the opposite iliac access sheath leaving the hydrophilic nitinol core through-and-through wire guide in place.
2. Advance a suitably sized up-and-over sheath over the hydrophilic nitinol core through-and-through wire guide until the tip of the sheath's dilator tip comes in contact with the tip of the pre-loaded catheter. (Figure 13)

NOTE: It may be necessary to apply gentle traction on both ends of the wire guide to aid in advancement.

CAUTION: If the pre-loaded catheter is not advanced over the aortic bifurcation use extra care during wire manipulation to avoid damage to the vasculature.

3. Clamp the hydrophilic through-and-through wire guide at both ends (where it emerges from the pre-loaded catheter and the up-and-over sheath).
4. Continue to advance the sheath over the secured wire guide while applying gentle, intermittent tension on the pre-loaded catheter so that the sheath and catheter will advance over the bifurcation, into the proximal opening of the device, and out through the sidebranch. (Figure 14)

NOTE: Dependent upon the length of the dilator tip of this sheath, it may be necessary to advance the sheath over its dilator tip, to ensure sheath placement within the sidebranch.

5. Remove the dilator of the up-and-over sheath.

10.7 Cannulation of Internal Iliac Artery

1. Puncture the hemostatic valve of the up-and-over sheath as far to one side as possible and insert a suitable wire guide and catheter combination. (Figure 15)
2. Advance the wire guide and catheter through the up-and-over sheath into the internal iliac artery. (Figure 16)
3. Replace the wire guide in the internal iliac artery with a support wire guide (e.g., Rosen or AES wire guide). Ensure that the wire guide placement is stabilized in the internal iliac artery.

10.8 Balloon-Expandable Covered Peripheral Bridging Stent Placement

1. Remove the catheter from within the up-and-over sheath, leaving the sheath and wire guides in place, and advance the pre-mounted balloon-expandable covered peripheral stent of suitable size over the support wire guide until it enters the internal iliac artery. (Figure 17)

NOTE: Pre-loading the pre-mounted balloon-expandable covered peripheral stent into an additional sheath may be desirable to provide protection as the stent is advanced through the up-and-over sheath and into the internal iliac artery. Suitable sheath/wire guides/stent combinations must be selected such that the inner sheath can pass through the up-and-over sheath in a co-axial manner while allowing the through-and-through wire to remain in place.

2. Check angiography as required to determine that the position of the balloon-expandable covered peripheral stent is satisfactory with respect to the sidebranch overlap and the intended internal iliac artery fixation site.

NOTE: When deployed, the balloon-expandable covered peripheral stent should overlap the Zenith Branch Endovascular Graft-Iliac Bifurcation sidebranch by 10-14 mm.

NOTE: Angiography should be possible through the up-and-over sheath during this stage of the procedure or through the catheter in the common iliac artery.

3. Withdraw the hydrophilic nitinol core through-and-through wire guide and the pre-loaded catheter.
4. Adjust the position of the Zenith Branch Endovascular Graft-Iliac Bifurcation, if necessary, so that the distal opening of the sidebranch is in a suitable location in relationship to the origin of the internal iliac artery.

10.9 Distal Zenith Branch Endovascular Graft-Iliac Bifurcation Deployment

1. Withdraw the introducer sheath of the Zenith Branch Endovascular Graft-Iliac Bifurcation until the external iliac component is fully deployed. (Figure 18)
2. Remove the safety lock from the black trigger-wire release mechanism. Under fluoroscopy withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 19)

NOTE: The distal stent is still secured by the trigger-wire.

10.10 Balloon-Expandable Covered Peripheral Stent Deployment

1. Withdraw the up-and-over sheath to a suitable position to allow deployment of the balloon-expandable covered peripheral stent.
2. Deploy the balloon-expandable covered peripheral stent and perform additional ballooning as required. (Figure 20)

NOTE: The balloon-expandable covered peripheral stent should be expanded to at least 8 mm in the sidebranch of the Zenith Branch Endovascular Graft-Iliac Bifurcation, and expanded to the appropriate size to achieve seal in the internal iliac artery.

NOTE: The balloon-expandable covered peripheral stent should overlap the Zenith Branch Endovascular Graft-Iliac Bifurcation sidebranch by a minimum of 10 mm, and a maximum of no further than the second marker on the distal end of the side branch, with a distal fixation site within the internal iliac artery of at least 10 mm (preferably 20-30 mm).

NOTE: Do not yet remove the balloon from the sidebranch segment of the Zenith Branch Endovascular Graft-Iliac Bifurcation.

3. Remove the safety lock from the white trigger-wire release mechanism. Under fluoroscopy withdraw and remove the trigger-wire by sliding the white trigger-wire release mechanism off the handle and then remove via its slot over the inner cannula. (Figure 21)
4. Re-inflate the balloon in the sidebranch segment to stabilize the device and withdraw the introducer system through the graft into the introducer sheath.

CAUTION: Care should be taken to avoid displacing or damaging the Zenith Branch Endovascular Graft-Iliac Bifurcation when removing the introducer system or when advancing other components through it.

NOTE: Ensure that the wire guide remains in a satisfactory position to cannulate the Zenith AAA graft.

5. Deflate and remove the balloon.
6. Withdraw the up-and-over wire guide and sheath into the opposite common iliac artery.

10.11 Zenith AAA Graft and Iliac Leg Graft Deployment

1. Introduce a catheter, and stiff wire guide (AUS or LES) standard 0.035 inch (0.89 mm), 260 cm long through the iliac artery access in the side opposite the Zenith Branch Endovascular Graft-Iliac Bifurcation.
2. Remove the existing sheaths and catheters, and introduce and deploy the main body of the Zenith bifurcated Graft according to standard Zenith protocol (see relevant Zenith Graft Instructions for Use).
3. Cannulate the short (contralateral) limb of the Zenith bifurcated graft through the Zenith Branch Endovascular Graft-Iliac Bifurcation.
4. Deploy an iliac leg graft (e.g., TFL-16-XX, ZSLE-16-XX-ZT, ZISL-16-XX) of suitable length through the Zenith Branch Endovascular Graft-Iliac Bifurcation such that the recommended overlap is achieved as shown below.

Recommended Overlap of Zenith Iliac Leg Graft connecting the Contralateral Limb of Zenith bifurcated Graft to the Zenith Branch Endovascular Graft-Iliac Bifurcation

Overlap into the contralateral limb of the Zenith bifurcated graft	Minimum 1 ¼ stents
	Maximum 1 ¾ stents
Overlap into the body of the Zenith Branch Endovascular Graft	Minimum 1 ½ stents
	Maximum of no further than the second marker on the proximal end of the device

5. Continue deployment of the Zenith bifurcated graft as outlined in the relevant Instructions for Use.

10.12 Molding Balloon Insertion

NOTE: Prior to molding balloon use, carefully remove any accessory devices.

1. For the Zenith bifurcated Graft use molding balloon as outlined in the relevant Instructions for Use.
2. For the Zenith Branch Endovascular Graft-Iliac Bifurcation prepare the molding balloon as follows:
 - Flush wire lumen with heparinized saline
 - Remove all air from balloon
3. In preparation of the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counter-clockwise.
4. Advance the molding balloon over the stiff wire guide and through the Captor Hemostatic Valve of the Zenith Branch Endovascular Graft-Iliac Bifurcation introduction system to the overlap between the common iliac segment of the Zenith Branch Endovascular Graft-Iliac Bifurcation and the Zenith® Iliac Leg grafts (e.g., TFLE, ZSLE, ZISL) and expand.

CAUTION: Confirm complete deflation of balloon prior to repositioning.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

NOTE: When ballooning in the vicinity of the bifurcation of the Zenith Branch Endovascular Graft-Iliac Bifurcation, ensure that the balloon does not cross the bifurcation.

11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

Please refer to **Table 11** for imaging and post-operative follow-up requirements.

Table 11 Imaging Schedule

	Pre-op	Intra-op	Post-Procedure	30 Day	6 Month	12 Month	24 Month
CT	X		X ¹	X ¹	X ¹	X ¹	X ¹
KUB (Supine radiograph, Kidneys, Ureters, and Bladder)			X	X	X	X	X
Angiography	X	X					

¹Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan.

NOTE: Do not balloon in the vicinity of the sidebranch to avoid disruption of the balloon expandable stent.

5. Withdraw the molding balloon to the Zenith Branch Endovascular Graft-Iliac Bifurcation distal fixation site and expand.

CAUTION: Do not inflate balloon in vessel outside of graft.

6. Remove molding balloon and replace it with an angiographic catheter to perform completion angiograms.
7. Remove or replace all stiff wire guides to allow iliac arteries to resume their natural position.

10.13 Final Angiogram

1. Position angiographic catheter just above the level of the aortic bifurcation. Perform angiographs to verify that the common iliac artery and corresponding internal iliac artery are patent and that there are no endoleaks.
2. Confirm there are no endoleaks or kinks. Remove the sheaths, wires and catheters.

NOTE: If endoleaks or other problems are observed and require intervention, consider use of additional Zenith AAA Ancillary Components.

3. Repair vessels and close in standard surgical fashion.

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© COOK 2020
WILLIAM A. COOK AUSTRALIA PTY. LTD.
COOK MEDICAL EUROPE LTD.