**EN**  
**2**

**Double Lumen Ovum Aspiration  
Needle (K-OPSD series)**

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**MANUFACTURER**

WILLIAM A. COOK AUSTRALIA PTY. LTD.  
95 Brandl Street, Brisbane Technology Park  
Eight Mile Plains, Brisbane  
QLD 4113 AUSTRALIA

## DOUBLE LUMEN OVUM ASPIRATION NEEDLE (K-OPSD SERIES)

### CONTAINS

A stainless steel double lumen needle with a bevelled tip which allows for puncture and an echogenic tip which enhances visualization under ultrasound.

**Needle Gauge:** 16 or 17G

**Needle Length:** 30, 33 or 35 cm

**Aspiration Tubing:** 75 cm, 1.32 mm OD

**Flushing Tubing:** 70 or 100 cm length, 1.83mm (16G) or 2.01mm (17G) OD

**Connector type** (e.g. for connection to pump): female luer fitting (6% taper)

### INDICATION FOR USE

The Double Lumen Ovum Aspiration Needles are used for laparoscopic or ultrasound guided transvaginal aspiration and flushing of oocytes from ovarian follicles for patients undergoing Assisted Reproductive procedures.

**NOTE:** Two cell mouse embryo tested and passed with 80% or greater Blastocyst rate at 72h. USP endotoxin (LAL) tested and passed with 20 EU's or less per device. Testing is conducted on a lot-to lot (batch) basis.

**CAUTION:** Sterile if the package is unopened and undamaged. Do not use if package is broken. Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

### CONTRAINDICATIONS

This needle should not be used on a patient with an active vaginal or intrauterine infection, a sexually transmitted disease, or who is currently pregnant.

### PRECAUTIONS

The assisted reproduction needle (and any other accessories used during these procedures having direct contact with oocytes) should be comprised of embryo compatible materials.

Pain may be associated with oocyte retrieval via percutaneous abdominal needle puncture or transvaginal puncture, but is usually perceived as tolerable when using smaller gage needles.

Local anaesthetic may be administered as necessary.

Hematuria may occur due to the aspiration needle penetrating a filled bladder during transvaginal ultrasound aspiration. This complication typically resolves spontaneously within a day. Extravasation of urine may occur within the abdominal cavity when the needle puncture traverses the bladder. Patients should be monitored for evidence of this known complication, however, there is typically no associated discomfort or adverse sequelae.

Infection may be introduced via needle puncture and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), uterine infection or cystitis.

Recommendations to minimise occurrence include use of only biocompatible materials, flushing the needle (and any other accessories used) with sterile, compatible culture media, and closely adhering to sterile techniques.

Vaginal bleeding has been reported to be associated with the transvaginal route for oocyte retrieval via needle aspiration. Bleeding is typically easily controlled with direct pressure.

### INSTRUCTIONS FOR USING OVUM ASPIRATION NEEDLE ULTRASOUND GUIDED PROCEDURE

1. Position the patient in the lithotomy or knee/chest position.
2. Place an ultrasound transducer into the vaginal fornix to visualize the ovary and follicles.
3. Insert the needle into the transducer needle guide.
4. Under ultrasound guidance, advance the needle into the ovarian follicle.
5. Using a vacuum unit or syringe, aspirate and/or flush the follicle to obtain the oocyte. Repeat steps 4 and 5 on remaining follicles. The vacuum pressure used with a specific gage and/or type of ovum pick-up needle is at the discretion of the clinician performing the procedure.
6. Remove needle and reposition the transducer to visualise the remaining ovary. Repeat steps 4 and 5 on remaining follicles.
7. When the desired number of oocytes have been aspirated, remove the needle and discard.

### LAPAROSCOPIC PROCEDURE

1. Position the patient in the lithotomy or knee/chest position.
2. Under laparoscopic guidance, nick the skin for ease of trocar/cannula placement.
3. Place trocar/needle assembly to the desired location. Remove trocar while leaving the needle cannula in place.
4. Under laparoscopic vision, place the needle through the abdominal cannula and advance into the ovarian follicle.
5. Using a vacuum unit or syringe, aspirate and/or flush the follicle to obtain the oocyte. Repeat steps 4 and 5 on remaining follicles and ovary.
6. When the desired number of oocytes have been aspirated, remove the needle and discard.

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