MINC+™

BENCHTOP INCUBATOR

K-MINC-2000

USER MANUAL
CE0123

English

Addendum to IFU-MINC3_EN-6, IFU-MINC3_EU-1, and IFU-MINC3_US-1

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IFU-MINC3_EN-6, IFU-MINC3_EU-1 ve IFU-MINC3_US-1'e ek

English

Section 5.6.2 of the IFU is no longer applicable.

This addendum is to replace section 5.6.2 of the IFU with the following:

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- This will trigger external alarm contacts to activate if the power is still off after two minutes.
- If the power is restored within two minutes, the external alarm will not activate and the MINC+ will recommence normal operation.
- If the power is restored after two minutes, the external alarm will reset and the MINC+ will recommence normal operation.
- If the gas flow was active before the power was lost, then the MINC+ will restart the flow sequence in gas purge mode when power is restored.

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Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

Email: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

Email: cau.custserv@cookmedical.com

WILLIAM A. COOK AUSTRALIA PTY. LTD.

95 Brandl Street, Eight Mile Plains, QLD 4113, AUSTRALIA

COOK MEDICAL EUROPE LTD.

O'Halloran Road, National Technology Park, Limerick, IRELAND

Language index

English	1
Español	2
Français	3
Português	4

General information

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

Please familiarize yourself with the safety instructions before using the device.

This device may only be used by physicians and medical assistants with the appropriate technical qualification.

1

Copyright

This manual contains information that is subject to copyright. All rights reserved. This manual should not be photocopied, duplicated on microfilm, or otherwise copied or distributed, completely or in part, without the approval of William A. Cook Australia Pty. Ltd.

Some of the parts and equipment referred to in this manual bear registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark indicates that any given designation is not subject to trademark protection.

Users of William A. Cook Australia Pty. Ltd. products should not hesitate to contact us if there are any unclear points or ambiguities in this manual.

© COOK 2022

2022-02

Document No: IFU-MINC3_US-1

Service address:

Please refer to your local Cook Medical distributor for details of your nearest authorized service agent.

Trademarks

MINC+ and DishTrace are trademarks of Cook Medical. Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and/or other countries. Bonjour is a registered trademark or trademark of Apple Inc. in the United States and/or other countries. All other trademarks and trade names are the property of their respective owners.

Contents

GENERAL INFORMATION	1-1	4.7 Gas supply	1-8
Copyright.....	1-1	4.7.1 Gas mixture required	1-8
Trademarks.....	1-1	4.7.2 Gas regulator selection.....	1-8
CONTENTS	1-1	4.7.3 Connection to gas cylinder	1-9
1. OVERVIEW OF MANUAL.....	1-3	4.7.4 Connection to single device	1-9
1.1 Purpose of this document	1-3	4.7.5 Connection for multiple devices.....	1-9
1.2 Explanation of pictograms.....	1-3	4.7.6 Testing for gas leaks	1-9
1.3 Intended audience	1-3	4.7.7 Other gas connections	1-10
1.4 Warnings and important notes	1-3	4.8 External alarm connection.....	1-10
1.5 Safety instructions.....	1-4	4.9 Ethernet connection.....	1-10
1.6 Table of definitions or actions.....	1-4	4.10 Humidification flasks	1-10
2. ABOUT THE MINC+ BENCHTOP INCUBATOR.....	1-4	4.10.1 Humidification flask components.....	1-10
2.1 Intended use	1-4	4.10.2 Filling the flasks	1-10
2.2 Indications for use	1-4	4.10.3 Flask installation	1-11
2.3 Intended users and use environment	1-4	4.10.4 Removing flasks for replacement	1-12
2.4 Contraindications	1-4	4.11 Activating the device	1-12
2.5 Device description	1-4	4.11.1 Unlock button function	1-12
2.5.1 Humidification flask	1-5	4.11.2 Entering standby and active modes	1-12
2.6 Precautions for device use	1-5	4.12 Initial configuration	1-12
3. PREREQUISITES	1-5	4.12.1 Assign a name to the incubator	1-12
3.1 MINC+ Benchtop Incubator	1-5	4.12.2 Assign a language	1-12
3.2 DishTrace PC	1-5	4.12.3 Network configuration	1-12
3.3 Network configuration	1-5	4.12.4 Verifying the date and time	1-12
3.4 Cybersecurity	1-5	4.13 Installing the DishTrace PC software	1-13
4. INSTALLATION AND SETUP.....	1-6	4.13.1 DishTrace PC software	1-13
4.1 MINC+ Benchtop Incubator	1-6	4.13.2 Before you begin	1-13
4.2 Unpacking and contents	1-6	4.13.3 Installing the DishTrace PC software	1-13
4.2.1 Supplied items	1-6	4.13.4 Add users to DishTrace PC users group	1-13
4.2.2 Must be supplied	1-6	4.13.5 Export DishTrace PC certificate	1-13
4.3 Front of the device	1-6	4.13.6 Import DishTrace PC certificate into MINC+	1-14
4.4 Rear of the device	1-7	4.13.7 Establishing a secure network	1-14
4.5 Electromagnetic compatibility	1-7	4.13.8 Run DishTrace PC	1-14
4.6 Incubator placement	1-7		

Contents (continued)

4.14 Chamber operation.....	1-14	5.5.12 Editing a dish record.....	1-23
4.14.1 Enable incubation chamber keypad.....	1-14	5.5.13 Dish records and dish reporting.....	1-23
4.14.2 Turning the incubation chamber heater on or off.....	1-15	5.5.14 Export a dish-specific incubator parameter report.....	1-24
4.14.3 Incubation chamber temperature selection.....	1-15	5.5.15 System tab	1-24
4.14.4 Lid openings	1-15	5.5.16 To view incubator alerts.....	1-24
4.14.5 Incubation chamber gas flow	1-15	5.5.17 DishTrace PC settings.....	1-24
4.14.6 Turning on water-level monitoring for the humidification flask.....	1-15	5.5.18 Back up the DishTrace PC database	1-24
4.14.7 Media warming chamber (MWC)	1-15	5.5.19 Restore DishTrace PC database.....	1-24
4.14.8 First time commissioning run	1-16	5.5.20 Purge the DishTrace PC database.....	1-24
4.15 Installation and setup checklist.....	1-16	5.6 Alarm conditions.....	1-25
5. DEVICE OPERATION 1-16		5.6.1 External alarm.....	1-25
5.1 Before use	1-16	5.6.2 Loss of mains power	1-25
5.2 Inserting culture dishes.....	1-16	5.6.3 Low inlet gas pressure	1-25
5.3 Patient identification plates	1-16	5.6.4 Low humidification flask water level	1-25
5.3.1 Placing test tubes inside the media warming chamber.....	1-17	5.6.5 No gas flow or gas flow out of range.....	1-25
5.4 DishTrace MINC+	1-17	5.6.6 Temperature out of range	1-25
5.4.1 Resetting the center console touchscreen.....	1-17	5.6.7 Chamber screen and center console error displays	1-25
5.4.2 DishTrace MINC+ setup.....	1-17	5.6.8 Internal functional error	1-25
5.4.3 Assigning dishes to a chamber	1-17	5.7 After use	1-25
5.4.4 Checking out a dish.....	1-17	6. SERVICE AND MAINTENANCE 1-26	
5.4.5 Check in a dish to its preassigned location.....	1-18	6.1 Changing the humidification flasks.....	1-26
5.4.6 Unassign a dish position	1-18	6.2 Cleaning the incubator.....	1-26
5.4.7 Edit dish record	1-18	6.2.1 Full cleaning	1-26
5.4.8 Lid open alert – Dish assignment	1-18	6.2.2 On-site cleaning	1-26
5.4.9 Records screen	1-18	6.3 Cleaning the media warming chamber	1-26
5.4.10 To search for a dish or patient record and generate a dish report..	1-19	6.4 Biannual functionality testing	1-27
5.4.11 To generate an incubator report	1-19	6.4.1 Temperature	1-27
5.4.12 Status screen	1-19	6.4.2 Gas flow rate	1-27
5.4.13 Settings screen – General tab	1-19	6.4.3 External alarm contacts	1-27
5.4.14 Configuration tab.....	1-20	6.5 Inspection by an authorized Cook Medical service agent	1-27
5.4.15 Network tab	1-20	6.6 Return procedure	1-27
5.4.16 Date and time settings	1-20	7. TROUBLESHOOTING 1-28	
5.4.17 Database and backup information	1-20	8. TECHNICAL DATA 1-30	
5.4.18 Purging data	1-20	9. LIMITED WARRANTY 1-33	
5.4.19 Backup data to an external USB drive	1-20	9.1 Liability	1-33
5.4.20 Restore data settings from an external USB drive	1-20	9.2 Life of the product	1-33
5.4.21 Upgrade DishTrace MINC+ software	1-21		
5.4.22 Restore factory settings	1-21		
5.4.23 Info – Device settings	1-21		
5.4.24 DishTrace MINC+ screen saver	1-21		
5.5 DishTrace PC software	1-21		
5.5.1 Using the DishTrace PC software	1-21		
5.5.2 Incubators screen	1-21		
5.5.3 Adding incubators	1-21		
5.5.4 Incubator display	1-21		
5.5.5 To edit incubator settings	1-22		
5.5.6 View Contents, History, Events and Information	1-22		
5.5.7 To generate a lid openings report	1-23		
5.5.8 To export a running log of the incubator	1-23		
5.5.9 Records	1-23		
5.5.10 Search for a dish record	1-23		
5.5.11 Adding a new dish record	1-23		

1. Overview of manual

1.1 Purpose of this document

This guide describes the operation of the MINC+ Benchtop Incubator (K-MINC-2000) system including the disposable humidification flask (K-MINC-2000-HF).

1.2 Explanation of pictograms

The following pictograms appear on the MINC+ Benchtop Incubator front and back panels, the device packaging, and/or the sterile humidification flask packaging.

	Caution: Read all warnings and precautions in the instructions for use before using the device
	Consult instructions for use
	Standby/on
	Unlock
	Enable/disable gas flow
	Increase/decrease temperature set point
	Enable/disable water level
	Enable/disable heater
	USB connection
30V ~ 2.0A	Normally open/normally closed alarm contact rating
	Gas cylinder
135 kPa - 165 kPa 600 mL/min	Inlet gas pressure tolerance/maximum flow rate
	Gas supply inlet
100-240 V~ 50-60 Hz 200 VA	Power supply/frequency/max power consumption
	Center console system reset button
	Ethernet connection
	Degree of enclosure protection from solid objects and liquids
	Dispose of in accordance with WEEE directive (2012/19/EU)
[SN]	Serial number
EC REP	Authorized Representative in the European Community/ European Union
REF	Catalogue number
	CE Mark approval
	UL Mark approval
[MD]	Medical device
	Temperature limit
	Humidity limitation

Two-person lift required

This way up

Fragile, handle package with care

Rx only

Prescription only

Keep dry

Manufacturer

Keep away from sunlight

Do not use if package is damaged and consult instructions for use

Date of manufacture

Lot Number

STERILE | EO

Sterilized using ethylene oxide (1)

Quantity per box (1)

Use-by date (1)

Do not reuse (1)

Do not resterilize (1)

Note (1): Symbols are on sterile disposable humidification flask packaging only.

1.3 Intended audience

This manual is intended for clinicians and authorized Cook Medical service agents involved with the operation or maintenance of the MINC+ Benchtop Incubator.

1.4 Warnings and important notes

Throughout these Instructions for Use, blocks of text may be accompanied by a pictogram and/or printed in bold type. These are warnings and important notes and are used as follows:

WARNING: The personal safety of the user may be involved. Disregarding this information could result in injury to the operator, device, or the contents!

WARNING: BIOLOGICAL HAZARD

WARNING: ELECTRIC SHOCK HAZARD

WARNING: EXPLOSION HAZARD

WARNING: RADIO FREQUENCY INTERFERENCE

IMPORTANT NOTE

This symbol indicates that this product may not be treated as municipal waste. Disposal of this product must be undertaken in accordance with the WEEE directive (2012/19/EU).

Please ensure that this product is properly disposed of. Inappropriate waste handling of this product may cause potential hazards to the environment and human health. For more detailed information about disposal of this product, please contact your local city office or Cook Medical service agent.

Cook fulfills its legal obligations with regards to WEEE and waste packaging compliance through its own takeback initiatives and through national takeback schemes.

Please see <https://www.cookmedical.com/about/sustainability-environmental-practices/waste-electrical-electronic-equipment-directive/> for details on how to properly recycle WEEE or waste packaging in your country.

1.5 Safety instructions

WARNING: It is essential to use this manual to familiarize yourself with the safety instructions, device functions, and operation before use. Not following these instructions can result in damage or breakdown of the device and device contents and/or user injury. This device should be used only by qualified personnel.

WARNING: ELECTRIC SHOCK HAZARD. The equipment is to be used only with electrical systems complying with all IEC, CEC, and NEC requirements.

Internal circuitry is energized whenever the MINC+ is connected to mains power, irrespective of the displays being illuminated. Always disconnect the device from mains power before cord replacement or cleaning. Should any power cord or plug associated with the MINC+ become cracked, frayed, broken, or damaged, it must be replaced immediately.

To reduce the risk of electric shock, do not remove covers. Please refer all servicing to the Cook Medical service agent.

WARNING: Do not sterilize the device.

WARNING: There are no user-serviceable parts inside the device. Any adjustment, modification, or repairs to the equipment should be carried out only by persons authorized to perform them.

WARNING: Protect the incubator from being splashed by liquid. Should any liquid enter the device, discontinue use immediately.

WARNING: EXPLOSION HAZARD. Do not use the MINC+ in an area where flammable gases are present.

The MINC+ can cause explosion in the presence of flammable gases.

WARNING: BIOLOGICAL HAZARD. Do not use a contaminated humidification flask in the MINC+.

If the humidification flasks are used continuously, they may be used until the water level reaches the MIN FILL line on the front of the flask.

If the humidification flasks are not used continuously, there may be an increased risk of biological contamination. If the device is not in use, remove the flasks and replace with new flasks upon startup.

Do not refill the humidification flasks.

WARNING: GAS UNDER PRESSURE. Use only the PTFE-lined gas connecting hose supplied with the MINC+ Benchtop Incubator to connect it to a gas supply. The use of another tubing type or hose may result in the loss of desired gas concentrations.

WARNING: For your own safety, use only Cook branded consumables and spare parts.

1.6 Table of definitions or actions

Term	Definition or action
DishTrace™ PC	Software used to manage dish records from the incubator via a PC.
DishTrace MINC+	Incubator-based dish-data management application that allows users to add dish information, assign a location, and review dish status and settings.
Keypad Enable	Unlock keypad by pressing the unlock button for 2 seconds; the keypad will beep on touch and the backlight will come on when the keypad is active.
MINC+	Abbreviated device name
MWC	Media warming chamber
Restart	When instructed to perform a restart, pull power cord out, wait at least 2 minutes, and plug back in. Incubator will reboot to its last known active state. In some cases the device may automatically perform a restart or a center console only restart to clear an internal error.

2. About the MINC+ Benchtop Incubator

2.1 Intended use

The MINC+ Benchtop Incubator (MINC+) is intended to store and preserve gametes and/or embryos close to body temperature.

Figure 1: MINC+ Benchtop Incubator

2.2 Indications for use

The MINC+ Benchtop Incubator is indicated to be used in the assisted reproduction technology as part of the process to treat infertile patients.

2.3 Intended users and use environment

The MINC+ Benchtop Incubator device is intended to be used by clinical embryologists. Service personnel may install the MINC+ Benchtop Incubator and will routinely service the device.

The MINC+ Benchtop Incubator is intended to be used within an IVF laboratory.

2.4 Contraindications

There are no known contraindications for this device.

2.5 Device description

The MINC+ includes the following features:

- Two incubation chambers (located on the left and right) to hold culture dishes.
- Two humidification flasks (located on the left and right) containing water for humidifying the gas supplied to each incubation chamber.
- DishTrace MINC+ – Dish data management via the center console.
- DishTrace PC – Software for dish record management on networked PC.
- Left and right keypads for controlling respective incubation chambers, parameters and humidification flasks' water level monitoring.
- Left and right incubator chamber status screens, which indicate current chamber operating environment.
- Media warming chamber (MWC) – Located between the incubation chambers.

The MINC+ contains two individual incubation chambers with independent temperature control that can be selected in the range of 35 °C to 40 °C. There is also gas flow provided independently to each chamber at a nominal fixed flow rate of 25 mL/min to optimize culture conditions inside the incubation chambers.

Premixed gas is supplied to each incubation chamber independently through a humidification flask. The gas is humidified as it passes through the water inside the flask. The use of humidified gas may reduce the evaporation rate of culture media during incubation.

Each incubation chamber has an individual control panel to turn the chamber on or off, set the desired temperature, activate the water level monitoring (in the humidification flasks), and control the gas flow status.

Each chamber has a display indicating the chamber status, the current temperature, temperature set point, gas flow status, water level sensing status, and any associated alarms.

The information displayed on each chamber screen can also be viewed on the center console. The center console does not control the chamber operational parameters. These operational parameters can be changed only by using the individual chamber keypads.

The MINC+ can accept Nunc® four well culture dishes, or Nunc and BD® Falcon® 35 mm and 60 mm single round culture dishes, or equivalent sized culture ware in two separate chambers. If using alternate culture ware, the user must ensure that dish fits in the chamber, sits flat on the base plate and locates correctly within the chamber base plate grooves.

The heater blocks under each culture dish are in direct contact with the bottom surface of the dish. Whenever the chamber lid is opened and closed, a rapid purge of the gas mix occurs to reestablish the required environment.

The MINC+ provides performance security by continuously monitoring critical functions. If a deviation from normal operation is detected, the user is alerted. The device can be connected to an external remote alarm to alert staff outside of normal working hours to any problems that the device is unable to self-correct in the allotted time.

The media warming chamber (MWC) can be used to maintain a stable temperature for test tubes containing IVF-related liquids such as media or overlay oil. The MWC is not gassed and only provides a constant temperature, which is manually set.

Dish data can be managed on the MINC+ via DishTrace MINC+ on the center console touchscreen on the front of the device. The touchscreen allows IVF clinicians to assign the locations of culture dishes within each chamber, check the history of dishes previously logged, and check the chamber temperature conditions.

The incubator may be connected to either an ethernet network or directly to a personal computer (PC) via LAN cable connection.

The DishTrace PC software shows the status of multiple MINC+ incubators connected on the PC network. It can connect up to 50 MINC+ Benchtop Incubators, including up to 10 (K-MINC-1000) MINC Benchtop Incubators (serial number A803560 onwards). DishTrace PC cannot control the MINC+ or the MINC incubator settings.

2.5.1 Humidification flask

The MINC+ uses a sterile, disposable humidification flask for each incubation chamber. The purpose of the flasks is to humidify and direct the gas into each chamber. The flask includes an internal filter to prevent backflow from filled flasks from entering the incubator's internal gas plumbing. The reorder code of the humidification flasks is K-MINC-2000-HF.

There are two humidification flasks included with the incubator. Each flask is supplied sterilized with EtO (ethylene oxide) in a sealed pouch.

Each pouch includes:

- One humidification flask.
- One disposable funnel.
- One flask T-piece.
- One temporary support clip to aid in filling the flask.

Refer to **Section 4.10** for more information

2.6 Precautions for device use

If any electrical or mechanical fault occurs while the MINC+ is in use, stop using the device until it has been checked by an authorized Cook Medical service agent. Complete the installation and setup checklist prior to operating the MINC+ (see **Section 4.15**).

3. Prerequisites

A detailed checklist of prerequisites is provided for each component further in the document. In summary, the following are needed:

3.1 MINC+ Benchtop Incubator (see Section 4 for details)

- Clean, secured laboratory environment.
- Stable mains power supply.
- Premixed gas supply.
- Sterilized distilled water for flasks.
- A 9/16" wrench (spanner) to connect the Braided Connecting Tube.
- Uninterruptible power supply (UPS) (recommended).

3.2 DishTrace PC (see Sections 4.13 and 5.5 for details)

- A PC running Microsoft Windows® 8.1 or 10 operating system.
- Mouse with center button scroll wheel.
- 4 GB of RAM.
- 10 GB of disc space for data storage.
- USB connection hub (up to 10 ports) for (K-MINC-1000) MINC Mini Incubators.
- Ethernet connection for MINC+ connectivity.
- USB drive for backups from MINC+ (FAT32 format).
- Antivirus software (e.g., Microsoft Windows Defender Antivirus or McAfee® for cybersecurity risk mitigation).
- All PC and software updates must be completed before DishTrace PC is installed. A PC restart may be required during an update process. Ensure that all Windows operating systems and application updates have been completed and installed successfully. Refer to **Troubleshooting (Section 7), DishTrace PC install and connectivity**, for instructions on how to verify that the updates have been successfully installed.

3.3 Network configuration

- DHCP server to assign dynamic IP addresses to each MINC+ incubator and PC.
- NTP server to provide date and time information to the incubators.
- Ethernet LAN connection for the MINC+ incubator and DishTrace PC.

3.4 Cybersecurity

- MINC+ incubators and the PC running DishTrace PC should be on a separate network segment that is firewalled from the business network and the Internet to ensure appropriate security.
- Use a dedicated PC for the DishTrace PC software and avoid the use of business productivity applications such as email, etc., on this PC.
- The PC used for DishTrace PC should be configured to require user logon (no automatic logon).
- The PC used for DishTrace PC should have an antivirus program and BitLocker® enabled.
- Ensure that physical access to the MINC+ incubators is restricted to only authorized users.
- MINC+ patching and security updates must be installed and kept up to date.

The MINC+ incubator has been designed to use only the following network port services:

- HTTPS (TCP 443, bidirectional).
- NTP client (random UDP port above 1023, bidirectional).
- Autodiscover (UDP 5353, bidirectional).
- DHCP (UDP 68, bidirectional).

The following network ports are used by DishTrace PC:

- HTTPS (TCP 443, bidirectional).
- SQL (TCP 1433, bidirectional, localhost only).
- Autodiscover (UDP 5353, bidirectional).

 **IMPORTANT NOTE:** In situations where it is necessary to retrieve system logs from DishTrace MINC+ (e.g., cybersecurity event investigation), contact your Cook Medical service agent for assistance.

4. Installation and setup

4.1 MINC+ Benchtop Incubator

An installation and setup checklist has been included in **Section 4.15**. This must be used to help ensure correct preparation before using the device.

IMPORTANT NOTE: To reduce the risk of electrical damage, the use of an uninterruptible power supply (UPS) with power conditioning capabilities is strongly recommended. To ensure uninterrupted device functionality, the use of the following is also strongly recommended:

- An automatic gas cylinder changeover unit that can maintain a minimum pressure of 135 kPa.
- A remote alarm response system to alert staff to any failures in the gas or power supply.

IMPORTANT NOTE: Use of cables other than those provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

IMPORTANT NOTE: This device has been certified by an independent testing authority for use with USB and with remote alarm cables less than 3 meters in length.

4.2 Unpacking and contents

IMPORTANT NOTE: It is important to retain packaging for future use. (Refer to **Section 6.6**.)

WARNING: The incubator is heavy – it weighs approximately 24 kg (53 lb). Use the handholds at the base of the incubator when lifting the device. To avoid injury, two persons should lift the device when it is being moved.

4.2.1 Supplied items

Figure 2: Supplied items

The following items are supplied inside the box:

- MINC+ Benchtop Incubator.
- MINC+ user manual.
- Disposable humidification flasks (two flasks – K-MINC-2000-HF).
- Braided Connecting Tube (3 m length) (K-MINC-BCT-10-300).
- Medical grade mains power cord.
- USB flash drive with DishTrace PC software.
- Ethernet cable for connection to a network, or PC.
- External alarm connector on rear of device.
- Magnetic patient identification plates (12) on each chamber lid for the recording of patient information.
- Gas inlet T-piece with thread protector and blanking cap (on rear of device).
- Test tube cradle in the media warming chamber.
- USB cover in the front panel USB port.

Check the MINC+ and all items immediately upon receipt to make sure that the contents are complete and that nothing is damaged. The manufacturer will only honor claims for compensation that are forwarded immediately to the sales representative or Cook Medical service agent.

IMPORTANT NOTE: Remove all items from packaging except the disposable humidification flasks, which need to be handled under sterile conditions.

4.2.2 Must be supplied

The following items are required in order to use the MINC+ but are not supplied:

- A clean source of mains power. The use of a UPS with power conditioning capabilities suitable for your local power supply system is strongly recommended.
- Premixed medical grade gas. All gas concentrations should be within a tolerance of $\pm 0.2\%$. Cylinders should be supplied with a certificate of analysis.
- A single-stage high-purity gas regulator that can supply the gas listed above at 150 kPa (1.5 bar, 22 psi) to the MINC+ gas inlet.
- An automatic gas cylinder changeover unit. To ensure the continuity of gas supply, an automatic changeover unit is strongly recommended.
- A source of sterile, distilled water.
- A pressure manometer to measure the input gas supply pressure. This item is not required if the regulator has gauges that display the gas output pressure.
- A 9/16 inch wrench (spanner) to tighten the connections of the Braided Connecting Tube.
- A PC/laptop with an ethernet connection for the DishTrace PC software.
- Ethernet network access if multiple incubators (up to a maximum of 50) will be used with the DishTrace PC software.

IMPORTANT NOTE: If a longer Braided Connecting Tube is required, contact your Cook Medical representative for further details.

4.3 Front of the device

Front controls

The front of the MINC+ incubator has two easy-to-use control panels with backlit buttons. Each function is displayed on the screen above the buttons so that what is happening in each chamber can be easily seen.

IMPORTANT NOTE: Front panel control panel buttons are touch-type proximity sensors and do not respond to increased finger pressure. The MWC temperature adjustment buttons are push-switch-type buttons and do not respond to increased finger pressure. Buttons do not auto-repeat if touch is sustained; move finger clear of pad and touch the button again to repeat input.


Figure 3: Front controls and features

The display icons, buttons, screens, and chambers are described below:

Buttons and features	Description
1	Standby/ON. Hold until a beep is heard. Keypad unlock is not required in order to switch from Standby to ON
2 Humidification flasks	A sterile, disposable humidification flask for each incubation chamber
3	Left-chamber status screen
4 Gas outlet port	Port that enables gas to exit the chamber
5 Chamber lid	Lid of the incubation chamber
6 Media warming chamber (MWC)	Chamber used to maintain a stable temperature for test tubes of IVF-related liquids such as media and overlay oil

1

Buttons and features	Description
7 Media warming chamber lights	MWC temperature status indicator
8	Keypad unlock (times out; will relock after 5 seconds of no keypad input)
9 USB port	USB port to download patient and dish information and for service maintenance updates
10	Center console touchscreen
11	Enable/disable gas flow
12	Enable/disable monitoring of the water level in the humidification flask
13	Increase/decrease chamber temperature set point
14	Enable/disable the incubation chamber heater
Display icons	Description
	Incubation chamber temperature set point
	Open incubation chamber lid
	Incubation chamber heater active
	Incubation chamber heater inactive
	Gas flow normal (flowing at a steady rate)
	Gas purge (actuates after the lid closes to reestablish environmental conditions)
	Gas flow disabled
	Water level monitoring disabled
	Water level empty
	Water level full

4.4 Rear of the device


Figure 4: Rear incubator controls

Description
1 Type B USB socket (for service only)
2 Center console reset (for DishTrace MINC+)
3 External alarm contacts (either normally open [NO] or normally closed [NC])
4 IEC power inlet socket
5 Gas inlet port
6 Ethernet connection socket

4.5 Electromagnetic compatibility

The MINC+ Benchtop Incubator complies with the electromagnetic compatibility (EMC) limits for medical devices as specified by IEC 60601-1-2: 2014. These limits provide reasonable protection against harmful interference in a typical medical environment.

Medical electrical equipment requires special precautions regarding EMC and must be installed and operated according to these instructions. It is possible that high levels of radiated or conducted radio-frequency electromagnetic interference (EMI) from portable and mobile RF communications equipment or other strong or nearby radio-frequency sources could result in performance disruption. Evidence of disruption may include erratic readings, equipment ceasing to operate, or other incorrect functioning. If any of these disruptions occur, cease using the incubator and contact your Cook Medical service agent.

For guidance and the manufacturer's declaration on electromagnetic emissions and immunity of the incubator, refer to **Section 8**.

4.6 Incubator placement

The MINC+ should be placed on a secure level surface away from heaters, coolers, air-conditioning outlets, mists, splashes, and exposure to direct sunlight. Position the device to ensure that it can be easily disconnected from the power supply. Once installed and placed into service, the MINC+ is not intended to be moved from one place to another. The lab environment should be clean and the device must be kept away from areas with dust.

Allow a reasonable amount of space around each MINC+ to ensure sufficient air circulation. At least 1 cm space is recommended around each device. Allow at least 60 cm of headroom, measured from the bottom of the incubator, to allow the lid to open.

WARNING: The MINC+ should not be used adjacent to or stacked with other equipment or be stacked during storage. If adjacent or stacked use is necessary, please ensure that appropriate shelving racks are used to allow the device to function normally.

WARNING: EXPLOSION HAZARD. Do not use the MINC+ in the presence of flammable gases.

WARNING: ELECTRIC SHOCK HAZARD. Do not immerse the MINC+.

WARNING: ELECTRIC SHOCK HAZARD. Determine if the available voltage corresponds with the MINC+. Connecting the device to the wrong voltage will cause the incubator to malfunction or may permanently damage the device. The power cord must be equipped with a safety plug. Use the enclosed power cord for the connection between the power point and the device socket.

WITHIN CANADA & USA – Use only a listed detachable power-supply cable, type SJT, minimum 18 AWG x 30, with 3 conductors, and one end configured for NEMA 5-15 and the other end for IEC 320/CEE22.

To avoid the risk of electric shock, this equipment must be connected to a supply main with protective earth.


Figure 5: Incubator placement and dimensions

4.7 Gas supply

4.7.1 Gas mixture required

⚠ WARNING: Select the appropriate high-purity gas mixture to suit the altitude and the culture media that is used.

👉 IMPORTANT NOTE: To maintain the correct operating pH level in bicarbonate-buffered media, the concentration of atmospheric CO₂ that is in contact with the media must be strictly controlled.

When the Cook Culture System is used at sea level, a 6% CO₂, 5% O₂, and 89% N₂ blend high-purity gas mixture is recommended.

Cook Medical recommends reducing the oxygen concentration from the normal atmospheric concentration to a value of 5%-8% in the atmosphere that is in contact with the media, because this reduced concentration may reduce the formation of reactive oxygen species.

It is recommended that the end user perform in-situ testing to determine the correct gas mix that is needed to attain the desired pH of the IVF media.

4.7.2 Gas regulator selection

The following information contains recommendations for the selection of the gas cylinder regulator that is to be used with the MINC+. Following these guidelines will provide a reliable gas connection between a gas cylinder and the MINC+.

👉 IMPORTANT NOTE: The MINC+ incubator must be connected to a regulated gas supply with a delivery pressure of 150 kPa ± 15 kPa at the MINC+ gas inlet. It cannot be directly connected to an unregulated gas supply. The gas regulator is an important part of the gas delivery system, and it must be designed and manufactured to supply medical-grade gases that are highly pure and accurately delivered to the device. It is also important that the regulator does not contaminate the gas stream.

Please request these features when ordering your gas regulator:

- Single-stage high-purity regulator.
- Metal-to-metal diaphragm seal.
- Stainless steel diaphragm, because stainless steel does not contaminate high-purity gas streams.

1

- Dual scale gauges (optional low-pressure gauge recommended, 200 kPa, full scale).
- Ability to connect to gas cylinders that contain a special blend of medical-grade gas.
- Delivery pressure of 150 kPa ± 15 kPa at the MINC+ gas inlet.
- Ability to maintain a 600 mL/min gas flow per MINC+ device during the gas purge cycle without undue supply-pressure sag.
- A Swagelok® SS-400-1-4RT outlet fitting to match the Braided Connecting Tube supplied with the MINC+ (assuming that the regulator outlet has a 1/4 inch ISO thread).


Figure 6: Gas regulator

4.7.3 Connection to gas cylinder

⚠️ WARNING: Use only a suitable medical gas-pressure regulator set to a nominal 150 kPa. Do not use flow restrictors or flow regulators in the gas line.

⚠️ WARNING: Do not use silicone tubing in the connection between the gas cylinder and the MINC+, or in series connections of multiple devices. Connect the regulator to the MINC+ with the supplied Braided Connecting Tube. It has a nontoxic, low-permeability PTFE lining that prevents the loss of CO₂ content associated with high-pressure use of more permeable materials such as silicone and PVC. Silicone tubing is relatively porous to pressurized CO₂ and must not be used anywhere in the connection from the gas cylinder, because the correct gas concentrations may not reach the MINC+ chambers.

The MINC+ and Braided Connecting Tube are fitted with Swagelok 1/4 inch series tube fittings ready for connection. The standard tube length is 3 meters. Alternative lengths are available in 60 cm, 100 cm, 6 m, 10 m, and 20 m. A union fitting is available to connect Braided Connecting Tubes. Refer to your local Cook Medical representative for more information.

The outlet fitting on the pressure regulator must be a Swagelok 1/4 inch series tube fitting. Any ferrule and backnut supplied with the regulator outlet fitting are surplus to requirements.

👉 IMPORTANT NOTE: All Swagelok 1/4 inch series tube fittings require a 9/16 inch wrench (spanner) for tightening the backnut. Ensure that all gas connections are tightly secured.

4.7.4 Connection to single device

- Connect the supplied Braided Connecting Tube to the regulator outlet fitting.
- Remove the thread protector cap from the gas inlet T-piece on the rear of the MINC+ device. (Retain the thread protector cap.)
- Connect the Braided Connecting Tube to the gas inlet T-piece, ensuring that it is tightly secured.
- Ensure that the blanking cap on the other side of the gas inlet T-piece is tightly secured also.


Figure 7: Back connection for a single device

4.7.5 Connection for multiple devices

If multiple MINC+ devices must be operated from a single gas supply point, the devices can be connected sequentially, as described below.

👉 IMPORTANT NOTE: Ensure that your pressure regulator can flow at a sufficient rate for the MINC+. Check that the regulated pressure is still a nominal 150 kPa at the last series connected MINC+ with all the devices purging.

A minimum flow capability (without causing undue supply pressure sag) of 600 mL/min per MINC+ device is required during the gas purge cycle. To ensure that the incubator at the end of the series receives sufficient gas pressure, it is recommended connecting no more than seven devices in series when the provided gas fittings are used. The Braided Connecting Tube can be ordered in different lengths to suit individual installation requirements. Contact your Cook Medical representative for further details.

⚠️ WARNING: Silicone tubing must not be used to connect multiple devices. For multiple connections, it is recommended to use only the Braided Connecting Tubes available from Cook.

- Remove the thread protector and blanking cap from the gas inlet T-piece connection points of all but the last device in the sequence. For the last device, remove only the thread protector.
- Connect the Braided Connecting Tube of the last device to the second gas connection point of the previous device.
- Repeat Step 2 until all devices are linked.
- Connect the supplied Braided Connecting Tube to the regulator outlet fitting.
- Connect the Braided Connecting Tube from the gas supply to the gas inlet T-piece of the first unit in the sequence.


Figure 8: Gas inlet T-piece connection for multiple devices


Figure 9: Gas connection to multiple devices

👉 IMPORTANT NOTE: Blanking caps must be tightened onto any unused gas connections. Store unused blanking caps in case the independent operation of the devices is required in the future.

4.7.6 Testing for gas leaks

The instructions below must be performed each time a Braided Connecting Tube connection is made to ensure that a reliable, leak-proof seal is established.

- Ensure that the fitting body, tube end, and ferrules are free of any foreign materials.
- Insert the tube adaptor with preswaged ferrules into the fitting body until the front ferrule contacts the fitting body.

3. Tighten the backnut firmly by hand. Use a 9/16 inch wrench (spanner) to securely tighten the backnut (usually 1/8 turn or less).
4. After applying $150 \text{ kPa} \pm 15 \text{ kPa}$ pressure to the line, check that the gas connections do not leak by brushing them with a brush that has been dipped in soapy water and then looking for the presence of bubbles. If bubbles are observed, turn off gas supply, disconnect the Braided Connecting Tube, and repeat **Step 1**, retightening the nut securely if no foreign material is found in the fitting.


Figure 10: Gas connection to incubator rear ports

4.7.7 Other gas connections

For other gas connections or existing gas delivery systems (e.g., gas cylinder changeover units or permanently installed gas supply systems), it is recommended that the user contact the relevant gas equipment supply center, e.g., BOC Gases or Air Liquide, for installation advice. For gas connection information, contact your Swagelok distributor (www.swagelok.com).

4.8 External alarm connection

The MINC+ can be connected to an external alarm system that will alert staff of alarms that occur outside of business hours. The MINC+ has two alarm contact configurations: normally open (NO) and normally closed (NC). Use the MINC+ only with alarm systems that can be activated by either contact closure or contact opening.

Refer to **Figure 4** in **Section 4.4** for connection points.

Battery-powered or safety extra-low voltage powered alarm equipment that does not exceed contact ratings can be connected to the external alarm terminals. Refer to **Section 8** for alarm contact rating.

4.9 Ethernet connection

If the DishTrace PC software will be used, connect the computer on which the DishTrace PC software is installed either directly to the device with the supplied ethernet cable or into a local area network connection port.

4.10 Humidification flasks

The MINC+ uses a sterile, disposable humidification flask for each incubation chamber. The purpose of the flasks is to humidify and direct gas into each chamber. The flask includes an internal inlet filter to prevent water from backflowing from the filled flask into the incubator's internal gas plumbing. The reference part number for reordering the humidification flasks is K-MINC-2000-HF.

IMPORTANT NOTE: Position the MINC+ in its final front-facing position before installing the humidification flasks.

WARNING: BIOLOGICAL HAZARD. Do not use a contaminated humidification flask in the MINC+. If the humidification flasks are used continuously, then they can be used until the water level reaches the MIN FILL line on the front of the flask. If the flasks are not used continuously, then the risk of biological contamination can increase. When the device is not in use, remove the flasks and replace them with new flasks at startup. Do not refill the humidification flasks.

Used humidification flasks are classed as infectious waste. Dispose of all infectious waste in a biohazard container or bag. Do not put sharps in biohazard bags.

IMPORTANT NOTE: Do not fill the flask above the max fill line. If overfilling occurs, the water from the flask may overflow into the chamber. If water overflows into the chamber, replace the whole humidification flask assembly (K-MINC-2000-HF).

IMPORTANT NOTE: To avoid water spillage inside the unit, remove the humidification flask assembly from the bay before moving the MINC+.

1

4.10.1 Humidification flask components

Two humidification flasks are included with the MINC+. They are supplied sterilized with ethylene oxide (EtO) in a sealed pouch. Each pouch contains the following items:


Figure 11: Humidification flask components


Description
1 Humidification flask reservoir
2 Humidification flask seal housing (including backflow prevention filter)
3 Humidification flask filter sleeve
4 Disposable funnel used when filling the flask
5 Humidification flask T-piece
6 Temporary support clip used when filling the flask

IMPORTANT NOTE: If the humidification flask packaging is damaged, the flask may be contaminated, so discard the flask and open a new one. To avoid contamination, do not open the humidification flask packaging until you are ready to install the flask.


4.10.2 Filling the flasks

WARNING: Failure to attach the flask T-piece correctly may result in a loss of gas flow to the chambers.


1. Under laminar flow conditions and using aseptic techniques, remove the components from the packaging.
2. Place the temporary support clip on the work surface and push the humidification flask down into the clip until engaged with the recess on the bottom edge of the reservoir. Slide the top latch back to the position shown below.


3. Insert the disposable funnel as shown below.


4. Fill the humidification flask reservoir to the MAX FILL line with sterile or filtered water (approximately 160 mL). Do not overfill.

1

5. Remove the funnel and then firmly push the flask T-piece onto the reservoir outlet spigot. Orient the outlet of the flask T-piece perpendicular to the flask reservoir to the left or right, to match the right or left humidification bay that it will be inserted into.

6. Place a finger on the support clip and lift the assembled humidification flask upward to detach it from the support clip. Check that no water is leaking from the flask.

4.10.3 Flask installation

1. Open the chamber lid. Ensure that the gas flow to the chamber is disabled before proceeding.
2. Hold the humidification flask upright with the flask T-piece oriented into the chamber in which the flask is to be installed. Hold the top latch back. While holding the end of the flask T-piece as high as possible to prevent water from spilling out, insert the flask into the humidification flask (HF) bay at a 10-15 degree angle, positioning the HF seal housing over the HF bay gas nozzle at the back of the HF bay (see **Figure 12**).
3. Rotate the flask, seating it into the HF bay until it is positioned so that the top latch is parallel with the surrounding surface of the chamber baseplate.
4. Insert the flask T-piece into the slot in the chamber baseplate. Ensure that the flask T-piece is straight and is pushed completely into the slot.
5. Push the top latch forward to lock the flask into position in the HF bay.
6. Gently push the rear edge of the HF seal housing into the HF bay gas nozzle to ensure that it is properly seated.
7. Gently push the top of the flask T-piece to ensure that it is seated.
8. Close the chamber lid.


Figure 12: Steps for humidification flask installation

4.10.4 Removing flasks for replacement

The humidification flasks must be replaced when the sterile water level reaches the MIN FILL line on the front of the flask. Or, if the MINC+ device is not in use, the flasks must be removed and replaced with new flasks at startup.

To remove the flask:

1. Ensure that the gas flow to the chamber is disabled before proceeding.
2. Open the chamber lid of the flask that needs to be replaced.
3. Slide the HF top latch back, and remove the flask T-piece from the slot in the chamber baseplate. The tab on the HF stand base may be used to assist in removing the T-piece tube from the chamber baseplate slot.
4. While holding the end of the flask T-piece as high as possible to prevent any water spillage, grasp the top latch or the front of the flask and rotate the flask out of the HF bay until the flask is at an angle of approximately 10-15 degrees.
5. Gently pull the flask back out of the HF bay to disengage the HF seal housing from the HF bay gas nozzle, ensuring that the HF filter sleeve remains with the flask to be disposed.
6. Dispose of the flask in a biohazard container or bag.

4.11 Activating the device

Turning on the MINC+

1. Connect the IEC plug on the power cord to the corresponding socket on the rear of the MINC+.
2. Plug the power cord into an earthed mains power supply socket.
3. Turn on the gas supply, ensuring that the gas pressure to the MINC+ is set to 150 kpa ± 15 kPa.
4. Switch on the power to the MINC+.

The MINC+ will perform a self-test and start the center console. The firmware version will briefly appear on each chamber display screen. The self-test sequence may take approximately one minute. After the self-test, upon first-time startup, the MINC+ should default to active mode. The incubator environment controls will show both heater chambers active and set to 37.0 °C, gas flow set to off, and the water level set to off. Refer to **Section 4.14**, if required, to set device parameters.

 **WARNING: ELECTRIC SHOCK HAZARD.** Internal circuitry is running whenever the incubator is connected to mains power, irrespective of whether the device is on or in standby mode.

 **IMPORTANT NOTE:** When the MINC+ has just been turned on, or if the temperature has been adjusted, the temperature alarm deactivates for 120 minutes while the device reaches stable conditions.

4.11.1 Unlock button function

The **Unlock** button sequence is as follows: beep on touch, and hold the button for two seconds until the keypad back-light illuminates to indicate that keypad is active. The keypad will remain active for 5 seconds after the last input. The **Unlock** button remains constantly backlit while the MINC+ is in active mode.

4.11.2 Entering standby and active modes

When the MINC+ is in standby mode, the center console and chamber display screens are off. When the MINC+ is in active mode, the center console touchscreen and chamber operational controls are active and may be used.

To enter standby mode, perform the following steps:

1. Touch the **Unlock** button until the keypad is active.
2. Touch and hold the **Power** button until the device shuts down and enters standby mode.

To enter active mode, perform the following steps:

1. Touch and hold the **Power** button until the MINC+ beeps and activates the chamber display screens, then release the **Power** button. The chamber status screen will start in approximately 6-10 seconds.
2. Wait approximately 40-60 seconds until the center console boot up sequence has completed.
3. The MINC+ is now in active mode.

4.12 Initial configuration

4.12.1 Assign a name to the incubator

By default, the incubator name of the MINC+ will be its serial number. It is important to give it a meaningful name that will identify it in a network of several MINC+ devices.

1

To assign a name, perform the following steps:

1. On the center console, tap the **Settings** tab.
2. Tap the **General** tab.
3. Tap the pencil icon in the text box adjacent to **Incubator Name** to add the name. Assign a name (maximum of 20 characters) that is meaningful to the lab or location of the device. The allowable characters for the name are: 0-9, a-z, A-Z, -, ., (,), #, @, &.
4. Tap **Save** to store the assigned name.

4.12.2 Assign a language

To assign a language, perform the following steps:

1. On the center console, tap the **Settings** tab.
2. In the **Settings** tab, tap **Language** and then slide up and down to scroll to highlight the desired language in the list.
3. Tap **Save** to keep the changes. The default language is English.

4.12.3 Network configuration (optional if connecting to DishTrace PC software)

Connectivity, either directly to the PC or through an ethernet switch, will require static IP addresses to be assigned. Alternatively, connection via an ethernet switch with a DHCP server will automatically assign IP addresses.

Dynamic allocation of the IP address

By default, the MINC+ device is automatically allocated an IP address when it is connected to a DHCP server. The allocated IP address can be seen in the **Settings** -> **Network** screen. DishTrace PC automatically finds each MINC+ device on the network.

If the default IP address 169.254.x.x is still seen on the **Network** screen, then the MINC+ device cannot access the DHCP server. Check the prerequisites **Section 3.3** or contact your network support service.

Static IP address allocation

This configuration must be selected if no DHCP server is available in the MINC+ network. All devices on the network should be on the same network subnet (the first three fields of the IP address). Most PCs support a subnet of the form 10.0.0.x, with x being any number from 2 to 253. The number in the x position is required to be unique to each PC and MINC+ device on the network.

Setting a static IP address on a PC (Sequence may vary depending on Windows configuration)

1. Via **Control Panel/Network and Sharing Center**, select **Change adapter settings**, then select **Ethernet**.
2. In the **Ethernet Properties** window, under the **Networking** tab, select **Internet Protocol Version 4 (TCP/IPv4)**.
3. Check **Use the following IP address** button.
4. In the **IP address** text box, enter the subnet 10.0.0.x, with x being a unique number for this device.
5. Enter the **Subnet Mask** 255.255.255.0.
6. Click **OK** to save and close the network sharing window.

Setting a static address on the MINC+

1. In the **Settings** tab, tap **Network**, and then tap the **Static IP** checkbox.
2. Edit the IP address and enter 10.0.0.x, with x being the a unique number for this MINC+. Tap **OK**.
3. Edit the **Subnet Mask** and enter 255.255.255.0. Tap **OK**.
4. Tap **Save**.

Adding MINC+ devices in DishTrace PC manually

1. In the **Incubators** tab, enter the static IP address assigned to the MINC+ device that is to be connected, and then tap **Add**. The incubator details should appear within 30 seconds.

4.12.4 Verifying the date and time

The date and time settings are important, because they are used when reports are generated by the MINC+. To verify that the date and local time are correct, perform the following steps:

1. In the **Settings**, tap the **Date & Time** tab, and then tap **Time Zone**. Verify that the local date and time that are shown are correct.
2. If the date and time are not correct, clear the **Network Time** checkbox and then edit the date and time as required.
3. Tap **Save** to save the changes.

4.13 Installing the DishTrace PC software (optional)

If the DishTrace PC software will be used, it can be installed as follows:

4.13.1 DishTrace PC software

DishTrace PC is dish-data software for use on a PC that has the Microsoft Windows operating system. It is used to view the status of the MINC+ incubators connected to the local PC network and to log information pertaining to the dishes in the incubators. See **Section 3.2** for details about the operating systems that can support it.

DishTrace PC can simultaneously support a recommended maximum of 50 MINC+ Benchtop Incubators, including 10 MINC Mini Incubators (K-MINC-1000). The DishTrace PC software functionality applies only to the MINC+ Benchtop Incubator. For connected MINC Mini Incubators, the software can only record their operating state. DishTrace PC can view and search all data records of registered dishes in any connected MINC+ incubators. The graphical history of any dish can be generated and viewed as well.

DishTrace PC can generate various reports for all MINC+ incubator activity. These reports are generated in an exportable PDF or CSV format and can be stored in other folders in the PC's file system or on external portable media such as a USB drive (FAT32 format).

DishTrace PC is described in detail in **Section 5.5** of this document.

4.13.2 Before you begin

IMPORTANT NOTE: Ensure that the installer has local administrative rights for the computer environment in which the software will be installed.

IMPORTANT NOTE: DishTrace PC is published by **William A. Cook Australia Pty. Ltd.** and Windows will automatically confirm the validity of the software during installation. If the publisher is different or cannot be verified by Windows, contact your service agent.

4.13.3 Installing the DishTrace PC software

The DishTrace PC software is on the USB drive provided. To install, perform the following steps:

1. Insert the Cook Medical USB drive into the computer's USB port and open the drive folder via Windows Explorer.
2. Open the **DishTrace_PC_xx.exe** file as an Administrator to start the installation process. The numbers in the xx position are the release version.
3. During first-time installation, the following software packages (high-level CBOM) will be automatically installed:

- USBXpress®
- Bonjour®
- .Net framework 4.5
- SQL Server®

These software packages are used by the installation application when it is executing. Subsequent reinstallations of DishTrace PC do not overwrite these files.


4.13.4 Add users to DishTrace PC users group

After DishTrace PC is installed, a new group called "DishTrace PC Users" is automatically created. All non-administrative users of DishTrace PC must be members of the "DishTrace PC Users" group to allow access to the DishTrace database. All users must have an existing valid log-in account.

To add entries to the "DishTrace PC Users" group, the administrator must perform the following steps:

1. Click on the Windows Start icon in the bottom left of the screen, type "**Computer Management**", and press **Enter**.
2. From the left-hand panel, expand the **Local Users and Groups** entry.

3. Click on the **Groups** folder to display all known groups.
4. Right-click on the **DishTrace PC Users** group and select **Add to Group...**
5. Click on the **Add** button at the bottom of the screen.
6. Enter the name of the new user into the panel that opens up similar to the display below:

7. Click on the **Check Names** button to select the correct user and then click on **OK**.
8. The original name entry will be updated to the correct valid user object name. Click **OK**.
9. The new user name will now be seen in the **Members** panel.
10. Repeat from step 5 for all users to be added to this group. When finished click **OK**.
11. Close the **Computer Management** utility.

4.13.5 Export DishTrace PC certificate

The network can only be secured once DishTrace PC has been installed and **Section 4.12.3** has been completed.

The security of the network is achieved using digital certificate technology. All connected PCs and MINC+ devices will generate their own public/private key pairs and create their own certificates. These certificates will be exported manually from each connected PC to the MINC+ devices. The corresponding export of the MINC+ device certificate to the PCs will be performed automatically.

Export DishTrace PC certificate to USB drive

After DishTrace PC is installed, a certificate for the PC will be automatically created and placed in the certificate repository. This certificate is a self-signed certificate that will need to be imported into every MINC+ connected to the PC.


The Windows operating system utility **mmc** will be used to export the PC certificate to a USB drive. Ensure that a USB drive is connected to the PC before the next steps.

To export the PC certificate with mmc, perform the following steps:


1. On the PC, click on the Windows Start icon, type "**mmc**" in the search dialog box, and press **Enter**. Click **Yes** to allow **mmc** to make changes to the device. The following screen will be displayed:

2. From the **File** menu select **Add or Remove Snap-ins**. From the **Available snap-ins:** panel, select **Certificates** and then click the **Add** button as shown below:

3. Select **Computer account** and click on the **Next** button.
4. With the default selection of **Local computer**, click on the **Finish** button.
5. Click on the **OK** button to complete to snap-in addition.
6. Expand **Certificates - Local Computer > Personal > Certificates** as shown below:

7. Right-click on the client certificate (PC name in **Issued To** and **Issued By** fields) and select **All Tasks > Export...**. This is a self-signed certificate created by the installation of DishTrace PC.
8. Press **Next** in Certificate Export Wizard to go to Export Private Key.
9. Select **No, do not export the private key** and click on **Next**.
10. Select **DER encoded binary X.509(.CER)** format and press the **Next** button.
11. Browse to the USB drive and provide the desired file name for the certificate (e.g., PC name) and click on **Save**.
12. Select **Next** to save the certificate file to the USB drive and press the **Finish** button to export the client certificate. The Certificate Export Wizard should indicate the export was successful. Click **OK**.
13. Exit the snap-in utility. Do not save the console settings.
14. Remove the USB drive containing the PC's certificate.

4.13.6 Import DishTrace PC certificate into MINC+

The PC certificate must be imported to each MINC+ to allow a secure communications link to be established.

Using the USB drive containing the PC's certificate created as per **Section 4.13.5**, perform the following steps:

1. Insert the USB drive into the front panel of the MINC+ device.
2. Navigate to the **Settings > Network** screen.
3. Ensure the IP address shown is valid as outlined in **Section 4.12.3**.
4. Press on the pencil icon on the right-hand side of the **Certificate** text box.
5. Press the **Import New Certificate** button and confirm the action by pressing the **Import** button.

If you press the magnifying glass symbol on the right-hand side of the certificate text box, a list of PC certificates already residing on the MINC+ device can be seen in a scroll box.

4.13.7 Establishing a secure network


Once the DishTrace PC certificate has been imported into all of the MINC+ devices, a secure network can be established.

The following steps must be completed when DishTrace PC is run by an administrator. Start DishTrace PC by right-clicking the DishTrace PC start icon or selecting DishTrace PC from the Cook Medical program menu, and then selecting **Run as administrator**.

At startup of the DishTrace PC program, a summary of the status of the connected MINC+ devices will be shown at the top on the screen as one of the following:

- **Unlocked Padlock N** – Where N is the number of detected MINC+ devices that do not have a secure connection.
- **Online: N** – Where N is the number of MINC+ devices which have a secure connection.
- **Offline: N** – Where N is the number of MINC+ devices which have a secure connection but are not currently operating.

To create a secure connection to the MINC+ devices, click the **Unlocked Padlock** symbol. This will display a screen as shown below:

The unsecured MINC+ devices are shown in the **List** box at the bottom of the screen. To create a secure connection to incubator COOK001, for example, perform the following steps:

1. Click on the (+) symbol in the **Connection** column. DishTrace will initiate a secure connection with the MINC+ device and show a dialog requesting that the user accept a certificate sent by the MINC+ device.
2. Accept the certificate. After you accept the certificate, the Windows operating system certificate utility will start and ask whether it should trust the certificate and store it in the PC's certificate repository.
3. Choose to trust the certificate. The certificate will be stored, and a secure communication link will be created.
4. Repeat the above steps for all MINC+ devices in the list. As a secure connection is created for each MINC+ device, the device will be moved to the incubator list.

4.13.8 Run DishTrace PC

To run the DishTrace PC software, do one of the following:

- Click on the DishTrace PC shortcut on the desktop.
- Click on the Windows operating system start button, scroll down to the Cook Medical folder, and click on the DishTrace PC entry.

The **Incubators** screen will appear, and after 30 seconds, it will detect any MINC+ devices connected to the network. The **Incubators** screen is shown by default.

4.14 Chamber operation

4.14.1 Enable incubation chamber keypad

Before any of the MINC+ incubation chamber settings can be altered, the incubation chamber keypad must be enabled.

To enable the incubation chamber keypad, perform the following steps:

1. Touch the **Unlock** button. A beep will sound to indicate that the MINC+ has recognized the button touch.
2. Continue to touch the **Unlock** button for two seconds and the keypad backlight will then illuminate.
3. Both left and right keypads are now active.

1

4.14.2 Turning the incubation chamber heater on or off

1. Touch the **Unlock** button for two seconds. The keypad will beep on touch and the backlight will come on when the keypad is active.
2. Touch and quickly release the **Enable/disable incubation chamber heater** button. This will toggle the chamber on or off, depending on the current state. The state will be displayed on the chamber's display screen.
3. The new changes to the heating state will be saved automatically.

IMPORTANT NOTE: When the chamber heater is turned on, the temperature display will be yellow as it warms up to within 0.2 °C of the temperature set point. When the temperature is within 0.2 °C of the temperature set point, the display will turn green.

4.14.3 Incubation chamber temperature selection

When the device is first turned on, it will default to a temperature of 37.0 °C. The incubation chamber displays will show both the temperature set points and actual temperatures of each chamber in degrees Celsius (°C).

Adjusting the incubation chamber temperature set point.

1. Touch the **Unlock** button for 2 seconds. The keypad will beep on touch and the backlight will come on when the keypad is active.
2. Repeatedly touch and release the **Increase/decrease chamber temperature set point** buttons for that chamber until the desired temperature set point is displayed on the chamber display screen. The set point is displayed on the chamber display screen in white text immediately above the increase/decrease temperature set point buttons. The temperature can be adjusted in 0.1 °C increments, within a range of 35.0 °C to 40.0 °C. The new changes to the state are saved automatically.

4.14.4 Lid openings

Opening the incubation chamber lid causes the open incubation chamber lid icon to display on the chamber display screen. A warning beep will emit from the MINC+ approximately every 15 seconds to alert the user that the lid is open. After the lid is closed and the latch is moved to the closed position, the chamber temperature will display on the chamber display screen, the warning beep will stop, and an entry will display on the center console to log the event.

Opening the lid

1. Unlock the latch to the fully open position.
2. Use the front lip to lift the lid open. The lid will maintain its opening position from approximately a 45° angle to the fully open position.

IMPORTANT NOTE: When a lid is opened and then closed, or when the gas flow is started, the chamber is automatically purged to quickly reestablish the appropriate gaseous environment. The gas-purge flow rate is preset. Each chamber performs the gas-purge cycle independently of the other incubation chamber.

Closing the lid

1. Ensure that the latch is in the **OPEN** position (pointing towards the front of the device).
2. Hold the front lip and lower the lid until fully closed.
3. After the lid is closed, move the latch from the **OPEN** to **CLOSE** position to complete the closure of the chamber.
4. The open incubation chamber lid icon will then disappear from the chamber display screen, the incubation chamber temperature will display, and the gas-purge cycle will begin for that chamber.


Figure 13: Closing the lid

4.14.5 Incubation chamber gas flow

The rate of gas flow into each incubation chamber is not adjustable and is calibrated at sea level to a normal flow rate of 25 mL/min and a gas purge flow rate of 300 mL/min. The gas supply pressure at the MINC+ gas inlet is 150 kPa ± 15 kPa.

IMPORTANT NOTE: The MINC+ measures the gas flow rate internally prior to the gas entering the humidification bay. There is no mechanism to verify that the rate of gas flow into each chamber is correct. Users must visually ensure that the humidification flask and Braided Connecting Tube are correctly installed so that the gas flows into each chamber. To verify that the flow is unobstructed, look for rising bubbles at the front of the humidification flasks. Each chamber's heater must be on to enable gas flow into the chamber.

Turning the gas flow ON

1. Touch the **Unlock** button for 2 seconds. The MINC+ will beep on touch and the backlight will come on when the keypad is active.
2. Touch and release the **Enable/disable gas flow** button. The chamber's status screen will display the gas purge icon, indicating that the gas flow is enabled. The MINC+ will begin a gas purge cycle for approximately 3 minutes and then go into normal flow operation. The state will be saved automatically.

Turn the gas flow OFF

1. Touch the **Unlock** button for two seconds. The MINC+ will beep on touch and the backlight will come on when the keypad is active.
2. Touch and release the **Enable/disable gas flow** button. The chamber status screen will display the **Gas flow disabled** icon indicating that the gas flow is disabled for that chamber. The state will be saved automatically.

4.14.6 Turning on water level monitoring for the humidification flask

1. Touch the **Unlock** button for two seconds. The MINC+ will beep on touch and the backlight will come on when the keypad is active.
2. Touch and release the **Enable/disable humidification flask water level monitoring** button. The water level indicator icon will show that the water monitoring is enabled and the amount of water in the humidification flask. The state will be saved automatically.

4.14.7 Media warming chamber (MWC)


Figure 14: Media warming chamber

Setting the media warming chamber temperature.

1. Touch the **Unlock** button for two seconds. The MINC+ will beep on touch and the backlight will come on when the keypad is active.
2. Lift the lid of the MWC.
3. Press the MWC temperature-adjust buttons until the desired temperature set point is displayed on the display screen in the MWC. The state will be saved automatically.

IMPORTANT NOTE: The MWC will stop regulating to its temperature set point if both left and right incubation chambers are turned OFF. If the media warming chamber temperature is out of range, open and close the MWC lid or toggle both incubation chambers OFF and then ON to reset it.

4.14.8 First time commissioning run

Ensure that the following actions are carried out before the first use:

- Both chamber heaters and gas flow to each chamber, should be turned on for 24 hours to ensure that any residual outgassing of components is complete.
- Install the humidification flasks in both incubation chambers and turn on the gas flow for at least 4 hours to saturate the water with CO₂.
- Use a culture media that contains phenol red indicator to test each chamber for pH maintenance. After overnight incubation, look to see if the phenol red indicator is the correct color (salmon pink).
- If the DishTrace PC software has been running during the commissioning run, download the PDF file for each device, and check that the performance and function of the device are satisfactory.

The MINC+ is now successfully installed and commissioned.

4.15 Installation and setup checklist

- All items have been supplied.
- The packaging has been safely stored for future use.
- All nonsterile items have been removed from packaging.
- The power cord is the correct format for the region.
- The MINC+ has been placed in a suitable location.
- The appropriate gas mixture has been determined.
- An appropriate gas regulator has been sourced.
- The gas regulator has been set to 150 ± 15 kPa.
- The gas connections have been made and checked.
- The humidification flask has been aseptically filled and fitted.
- External alarms have been connected.
- The MINC+ has been switched on.
- The MINC+ has been configured.
- The gas flow to the chambers has been activated.
- The MINC+ has been left to operate for 24 hours to ensure that any residual outgassing of components is complete.

Check the following for the DishTrace PC software (optional)

- The installation USB flash drive has been supplied.
- A computer with the appropriate specification is available for using the DishTrace PC software.
- The DishTrace PC software has been installed on the desired computer.
- A direct or LAN network connection between the incubator and PC exists.
- The MINC+ devices have been activated.
- Users group has been populated and secure connectivity established between each PC with DishTrace PC and each MINC+ on the closed network.

5. Device operation

1

This section describes how to use the MINC+. When operating the MINC+ for the first time, follow the instructions below to ensure that the device performs reliably.

Only the clinician or embryologist can evaluate the clinical factors for each patient and determine if the use of this device is indicated. The clinician or embryologist must decide on the specific technique and procedure that will accomplish the desired outcome.

⚠️ WARNING: To guarantee safe operation, properly care for and maintain the device and perform regular checks to confirm that the device is functioning correctly. New and repaired products must be prepared and tested according to the MINC+ instructions before use.

👉 IMPORTANT NOTE: Ensure that gas supply is available to the MINC+. Each chamber's heater must be on to enable gas flow into the chamber.

5.1 Before use

Prior to commencing the culturing procedure, perform the following steps:

- Use the temperature-selection touchpads to select the desired chamber temperature (refer to **Section 4.14.3**).
- Fill the humidification flasks with sterile water and install them as detailed in **Section 4.10**.
- Confirm that the humidification flasks and gas tubes are correctly installed.
- Activate the required incubation chambers.
- Activate the gas flow (see **Section 4.14.5**).
- Ensure that gas bubbles can be seen in the humidification flask during both the gas purge and normal flow cycles.
- Wait a minimum of four hours before using the incubator to allow for equilibration.

5.2 Inserting culture dishes

The MINC+ has been designed to be compatible with Nunc four-well dishes and Nunc and BD Falcon 35 mm and 60 mm single round culture dishes. Equivalent dishes maybe used, however the user must ensure correct fitment as described below.

Four-well dishes or single round dishes can be placed on the chamber base. Ensure that they sit securely in the grooves designed to accept them. There should be direct contact between the base of the dish and the surface of the chamber base. The height of the dish should not interfere with the lid closure.

5.3 Patient identification plates

Twelve magnetic patient identification plates are provided with the MINC+. Use the magnetic patient identification plates on the chamber lids to record culture dish contents with a marking pen. Pen markings can be removed with an alcohol solution.


Figure 15: Patient identification plates

5.3.1 Placing test tubes inside the media warming chamber (MWC)

IMPORTANT NOTE: The MWC is not a gassed environment and is heated independently of the incubation chambers.

1

The removable MWC cradle can hold six (6) test tubes. The MWC cradle is compatible with BD Falcon tubes (17 mm diameter, 100 mm long) with the test tube caps in either the raised or lowered position. Equivalent test tubes can be used, however the user is responsible to ensure the tube fitment within the MWC cradle and the MWC lid can fully close.

The MWC lid must be kept closed when warming media. The MWC temperature alarm is deactivated when the lid is opened, and remains so for 120 minutes on lid closure to allow media to equilibrate without causing a temperature out of range alarm.

The colors of the three lights on the front panel below the MWC lid indicate the chamber temperature status as follows:

- Green indicates that the MWC is within ± 0.5 °C of the set point temperature.
- Amber indicates that the MWC is equilibrating to temperature when the alarm is deactivated or within ± 0.5 °C to ± 1 °C of the set point temperature during normal operation.
- Red indicates that the MWC is out of range of the set point temperature by more than ± 1 °C during normal operation. Refer to **Section 5.6.6** for alarm condition details.

DishTrace

The DishTrace software can be used to assign logged data to identified dishes in the MINC+ incubation chambers. DishTrace has two components:

- DishTrace MINC+ – Installed on the incubator and accessed via the center console.
- DishTrace PC – Software that can be installed on a PC to access and make available dish information to connected MINC+ devices.

5.4 DishTrace MINC+

The DishTrace MINC+ enables environmental and placement information to be logged via the center console touchscreen to a dish record and to record when dishes are placed in and taken out of the incubation chambers.

Dish logging actions in DishTrace MINC+ are performed under the **Contents** menu option on the center console. This menu allows users to add new dishes, edit dish records, assign dishes to a location, and review dish status and settings.

Dish data entered via the DishTrace MINC+ can also be viewed with the DishTrace PC software described in **Section 5.5**.

WARNING: DishTrace MINC+ is only an aid to logging and viewing the location of dishes within the MINC+. DO NOT use DishTrace to replace laboratory identification protocols for dishes and their contents. Always use dish identity checking processes such as standard laboratory protocols (including, for example, individual dish labeling and identification verification processes) to minimize the risk of incorrect dish identification.

IMPORTANT NOTE: Dishes can be assigned to a location in the incubation chambers only via the center console of the MINC+ in which the dish is being placed. They cannot be assigned to a specific location in a MINC+ via the DishTrace PC software (see **Section 5.5**).

5.4.1 Resetting the center console touchscreen (DishTrace MINC+)

If the center console touchscreen does not respond at any time, it can be reset without interrupting the incubation chambers.


- On the back of the MINC+, locate the DishTrace MINC+ reset button. Refer to the diagram under **Section 4.4** for more information.
- Press and release the button once and wait for the center console to reset.

5.4.2 DishTrace MINC+ setup

IMPORTANT NOTE: It is necessary to set the **Time and Date**, because event logs are timestamped with the time settings in the MINC+. Refer to **Section 4.12.4** for more information.

The **Contents** tab displays the dishes assigned to locations in each incubation chamber. The following describes dish management actions that can be performed on the center console. Refer to **Section 5.4.8**.

to automatically display dish assignment options on lid opening. Dish locations can have three different statuses: **Unassigned**, **Assigned/Checked In**, and **Assigned/Checked Out**:


Dish status	Description	Icon
Unassigned	Empty black square or circle (varies according to the type of dish configuration selected under the Settings screen) – No dish is assigned to the location.	
Assigned/Checked In	Solid light green square or circle – A dish has been added to the incubator chamber at the highlighted location.	
Assigned/Checked Out	Outlined circle or square – This space is assigned to a dish that has been temporarily removed from the chamber.	

5.4.3 Assigning dishes to a chamber

To assign a new dish to a chamber

- Tap an unassigned dish location on the **Contents** screen. The following menu appears:

- Tap **Assign New Dish** to allocate a new dish to a position in the chamber. A new dish is a dish that has not had its details entered into the MINC+.
- Enter dish identification: **First Name**, **Last Name**, **Patient ID**, and/or a **Dish ID**. The chosen **Primary Dish ID** must be entered as minimum to correctly display the dish identification. Refer to **Section 5.4.13** to set **Primary Dish ID**. Other fields are optional or maybe left blank.
- Press **Save** to continue. The new dish has now been assigned and checked into the incubation chamber.

To assign an existing dish


An existing dish is a dish that has already had its details entered using the DishTrace PC software (see **Section 5.5.11**).

- Tap an unassigned dish location on the center console.
- Select **Assign Existing Dish**.
- Select the desired dish from the list of existing dishes, and assign it to the new location.

5.4.4 Checking out a dish

Use this function to check out a dish from its assigned location when a dish is physically removed from the incubator chamber.

- Tap the dish location of a dish that is checked in. The following menu will appear:


- Select **Check Out Dish**.

3. The dish is checked out from the current position. The dish status displayed will change to an outline of a circle or square, depending on the dish type that was selected.

5.4.5 Check in a dish to its preassigned location

1. Tap on the empty dish location where the dish was previously assigned. The following menu will appear:

2. Select the **Check In Dish** option from the drop-down list. The dish is now checked in to its previous location. The dish status displayed will change to a solid light green filled circle or square, depending on the dish type that was selected.

5.4.6 Unassign a dish position

1. Tap an **Assigned/Checked In** or an **Assigned/Checked Out** dish location.
2. Select **Unassign Dish Position** from the drop-down menu. A second screen will display with the message '**Are you sure you wish to unassign dish: x from the selected dish position?**'.
3. Select **Unassign Dish Position** to confirm. The dish status displayed will change to a empty black circle or square, depending on the dish type that was selected.

5.4.7 Edit dish record

1. Tap a dish location that is currently assigned to a dish, even if the dish is checked out.
2. Select **Edit Record**. The **Edit Record** screen will appear:

3. Update the required details.
4. Select **Save** when the details are correct.

View and export a patient record

1. Tap a dish location that is currently assigned to a dish, even if the dish is checked out.
2. Select **View Record**. The **Dish Record** screen appears. The dish temperature history is displayed graphically in the **Chart** tab. The **Alerts/Events** tab displays a list of alerts and events for the dish:


3. To save a copy of the patient record as a report, insert a USB drive into the USB port in the front of the MINC+ to save the report. Select **Export Report**, and choose the option to export the report as a PDF or CSV file.

IMPORTANT NOTE: The USB drive that is inserted into the front-panel USB port for reports and backup must be formatted to FAT32.

1

5.4.8 Lid open alert – Dish assignment

If **Lid open alert** in the **Settings - General** tab is checked, the dish assignment menu options to check dishes in and out of each chamber will be enabled, and is displayed as shown in the figure below. The dish check in/out menu will be displayed on the opposite side of the screen relative to the chamber that was opened, allowing the chamber contents to remain displayed.

Assigning a new dish

1. Tap an unassigned dish location. A check mark will appear on the selected location and the dish assignment options will appear.
2. Select **Assign New Dish**.
3. Enter the dish identification and/or patient details: **First Name**, **Last Name**, **Patient ID**, and/or **Dish ID**. Any field can be left blank; however, at least one field must be filled in to create the record.
4. Press **Save** to continue. The new dish has been assigned and checked into the incubation chamber.

Checking in one or more dishes

1. Tap any checked out dish locations. A check mark will appear on the selected locations.
2. Select **Check Dish In**.

Checking out one or more dishes

1. Tap any assigned dish locations. A check mark will appear on the selected locations.
2. Select **Check Dish Out**.

5.4.9 Records screen

The **Records** screen enables users to search, view, and export records and reports. It also provides a scalable graphical representation of the temperature logs for each of the incubation chambers.

To view a graph of the incubation chamber's history.

1. Tap the **Chambers** icon. The following screen will appear:

2. Tap the plus (+) icon on the top right of each incubation chamber graph to show only the selected chamber. Tapping the magnify +/- will zoom/expand the time scale.
3. Use the < and > icons on the top right-hand corner to scroll through the recorded history.

To view logged alerts and other events

1. Tap the **Alerts/Events** icon. The following screen will appear showing all events and alerts for the MINC+:

Alerts/Events			<input type="checkbox"/> Display Alerts Only
Description	Occurrence	Acknowledged	
Right Chamber lid closed	16-Dec-2021 20:42:46		
Right Lid Open Alert	16-Dec-2021 20:42:28	16-Dec-2021 20:42:47	
Right Chamber lid open	16-Dec-2021 20:42:12		
Left Chamber lid closed	16-Dec-2021 20:42:05		
Left Chamber lid open	16-Dec-2021 20:42:00		
Right Chamber lid closed	16-Dec-2021 20:41:48		
Right Chamber lid open	16-Dec-2021 20:41:44		
Right Chamber lid closed	16-Dec-2021 19:57:40		
Right Chamber lid open	16-Dec-2021 19:57:36		
System Restart	16-Dec-2021 11:23:18		

1

2. Tap the up and down arrows to scroll through the occurrences.
3. To display alerts only (i.e., omitting events such as lid openings) tap the **Display Alerts Only** checkbox.

5.4.10 To search for a dish or patient record and generate a dish report

1. Ensure that the USB drive is inserted into the front USB port below the center console touchscreen.
2. Tap the **Reports** icon. The following screen will appear:

3. Tap **Search for Dish/Patient Report**.
4. Enter a search term. Records can be searched by **Last Name, First Name, Dish ID, and Patient ID**.
5. Select **Export Report**.
6. Select **PDF** or **CSV** format. A report will be generated for the specific dish and saved to the external USB drive.

5.4.11 To generate an incubator report

1. Ensure that a USB drive is inserted into the front USB port below the center console touchscreen.
2. Tap the **Reports** icon. The following screen will appear:

3. Tap **Generate Incubator Report**.
4. Select the date range for the report by tapping the calendar icon. The maximum date range is three months.

5. Select **Export Report**.
6. Select **PDF** or **CSV** format. A report of all the MINC+ logs, events, and alerts within the specified date range will be generated and saved to the external USB drive.

5.4.12 Status screen

The status screen provides a current snapshot of the left and right incubation chambers, the media warming chamber (MWC), the gas flow status, and the water level of each humidification flask. There are no available actions on this screen to change the operational settings of the MINC+ via the center console.

5.4.13 Settings screen – General tab

The **Settings** screen enables the user to change general settings, dish configuration, network settings, date and time settings, manage the database and display incubator information.


General tab

This screen displays the incubator name and dish settings.

- Device Settings** – This allows the MINC+ to be assigned a name. Select the pencil icon to add the incubator name. Assign a name (maximum of 20 characters) that is meaningful to the lab or describes the location of the device. The characters allowed for the incubator name are: 0-9, a-z, A-Z, -, ., (,), #, @, &.
- Lid Open Alert** – If this option is selected, the **Dish Assignment** menu option is enabled. Refer to **Section 5.4.8**.
- Beep on Touch** – If this option is enabled, then when the center console is touched, a beep will confirm user actions.
- Language** – Select language from the drop-down list. The default is English.
- Primary Dish ID** – This option determines the text that is displayed for the dish on the **Contents** screen. This must be defined to correctly identify dishes.
- Secondary Dish ID** – This option sets the text for secondary information that is displayed for each dish on the **Contents** screen. This information is not displayed for small (35 mm) dishes.
- Name Format** – Sets the order in which the patient's first and last names are displayed.

5.4.14 Configuration tab

This screen enables the dish configuration to be set for each incubation chamber. Three options are available: 60 mm, 35 mm, and square four-well dishes.

To set the dish configuration:

1. Press **Edit** on the top right-hand side. The **Chamber Configuration** screen will open.
2. Press the left and right arrows < > on each row to change the dish type for each chamber.
3. Press **Save** to save the changes. A second screen will display to confirm the changes. Press **OK**.

5.4.15 Network tab

This screen displays the network settings for the MINC+. The network settings may be changed from this screen. It is recommended that you consult your local network administrator for help with these settings.


The small icon under the incubator name indicates that there has been, or is, connectivity to the DishTrace PC software. An orange display indicates a lost connection to the DishTrace PC software. A green display indicates the MINC+ is connected to the DishTrace PC software.

- **Static IP** – Allows for the network information to be set manually if no DHCP server is present.
- **Dynamic IP** – Set by default to allow a DHCP server to allocate an IP address to the MINC+.
- **IP Address** – The incubator's own IP address, which is used to connect to the DishTrace PC software.
- **Subnet Mask** – Restricts settings in the IP address.
- **Certificate** – Displays and manages imported PC certificates residing on the MINC+

5.4.16 Date and time settings

This screen displays the date and time settings for the MINC+:

- **Network Time** – Select this option to synchronize the incubator's time to the network's NTP Server.
- **12 Hour Format** – Select this option to change the time from 24-hour to 12-hour format.

- **Date, Time, and Time Zone** – Tap the pencil icon for each option to change the date, time, and time zone. Press Save to confirm the changes.
- **Date & Time Format** – Tap the drop-down boxes to change the date and time formats.

5.4.17 Database and backup information

This screen displays the database operations and management options. The MINC+ database can hold up to 2 GB of data, which under typical operating conditions is approximately 12 months of information. Data can be deleted (purged), backed up, or restored.

IMPORTANT NOTE: The MINC+ will display a warning when the database is approximately 80% full. If the database becomes full, any new data will automatically overwrite the oldest data in the database. Regular backups of the database are recommended. It will take approximately 3.5 minutes to backup a full 12 months of data.

5.4.18 Purging data

This function deletes selected data in the MINC+. This operation should be used if the database is reaching its capacity.

1. To purge data, tap the **Settings** screen and then tap the **Database** tab.
2. Tap **Purge**.
3. Select the date range to delete:
 - Oldest 3 months
 - Oldest 6 months or
 - All data
4. A confirmation message will appear: '**What would you like to do?**' **Backup & Purge**, **Purge** or **Cancel**. Select **Purge** or **Backup & Purge** to continue.
5. Another confirmation message will appear: '**Are you sure you want to delete the stored data for the selected duration?**'
6. Press **OK** to confirm the data purge and return to the **Database** tab. The USB drive can now be removed.

5.4.19 Backup data to an external USB drive

This function backs up the database (dish and patient records) and MINC+ configuration settings.

1. Insert a USB drive into the USB port at the front of the MINC+.
2. Tap the **Settings** screen and then tap the **Database** tab.
3. Tap **Backup**. The device will begin to save the data to the USB drive.
4. When the report has finished saving, a message will appear: '**Database successfully saved to USB drive**'.
5. Press **OK** to return to the **Database** tab. The USB drive can now be removed.

IMPORTANT NOTE: The database backup file can be restored only to the same MINC+ incubator. It cannot be opened on a computer or another MINC+ incubator.

5.4.20 Restore data settings from an external USB drive

Use this function to restore a database from a saved backup on an external USB drive.

Data cannot be transferred to another incubator and can be restored only to the original device it was backed up from.

1. Insert a USB drive into the USB port at the front of the incubator.
2. Tap the **Settings** screen and then tap the **Database** tab.
3. Tap **Restore**. A message will appear: '**Restore operation will overwrite the current database. Do you wish to proceed anyway?**'
4. Tap **Yes**.

5. A message will appear: **'Database restore has been completed successfully'**.

6. Press **OK** to return to the **Database** tab. The USB drive can be removed from the port.

5.4.21 Upgrade DishTrace MINC+ software

The DishTrace MINC+ software may occasionally be required to be updated. DishTrace MINC+ is published by **William A. Cook Australia Pty. Ltd.** and MINC+ will automatically confirm the validity of the software during installation. If the publisher is different or cannot be verified by the device, contact your service agent. An option to downgrade the DishTrace MINC+ software to the previous version is available, however should only be used under instruction from Cook Medical or an authorised service agent.

To update the DishTrace MINC+ software:

1. Insert the USB drive containing the software upgrade files into the USB port on the front of the MINC+.
2. Tap the **Settings** screen and then tap the **Database** tab.
3. Next to the **Upgrade Software**, tap **Upgrade**. A message will appear asking for conformation to upgrade the software.
4. Tap **Yes** to continue with the upgrade. The MINC+ will reboot once the software update is complete.

5.4.22 Restore factory settings

 **WARNING:** Selecting **Restore Factory Settings** will delete all logged data and restore default settings.

The **Restore Factory Settings** function deletes all dish and patient records. Performing this function before requesting a service operation is recommended. This function also resets the DishTrace MINC+ function settings to the default.

1. Tap the **Settings** screen and then tap the **Database** tab.
2. Tap **Reset**, which is next to **Restore Factory Settings**. A message will appear: **'All data including user settings, database will be deleted. Do you want to proceed anyway?'**
3. Tap **Yes** to continue. The device settings will be restored, and the center console will restart automatically.
4. After the device restarts, the date and time and possibly network settings will need to be reconfigured. See **Section 4.12** for details.

5.4.23 Information – Device settings

This screen shows information about the incubator. There are no user actions on this screen.

- **Incubator Name** – Incubator name.
- **Serial No.** – Incubator serial number.
- **BIOS Version** – Operating system BIOS version for center console.
- **OS Version** – Operating system version for the center console.
- **DishTrace MINC+ Version** – DishTrace MINC+ software version number.
- **FW Version** – The firmware version for the MINC+ core control system.
- **Manufacture Date** – The date on which the device was manufactured.
- **Last Serviced Date** – The date on which the device was last serviced.

Contents	Records	Status	Settings	90987680014
General	Device Settings			
Configuration	Incubator Name	90987680014		
Network	Serial No.	90987680014		
Date & Time	BIOS Version	BTI_3.1.0.314_X64		
Database	OS Version	4.19.80-intel-pk-standard		
	DishTrace MINC+ Version	1.1.0.2		
Information	FW Version	1.1.0.1		
	Manufacture Date	24-Apr-2019		
	Last Serviced Date	06-Dec-2021		

5.4.24 DishTrace MINC+ screen saver

The screensaver is a Cook Medical image that appears on the center console touchscreen every 6 hours when the touchscreen is not used. It appears on the center console touchscreen for only a few seconds to reduce image burn-in. The screensaver cannot be modified or disabled. The two outer incubator chamber screens also display a black screen as a screensaver for a few seconds, every 6 hours, if the chamber controls have not been used. The screensavers will not affect the operational aspects of the MINC+.

5.5 DishTrace PC software

The DishTrace PC software is a comprehensive software tool that shows the status of multiple MINC+ incubators connected to a local network. It shows a list and the status of all the active MINC+ devices currently connected as well as any previously connected incubators that may be currently offline. Offline status may indicate that the device is turned off or has been sent away for servicing.

The DishTrace PC software must be installed in accordance with **Section 4.13** to establish a secure connection between the PC with DishTrace PC installed and each MINC+ on the network. Once the connection process is completed, each device will be automatically detected when turned on and added to the list of active incubators. If a MINC+ is later turned off, it will still appear on the list but will be shown as inactive on the network.

For each incubator, DishTrace PC displays the status of the following:

- Incubator chambers
- Humidification flasks
- Gas
- Dish records
- MINC+ administration information

5.5.1 Using the DishTrace PC software

To run DishTrace PC, do one of the following:

- Click on the **DishTrace PC** shortcut on the desktop.
- Click on the Windows operating system start button, scroll down to the Cook Medical folder, and click on the **DishTrace PC** file.

The main screen will appear, and the **Incubators** screen will be shown by default.

5.5.2 Incubators screen

The **Incubators** screen is shown by default when DishTrace PC is launched. Any active incubator on the network will be automatically connected and shown.

5.5.3 Adding incubators

If an incubator fails to appear on the list after several minutes have elapsed, it can be manually added:

1. On the top left-hand side, in the MINC+ address bar (as shown in the screenshot above), enter the IP address of the incubator you want to add. (To find the IP address on the MINC+ incubator, go to **Settings > Network**.)
2. Press the **Add** button. After the incubator is found, a message will say that the incubator has been added successfully. Click **OK**.

5.5.4 Incubator display


The total number of online and offline incubators is shown on the top right-hand side. The incubator list can be shown as a graphical representation or a list view, either of which can be selected by clicking any of the icons next to the **VIEW** label.

To view further details of an incubator, double-click on any incubator entry. The following screen will appear:

5.5.5 To edit incubator settings

- Double-click on the incubator to edit.
- Click on the edit button on the top right-hand side (shown in the screenshot above).
- A new screen will open on which the following settings can be edited:
 - Edit Incubator Name:** Edit the incubator name in this textbox, and click **Save** when finished. The characters allowed for the incubator name are: 0-9, a-z, A-Z, -, ., #, @, &. (Maximum of 20 characters.)
 - Remove Incubator:** This option will remove the incubator from the incubator list and cause the incubator to no longer be autodetected by DishTrace PC. To re-list an incubator that has been removed, manually enter its IP address as described in **Section 5.5.3**.
 - Resync Incubator Data:** Any major change to the MINC+ database, such as restore data from back-up, purge, or restore factory settings, will require the DishTrace PC database to be resynchronized with the MINC+ database. Clicking **Resync** will erase the existing data from the selected MINC+ on the DishTrace PC database. The incubator will require the secure connection to be established as per **Section 4.13.7**. Upon resynchronizing, data from the selected MINC+ device's database will transfer to the DishTrace PC database. This data transfer may take some time and will happen in the background.
- To go back to the main **Incubators** screen, click the back arrow on the top left-hand side.

5.5.6 View Contents, History, Events and Information

To view the contents, history, events and information of an incubator, double-click on the incubator. The following four tabs will appear:

Contents, History, Events, and Information:

Contents

The **Contents** tab shows the dish layout for each incubation chamber and its contents.

- Clicking on the green **List View** icon (top right-hand side of each chamber field) expands the opposite column to show more detail on each dish, such as **Patient Name**, **Patient ID**, **Dish ID**, **Position**, **First Check In**, and the **Stage** of each dish.

History


1

The **History** tab shows the history of incubator temperature and alerts:

- The graphs show the events that have been logged for each incubator chamber.
- The **Date Range** button allows the displayed date range to be adjusted as desired. After the date range is adjusted, the **Date Range** button will be grayed out and only the **Reset** button will be active. Selecting the **Reset** button will restore the original date range and display the **Report** and **Export** buttons.
- The **Zooming** and **Panning** functions described below also may be used in combination with the **Date Range** button.
- Zooming:** For viewing the data values by different timespans along the X-axis, scroll the mouse wheel up or down over the graph that is to be adjusted. Select **Show All** to reset the zoomed graph to normal.
- Panning:** Each temperature graph provides an option to pan left or right by holding the left mouse button down and moving the mouse left to right, but this functionality is available only when the graph is zoomed. The panning of the temperature graph will stop if there is no further data beyond the time limit.
- The **Report** button generates a report of lid openings. Refer to **Section 5.5.7**.
- The **Export** button generates a CSV file with a running log of detailed information. Refer to **Section 5.5.8**.

Events

The **Events** tab shows the events or alerts that have been logged for each incubation chamber and any general alerts that have been logged for the incubator:

- Event/Alert Details** – What events occurred and where.
- Timestamp** – When the event occurred.
- Time Acknowledged** – When the event was acknowledged by the incubator control system.

Information

The **Information** tab shows the MINC+ details and information about the network connection.

Under the **Device** heading:

- **Incubator Name** – Incubator name assigned by the clinic.
- **Serial Number** – Incubator serial number.
- **FW Version** – The firmware version for the MINC+ core control system.
- **DishTrace MINC+ Version** – DishTrace MINC+ software version number.
- **Last Service Date** – Date of the last service.
- **Manufacture Date** – Date of device manufacture.
- **First Connected on** – Date of first connection to network.

Under the **Network Connection** heading of the incubator:

- **IP Address** – Network IP address of the incubator.
- **Connected** – Date and time of connection of selected incubator.
- **Duration** – Duration of connection of selected incubator.

5.5.7 To generate a lid openings report

1. Under the **Incubators** screen, select the **History** tab.
2. Click the **Report** button. If the Report button is not displayed, click the **Reset** button.
3. Select the start date for the report. The end date of the report must be within three months of the start date. Click **OK**.
4. Select the folder or a drive location other than the Desktop or My Documents folder to which you want to save the report.
5. Click **OK**. After the report has been successfully saved, a message will appear. The report will be in PDF format.

5.5.8 To export a running log of the incubator

1. Under the **Incubators** screen, select the **History** tab.
2. Click the **Export** button. If the Export button is not displayed, click the **Reset** button.
3. Select the start date for the report. The end date of the report must be within three months of the start date. Click **OK**.
4. Select the folder or a drive location other than the Desktop or My Documents folder to which you want to save the report.
5. Click **OK**. After the report has been successfully saved, a message will appear. The report will be in CSV format.

5.5.9 Records

Dish records for all the MINC+ devices that are connected on the network can be searched and viewed under the **Records** tab.

5.5.10 Search for a dish record

Dish records are displayed under the **Recent Records** heading, which shows the **Dish ID**, **Patient ID**, **Patient Name**, the **First Check In** for each dish, **Status**, **Location**, and the **View** column. To view individual records, click the arrow under **View** for each row.

1. Select one of the following fields from the drop-down menu in the top left-hand corner of the screen:
 - **Dish ID**
 - **Patient ID**
 - **Patient Name** (first or last name only)
 - **Incubator Name**
2. Type information into its corresponding field (e.g., Dish ID, Patient ID), or type the beginning of the information to show all exact matches in the dish records.
3. Select the required **Start** and **Finish** dates for the search.

4. Click the **Search** button to display the results. Click **Clear** to redisplay the dish records.

5.5.11 Adding a new dish record

This option can be used to add multiple dish records. After the records are saved, they are broadcast to all networked incubators.

1. Under the **Records** screen, select the **Add** button on the bottom right-hand corner. The **Add Record** screen will appear.
2. Enter dish information into at least one of the following fields: **First Name**, **Last Name**, **Patient ID**, and/or **Dish ID**. After the required fields have been filled, click **Save**.
3. A dialog box will appear: '**Are you sure you wish to save the data?**' Click **Yes**.
4. A dialog box will appear to confirm that the dish record has been created and broadcast to all connected incubators. Click **OK**. The new dish record will be added to the **Recent Records** table with an unassigned status.

5.5.12 Editing a dish record

1. On the **Records** screen, either double-click the row, or under **View**, select the '**Click to view details**' icon of the record row to be edited.
2. The **Dish Record** screen will open as shown in **Section 5.5.13**.
3. Select the edit icon on the top right of the screen in the dish information bar.
4. The **Edit Record** screen will appear. Place the cursor in the field(s) to be edited and amend the information as required.
5. Click **Save**. A dialog box will appear: '**Are you sure you wish to save the data?**' Click **Yes**.
6. A dialog box will confirm that the dish record has been successfully modified and will be broadcast to all connected incubators. Click **OK**.
7. The **Dish Record** screen will reopen, showing the amended information. Click the **Return** icon in the top left of the information bar to return to the **Recent Records** screen.

5.5.13 Dish records and dish reporting

On the **Records** screen, either double-click the row, or under **View**, select the '**Click to view details**' icon of the dish record row to be viewed or reported on. The **Dish Record** screen shown below will be displayed. The **Dish Record** information bar shows incubation information associated with the dish selected. The **History** tab shows a graphical temperature profile with alerts for the incubation period of the dish. Hovering the cursor over an individual alert marker will provide additional information. The **Events** tab shows a list of events associated with the selected dish, such as the dish being checked in or out, lid openings, and any alerts that occurred during the incubation period.

A dish report can be generated by clicking the **Report** button on the **History** tab. Two options can be selected:

- **Total Duration Checked Out** – Total time the dish was checked out over the assigned period.
- **Average Checked Out Duration** – The total time that the dish was checked out divided by the number of checkouts.

Reports are generated in PDF format. A detailed log of all the MINC+ operating parameters from the selected period for the assigned dish can be generated by using the **Export** button. Refer to **Section 5.5.14**.

To generate a dish report

1. Under the **Records** screen, double-click a dish record.
2. Select the **History** tab.
3. Click the **Report** button. A dialog box will appear. Choose the options to include in the report by clicking on the boxes, and then add any comments to be included in the report. Click **OK** to continue.
4. Select a location to save the **Dish Report**. Select the folder or a drive location other than the Desktop or My Documents folder to which you want to save the report. Click **OK** to continue.
5. Once the dish report is saved, a dialog box will appear: '**Report Successfully Exported.**' Click **OK**.

5.5.14 Export a dish-specific incubator parameter report

This function generates a CSV file with a detailed log of the MINC+ operating parameters over the assigned period for the selected dish. The file can be used for detailed analysis, if required.

1. Under the **Records** screen, double-click a dish record.
2. Select the **History** tab.
3. Click the **Export** button. A dialog box will appear. Select the folder or a drive location other than the Desktop or My Documents folder to which you want to save the report. Click **OK** to continue.
4. After the report is saved, a dialog box will appear: '**Report Successfully Exported.**' Click **OK**.

5.5.15 System tab

The system screen provides the following information for all the active networked incubators.

Temperature – Graphs the number of left and right chamber temperature set points.

Incubation Stage – Graphs the number of active dishes and their incubation time in days.

Database Management – Shows the amount of storage space used on the Dishtrace PC database.

Recent Dishes – Shows when dishes were added to the list of available dishes.

Recent Alerts & Events – Shows all alerts and events of all active incubators.


5.5.16 To view incubator alerts

1. Under the **System** tab, click the **View Alerts** button in the middle of the right-hand side of the screen. A new screen will appear that shows the incubator alerts:

Incubator Name	Description	Time	Status	View
COOK MEDICAL-E1	CO2 Warning Low Pressure Low	2021-11-05 18:19:30	inactive	<input checked="" type="checkbox"/>
COOK MEDICAL-E1	CO2 Out of Range High Pressure Low	2021-11-05 18:19:47	inactive	<input checked="" type="checkbox"/>
	Right Lai Cover Alert	2021-11-05 18:20:04	inactive	<input checked="" type="checkbox"/>

2. Double-click on each row to view the alert details.

5.5.17 DishTrace PC settings

Changing general settings

The **Settings** tab allows you to customize the way that the DishTrace PC software presents information.

1. To change the **Language**, **Date**, **Time**, **Name**, or **Record Identification**, use the drop-down menus for each field.
2. After the changes are complete, click **Save**.

5.5.18 Back up the DishTrace PC database

1. In the **Settings** tab, click the **Backup** button.
2. Choose a location to save the database file. This file cannot be saved on the Desktop or My Documents folder. It must be saved on an external drive or within the computer's C drive (e.g. C:\MINC3 data).
3. Press **Save**.

IMPORTANT NOTE: This backup will include all the incubators on the network that are currently detected and used by DishTrace PC.

IMPORTANT NOTE: The database backup file can be used only in DishTrace PC. It cannot be opened on a PC or Mac.

5.5.19 Restore DishTrace PC database (from a saved backup file)

1. In the **Settings** tab, click the **Restore** button.
2. Browse to the backup file from its current location. A dialog box will ask for confirmation. Click **Yes** to continue.
3. The database will start to restore. When it is complete, a message will appear: '**Database restore has been completed successfully.**' Click **OK** to continue.

5.5.20 Purge the DishTrace PC database

1. In the **Settings** tab, click the **Purge** button. A pop-up window will provide a choice of three time spans of data to be purged. **Note:** For oldest 3 and 6 month purges, the last 14 days data will be retained if the time period includes the last 14 days. Unassigned dishes in the selected time period will be purged. Purge **All Data** will delete all data except for assigned checked in/out dishes.
 - **Oldest 3 months**
 - **Oldest 6 months**
 - **All Data**
2. Select any option to continue.
3. A pop-up window will ask you if you want to delete the incubator data on the PC for the selected duration. Click **Yes** to continue.
4. The purging process will commence. When it is complete, a message will appear: '**Database Purge has been completed successfully.**' Click **OK** to continue.
5. After a purge **All Data**, the MINC+ will need to be reconnected to DishTrace PC as per **Section 4.13.7**. The certificate will already be in the PC's repository and only requires acknowledgment that the incubator has been added.

5.6 Alarm conditions

5.6.1 External alarm

The MINC+ can connect to an external alarm response system that will alert staff to alarms that activate outside of business hours. All alarm conditions are recorded in the log files.

This external alarm connection can be one of two types of circuits:

- Normally open circuit (closes when a fault is detected).
- Normally closed circuit (opens when a fault is detected).

The external alarm connection will activate under the following conditions:

- Loss of mains power (after two minutes).
- Low inlet gas pressure (less than 50 kPa for more than 15 minutes).
- Low water level in the humidification flask (if the water level is less than the minimum for more than 75 minutes).
- No gas flow or gas flow out of range during normal flow (< 21 mL/min or > 35 mL/min) for more than 15 minutes.
- Incubation chamber temperature out of range from set point by 0.4 °C for more than two minutes.
- Media warming chamber temperature out of range from set point by 1 °C for more than two minutes.
- Internal function error is detected.

Customer requirements for the external alarm response system should be referred to a company that specializes in this type of equipment. Refer to **Section 8** for the alarm contact rating.

If any alarm condition is unable to be rectified, contact your local Cook Medical representative or Cook Medical service agent.

5.6.2 Loss of mains power

If a loss of mains power is encountered, the following sequences will occur:

- If the MINC+ is switched on and one or both incubation chambers are active, loss of mains power will cause the external alarm contacts to activate if the power is still off after two minutes.
- If mains power is restored within two minutes, the external alarm will not activate and the MINC+ will recommence normal operation.
- If mains power is restored after two minutes, the external alarm will reset and the MINC+ will recommence normal operation.
- If the gas flow was active before the power was lost, then the MINC+ will restart the flow sequence in gas purge mode when power is restored.

5.6.3 Low inlet gas pressure

When the inlet gas pressure goes below 50 kPa the gas flow status will display a **Low inlet gas pressure** icon and emit an audible alarm every 30 seconds.

The external alarm contacts will activate 15 minutes later if the correct inlet gas pressure is not restored. The alarm will reset when an inlet gas pressure above 60 kPa is restored.

A warning icon will also be displayed on the center console while the alert is active.

5.6.4 Low humidification flask water level

If a low water level is detected in the humidification flasks, an immediate alert on the center console screen is displayed. The **Low humidification flask water level** icon will start flashing in red and emit an audible alarm after 15 minutes where the low water-level alert is detected.

If the water level remains low for longer than 75 minutes, the external alarm contacts will activate and the **Low humidification flask water level** icon will continue to be displayed and the audible alarm emitted.

Alarm contacts will reset if the humidification flask is replaced with a full flask or if the water level alert is disabled.

5.6.5 No gas flow or gas flow out of range

If the normal gas flow is less than 21 mL/min or greater than 35 mL/min (including no gas flow) for longer than 10 minutes, the normal gas flow status will display the **Normal gas flow out of range** icon and emit an audible alarm every 30 seconds. The external alarm contacts will activate 5 minutes later.

To reset the alarm after rectification of the cause, press the **Unlock** button and then toggle the **Enable/disable gas flow** button. The normal flow out of range alarms will remain active until they have been manually reset after rectification of the cause.

If the gas-purge flow rate is less than 259 mL/min or greater than 450 mL/min for more than 2 minutes during the 3 minute purge cycle, a warning icon will also be displayed on the center console and logged as a purge flow out of range notification.

To conserve gas, the flow will revert to normal gas flow status after the elapsed 3 minute purge cycle regardless of the attained purge flow. Toggle the associated chamber **Enable/disable gas flow** button to acknowledge the notification. If the problem repeats during the subsequent purge cycle, contact your local Cook Medical service agent.

5.6.6 Temperature out of range

When the MINC+ has just been turned on or the set temperature has been adjusted, the temperature alarm is deactivated for 120 minutes to allow the device to reach stable conditions without the alarm activating.

If after 120 minutes the incubation chamber temperature differs from the set point temperature by greater than ± 0.4 °C, for longer than 2 minutes, the chamber display will show a **Temperature out of range** icon, the temperature value will be red, an audible alarm will sound, and the external alarm contacts will be activated. To reset the alarm, touch the **Unlock** button and then toggle the **Enable/ disable heater** button off and on.

If the media warming chamber temperature is out of range by more than +/- 1 °C for longer than 2 minutes, when the MWC temperature alarm is activated, the three LEDs in the front panel below the MWC lid will flash red, the MWC temperature value will be red, an audible alarm will sound, and the external alarm contacts will be activated. To reset the MWC temperature alarm, open and close the MWC lid or toggle both incubation chambers **Enable/disable heater** off and on.

A warning icon will also be displayed on the center console while the alert is active.

 **IMPORTANT NOTE:** Resetting an incubation chamber or MWC temperature alarm will deactivate that alarm for a further 120 minutes.

5.6.7 Chamber screen and center console error displays

The table below shows the icons that will display on the incubation chamber screens. A warning icon will also simultaneously display on the center console.

Icon name	Icon	Display color and state
Low inlet gas pressure		Flashing red icon with audible alarm
Normal gas flow rate out of range		Flashing red icon with audible alarm
Temperature out of range		Flashing red icon with audible alarm
Low humidification flask water level		Flashing red icon with audible alarm

5.6.8 Internal functional error

If an internal error is detected, a five-digit number will be displayed in the top left-hand corner of either or both chamber display screens. The alert symbol will also be displayed on the center console.

Should an internal functional error be detected, the MINC+ will try to rectify the error, in some cases performing a restart. If the error persists after the restart, the external contacts will activate to alert the user of an abnormal condition. If the alarm condition is unable to be rectified, contact your local Cook Medical service agent.

5.7 After use

If the MINC+ will not be used for an extended period of time, follow the steps below.

1. Place the MINC+ in standby mode.
2. Turn off power to the device.
3. Remove the humidification flasks.
4. Disconnect the power cord.
5. Disconnect the gas supply.
6. Clean the device (refer to **Section 6.2**).

6. Service and maintenance

⚠️ WARNING: To guarantee safe operation, it is necessary to carry out proper care and maintenance of the MINC+ and disposables. Regular checks are recommended to confirm that the device is functioning correctly.

⚠️ WARNING: New and repaired products must be prepared and tested according to the user manual before use.

6.1 Changing the humidification flasks

⚠️ WARNING: BIOLOGICAL HAZARD. Do not use a contaminated humidification flask in the MINC+.

If the flask is used continuously, it can be used until the water level reaches the MIN FILL line on the front of the flask.

If the flasks are not used continuously, there may be an increased risk of biological contamination. If the device is not in use, remove the flasks and replace with new flasks upon startup.

Do not refill the flasks.

To change the humidification flasks:

1. Ensure that the gas flow to each chamber is turned off. It is also recommended that the power cord to the MINC+ be removed from the mains supply wall socket.
2. Lift both heater chamber lids and slide out the humidification flasks. For detailed instructions refer to **Section 4.10.4**.
3. It is recommended that the entire device be cleaned prior to the installation of a new humidification flask. Refer to **Section 6.2**.
4. Prepare and insert a new flask into each HF bay. Refer to **Section 4.10.2 and 4.10.3**.
5. To recommence operation of the MINC+, turn on the gas flow to begin the gas purge process.
6. Allow 4 hours for the water to heat and saturate with CO₂.
7. Accelerate the equilibration process, if desired, by following these steps:
 - a. Preheat the aseptically prepared humidification flask to 35 °C in **Step 4**.
 - b. Purge the system, as described in **Step 5**, at least four times in quick succession.

Performing these steps will allow immediate use of the device.


Figure 16: Changing humidification flasks

6.2 Cleaning the incubator

6.2.1 Full cleaning

⚠️ WARNING: ELECTRIC SHOCK HAZARD. Do not immerse the MINC+ in any type of liquid.

Before periodic cleaning, or if the MINC+ is not in use for extended periods of time, follow the steps below.

1. Remove the contents of the incubation chambers.
2. Place the MINC+ in standby mode.
3. Disconnect the power cord from the mains supply outlet, and remove the power cord from the power inlet at the rear of the device.
4. Remove the humidification flasks.
5. Remove the gas supply.
6. Clean the device per the instructions below.

It is recommended that the MINC+ be cleaned with aqueous 70% isopropyl alcohol, purified water, or 3% hydrogen peroxide solution. Moisten a cloth and wipe all internal and external surfaces of the device. Do not immerse the device in the cleaning solution.

Clean the gas vent tubes located in the center of the front of each chamber base with a small pipe cleaner wetted with the aqueous 70% alcohol solution.

After cleaning, leave the lids of the unit open for 30 minutes to ensure that all alcohol fumes have dissipated.

👉 IMPORTANT NOTE: Purified water may be used to wipe device surfaces when an alcohol solution is not available or suitable.

👉 IMPORTANT NOTE: 70% aqueous ethanol maybe used only to clean the aluminum anodized surfaces of the incubation chambers, the MWC chamber, and the front decal, using a moistened cloth. Avoid contacting all other plastic surfaces.

6.2.2 On-site cleaning

To perform an on-site cleaning of the incubator in the laboratory when alcohol cannot be used, use a cloth moistened with purified water to wipe the incubator surfaces. Removing the flasks is optional. It is also recommended that the power to MINC+ be turned off by removing the power cord from the power inlet at the rear of the device.

6.3 Cleaning the media warming chamber (MWC)

⚠️ WARNING: ELECTRIC SHOCK HAZARD. Ensure that the MWC cradle and bay are dry before use. Before periodic cleaning of the MWC, remove its contents. Shut down the device by removing the power cord. It is recommended that the MWC be cleaned with aqueous 70% isopropyl alcohol, purified water, or 3% hydrogen peroxide solution if alcohol is not available. 70% aqueous ethanol maybe used only to clean the aluminum anodized surfaces of the MWC.

To clean the media warming chamber:

1. Lift the MWC cradle from its bay.
2. Moisten a cloth and wipe all internal and external surfaces of the MWC cradle and its bay. The MWC cradle can be immersed in purified water for cleaning.
3. After cleaning, ensure that the bay and the MWC cradle are dry. Replace the MWC cradle in the bay.


Figure 17: Removing the MWC cradle

1**6.4 Six monthly functionality testing**

To preserve the MINC+ and maintain its performance, regular inspections are necessary for early detection of possible device malfunctions. The user or a qualified technician must regularly test the device to assess its functionality.

The following must be tested every six months:


- Temperature
- Gas flow rate
- External alarm contacts. The functionality of the external alarm contacts is also checked internally each time the MINC+ is powered on. It is therefore recommended that the power to the device be turned off at regular intervals, such as during cleaning and flask replacement as a suggested minimum.

6.4.1 Temperature

Set each chamber temperature to 37.0 °C. Insert a thin wire thermocouple through the gas outlet tube into the chamber, as shown in the diagram below. Ensure that the sensor wire is thin enough that it does not occlude the gas flow. Place a Nunc four-well dish into the front outer dish location of each chamber so that the well with the thermocouple is in the position shown below. Fill the well that has the thermocouple with 1 mL of water, close the chamber lid, and allow enough time for the thermocouples to thermally equilibrate. The temperature reading should be 37.0 ± 0.2 °C.

If the measurement lies outside the acceptable range, contact your Cook Medical service agent.

To measure the temperature, the thermocouple must be secured so that the tip is immersed in the fluid and touching the bottom of the dish and remains in that position when the dish lid is fitted and the chamber lid is closed. It is recommended that the thermocouple and thermometer readout be calibrated at 37.0 °C to compensate for any offset between the actual temperature and the indicated value, if required.

The temperatures measured will be the temperature set point ± 0.2 °C.

IMPORTANT NOTE: To avoid damaging the lid seal, do not pass the thermocouple wire through the lid seal.

6.4.2 Gas flow rate

To test the gas flow rate:

1. Insert a fully filled flask into each HF bay.
2. Wait until a normal flow cycle is initiated and the rate of bubbling appears stable.
3. Count the number of bubbles formed against each flask front face over a one-minute period. The number of bubbles formed in each flask should be at least 9 per minute. The bubble count per minute can vary significantly between flasks; it may be in excess of 25 bubbles per minute.

6.4.3 External alarm contacts

To test the external alarm contacts, ensure that at least one heater chamber is active and then disconnect the MINC+ from mains power. The contacts should open or close (depending on which one is the normal position) approximately 2 minutes later. The contacts should actuate to the opposite of their normal state. After you confirm that the contacts have actuated, reconnect the device to the mains power. The contacts should then reset to their normal position (i.e., normally open or normally closed).

If the contacts fail to actuate, contact your Cook representative.

6.5 Inspection by an authorized Cook Medical service agent

WARNING: No user serviceable parts inside!

Inspections at least once a year

For ongoing operational performance of the MINC+, an authorized Cook Medical service agent must maintain the device annually. The service agent will assess the operational functionality of the water level detector, gas flow and temperature control systems, and external alarm contacts.

Authorized service agents

All services such as alterations, repairs, calibrations, etc., may be performed only by the manufacturer or by service agents who are authorized by the manufacturer.

Liability

The manufacturer is free from all liability for the operational safety of the MINC+ if the device is willfully opened and unauthorized persons perform repairs or alterations on it during the warranty period.

Certification

The MINC+ owner will receive a signed certificate from the service agent for all inspections or repairs. This certificate states the type and scope of the services rendered, the service date, and the name of the service company.

Technical documentation

If the manufacturer provides technical documentation, this does not authorize the user to perform repairs, adjustments, or alterations to the MINC+ or disposables.

6.6 Return procedure

WARNING: BIOLOGICAL HAZARD. Before you return a product, you must clearly mark it with a contamination warning and seal it in a plastic bag that is sealed within a second plastic bag. When you ship the MINC+, remove the humidification flasks and all chamber contents prior to transport.

IMPORTANT NOTE: When returning goods, use the original packaging. The manufacturer does not take responsibility for damage that occurs during transportation if the damage is caused by inadequate transport packaging.

All devices or disposables that are to be returned must be prepared as described below for the protection of the service agent and for safety during transportation.

To return the MINC+

1. Reset the MINC+ to its factory settings to remove patient and dish information. Refer to **Section 5.4.22**.
2. Clean the unit as detailed in the cleaning **Section 6.2**.
3. Place in the original packaging. Ensure that the device is in the original bag and that the bag is sealed with tape.
4. Enclose the following information:
 - Owner's name.
 - Owner's address.
 - Model type.
 - Serial number of the equipment (located on the rear of the MINC+, or found by tapping **Settings** and then tapping **Information** on the center console).
 - Description of the damage or fault.

The manufacturer has the right to refuse to carry out repairs if the products received are contaminated.

7. Troubleshooting

The following table displays error conditions and alarms and provides possible solutions. If solutions provided fails to remedy the problem, contact a Cook Medical service agent.

Error and alarm indicator	Source of error	Elimination of error
MINC+ will not turn on.	Power supply cord is not connected.	Check the power supply connection.
Displays do not illuminate.	Mains power is not turned on. The device is in standby.	Turn on the mains power. Check that the standby indicator is illuminated. Press the standby touchpad to activate the MINC+.
Numerical error displayed (e.g., 10101).	Device has internal fault.	Allow the MINC+ to run for five minutes, if the error remains, restart the incubator. If the problem persists, contact a Cook Medical service agent.
MINC+ will not heat to correct temperature.	Temperature settings not entered correctly.	Enter correct temperature settings. If problem persists, contact a Cook Medical service agent.
'Temperature out of range' alarm. Chamber temperature value displayed in red.	Set temperature not achieved in time. Lids open for extended periods. Device is not operating in functional environmental conditions.	Turn heater bay off and on again Do not leave lids open for extended periods. Turn heater bay off and on again. Evaluate device placement.
'Low inlet gas pressure' alarm. 'CO ₂ ' displayed.	Gas supply not connected. Inlet gas pressure below 50 kPa. Gas lines are kinked.	Check gas supply connection. The gas flow should recommence in gas purge mode when the gas supply is reconnected. Check gas supply pressure is above 50 kPa. The low inlet gas pressure alarm icon will disappear when the inlet gas pressure is reinstated above 60 kPa. The gas flow will recommence in gas purge mode when the gas supply is reconnected. Check that gas lines are not obstructed. The gas flow should recommence in purge mode when the gas supply is reconnected.
	Gas cylinder changeover unit incompatible with the MINC+.	Ensure that the inlet gas pressure to the MINC+ does not fall below 135 kPa. The gas flow should recommence in gas purge mode when the gas supply is reconnected.
	Resistance to gas flow from series connection of multiple devices is too high to allow multi-unit purging.	Check the regulated pressure is still a nominal 150 kPa ±15 kPa at the last series connected MINC+, with all series-connected devices purging.
No gas flow or 'Gas flow out of range' alarm.	Flask is incorrectly installed and is occluding the flow.	Check the installation of the humidification flask and that the flask T-piece is correctly installed. Ensure that the HF filter sleeve from previous flask has been removed from the HF bay gas nozzle.
No bubbles visible in the front of the humidification flasks.	Flasks not installed correctly.	Turn flow off and on again. Check the installation of the humidification flasks. Ensure that the HF filter sleeve from previous flask has been removed from the HF bay gas nozzle.
Excessive gas consumption.	High supply pressure. Gas connections not secure or damaged. Gas line damaged.	Ensure that the inlet gas pressure to the MINC+ does not rise above 165 kPa. Ensure all gas fittings from the cylinder to the MINC+ gas inlet are tight. If error persists, inspect all fittings for damage to seal surfaces, and replace as necessary. Replace the gas line.
Low humidification flask water level alarm.	Not enough water in the humidification flask. Wrong water levels are displayed.	Replace the humidification flask with a new filled flask. Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Lid open icon displayed and device emitting beep every 15 seconds.	A lid latch is not in the closed position. Lid latch is closed but not engaged with the bayonet (the lid is propped open slightly).	Close the lid, and move the latch to the closed position. Open the latch and ensure that the lid is fully closed, and then move latch to the closed position.

Error and alarm indicator	Source of error	Elimination of error
DishTrace MINC+		
Network error (orange icon under incubator name).	Network cable connections are not correctly connected.	Check the network connections to the MINC+ and that DishTrace PC is running.
MINC+ does not recognize the USB drive.	USB drive is not correctly formatted.	Format the USB drive to FAT32.
Center console		
Center console fails to respond.	System has frozen.	Press the Reset button on the rear panel of the MINC+. If the problem persists, contact a Cook Medical service agent.
Center console does not display icons.	Device has an internal fault.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Alarm		
Alarm fails to actuate: no sound emitted, no alarm icon is displayed, alarm contacts do not actuate.	Internal error.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Buttons		
Power button not turning the incubator on or putting it into standby mode.	Keypad functions not working properly, system has frozen.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Power button flashing.	Internal system error.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Lock/Unlock button not working, keypad cannot be locked or unlocked for user input.	Lock button failure due to internal malfunction.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Heater Enable button not responding – heater will not turn on or off.	Heater button failure due to internal malfunction.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Gas Flow Enable button not responding – gas flow inactive.	Gas Flow Enable button failure due to internal malfunction. Keypad sensor sensitivity is wrongly configured.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Water Level Enable button not responding – water level cannot be enabled or disabled.	Water level button failure due to internal malfunction.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Temperature Up and Down buttons not responding.	Buttons not responding due to internal malfunction.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Temperature Up and Down buttons not functioning correctly. Up decreases the temperature and Down increases it.	Up function is detected when the down button is pressed.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Media warming chamber (MWC) Up and Down buttons not responding.	Buttons not responding due to internal malfunction.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
MWC temperature Up and Down buttons not functioning correctly. Up decreases the temperature and Down increases it.	Up function is detected when the down button is pressed.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
MWC temperature not displaying at all, or not displaying correctly while the chamber is on and warming test tubes.	Temperature not displaying due to internal display corruption.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
Chamber display screens		
Set Temperature display on screen is wrong.	Device has internal malfunction.	Restart the incubator. If the problem persists, contact a Cook Medical service agent.
DishTrace PC install and connectivity		
Windows 8.1/10 "Application Installation was unsuccessful." DishTrace PC software won't install or the installation is unsuccessful.	Windows has not completed the installation of Windows updates in the background. Ensure that the updates are completed, restart the PC, and then begin the DishTrace PC installation again.	Refer to Section 3.2 , DishTrace PC. To view the update history on your Windows operating system, do the following: 1. From the Start menu, open Settings > Update & Security > Windows Update. 2. Click the View Update History link to open the Update History window. If the above instructions do not work, contact your Cook Medical service agent.
MINC+ fails to automatically connect to the DishTrace PC software.	DHCP server is not present or is not configured. Ethernet cables are not connected correctly.	If the DHCP server should be available, contact your network provider for advice on DHCP configuration. If DHCP is not available, assign static IP addresses with a subnet of 10.0.0.x, where x is a number from 1 to 254 that is unique to each connected device. If the above instructions do not work, contact your Cook Medical service agent. Reseat the ethernet cable connections or replace the cable.
	DishTrace PC certificate exchanges not completed or user group incomplete.	Complete DishTrace PC installation process in accordance with Section 4.13 . If problem persists contact a Cook Medical Service agent.
Events and dish information not displayed on DishTrace PC screen.	MINC+ database not synchronized with DishTrace PC database.	Refer to Section 5.5.5 to resynchronize the DishTrace PC database with connected MINC+ device databases.

8. Technical data

1

Classification according to IEC 60601-1

Type of protection against electric shock	Class I equipment
Degree of protection against harmful ingress of solids and water	IP31

General specifications

Power supply	100 - 240 VAC
Frequency	50 - 60 Hz
Oversupply category	II
Pollution degree	2
Mains supply voltage fluctuations	+/-10%
Maximum power consumption	200 VA
Alarm contact rating	30V \sim 2.0A
Environmental operating conditions	Indoor use only: + 18 °C to + 32 °C at temperature set point of 37.0 °C + 20 °C to + 28 °C for all other temperature settings 10% to 85% RH 700 hPa to 1,060 hPa (Altitude up to 3000m)
Storage and transport directions	+ 5 °C to + 40 °C 10% to 85% RH
Manufactured and tested to the following standards	IEC 60601-1: 2005 + AMD1:2012 IEC 60601-1-2: 2014 IEC 61010-1: 2010 + AMD1:2016, IEC 61010-2: 010: 2019 UL 61010-1 (3rd Ed.): Am.1 CSA C22.2 No. 61010-1 3ED + Amd 1 CSA C22.2 No. 61010-2-010:15
Dimensions	Door closed: 556 mm wide x 197 mm high x 449 mm deep. Door opened: 556 mm wide x 528 mm high x 449 mm deep.
Weight (device only)	23.4 kg (51.6 lb)
Gas supply type	Blend of 6% CO ₂ , 5% O ₂ , 89% N ₂ (Cook Culture System at sea-level) or high purity 6% CO ₂ in air (recommended tolerances ± 0.2%).
Gas supply pressure	150 kPa ± 15 kPa (21.8 psi ± 2.2 psi) (1,500 mbar ± 150 mbar)
Gas flow rate, nominal	25-29 mL/min per chamber (both chambers in normal flow status). Gas-purge flow rate at 300 mL/min per chamber for 3 minutes (both chamber in purge flow status).
Gas flow rate, accuracy	4 mL/min per chamber (normal flow). ± 30 mL/min per chamber (purge).
Chamber temperature capability	35.0 °C to 40.0 °C in 0.1 °C increments in an ambient temperature range of + 20 °C to + 28 °C. At set point of 37.0 °C, the ambient temperature range is extended to + 18 °C to + 32 °C.
Incubation chamber temperature accuracy	± 0.2 °C at calibration point.
Media warming chamber	Temperature set point range 35.0 °C to 40.0 °C. Accuracy +/- 1 °C.

Computer requirements for the DishTrace software

Operating system	Windows 8.1, Windows 10.
Processor and RAM	Minimum configuration as specified by the host OS.
Free hard disc space	10 GB minimum.
Processor	Pentium 4, 1.6 GHz or better.
Display	1024 x 768, 256 color or better.
I/O ports	Ethernet port for connecting to networked MINC+ devices. USB port for connecting up to ten K-MINC-1000 incubators via hub. USB port for installation of software.

Specifications

Maximum number of connected MINC+ devices	50 (including ten K-MINC-1000)
Maximum logging time period	Limited only by free disc space.
Log file format	Comma-separated values (CSV)
Report format	PDF
USB drive format for backup and reports	FAT32
Graph update interval	25 seconds
Maximum displayed time period on screen	5 days
Temperature resolution	0.1 °C

Guidance and manufacturer's declaration – Electromagnetic immunity

The MINC+ Benchtop Incubator is intended for use in the electromagnetic environment specified below. The customer or the end user of the MINC+ Benchtop Incubator should ensure that it is used in such an environment.

1

<i>Immunity test</i>	<i>IEC 60601 test level</i>	<i>Compliance Level</i>	<i>Electromagnetic environment guidance</i>
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient (EFT) IEC 61000-4-4	± 2 kV at 100 kHz for AC mains ± 1 kV at 100 kHz for I/O ports	± 2 kV at 100 kHz for AC mains ± 1 kV at 100 kHz for I/O ports	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line-to-line ± 2 kV line-to-ground	± 1 kV line-to-line ± 2 kV line-to-ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% Ut for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% Ut for 1 cycle and 70% Ut for 25 cycles at 50 Hz or 30 cycles at 60 Hz 0% Ut for 250 cycles at 50 Hz or 300 cycles at 60 Hz	0% Ut for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% Ut for 1 cycle and 70% Ut for 0.5 seconds 0% Ut for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MINC+ Benchtop Incubator requires continued operation during power mains interruptions, it is recommended that the MINC+ Benchtop Incubator be powered from an uninterruptible power supply, institutional standby generator, or an external battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms for 0.15 to 80 MHz 6 Vrms in ISM bands between 0.15 and 80 MHz 80% amplitude modulation at 1 kHz	3 Vrms for 0.15 to 80 MHz 6 Vrms in ISM bands between 0.15 and 80 MHz 80% amplitude modulation at 1 kHz	Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the MINC+ Benchtop Incubator, including cables. ISM bands between 0.15 and 80 MHz are considered to be: 6.765 to 6.795 MHz 13.553 to 13.567 MHz 26.957 to 27.283 MHz 40.66 to 40.7 MHz Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ¹ should be less than the compliance level in each frequency range, ² Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80% amplitude modulation at 1 kHz	3 V/m 80 MHz to 2.7 GHz 80% amplitude modulation at 1 kHz	 Refer to table below on proximity fields for minimum separation distance.
Proximity field from wireless transmitters IEC 61000-4-3	Refer to table	Refer to table	Refer to table below on proximity fields for minimum separation distance.

At 80 MHz and 800 MHz, the higher frequency range applies.

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

1. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MINC+ Benchtop Incubator is used exceeds the applicable RF compliance level above, the MINC+ Benchtop Incubator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the MINC+ Benchtop Incubator.
2. Over the frequency range of 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Guidance and manufacturer's declaration – Electromagnetic emissions

The MINC+ Benchtop Incubator is intended for use in the electromagnetic environment specified below. The customer or the end user of the MINC+ Benchtop Incubator should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The MINC+ Benchtop Incubator uses RF energy only for its internal function. Therefore, its RF emissions are very low and unlikely to interfere with nearby electronic equipment. The MINC+ Benchtop Incubator is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 32	Class B	
Harmonic emissions IEC 61000-3-2	Class A, Class D	
Conducted emissions CISPR 11	Group 1 Class B limits	
Conducted emissions CISPR 32	Class B limits	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Proximity fields from RF wireless communication equipment

Test frequency (MHz)	Band (MHz)	Service ¹	Modulation ²	Maximum power (W)	Distance (m)	Immunity test level (V/m)	Minimum separation distance (m)
385	360-390	TETRA 400	Pulse modulation ² 18 Hz	1.8	0.3	27	0.3
450	430-470	GMRS 460 FRS 460	FM ³ ± 5 kHz deviation 1 kHz sine	2	0.3	28	0.3
710	704-787	LTE band 13,17	Pulse modulation ² 217 Hz	0.2	0.3	9	0.3
745							
780							
810	800-960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE band 5	Pulse modulation ² 18 Hz	2	0.3	28	0.3
870							
930							
1,720	1,700-1,990	GSM 1800 CDMA 1900 DECT LTE band 1, 3, 4, 25; UMTS	Pulse modulation ² 217 Hz	2	0.3	28	0.3
1,845							
1,970							
2,450	2,400-2,570	Bluetooth WLAN 802.11 b/g/n WiFi 2450 LTE band 7	Pulse modulation ² 217 Hz	2	0.3	28	0.3
5,240	5,100-5,800	WLAN 802.11 a/n	Pulse modulation ² 17 Hz	0.2	0.3	9	0.3
5,500							
5,785							

Note: If necessary to achieve the IMMUNITY test level, the distance between the transiting antenna and the ME equipment or ME system can be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

1. For some services, only the uplink frequencies are included.
2. The carrier shall be modulated using a 50% duty cycle square wave signal.
3. As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

⚠ WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

⚠ WARNING: Use of cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

⚠ WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MINC+ device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment will result.

9. Limited warranty

Oceania regional only warranty

1

This warranty is effective only in the Oceania region, including Australia, New Zealand, and the Pacific islands. Products purchased in the Oceania region shall be subject to the terms of this warranty and not 'Company' worldwide warranty.

William A. Cook Australia's goods and services come with guarantees that cannot be excluded under the Australian Consumer Law. For major failures with the service, you are entitled:

- to cancel your service contract with us; and
- to a refund for the unused portion, or to compensation for its reduced value.

You are also entitled to choose a refund or replacement for major failures with goods. If a failure with the goods or a service does not amount to a major failure, you are entitled to have the failure rectified in a reasonable time. If this is not done you are entitled to a refund for the goods and to cancel the contract for the service and obtain a refund of any unused portion. You are also entitled to be compensated for any other reasonably foreseeable loss or damage from a failure in the goods or service.

Worldwide limited warranty

William A. Cook Australia warrants to the purchasers of this device that at time of manufacture, the product was prepared and tested in accordance with good manufacturing practices and guidelines specified by the Australian Therapeutic Goods Administration or relevant competent authority.

In the event of product failure under normal use, due to defects in material or workmanship, within a period of one (1) year from the date of purchase, the product will be repaired, or at Cook's option, replaced, at no charge. This limited warranty does not apply to products subjected to abnormal use or conditions, improper storage, damaged by accident, misuse or abuse, improper line voltage or to products altered or serviced by anyone other than William A. Cook Australia or its authorized agent.

The foregoing limited warranty is exclusive and in lieu of all other warranties whether written, oral, expressed, or implied. In particular, William A. Cook Australia does not warrant that the product is suitable for the needs of the purchaser and there are no warranties given as to merchantability or fitness for a particular purpose. William A. Cook Australia's representations concerning fitness for purpose or suitability for use by any purchaser do not extend beyond those representations set out in the William A. Cook Australia literature that accompanies the product. William A. Cook Australia assumes that the purchaser is experienced in the use of this device and is able to judge from his/her own expertise the suitability or otherwise of the product for the intended use. William A. Cook Australia conducts a technical advisory service, which can be consulted by a purchaser or intended purchaser on an advisory basis.

After one (1) year from the date of purchase, this device will be repaired for a repair charge equal to the cost of parts, labor and transport. Before returning a product for any reason, please contact your nearest Cook distributor for assistance and instructions.

William A. Cook Australia reserves the right to change or discontinue this product without notice.

9.1 Liability

Because William A. Cook Australia has no control or influence over the conditions under which this device is used, over its method of use or administration, or on handling of the product after it leaves its possession, William A. Cook Australia takes no responsibility for the results, use, and/or performance of the product. William A. Cook Australia expects that use of the product will be confined to trained and expert users.

In no event will William A. Cook Australia be liable for any direct or indirect damages, including incidental, consequential, or special damages, arising out of or in connection with the use or performance of the product.

If the manufacturer provides you with technical documentation, this does not authorize you to perform repairs, adjustments, or alterations on the device or disposables.

No representative of William A. Cook Australia and no vendor or lessor of the product is authorized to change any of the foregoing terms and conditions, and the purchaser accepts the product subject to all terms and conditions herein, subject always to any contrary provisions which are necessarily implied by statute or law notwithstanding the within terms and conditions.

9.2 Life of the product

The expected service life of this product is deemed to be twelve (12) years, subject to the service schedule being maintained throughout the life of the device. After this time William A. Cook Australia will no longer be responsible for this product.

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The screenshot shows a table with columns: Contenu, Enreg., Statut, Paramètres, and a serial number 90987680014. The 'Alertes/Evénements' tab is selected. A checkbox labeled 'Afficher les alertes uniquement' is checked. The table lists various events and alerts with their descriptions, dates, and times.

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The screenshot shows a reporting interface with tabs: Annuler, Rapport incubateur, and Exporter rapport. It includes date selection fields 'De' (03-Nov-2021) and 'À' (03-Nov-2021).

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