
EN 3	Murdoch Mechanical Arm Instructions for Use
CS 5	Mechanické rameno Murdoch Návod k použití
DA 7	Murdoch mekanisk arm Brugsanvisning
DE 9	Murdoch mechanischer Arm Gebrauchsanweisung
EL 11	Μηχανικός βραχίονας Murdoch Οδηγίες χρήσης
ES 13	Brazo mecánico Murdoch Instrucciones de uso
ET 15	Murdochi mehaaniline statiiv Kasutusjuhised
FI 17	Murdoch- mekaaninen varsi Käyttöohjeet
FR 19	Bras mécanique Murdoch Mode d'emploi
IT 21	Braccio meccanico Murdoch Istruzioni per l'uso
LV 23	Murdoch mehāniskā svira Lietošanas norādījumi
NL 25	Murdoch mechanische arm Gebruiksaanwijzing
NO 27	Murdoch mekanisk arm Bruksanvisning
PL 29	Wysięgnik mechaniczny Murdoch Instrukcja użycia
PT 31	Braço mecânico Murdoch Instruções de utilização
SK 33	Mechanické rameno Murdoch Návod na použitie
SL 35	Mehanska roka po Murdochu Navodila za uporabo
SV 37	Murdoch mekanisk arm Bruksanvisning

MURDOCH MECHANICAL ARM

DEVICE DESCRIPTION

The device is comprised of the following components:

The device shall be used for holding retractors in patients undergoing a laparoscopic surgery.

The intended user of the device shall be an experienced/trained physician in laparoscopic surgery.

INTENDED USE

The Murdoch Mechanical Arm (MMA) is intended to hold the Nathanson Liver Retractor (NLRS) or Laparoscopic Pelvic Retractor (LPR) in position during Laparoscopic Surgery.

DEVICE AND CLINICAL BENEFIT

The Murdoch Mechanical Arm has been developed to minimize the risk of injury to the patient by holding the retractor stable for the duration of the procedure.

INDICATIONS

The Murdoch Mechanical Arm is indicated for use during laparoscopic surgery where retraction is required.

CONTRAINDICATIONS

There are no known contraindications for this device.

WARNINGS

Over-tightening the locking handles and bolts or reassembling components incorrectly will lead to irreparable damage to the unit.

PRECAUTIONS

- Ensure that the device is operated and used only by the persons with the requisite training, knowledge and experience.
- Read, follow, and keep these instructions for use.
- Use the device only in accordance with its intended use.
- Remove the packaging and clean the MMA with the cleaning instructions prior to its initial sterilization and usage.
- Do not use if MMA is damaged or defective.

CLEANING AND STERILIZING INSTRUCTIONS

Warnings

- Do not use acidic detergents in automated or manual cleaning as this will cause corrosion of the MMA.
- It is recommended that the MMA is reprocessed as soon as reasonably practical, or within 30 minutes following use.

Initial treatment at the point of use

- Keep the device moist immediately after use, or as soon as practical.
- Remove excess soil from MMA with water and wipe with a non-linting cloth.

Preparation before cleaning

Equipment: non-abrasive brushes, detergent (enzymatic detergent pH 7.0-10.0).

- Disassemble the components of MMA.
- Rinse the components under ambient running tap water for 3 minutes.
- Remove excess soil while rinsing using non-abrasive brush or non-linting cloth.
- Soak the device in prepared enzymatic detergent solution (detergent concentration: 8 mL/Liter (1 oz per gallon))² for 20 minutes.
- Scrub the components with a non-abrasive brush while they are immersed in the detergent solution.
 - Scrub external and internal surfaces for a minimum of 3 minutes each to thoroughly clean the components. Pay attention to threaded areas, hinges, curves and internal areas.
 - Using a nylon channel brush (non-abrasive style) push the brush back and forth in twisting motion through all the individual components and repeat the process 5 times.
- Rinse the components under tap water for 3 minutes and aggressively flush threaded parts and internal/difficult to reach areas.

Cleaning (Automated)

Equipment: Washer/disinfector, detergent (enzymatic detergent pH 7.0-10.0).

- Load products into washer/disinfector.
- Complete a cleaning cycle using prepared enzymatic detergent solution (Detergent concentration: 4 mL/Liter (0.5 oz per gallon))².
- Minimum cycle parameter set points listed in **Table 1**, or alternatively use a cleaning cycle validated for non-complex re-useable medical devices (RMDs).
- When unloading, check for complete removal of visible contaminants. If necessary, repeat cycles or use manual cleaning.
- Dry the components using a non-linting cloth and compressed air.

Table 1: Process Parameters – Automated Cleaning

Phase	Time (minutes)	Temperature – Treated water
Pre-wash	2.2	Cold water
Wash	4.5	Set Point: 72 °C (162 °F)
Rinse	15 seconds	Set Point: 91 °C (196 °F)
Thermal rinse	1	Set Point: 82.2 °C to 93 °C (180 °F to 199 °F)
Drying	7	Set Point: 98.9 °C (210 °F) High

Cleaning (Manual)

Equipment: Sonicator, detergent (enzymatic detergent pH 7.0-10.0).

- Sonicate components in prepared detergent solution (Detergent concentration: 4 mL/Liter (0.5 oz per gallon))² in ultrasonic unit at 45°C (113°F) for 10 minutes with ultrasonic frequency at 40 kHz.
- Rinse the components using treated water¹ for 2 minutes and aggressively flush threaded parts and internal/difficult to reach areas.
- Dry the components using a non-linting cloth and compressed air.

¹ Treated water as defined in AAMI TIR34 Water for Reprocessing of Medical Devices.

² Detergent concentration validation was performed using Cidezyme Xtra Multienzyme Detergent. Detergent concentration used as recommended by manufacturer. Manual: 8 mL/Liter (1 oz per gallon); Automated: 4 mL/Liter (0.5 oz per gallon). Any equivalent detergent can be used for cleaning. Detergent concentration to be as per that manufacturer's recommendation for cleaning processes.

Drying

Ensure the MMA components are dry before packaging.

Sterilization

Moist heat sterilize as per ISO/TS 17665 Part 1 and Part 2.

- Disassemble for sterilization.
- Prepare the components by packing them in a sterilization pouch or place them in a tray and ensure that steam can penetrate through all surfaces.
- Run sterilization cycles in pre-vacuum mode with either of the parameters described in **Table 2**.

Table 2: Moist Heat Sterilization Parameters

Process steps	Parameters	
Method	Moist heat sterilization	
Pre-vacuum phases	3	
Temperature	132 °C (270 °F)	134 °C (273 °F)
Exposure time	4 minutes	3 minutes
Drying time	30 minutes	

Packaging – Post-Sterilization

Carefully examine the sterile instrument packaging prior to use, ensuring package integrity is maintained.

Do not use if package is opened or damaged. If so, repeat the reprocessing steps. The sterile barrier system should comply with ISO/TS 16775 and ISO 11607 Part 1 and 2.

HOW SUPPLIED

The device is supplied non-sterile and is reusable.
Store individual component sets together. Store in a dark, dry, cool place.
Avoid extended exposure to light.

INSPECTION OF DEVICE

- Visually inspect for damage and wear.
- Identify and keep individual sets together.
- Assemble to ensure all parts of the Murdoch Mechanical Arm are functioning and complete.

Component Identification

Current sets:

- Both bolts have engraving on both sides of the head.
- Locking Handle has an "N" engraved on the hub.


Original Sets:

- Both bolts have single engraving on the head.
- Locking Handle is not engraved.

DEVICE PREPARATION

Assembly Instructions

To ensure correct operation it is necessary to correctly assemble the device, as follows:

Retractor Insertion Point Assembly

1. Take the Retractor Clamp Bolt (1), the shorter of the two bolts, and pass it through the "eye" feature at the end of the Positioning arm. (2)
2. Loosely screw a Locking Handle to the bolt to retain it within the positioning arm.
3. Position the Nathanson Liver Retractor within the patient. The retractor handle is clamped against the face of the Positioning Arm by tightening the Locking Handle.

Main Body Clamp Assembly

1. Note that on the Main Clamp (C) there is a tapered face. Place the tapered end of the Main Collar Washer (B) into this hole.
2. Take the Main Clamp bolt (A), which is the longer of the two bolts and pass this through the assembled clamp and washer. Loosely screw the remaining Locking Handle to the bolt to retain the assembly.
3. Align the hole through the Main Clamp Bolt (A) with the semi-circular notch in the Main Collar Washer. Pass the end of the positioning arm through this hole.
4. Pass the assembled clamp over the support pole and gently tighten the Locking Handle to grip the pole.
5. By loosening the Main Clamp Locking Handle sufficiently, adjustments to the placement of the Positioning Arm may be made in multiple directions simultaneously.

NOTE: There is a difference between the bolts and the Locking Handle that excludes the possibility of mixing components between product generations. Please note the numbers on the poles and keep parts from the same set together.

INSTRUCTIONS FOR USE

The MMA consists of eight components that are to be cleaned and sterilized prior to assembly (see **Cleaning and Sterilization Instructions**).

The components are to be assembled into an adjustable arm that holds a retractor, such as the NLR5 or LPR (see **Assembly Instructions**).

1. Insert the thicker 1/2" support pole into the OT table rail clamp.
2. Loosen the locking handle and swing the positioning arm at right angles to the support pole.
3. Adjust and gently tighten the main clamp to secure the positioning arm.
4. Insert the handle from the Nathanson Liver Retractor or Laparoscopic Pelvic Retractor, position as desired and tighten the locking handle of the positioning arm.
5. Adjust the positioning arm to the desired position and tighten the locking handle.
6. After final positioning fully tighten both locking handles.

Limitations on Reprocessing

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use. The expected lifetime of the MMA and spare parts is 10 years.

Reorder codes

MMA-1000 Murdoch Mechanical Arm

Spare parts

- Main Body Clamp: MMA-1000-MBC
- Main Collar Washer: MMA-1000-MCW
- Bolt and Handle Set: MMA-1000-BHS

The instructions provided above have been identified by the manufacturer of the medical device as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This normally requires validation and routine monitoring of the process.

MAINTENANCE

Apply surgical lubricants to threads, bolts and moving parts. Lubricant must be applied if components have been cleaned in an ultrasonic bath.

DISPOSAL OF DEVICES

Damaged components must be sterilized before disposal. This device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

The Cidezyme Xtra Multienzyme is an Advanced Sterilization Products product.

SERIOUS INCIDENT REPORTING

Any adverse event (clinical incident) involving this device should be reported to Cook Australia and the competent authority of the Member State immediately. To report an incident, call your local Cook representative or contact the Customer Relations Department on 1800 777 222 (toll free within Australia) or +61 7 3434 6777.

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© COOK 2021
COOK MEDICAL EUROPE LTD.