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WILLIAM A. COOK AUSTRALIA PTY. LTD.
95 Brandl Street, Brisbane Technology Park
Eight Mile Plains, Brisbane
QLD 4113 AUSTRALIA

www.cookmedical.com
© COOK 2014

IFU-TIVM/1

IFU-TIVM/1

ENGLISH

OTRIEVA TAPERED OVUM ASPIRATION NEEDLE

CONTAINS

A stainless steel tapered single lumen needle where the cannula is 17 G tapering down to 20 G at the distal tip. The needle has a bevelled tip which allows for puncture and an echogenic tip which enhances visualization under ultrasound.

Needle Gauge: 17 G tapering to 20 G at distal tip. The 17 G section is approximately 80% of the needle cannula length.

Needle Length: 30 or 35 cm

Aspiration Tubing: 90 cm length, 1.32 mm OD

Vacuum Tubing: 50 cm length, 2.06 mm OD

Connector type (e.g. for connection to pump): female luer fitting (6% taper)

INDICATIONS FOR USE

Used for laparoscopic or ultrasound guided transvaginal aspiration and flushing of oocytes from ovarian follicles.

CONTRAINDICATIONS

Not for use on a patient:

- With an active vaginal or intrauterine infection
- With a sexually transmitted disease
- With a recent uterine perforation
- With a recent caesarean section
- After recent pelvic inflammatory disease
- Who is currently pregnant

Not for intrafallopian tube procedures.

ADDITIONAL NOTES

Intended for one time use only. This is a sterile device and should be stored at room temperature away from direct sunlight. If the product package is open or damaged when received, do not use this device. The shelf life of the product is 3 years from the date of manufacture. The device has been two-cell MEA tested and passed with 80% or greater Blastocyst rate in 72 hrs. USP Endotoxin (LAL) tested and passed with 20 EU or less per device.

PRECAUTIONS

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Where possible, the needle tip should be kept within the stroma or follicles to prevent the aspiration of air into the needle. This minimises the potential for oocyte damage and 'frothing' in the test tube.

Hematuria may occur due to the aspiration needle penetrating a filled bladder during transvaginal ultrasound aspiration. This complication typically resolves spontaneously within a day. Extravasation of urine may occur within the abdominal cavity if a needle puncture traverses the bladder. Patients should be monitored for evidence of this known complication; however, infection may be introduced via needle puncture and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), uterine infection or cystitis. Vaginal bleeding has been reported to be associated with the transvaginal route for oocyte retrieval via needle aspiration. Bleeding is typically easily controlled with direct pressure.

INSTRUCTIONS FOR USE – ULTRASOUND GUIDED PROCEDURE

1. Position the patient in the lithotomy position on operating table. Local or general anaesthetic may be administered as necessary.
2. Carefully remove the needle from the packaging maintaining the sterility of the product.
3. The sterile needle should be inspected for tip sharpness and kinking of any supplied tubing.
4. Connect the needle's vacuum tubing to a vacuum pump.
5. Fit the silicone stopper onto the collection tube (designed to fit 15 mL Falcon tubes).
6. The aspiration system should be tested for patency by placing the tip in a spare test tube containing approximately 5 mL of culture medium and applying vacuum. If system is blocked or leaks, do not use. Replace device with another. Before proceeding, change the collection test tube.
7. Introduce an ultrasound transducer into the vaginal fornix to visualise the ovary and follicles. Identify the follicles to be aspirated. Check for the presence of blood vessels in and

around the ovary and determine a direct path into the ovarian follicles to be aspirated.

8. Insert the aspiration needle into the transducer needle guide and puncture vaginal wall. Ensure the tubing does not become kinked during use.
 9. Following visualisation of the follicle to be aspirated, line up the follicle using the needle guide on the ultrasound monitor and advance the needle tip into the centre of an ovarian follicle via a rapid, stabbing motion. The combination of the needle bevel and Echotip® enhances visualisation of the position of the needle tip. The handle indent indicates bevel orientation as well as facilitating grip.
 10. Apply vacuum to aspirate the follicular contents into the test tube. As the follicle collapses, rotate the needle tip within the follicle to curette the follicular walls. If required the follicle can be flushed, as described below. Repeat **steps 9 and 10** for remaining follicles.
- Follicle Flushing:
- (a) Use a non-toxic syringe filled with a follicle-flushing buffer (fluid volume of the needle and aspiration tube is approximately 5 mL). Remove the stopper from the test tube and insert the syringe into the aspiration line which can be identified as the protruding section from the base of the stopper.
 - (b) With the needle tip in the collapsed follicle, slowly inject (1-2 mL per second) the flushing medium to refill the follicle.
 - (c) Replace the stopper and aspirate the follicular contents.
11. Withdraw the needle from the patient and reposition the transducer to visualise the remaining ovary. Repeat **steps 8 to 10** then withdraw the needle.
 12. On completion of the aspiration procedure, remove the needle from the ultrasound guide, rinse with flushing buffer and then discard in an appropriate sharps container.

LAPAROSCOPIC PROCEDURE

1. Position the patient in the lithotomy position.
2. Follow **steps 2** through **6** as previously described in the **Ultrasound Guided Procedure Section**.
3. Under laparoscopic guidance, make a small skin nick to assist placement of a trocar/cannula.
4. Place trocar/needle assembly through the abdominal wall in the desired location. Remove trocar if used while leaving the cannula in place.
5. Under laparoscopic vision, place the needle through the abdominal cannula and advance the needle tip into the ovarian follicle.
6. Apply vacuum to aspirate the follicular contents into the test tube. If required, the follicle can be flushed, as described in **Step 10 Follicle Flushing** in the **Ultrasound Guided Procedure Section**.
7. When the desired number of oocytes has been aspirated, remove the needle, rinse, and discard.

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