| English                                               | IFU-TOWAKO/1 |
|-------------------------------------------------------|--------------|
| <b>THE TOWAKO TRANSMYOMETRIAL EMBRYO TRANSFER SET</b> |              |
| <span></span>                                         |              |

#### THE TECHNIQUE OF TRANSMYOMETRIAL EMBRYO TRANSFER INSTRUCTIONS FOR USE

The technique described here is that used by the groups from Birmingham, England, and Komatsu, Japan. Between them they have reported more than 1000 cycles using transmyometrial embryo transfer, with a clinical pregnancy rate of around 25 per cent (Kato et al., 1993; Kato, 1995; Sharif et al., 1996; Sharif and Kato, 1998).

#### SEDATION AND ANAESTHESIA

- No sedation or anaesthesia is used as the discomfort associated with the technique is minimal.
- The duration of the procedure is brief, with a mean "needle-in-uterus" time of 35 seconds (range 31 to 47 seconds).

#### PATIENT'S POSITION

- The lithotomy position, as used for transcervical transfer.
- The patient should empty her bladder immediately before the transfer, as this will allow optimal uterine visualisation using a transvaginal (TV) ultrasound scan.

#### EQUIPMENT USED

##### Ultrasound Scan Equipment:-

The same equipment used for TV ultrasound-guided oocyte recovery can be used for transmyometrial transfer. The availability of a dual-frequency transducer with zoom facility will enable the operator to concentrate on the area of interest (uterus) and reduce the "background noise" in the picture.

##### The Towako Needle Set contains:-

An echotipped needle, a matching stylet and a transfer catheter. The needle comes in either 18 gauge (1.25 mm) or 19 gauge (1.06 mm) and is very sharp to make insertion more smooth and hence less painful. The EchoTip® needle enhances ultrasound visualisation (**Figure 1**). The matching stylet fits into the needle, with both bevels flush. The needle is inserted with the stylet in situ to prevent blockage of the needle lumen with tissue debris during insertion. The 2 French polyethylene transfer catheter protrudes 1 mm beyond the needle tip when its hub is pushed firmly into the Luer lock fitting of the needle.

**USAGE:** Used in conjunction with a transvaginal ultrasound transducer to place embryos into the endometrium or uterine cavity.

**CONTRAINDICATIONS:** This device should not be used on a patient with an active vaginal or intrauterine infection, a sexually transmitted disease, a recent uterine perforation, a recent caesarean section, or who is currently pregnant.

##### ADDITIONAL NOTES:

- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- If the product package is open or damaged when received, do not use this device.
- The shelf life of the product is 3 years from the date of manufacture.

#### EMBRYO TRANSFER


- Position the patient carefully on the operating table or chair. The patient should empty her bladder immediately before the transfer, as this will allow optimal uterine visualisation using a transvaginal ultrasound scan.
- The vagina is cleaned and the TV transducer is prepared (sterile condom and needle guide) exactly as for TV ultrasound-guided oocyte collection.
- A TV scan is performed to visualize the uterus and locate the endometrial cavity. Characteristically, the endometrial ultrasound pattern is usually triple-lined (see **Figure 2**). The outer lines indicate the junction between the myometrium and endometrium, while the middle line indicates the endometrial cavity.
- The embryos are loaded into the tip of the catheter in 10 µl to 15 µl of culture medium.
- The needle fitted with the stylet, is then inserted transvaginally under ultrasound guidance aiming for the middle of the uterine cavity. This insertion should be quick and deliberate in one move.

**NOTE:** In an anteverted uterus the needle is inserted quickly (but precisely) through the anterior fornix, anterior uterine wall, anterior myometrium, anterior endometrial layer, cavity, posterior endometrial layer and then posterior myometrium. The needle is then pulled back till its tip is in the cavity (middle endometrial line, **Figure 1 & 3**). This is the pull-back technique, and it allows more accurate placement of the tip into the cavity than if it were attempted while the needle was first being pushed into the uterus. For a retroverted uterus the process is reversed, i.e. posterior fornix, posterior wall, and so on.


- When the needle is in place, the stylet is removed and the loaded catheter is inserted until its hub is pushed firmly against the Luer lock fitting of the needle.

**REMEMBER:** The catheter tip extends 1 mm beyond that of the needle (which is more echogenic) and a final adjustment of about a 1 mm pull of the needle is necessary before injecting the embryos.


- The transfer is evident on the ultrasound monitor by the brightness caused by the culture medium in which the embryos are contained — the flashing-light sign.
- The needle and catheter are then removed in one quick move and checked for inadvertent retention of embryos.

**Figure 1:** Transvaginal ultrasound scan showing needle echogenicity.

**Figure 2:** Triple line endometrial pattern on ultrasound.

**Figure 3:** Diagram showing placement of needle tip.

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