EN 17	Zenith® Universal Distal Body Endovascular Graft Instructions for Use
CS 22	Endovaskulární graft Zenith® Universal Distal Body Návod k použití
DA 27	Zenith® Universal Distal Body endovaskulær protese Brugsanvisning
DE 32	Zenith® Universal Distal Body endovaskuläre Prothese Gebrauchsanweisung
EL 38	Ενδαγγειακό μόσχευμα Zenith® Universal Distal Body Οδηγίες χρήσης
ES 44	Endoprótesis vascular Zenith® Universal Distal Body Instrucciones de uso
ET 50	Zenith® Universal Distal Body endovaskulaarne protees Kasutusjuhised
FI 55	Zenith® Universal Distal Body -endovaskulaarinen grafti Käyttöohjeet
FR 60	Endoprothèse vasculaire Zenith® Universal Distal Body Mode d'emploi
IT 65	Protesi endovascolari Zenith® Universal Distal Body Istruzioni per l'uso
LT 71	„Zenith® Universal Distal Body“ endovaskulinis protezas Naudojimo nurodymai
LV 76	Zenith® Universal Distal Body endovaskulārā protēze Lietošanas instrukcija
NL 81	Zenith® Universal Distal Body endovasculaire prothese Gebruiksaanwijzing
NO 87	Zenith® Universal Distal Body endovaskulært implantat Bruksanvisning
PL 92	Stent-graft wewnątrznaczyniowy Zenith® Universal Distal Body Instrukcja użycia
PT 98	Prótese endovascular Zenith® Universal Distal Body Instruções de utilização
RO 104	Grefă endovasculară Zenith® Universal Distal Body Instrucțiuni de utilizare
SK 109	Endovaskulární štep Zenith® Universal Distal Body Návod na použitie
SL 114	Endovaskularni vsadek Zenith® Universal Distal Body Navodila za uporabo
SV 119	Zenith® Universal Distal Body endovaskulärt graft Bruksanvisning

ENGLISH

TABLE OF CONTENTS

Illustrations	13-16
1 DEVICE DESCRIPTION	17
1.1 Overall Device Description	17
1.2 Implant	17
1.2.1 Implant Sizes and Configurations	17
1.2.2 Gold Marker Placement	17
1.3 Delivery System	17
1.4 Devices used In Combination with the Zenith Universal Distal Body	17
2 INDICATIONS FOR USE	17
3 CONTRAINDICATIONS	17
4 WARNINGS AND PRECAUTIONS	17
4.1 General	17
4.2 Patient Selection, Treatment and Follow-Up	17
4.3 Implant Procedure Requirements	18
4.4 MRI Safety and Compatibility	18
5 POTENTIAL ADVERSE EVENTS	18
5.1 Device Related Adverse Event Reporting	18
6 PATIENT SELECTION AND TREATMENT	18
7 PATIENT COUNSELING INFORMATION	18
8 HOW SUPPLIED	19
8.1 General	19
8.2 Device Sizes and Sizing Guidelines	19
8.2.1 Device Options	19
8.2.2 Sizing Guidelines	19
9 CLINICAL USE INFORMATION	19
9.1 Physician Training	19
9.2 Inspection Prior to Use	19
9.3 Materials Required	19
9.4 Recommended Devices	19
10 DIRECTIONS FOR USE	19
10.1 General Use Information	19
10.2 Pre-Implant Review	19
10.3 Patient Preparation	19
10.4 Device Preparation/Flush	20
10.5 Vascular Access and Angiography	20
10.6 Placement of Proximal Devices	20
10.7 Zenith Universal Distal Body Placement	20
10.8 Iliac Leg Graft or Iliac device (Contralateral) Placement	20
10.9 Rest of Zenith Universal Distal Body Deployment	20
10.10 Iliac Leg Graft or Iliac Device (Ipsilateral) Placement	20
10.11 Final Angiogram	20
10.12 Zenith Universal Distal Body Moulding (Optional)	20
11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP	20
11.1 General	20
11.2 CT Recommendations	21
11.3 Abdominal Radiographs	21
11.4 Ultrasound	21
11.5 MRI Safety and Compatibility	21
11.6 Additional Surveillance and Treatment	21

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
Figure 8

ZENITH® UNIVERSAL DISTAL BODY ENDOVASCULAR GRAFT

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

For the Zenith product line, there are several applicable Instructions for Use (IFU). This IFU describes the Zenith Universal Distal Body Endovascular Graft. For information on other Zenith devices, please refer to the appropriate Instructions for Use.

1 DEVICE DESCRIPTION

1.1 Overall Device Description

The Zenith Universal Distal Body Endovascular Graft (Zenith Universal Distal Body) consists of an implant loaded into a delivery system. The device is designed to be used in combination with other Zenith proximal grafts (for example, Zenith® t-Branch Thoracoabdominal Endovascular Graft (t-Branch) or Zenith® Fenestrated AAA Endovascular Graft (ZFEN-P)) and Iliac Leg Grafts (e.g., Zenith® Flex AAA Endovascular Graft Iliac Leg or Zenith Spiral-Z AAA Iliac Leg). Refer to **Figure 1**. The device has been tested to exhibit a functional design lifetime of ten years.

1.2 Implant

The Zenith Universal Distal Body implant is a bifurcated endovascular graft with one long ipsilateral iliac limb and one short contralateral iliac limb. The graft is constructed of woven polyester fabric sewn to self-expanding stainless steel Cook Z-Stents® with braided polyester and monofilament polypropylene suture, providing a conduit that is intended to exclude the aneurysm from blood flow. The graft is fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Refer to **Figure 2**.

The Zenith Universal Distal Body graft is intended to work as a modular system in combination with other available Zenith proximal devices (e.g. Zenith t-Branch) and iliac leg grafts (e.g. Zenith Flex AAA Endovascular Graft Iliac Leg, Zenith Spiral-Z AAA Iliac Leg).

Devices such as main body extensions, iliac leg extensions, converters, and iliac plugs may also be required during an aortic endovascular procedure. Each individual device has its own delivery system and Instructions for Use booklet.

1.2.1 Implant Sizes and Configurations

A range of proximal graft diameters and graft lengths are available to accommodate a range of Zenith proximal devices and to suit individual patient anatomies. Available sizes include a proximal graft diameter of 22 and 24 mm each with lengths of 81, 98, 115 and 132 mm.

1.2.2 Gold Marker Placement

To facilitate fluoroscopic visualisation, gold radiopaque markers are positioned on the graft. Four gold markers are positioned in a circumferential orientation within 2 mm of the most superior aspect of the graft material (proximal markers) and one gold marker is positioned at the end of the second proximal stent to aid with overlap positioning. There is a radiopaque marker at the graft bifurcation and gold markers are positioned in the shape of a tick (check) marker on the distal end of the contralateral limb. Refer to **Figure 3**.

1.3 Delivery System

The Zenith Universal Distal Body graft is shipped loaded onto the H&L-B One-Shot Introduction System. It has a sequential deployment method with built-in features to provide continuous control of the graft throughout the deployment procedure.

The proximal end of the graft is attached to the delivery system by three nitinol trigger wires. The distal end of the graft is also attached to the delivery system and held by an independent stainless steel trigger-wire. The trigger wires enable precise positioning and allow readjustment of the final graft position before full deployment.

The delivery system uses a 7.7 mm O.D. (20 Fr) H&L-B One-Shot Introduction System and is compatible with a 0.89 mm (0.035 inch) wire guide.

For added haemostasis, the Captor Haemostatic Valve can be opened or closed for the introduction and/or removal of other devices into and out of the sheath. The delivery system features a Flexor® introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta. Refer to **Figure 4**.

1.4 Devices used in Combination with the Zenith Universal Distal Body

The Zenith Universal Distal Body can be used in conjunction with other approved Zenith devices. Examples include the following:

- Zenith t-Branch Thoracoabdominal Endovascular Graft
- Zenith Fenestrated AAA Endovascular Graft
- Zenith Flex AAA Endovascular Graft Iliac Leg
- Zenith Spiral-Z AAA Iliac Leg.

Refer to the respective Instructions for Use for any of these devices.

2 INDICATIONS FOR USE

The Zenith Universal Distal Body Endovascular Graft is indicated for the endovascular treatment of patients with abdominal or thoracoabdominal aortic aneurysms, when used with the appropriate Zenith proximal device and Zenith iliac leg grafts.

The patient must have anatomy suitable for endovascular repair, including:

- Adequate femoral/iliac access compatible with the required delivery systems;
- An appropriate length of abdominal aorta, distal from the most proximal point of the overlap with the proximal piece to the aortic bifurcation.

3 CONTRAINDICATIONS

The Zenith Universal Distal Body Endovascular Graft is contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, nitinol, polyester, solder (tin, silver), polypropylene, urethane or gold
- Patients with a systemic infection may be at increased risk of endovascular graft infection.
- Patients who cannot tolerate contrast agents necessary for intra-operative and post-operative follow-up imaging.

- Patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- The Zenith Universal Distal Body should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. Specific training expectations are described in **Section 9.1, Physician Training**.
- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenosis from CT.
- The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up**.
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including:
 - 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and;
 - 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease.
 If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.
- The Zenith Universal Distal Body is not recommended in patients unable to undergo, or who will not be compliant with the necessary preoperative and post-operative imaging and implantation studies as described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up**.
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- A vascular surgery team should always be made available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

4.2 Patient Selection, Treatment and Follow-Up

The Zenith Universal Distal Body has not been evaluated in the following patient populations:

- Traumatic aortic injury
- Leaking, pending rupture or ruptured aneurysms
- Mycotic aneurysms
- Pseudoaneurysms resulting from previous graft placement
- Revision of previously placed endovascular grafts
- Uncorrectable coagulopathy
- Indispensable inferior mesenteric artery
- Genetic connective tissue disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
- Concomitant thoracic aortic aneurysms
- Patients with active systemic infections
- Pregnant or nursing females
- Morbidly obese patients
- Less than 18 years of age.

Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 7.7 mm OD (20 Fr) to 8.5 mm OD (22 Fr) vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolisation.

Key anatomic elements that may adversely affect successful exclusion of the aneurysm include severe aortic angulation, tortuosity and calcification.

The Zenith Universal Distal Body is not recommended in patients who cannot tolerate contrast agents necessary for intra-operative and post-operative follow-up imaging.

The Zenith Universal Distal Body is not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.

The Zenith Universal Distal Body is not recommended in patients with known sensitivities or allergies to stainless steel, polyester, solder (tin, silver), polypropylene, urethane or gold.

Patients with a systemic infection may be at increased risk of endovascular graft infection.

Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.

Patients with uncorrectable coagulopathy may have an increased risk of endoleak or bleeding complications.

4.3 Implant Procedure Requirements

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Universal Distal Body graft.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Universal Distal Body requires administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure post-operatively. Care should be taken to limit the amount of contrast media used during the procedure.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
- Inaccurate placement, incomplete sealing, or inadequate fixation of the Zenith Universal Distal Body graft within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the visceral vessels. Visceral vessel patency must be maintained to prevent/reduce the risk of renal failure, mesenteric ischemia, and subsequent complications.
- Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Inadequate overlap of the Zenith Universal Distal Body graft with a Zenith proximal device may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance. Vessel or catheter damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.

Take care during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus which can cause distal embolisation.

Before deployment of the proximal trigger wires, verify that the position of the access wire guide extends sufficiently into the aortic arch.

Care should be taken not to damage or disturb positioning of the graft after placement in the event recannulation of the graft is necessary.

Moulding Balloon Use (Optional)

Confirm complete deflation of balloon prior to repositioning.

Do not inflate balloon in any vessel outside of graft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labelling.

Use care when inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.

The Captor Haemostatic Valve should be opened prior to insertion and withdrawal of a moulding balloon.

4.4 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the Zenith Universal Distal Body graft is MR Conditional according to ASTM F2503. A patient with this device may be scanned safely after placement under the following conditions:

- Static magnetic field of 3.0 Tesla or 1.5 Tesla
- Maximum spatial magnetic gradient of 720 Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning or less (i.e., per scanning sequence)
- Normal operating mode.

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Systems:

In non-clinical testing, the Zenith Universal Distal Body graft produced a maximum temperature rise of 1.8°C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Magnetom, Software Numaris/4.) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.1 W/kg).

3.0 Tesla Systems:

In non-clinical testing, the Zenith Universal Distal Body graft produced a maximum temperature rise of 1.9°C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (General Electric Excite, HDx, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 50 mm of the position of Zenith Universal Distal Body graft as found during non-clinical testing using T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this stent.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

NOTE: For the Zenith Universal Distal Body graft, the clinical benefit of an MRI scan should be balanced with the potential risk of the procedure.

5 POTENTIAL ADVERSE EVENTS

Adverse events associated either with the Zenith Universal Distal Body or the implantation procedures that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolisation (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Extremity ischemia or neurologic complications (e.g., brachial plexus injury, claudication)
- Fever and localised inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, haematuria, infection)
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Occlusion of device or native vessel
- Open surgical conversion
- Radiation exposure
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Splenic injury (e.g., infarction, ischemia)
- Vascular access site complications, including infection, pain, haematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage
- Visceral vessel occlusion
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death).

5.1 Device Related Adverse Event Reporting

Any adverse event (clinical incident) involving the Zenith Universal Distal Body should be reported to **Cook Medical Australia** immediately. To report an incident, call your local Cook Aortic Intervention representative or contact the Quality Assurance Department on 1 800 777 222 (toll free within Australia) or +61 7 3841 1188.

6 PATIENT SELECTION AND TREATMENT

(Refer to **Section 4, Warnings and Precautions**)

The Zenith Universal Distal Body sizes are to be selected as described in this document and all warnings and precautions outlined in **Section 4** are to be taken into consideration.

All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. Additional considerations for patient selection include but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith Universal Distal Body
- Ability to tolerate general, regional or local anaesthesia
- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 7.7 mm OD (20 Fr) and 8.5 mm OD (22 Fr) vascular introducer sheath
- Aneurysm anatomy suitable for endovascular repair as specified in **Section 2** above
- The final treatment decision is at the discretion of the physician and patient.

7 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure including:

- Risks and differences between endovascular repair and surgical repair
- Potential advantages of traditional open surgical repair

- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, Imaging Guidelines and Post-Operative Follow-Up.**
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of abdominal aortic aneurysms (AAA) or thoracoabdominal aortic aneurysms (TAAA). At a minimum, annual imaging and adherence to routine post-operative follow-up requirements is required and should be considered a life-long commitment to the patient's health and well-being.
- Physicians must advise all patients that it is important to seek prompt medical attention if he/she experiences signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.

Product Code	Proximal Diameter	Length	Distal Limb Diameter	French Size	Introduction Sheath (I.D. / O.D.)	Length
UNIBODY-22-81	22 mm	81 mm	11 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
UNIBODY-22-98	22 mm	98 mm	11 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
UNIBODY-22-115	22 mm	115 mm	11 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
UNIBODY-22-132	22 mm	132 mm	11 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
UNIBODY-24-81	24 mm	81 mm	12 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
UNIBODY-24-98	24 mm	98 mm	12 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
UNIBODY-24-115	24 mm	115 mm	12 mm	20 Fr	6.7 mm / 7.7 mm	40 cm
UNIBODY-24-132	24 mm	132 mm	12 mm	20 Fr	6.7 mm / 7.7 mm	40 cm

8.2.2 Sizing Guidelines

A 2 mm–4 mm interference fit, combined with a minimum 2 stent overlap, is typical in the use of a distal bifurcated endovascular graft with a proximal component to reduce the risk of component separation.

9 CLINICAL USE INFORMATION

9.1 Physician Training

CAUTION: Always have a vascular surgery team available during implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Universal Distal Body should only be used by physicians and teams trained in aortic endovascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Universal Distal Body are outlined below:

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) or thoracoabdominal aortic aneurysms (TAAA) and co-morbidities associated with endovascular aneurysm repair
- Knowledge of radiographic image interpretation, device selection and sizing
- A multi-disciplinary team that has combined procedural experience with:
 - Femoral cutdown, arteriotomy and repair
 - Percutaneous access and closure techniques
 - Non-selective and selective wire guide and catheter techniques
 - Fluoroscopic and angiographic image interpretation
 - Embolisation
 - Angioplasty
 - Endovascular stent placement
 - Appropriate use of radiographic contrast material
 - Techniques to minimise radiation exposure
 - Expertise in necessary patient follow-up modalities

9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilisation barrier has been damaged or broken. If damage has occurred, do not use the product and return to **Cook Medical Australia**. Prior to use, verify correct device (size) has been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient. Also ensure the expiration date has not been surpassed and the device has been stored as per storage conditions printed on the device label.

The device is to be used in conjunction with materials and devices outlined in the following sections.

9.3 Materials Required

(Not included in system)

The following products are recommended to assist in the implantation of the Zenith Universal Distal Body. For information on the use of these products, refer to the individual product's Instructions for Use:

- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinised saline solution
- Sterile gauze pads

8 HOW SUPPLIED

8.1 General

The Zenith Universal Distal Body is supplied STERILE (100% Ethylene Oxide) in peel-open pouches (sterile inner pouch).

This device is intended for single patient use only. Do not reuse, reprocess or resterilise any part of this device. Reuse, reprocessing, or resterilisation may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death.

Reuse, reprocessing, or resterilisation may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilisation barrier has been damaged or broken. If this occurs, return the device to William A. Cook Australia Pty Ltd.

The label includes the device code, components supplied, expiry date and storage requirements to be followed for this device.

- The device is loaded into a 7.7 mm OD (20 Fr) introducer sheath.
- The sheath's surface is treated with a hydrophilic coating that, when hydrated, enhances trackability.
- The device is supplied with a tip protector (tip end), needle hub protector (handle end) and shipping stylet (handle end) that are to be removed prior to use.

8.2 Device Sizes and Sizing Guidelines

8.2.1 Device Options

The product code for the Zenith Universal Distal Body Endovascular Graft is UNIBODY-[D1]-[L1]. All graft options for the Zenith Universal Distal Body are listed here:

9.4 Recommended Devices

(Not included in system)

- 0.89 mm (0.035") extra stiff wire guide, 260 cm; for example:
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- 0.89 mm (0.035") standard wire guide; for example:
 - Cook 0.035 inch wire guides
 - Cook Nimble™ Wire Guides
- Moulding Balloons; for example:
 - Cook Coda® Balloon Catheter
- Introducer sets; for example:
 - Cook Check-Flo® Introducer Sets
 - Cook Extra Large Check-Flo® Introducer Sets
- Sizing catheter; for example:
 - Cook Auros® Centimetre Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
 - Cook Beacon® Tip Angiographic Catheters
 - Cook Beacon® Tip Royal Flush Catheters
- Entry needles; for example:
 - Cook single wall entry needles.

10 DIRECTIONS FOR USE

Prior to use of the Zenith Universal Distal Body, review this Instructions for Use booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

The suggested general order of implantation is:

1. Zenith Proximal Body Graft
2. Covered Visceral Stents
3. Universal Distal Body Graft
4. Contralateral Iliac Leg or Iliac Device
5. Ipsilateral Iliac Leg or Iliac Device.

10.1 General Use Information

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Universal Distal Body. The Zenith Universal Distal Body is compatible with 0.89 mm (0.035") diameter wire guides.

Endovascular stent grafting is a surgical procedure and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the haemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the grey positioner. After the grey positioner has been removed, if blood loss is excessive consider placing an uninflated moulding balloon or an introduction system dilator within the valve, restricting flow.

10.2 Pre-Implant Review

Verify from planning that the correct device has been selected.

10.3 Patient Preparation

1. Refer to institutional protocols relating to anaesthesia, anticoagulation and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualisation from the aortic arch to the femoral bifurcations.
3. Gain access to both femoral arteries using standard surgical technique or percutaneous technique.
4. Establish adequate proximal and distal vascular control of both femoral arteries.

10.4 Device Preparation/Flush

1. Remove black-hubbed shipping stylet (from the inner cannula), cannula protector tube (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away® sheath from back of the haemostatic valve. See **Figure 5**.
2. Elevate tip of system and flush through the stopcock on the haemostatic valve with heparinised saline until fluid emerges from the flushing groove in the tip. Refer to **Figure 6**. Continue to inject a full 20 ml of flushing solution through the device. Discontinue injection and close stopcock on connecting tube.
3. Attach syringe with heparinised saline to the hub on the inner cannula. Elevate the tip of the system and flush until fluid exits the dilator tip. Refer to **Figure 7**.
4. Soak a gauze pad in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate sheath liberally.

10.5 Vascular Access and Angiography

1. Puncture the selected common femoral arteries using standard technique with an 18 or 19UT gauge arterial needle. Upon vessel entry, insert:
 - Wire guides – standard 0.89 mm (0.035") diameter, 260 cm long, soft tip Wire Guide
 - Appropriate size sheaths (e.g., 2.0 mm I.D. (6 Fr) or 2.7 mm I.D. (8 Fr))
 - Flush catheter (often radiopaque sizing catheters – e.g., Centimetre Sizing Catheter or straight flush catheter)
2. Perform angiography as needed to identify positions of the coeliac, SMA, renal arteries, and the aortic bifurcation.

NOTE: If fluoroscope angulation is used with an angulated neck it may be necessary to perform angiograms using various projections.

10.6 Placement of Proximal Devices

1. Implant any proximal grafts (example Zenith t-Branch or Zenith Fenestrated AAA Endovascular Grafts) as required. Refer to the Instructions for Use of these devices for preparation and deployment instructions.

10.7 Zenith Universal Distal Body Placement

1. Ensure the delivery system has been flushed with heparinised saline and that all air is removed from the system.
2. Give systemic heparin if required. Flush each catheter and/or wipe the wire guide after any exchange.

NOTE: Monitor the patient's coagulation status throughout the procedure.

3. A 0.89 mm (0.035") stiff wire guide, at least 260 cm long, should be in place on the ipsilateral side to guide the delivery system into the aorta. Maintain wire guide position.
4. Before insertion, position delivery system on patient's abdomen under fluoroscopy to assist with orientation and positioning. The side arm of the haemostatic valve may serve as an external reference to the contralateral limb radiopaque marker.
5. Advance delivery system over the wire guide until the contralateral limb is positioned above and slightly anterior to the origin of the contralateral common iliac. If the contralateral limb radiopaque marker is not properly aligned, rotate the entire system until it is correctly positioned half way between a lateral and an anterior position on the contralateral side.

NOTE: Zenith Universal Distal Body delivery system will not pass through the sheath used to deliver the proximal body. The proximal body delivery sheath must be removed prior to insertion of the Zenith Universal Distal Body delivery system.

CAUTION: Maintain wire guide position during delivery system insertion.

CAUTION: To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).

6. Verify position of wire guide in the thoracic aorta.
7. Perform an angiogram to verify the position and orientation, including:
 - The degree of overlap with proximal body (no less than 2 stents)
 - The position of the contralateral limb (on contralateral side and slightly anterior) and proximal to the aortic bifurcation.

8. Ensure the Captor Haemostatic Valve is turned to the open position.

NOTE: The Captor Haemostatic Valve must only be turned until it stops and the open/close marks must be correctly aligned. Do not overturn the Captor Haemostatic Valve as overturning will damage the valve which may result in significant blood loss.

9. Stabilize the gray positioner (the shaft of the delivery system) and begin withdrawal of the outer sheath. Deploy the first two stents by withdrawing the sheath while monitoring device location. Proceed with deployment until contralateral limb is fully deployed.

NOTE: Tick marker on the contralateral limb of the distal bifurcated body is used to determine anterior/posterior orientation of the contralateral limb. It is not intended to line up with the anterior tick mark on the proximal body.

10. Advance a contralateral catheter and wire guide into the common iliac artery to a level below the short contralateral limb and then manipulate the wire guide into the contralateral limb and into the Zenith Universal Distal Body. AP and oblique fluoroscopic views can aid in verification of device cannulation.
11. Advance an angiographic catheter into the body of the graft. Perform angiography to confirm correct position inside the Zenith Universal Distal Body. Advance the catheter to where the proximal end of the Zenith Universal Distal Body graft is attached to the introducer.
12. Remove the safety lock from the black trigger-wire release mechanism. Under fluoroscopy, withdraw the release mechanism off the handle to remove the trigger-wire and release the graft's proximal attachment. Remove via its slot over the inner cannula. Refer to **Figure 8**.

NOTE: The distal stent is still secured by the distal trigger-wire.

NOTE: If resistance is felt or system bowing is noticed, the trigger wire is under tension. Excessive force may cause the graft position to be altered. If excessive resistance or delivery system movement is noted, stop and assess the situation.

NOTE: Technical Assistance from a Cook product specialist may be obtained by contacting your local Cook representative.

CAUTION: During trigger-wire removal and subsequent deployment, verify that the position of the access wire guide extends sufficiently into the aortic arch.

10.8 Iliac Leg Graft or iliac device (Contralateral) Placement

1. See IFU of the iliac leg graft or iliac device that has been selected to be placed in the contralateral side of the Zenith Universal Distal Body graft.

10.9 Rest of Zenith Universal Distal Body Deployment

1. Withdraw the rest of the Zenith Universal Distal Body sheath to expose the ipsilateral limb.
2. Remove the safety lock from the white trigger-wire release mechanism. Withdraw the release mechanism off the handle to remove the trigger-wire and detach the distal end of the endovascular graft from the delivery system. Remove via its slot over the device inner cannula.
3. Remove the grey positioner and dilator tip as one unit by pulling on the inner cannula. Keep the wire guide and sheath immobilized.
4. Close the Captor valve immediately after removal of the grey positioner and dilator tip.

NOTE: The Captor Haemostatic Valve must only be turned until it stops and the open/close marks must be correctly aligned. Do not overturn the Captor Haemostatic Valve as overturning will damage the valve which may result in significant blood loss.

10.10 Iliac Leg Graft or Iliac Device (Ipsilateral) Placement

See IFU of the iliac leg graft or device that has been selected to be placed in the ipsilateral side of the Zenith Universal Distal Body graft.

10.11 Final Angiogram

Position an angiographic catheter just above the level of the visceral arteries. Perform angiography to verify that all the arteries are patent and that there are no endoleaks. If endoleaks or other problems are observed, address as appropriate.

10.12 Zenith Universal Distal Body Moulding (Optional)

1. Prepare moulding balloon as follows:
 - Flush wire lumen with heparinised saline.
 - Remove all air from balloon.
 2. In preparation for the insertion of the moulding balloon, open the Captor Haemostatic Valve by turning it counter-clockwise.
 3. Advance the moulding balloon over the wire guide and through the haemostatic valve of the delivery system to the level of the overlap. Maintain proper sheath positioning.
- NOTE:** Captor Haemostatic Valve may be utilised to assist with haemostasis by turning in a clockwise rotation to the "close" position.
- NOTE:** Captor Haemostatic Valve should always be in the "open" position when repositioning of moulding balloon.
4. Expand the moulding balloon with diluted contrast media (as directed by the manufacturer) in the area of the overlap between the proximal device and Zenith Universal Distal Body, starting proximally and working in the distal direction.

NOTE: Do not inflate the moulding balloon outside the graft.

NOTE: Care should be taken when manipulating a moulding balloon in the vicinity of the visceral stents (if applicable).

NOTE: Confirm complete deflation of balloon prior to repositioning.

5. Withdraw the moulding balloon to the ipsilateral limb device overlap site and expand. Withdraw to the distal fixation site and expand.

CAUTION: Do not inflate balloon in iliac vessel outside of graft.

6. Deflate and remove moulding balloon. Transfer the moulding balloon onto the contralateral wire guide and into the contralateral iliac leg delivery system. Advance moulding balloon to the contralateral limb overlap and expand.

CAUTION: Confirm complete deflation of balloon prior to repositioning.

7. Withdraw the moulding balloon to the contralateral iliac leg/vessel distal fixation and expand.

CAUTION: Do not inflate balloon in iliac vessel outside of graft.

8. Remove moulding balloon and replace it with an angiographic catheter to perform completion angiograms.
9. Repair vessels and close access sites in standard surgical fashion.

11 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

11.1 General

The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and performance of their endovascular graft. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up. Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAA and TAAA.

Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in **Table 11.1**.

This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.

Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs, non-contrast CT and duplex ultrasound may be used.

- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (separation between components, stent fracture and barb separation).

• Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound. Ultrasound may be a less reliable and sensitive diagnostic method compared to CT.

Table 11.1 lists the minimum imaging follow-up for patients with the Zenith Universal Distal Body. Patients requiring enhanced follow-up should have interim evaluations.

Table 11.1 Recommended Imaging Schedule for Endograft Patients

	Angiogram	CT	Abdominal Radiographs
Pre-procedure	X ¹	X ¹	
Procedural	X		
Pre-discharge (within 7 days)		X ^{2,3,4}	X
1 month		X ^{2,3,4}	X
6 month		X ^{2,4}	X
12 month (annually thereafter)		X ^{2,4}	X

¹Imaging should be performed within 6 months before the procedure.

²Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan. With ultrasound, non-contrast CT is still recommended.

³Either pre-discharge or 1 month CT recommended.

⁴If Type I or III endoleak, prompt intervention and additional follow-up post-intervention recommended. See **Section 11.6, Additional Surveillance and Treatment**.

11.2 CT Recommendations

Image sets should include all sequential images at lowest possible slice thickness (≤ 3 mm). Do NOT perform large slice thickness (> 3 mm) and/or omit consecutive CT images sets, as it prevents precise anatomical and device comparisons over time. When using a multiphase scan, the table positions must match.

Non-contrast and contrast enhanced baseline and follow-up imaging are important for optimal patient surveillance. It is important to follow acceptable imaging protocols during the CT exam. **Table 11.2** lists examples of acceptable imaging protocols.

Table 11.2 Acceptable Imaging Protocols

	Non-Contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral CT or high performance MDCT capable of > 40 seconds	Spiral CT or high performance MDCT capable of > 40 seconds
Injection volume	n/a	150 cc
Injection rate	n/a	> 2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E. or equivalent
Coverage - start	Humeral Head	Humeral Head
Coverage - finish	Proximal femur	Proximal femur
Collimation	< 3 mm	< 3 mm
Reconstruction	2.5 mm throughout – soft algorithm	2.5 mm throughout – soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

11.3 Abdominal Radiographs

The following views are required:

- Four images: supine-frontal (AP), cross-table lateral, 30 degree LPO and 30 degree RPO views centred on umbilicus
- Record the table-to-detector distance and use the same distance at each subsequent examination.

Ensure entire device is captured on each single image format lengthwise.

If there is any concern about the device integrity (e.g., kinking, stent breaks, relative component migration), it is recommended to use magnified views. The attending physician should evaluate images for device integrity (entire device length including components) using 2-4X magnification visual aid.

11.4 Ultrasound

Ultrasound imaging may be performed in place of contrast CT when patient factors preclude the use of image contrast media (**NOTE:** Imaging is limited to the abdominal aorta). Ultrasound may be paired with non-contrast CT. A complete aortic duplex is to be videotaped for maximum aneurysm diameter, endoleaks, stent patency and stenosis. Included on the videotape should be the following information as outlined below:

- Transverse and longitudinal imaging should be obtained of the abdominal aorta demonstrating the celiac, mesenteric and renal arteries to the iliac bifurcations to determine if endoleaks are present utilizing colour flow and colour power Doppler (if accessible)
- Spectral analysis confirmation should be performed for any suspected endoleaks
- Transverse and longitudinal imaging of the maximum aneurysm should be obtained.

11.5 MRI Safety and Compatibility

For MRI safety and compatibility information refer to **Section 4.4**.

11.6 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, ≥ 5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy and the patient's personal choices. Patients should be counselled that subsequent re-interventions including catheter-based and open surgical conversion are possible following endograft placement.

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© COOK 2020
WILLIAM A. COOK AUSTRALIA PTY. LTD.
COOK MEDICAL EUROPE LTD.