COOK
MEDICAL


CE
0123

EN 4	Cystotome Instructions for Use
BG 6	Цистотом Инструкции за употреба
HR 9	Cistotom Upute za uporabu
CS 12	Cystotom Návod k použití
DA 15	Cystotom Brugsanvisning
NL 18	Cystotoom Gebruiksaanwijzing
ET 21	Tsüstotoom Kasutusjuhised
FR 24	Cystotome Mode d'emploi
DE 27	Zystotom Gebrauchsanweisung
EL 30	Κυστεοτόμος Οδηγίες χρήσης
HU 34	Hólyagmetsző Használati utasítás
IT 37	Cistotomo Istruzioni per l'uso
LV 40	Cistotomijas ierīce Lietošanas instrukcija
LT 43	Cistotomas Naudojimo instrukcija
NO 45	Cystotomi-instrument Bruksanvisning
PL 48	Cystotom Instrukcja użycia
PT 51	Cistótomo Instruções de utilização
RO 54	Cistotom Instrucţiuni de utilizare
SK 58	Cystotóm Návod na použitie
ES 60	Cistótomo Instrucciones de uso
SV 64	Cystotom Bruksanvisning

I F U 0 0 0 5 - 1 2

Fig. 1

Proximal Electrode (Part A)
Проксимален електрод (Част А)
Proksimalna elektroda (dio A)
Proximální elektroda (část A)
Proksimal elektrode (del A)
Proximale elektrode (deel A)
Proksimaalne elektrood (osa A)
Électrode proximale (partie A)
Proximale Elektrode (Teil A)
Εγγύς ηλεκτρόδιο (μέρος Α)
Proximális elektróda („A” rész)
Elettrodo prossimale (A)
Proksimālais elektrods (A daļa)
Proksimalinis elektrods (A dalis)
Proksimal elektrode (del A)
Elektroda proksymalna (część A)
Eléctrodo proximal (componente A)
Electrod proximal (Partea A)
Proximálna elektróda (časť A)
Electrodo proximal (parte A)
Proximal elektrod (del A)

Distal Electrode (Part B)
Дистален електрод (Част Б)
Distalna elektroda (dio B)
Distální elektroda (část B)
Distal elektrode (del B)
Distale elektrode (deel B)
Distaalne elektrood (osa B)
Électrode distale (partie B)
Distale Elektrode (Teil B)
Περιφερικό ηλεκτρόδιο (μέρος Β)
Disztális elektróda („B” rész)
Elettrodo distale (B)
Distālais elektrods (B daļa)
Distalinis elektrods (B dalis)
Distal elektrode (del B)
Elektroda dystalna (część B)
Eléctrodo distal (componente B)
Electrod distal (Partea B)
Distálna elektróda (časť B)
Electrodo distal (parte B)
Distal elektrod (del B)

Diathermic ring	.038 inch Needle knife
Диатермичен пръстен	0,038 inch иглен нож
Dijatermički prsten	Igličasti skalpel od 0,038 inch
Diatermický kroužek	Jehlový papilotom 0,038 inch
Diatermiring	0,038 inch kanylekniv
Diathermische ring	0,038 inch naaldmes
Diatermiline rōngas	0,038-tolline (inch) punkteerimisnuga
Anneau diathermique	Papillotome-couteau-aiguille de 0,038 inch
Diathermie-Ring	Nadelmesser 0,038 Inch
Διαθερμικός δακτύλιος	Βελονοτόμος 0,038 inch
Diatermiküs gyűrű	0,038 inch-es tűkés
Anello diatermico	Ago di taglio da 0,038 inch
Diatermijas gredzens	0,038 inch adatas asmens
Diaterminis žiedas	0,038 inch adatinis peilis
Diatermisk ring	0,038 inch nålekniv
Pierścień diatermiczny	Nóż igłowy o średnicy 0,038 inch
Anel diatérmico	Agulha cortante de 0,038 inch
Inel diathermic	Partea activă a acului, de 0,038 inch
Diatermický krúžok	Nôž s ihlou veľkosti 0,038 inch
Anillo diatérmico	Papilótomo de aguja de 0,038 inch
Diatermiring	0,038 inch nålkniv

10 Fr Outer catheter	5 Fr Inner catheter
Външен катетър 10 Fr	Вътрешен катетър 5 Fr
Vanjski kateter veličine 10 Fr	Unutarnji kateter veličine 5 Fr
Vnější katetr 10 Fr	Vnitřní katetr 5 Fr
10 Fr ydre kateter	5 Fr indre kateter
10 Fr buitenste katheter	5 Fr binnenste katheter
10 Fr väline kateeter	5 Fr sisemine kateeter
Cathéter externe 10 Fr	Cathéter interne 5 Fr
Außenkatheter (10 Fr)	Innenkatheter (5 Fr)
Εξωτερικός καθετήρας 10 Fr	Εσωτερικός καθετήρας 5 Fr
10 Fr méretű külső katéter	5 Fr méretű belső katéter
Catetere esterno da 10 Fr	Catetere interno da 5 Fr
10 Fr ārējais katetrs	5 Fr iekšējais katetrs
10 Fr išorinis kateteris	5 Fr vidinis kateteris
10 Fr ytre kateter	5 Fr indre kateter
Cewnik zewnętrzny 10 Fr	Cewnik wewnętrzny 5 Fr
Cateter externo de 10 Fr	Cateter interno de 5 Fr
Cateter extern de 10 Fr	Cateter intern de 5 Fr
Vonkajší katéter veľkosti 10 Fr	Vnútorný katéter veľkosti 5 Fr
Catéter exterior de 10 Fr	Catéter interior de 5 Fr
10 Fr ytterkateter	5 Fr innerkateter

INTENDED USE

This device is designed to electrosurgically puncture a hole in the transgastric or transduodenal wall and into a pancreatic pseudocyst, when it is visibly bulging into the gastrointestinal tract.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than the stated intended use.

If the package is open or damaged when received, do not use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

This device must only be used with an active cord compatible with a 3mm diameter connector. This device is compatible with the following Cook active cords: ACU-1 and ACU-1-VL.

CAUTION

U.S. Federal Law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

Contraindications include those specific to blood coagulation disease, interposing vessels between the pseudocyst wall and that of the stomach or the duodenum. If the pseudocyst is < 4cm in diameter do not proceed.

WARNINGS

This device contains nickel, which may cause allergic reaction in individuals with nickel sensitivity.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with GI endoscopy include, but are not limited to: allergic reaction to medication • allergic reaction to nickel • aspiration • burn • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest • sepsis.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Maximum rated input voltage for this device is 2.0 kVp-p.

Any electrosurgical accessory constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration • burns • nerve and/or muscle stimulation • cardiac arrhythmia.

Before using this device, follow recommendations provided by electrosurgical unit manufacturer to ensure patient safety through proper placement and utilisation of patient return electrode. Ensure a proper path from patient return electrode to electrosurgical unit is maintained throughout procedure.

Examine all components to be used during procedure. Do not use a device that has been cut, burned or damaged. Damaged device insulation may cause unsafe currents in either patient or operator.

Switch electrosurgical unit to "off" position when not in use.

When applying current, ensure needle knife is completely out of endoscope. Contact of needle knife with endoscope may cause grounding, which can result in patient injury, operator injury, a broken needle knife and/or damage to endoscope.

INSTRUCTIONS FOR USE

Illustrations

1. Upon removing from package, uncoil device. **Note:** Do not extend or retract the needle knife while the device is coiled as this may cause damage to the device. Extend and retract the needle when the device is straight.
2. Examine features of device. It consists of a 10 French outer catheter with a diathermic ring, a 5 French inner catheter with a removable .038 inch needle knife and electrodes located distal and proximal to the handle. Proximal electrode (Part A) is indicated for needle knife, the distal electrode (Part B) is indicated for the diathermic ring (see Fig. 1).
3. With the electrosurgical unit off, prepare equipment. Active cord fittings should fit snugly into both device handle and electrosurgical unit. **Note:** Set the unit to pure cutting current (80-120 watts).
4. With the needle knife fully retracted into sheath, introduce device into endoscope accessory channel. Advance device in short increments until it is endoscopically visualized exiting the scope.
5. Once tip is visualized exiting scope, position catheter and extend needle knife to desired length. Identify position of cyst bulging into stomach and/or duodenum. With electrosurgical unit off determine position of puncture site.
6. Following electrosurgical unit manufacturer's instructions, verify desired settings and activate electrosurgical unit.
7. Puncture pseudocyst with needle knife using the pure cutting current.
8. Once puncture is complete, turn electrosurgical unit off and disconnect active cord from device handle.

9. With needle in the cyst, inject contrast to fill pseudocyst under fluroscopic visualisation.
10. Remove needle knife by unlocking proximal electrode from handle. **Note:** Before each introduction of a wire guide into a device, wet wire guide in a sterile water or saline bath. A contrast filled catheter may make wire guide difficult to advance. Flush contrast from catheter with sterile water or saline before each introduction of wire guide. Advance wire guide into pancreatic pseudocyst to facilitate the introduction of a stent or drainage set.
11. Push 10 French outer catheter while slightly pulling back inner catheter.
12. When diathermic ring comes in contact with digestive wall, place active cord on distal electrode (Part B) and advance into cyst, using pure cutting current (80-120 watts). **Caution: To avoid electric shock, ensure that ring does not come in contact with wire guide. When electrosurgical current is being applied, do not advance the diathermic ring beyond the proximal catheter marking.**
13. Turn off electrosurgical unit and disconnect active cord from side arm and from electrosurgical unit. Wipe active cord with damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping active cord tightly may damage device.
14. Fluid from cyst can be collected for analysis. Device can be removed leaving wire guide in place to facilitate introduction of stents or Drainage set.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© COPYRIGHT COOK 2019 ©
COOK IRELAND LTD.