Cystotome

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Cystotom

Cystotom

Cystotoom

Tsüstotoom

Cystotome

Cystotome

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English

INTENDED USE

This device is designed to electrosurgically puncture a hole in the transgastric or transduodenal wall and into a pancreatic pseudocyst, when it is visibly bulging into the gastrointestinal tract.

NOTES

Do not use this device for any purpose other than the stated intended use.

If the package is open or damaged when received, do not use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CAUTION

U.S. Federal Law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

Contraindications include those specific to blood coagulation disease, interposing vessels between the pseudocyst wall and that of the stomach or the duodenum. If the pseudocyst is < 4cm in diameter do not proceed.

POTENTIAL COMPLICATIONS

Potential complications associated with GI endoscopy include, but are not limited to: sepsis, perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Maximum rated input voltage for this device is 2.0 kVp-p.

Any electrosurgical accessory constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

Before using this device, follow recommendations provided by electrosurgical unit manufacturer to ensure patient safety through proper placement and utilisation of patient return electrode. Ensure a proper path from patient return electrode to electrosurgical unit is maintained throughout procedure.

Examine all components to be used during procedure. Do not use a device that has been cut, burned or damaged. Damaged device insulation may cause unsafe currents in either patient or operator.

Switch electrosurgical unit to "off" position when not in use.

When applying current, ensure needle knife is completely out of endoscope. Contact of needle knife with endoscope may cause grounding, which can result in patient injury, operator injury, a broken needle knife and/or damage to endoscope.

INSTRUCTIONS FOR USE

Illustrations

1. Upon removing from package, uncoil device. **Note:** Do not extend or retract the needle knife while the device is coiled as this may cause damage to the device. Extend and retract the needle when the device is straight.
2. Examine features of device. It consists of a 10 French outer catheter with a diathermic ring, a 5 French inner catheter with a removable .038 inch needle knife and electrodes located distal and proximal to the handle. Proximal electrode (Part A) is indicated for needle knife, the distal electrode (Part B) is indicated for the diathermic ring (see Fig. 1).
3. With the electrosurgical unit off, prepare equipment. Active cord fittings should fit snugly into both device handle and electrosurgical unit. **Note:** Set the unit to pure cutting current (80-120 watts).
4. With the needle knife fully retracted into sheath, introduce device into endoscope accessory channel. Advance device in short increments until it is endoscopically visualized exiting the scope.
5. Once tip is visualized exiting scope, position catheter and extend needle knife to desired length. Identify position of cyst bulging into stomach and/or duodenum. With electrosurgical unit off determine position of puncture site.
6. Following electrosurgical unit manufacturer's instructions, verify desired settings and activate electrosurgical unit.
7. Puncture pseudocyst with needle knife using the pure cutting current.
8. Once puncture is complete, turn electrosurgical unit off and disconnect active cord from device handle.

9. With needle in the cyst, inject contrast to fill pseudocyst under fluroscopic visualisation.

10. Remove needle knife by unlocking proximal electrode from handle.

Note: Before each introduction of a wire guide into a device, wet wire guide in a sterile water or saline bath. A contrast filled catheter may make wire guide difficult to advance. Flush contrast from catheter with sterile water or saline before each introduction of wire guide. Advance wire guide into pancreatic pseudocyst to facilitate the introduction of a stent or drainage set.

11. Push 10 French outer catheter while slightly pulling back inner catheter.

12. When diathermic ring comes in contact with digestive wall, place active cord on distal electrode (Part B) and advanced into cyst, using pure cutting current (80-120 watts). **Caution: To avoid electric shock, ensure that ring does not come in contact with wire guide. When electrosurgical current is being applied, do not advance the diathermic ring beyond the proximal catheter marking.**

13. Turn off electrosurgical unit and disconnect active cord from side arm and from electrosurgical unit. Wipe active cord with damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping active cord tightly may damage device.

14. Fluid from cyst can be collected for analysis. Device can be removed leaving wire guide in place to facilitate introduction of stents or Drainage set.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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Cook is a registered trademark of Cook Incorporated.
Cook Ireland Ltd.