COOK®

MEDICAL

IFU0019-8

Universa™ Stent Set
Комплект стент Universa™
Souprava stentu Universa™
Universa™ stentsæt
Universa™ stentset
Stentide komplett Universa™
Set d'endoprothèse Universa™
Universa™-Stent-Set
Σετ ενδοπροσθέσεων Universa™
Universa™ sztentkészlet
Set con stent Universa™
Universa™ stenta kompletst
Universa™ stento rinkinys
Universa™ stentsett
Zestaw stentu Universa™
Conjunto de stent Universa™
Setul de stenturi Universa™
Súprava stentov Universa™
Equipo de stents Universa™
Universa™ stentset

English	1
Bulgarian/Български	5
Czech/Česky	10
Danish/Dansk	14
Dutch/Nederlands	18
Estonian/Eesti	22
French/Français	26
German/Deutsch	30
Greek/Ελληνικά	35
Hungarian/Magyar	40
Italian/Italiano	44
Latvian/Latviski	48
Lithuanian/Lietuviškai	52
Norwegian/Norsk	56
Polish/Polski	60
Portuguese/Português	64
Romanian/Română	68
Slovak/Slovenčina	72
Spanish/Español	76
Swedish/Svenska	80

English

INTENDED USE

Used for temporary internal drainage from the ureteropelvic junction to the bladder. Ureteral stents have been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic, percutaneous, or open surgical techniques. The multi-length stents will accommodate ureters from approximately 22-32cm in length.

NOTES

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Inventory rotation of sterile products is essential. Verify the expiration date of the package label prior to using the product. If the expiration date has lapsed, do not use or resterilize the device. Store in a dark, dry location, away from temperature extremes.

SET CONTAINS *(set components may vary)*

- Radiopaque double pigtail/multi-length stent
- Stent Positioner
- Roadrunner hydrophilic wire guide (-RPC) or stainless steel wire guide
- Lockable Pusher (-LP stents)
- Closed End Stent (-CE stents)
- Pigtail Sideports only (-P stents)
- Distal Pigtail Sideports only (-PD stents)
- Tether (-T)
- Stent only (-B)

Note:

If no wire guide is provided with this set the following is recommended;

- 5.0 French stents accept .035" diameter
- 6.0, 7.0, 8.0 French stents accept .038" diameter

CONTRAINDICATIONS

There are no known contraindications.

CAUTIONS

Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to individual patients. Informed consent should be obtained to maximize patient compliance with follow-up procedures.

U.S. Federal Law restricts this device to sale by or on the order of a physician.

SUGGESTED INSTRUCTIONS FOR ACTIVATING HYDROPHILIC COATING

The hydrophilic coating on this wire guide is activated by immersion in sterile water or sterile saline solution. **Note:** Hydrophilic coated wire guides are intended for one-time use.

SUGGESTED INSTRUCTIONS FOR USE

1. Prior to using the wire guide, fill a 10 ml syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder.
 2. Inject sterile liquid into the holder until it is full and liquid can be seen exiting around the tip of the wire guide. **Note:** Hydrophilic coating is not permanent. If after extended use the coated wire guide does not perform smoothly, replace the wire guide with a new coated guide.
- 

ASSEMBLY FOR STENT INSERTER AND RELEASE SLEEVE

CAUTION: It is imperative the filiform flexible end of the wire guide is introduced into the ureter first.

1. Using a baseline pyelogram, estimate the proper stent length; add 1 cm to that estimated ureteral measurement. Accurate measurement enhances drainage efficiency and patient comfort.
2. Place the stent inserter inside the release sleeve and secure the Luer lock fitting.
3. Introduce the flexible end of the wire guide into the fitting end of the stent inserter and advance it until the flexible end of the wire guide extends half its length beyond the end of the inserter/release sleeve assembly. **Note:** The lockable fitting on the distal end of the stent inserter must be loosened to allow the wire guide to be introduced.

- 
- 
4. Carefully introduce the flexible end of the wire guide into the blunt end of the stent and advance it until the stent is straightened and the flexible tip of the wire guide extends beyond the tapered end of the stent.
 5. While maintaining the position of the stent and the wire guide, advance the inserter/release sleeve assembly to the blunt end of the stent. Push the tip of the stent inserter into the end of the stent. **Note:** The stent is pushed until it reaches the ink mark on the release sleeve.
 6. Check to ensure the flexible tip of the wire guide extends beyond the tip of the straightened stent and then tighten the lockable fitting to secure the position of the stent inserter on the wire guide.

SUGGESTED INSTRUCTIONS FOR USING UNIVERSA™ STENT SETS ENDOSCOPIC PLACEMENT

- 
1. Pass a flexible wire guide tip to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using a wire guide and an open-ended ureteral catheter in combination.
 2. Using a baseline pyelogram, estimate the proper stent length; add 1 cm to that estimated ureteral measurement. Accurate measurement enhances drainage efficiency and patient comfort.
 3. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Care should be taken not to allow the wire guide to advance into the renal parenchyma.
 4. When the distal end of the stent is observed at the ureterovesical junction, halt advancement of the stent. While stabilizing the stent with the positioner, remove the wire guide. The stent pigtail will form spontaneously. Carefully remove the positioner from the cystoscope. **Note:** If necessary, final adjustment can be performed with endoscopic forceps. The stent may be removed easily by gentle withdrawal traction using endoscopic forceps. Fluoroscopy facilitates stent placement; however, standard radiography may be used.
- 

CAUTIONS

Intended for one-time use only. The Universa™ Firm stent must not remain indwelling more than twelve (12) months. The Universa™ Soft stent must not remain indwelling more than six (6) months.

These stents are not intended as permanent indwelling devices. Periodic evaluation via cystoscopic, radiographic or ultrasonic means is suggested. Individual variations of interaction between

stents and the urinary system are unpredictable. The stent must be removed if encrustation hampers drainage. The stent may be replaced with a new stent, if the patient's status permits.

Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered. A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.

Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended. Scopes larger than 21.0 French are suggested.

REFERENCES

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," *The Journal of Urology*, 119 (1978), 731.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," *The Journal of Urology*, 121 (1979), 17.

M. F. Camacho, R. Pereiras, H. Carrion, M. Bondhus, V. A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single End Ureteral Stent," *Urology*, 13 (1979), 516-520.

M. A. Koyle, A. D. Smith: "Ureteral Stents," Chapter 13 in A. D. Smith, W. R. Castañeda-Zuñiga, J. G. Bronson: *Endourology: Principles and Practice*. Thieme Inc., New York, New York, 1986.

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© 2012 Cook Medical
Cook is a registered trademark of Cook Incorporated.
Cook Ireland Ltd.