EN 3	Resonance® Metallic Ureteral Stent Set Instructions for Use
BG 7	Комплект метален уретрален стент Resonance® Инструкции за употреба
CS 11	Souprava kovového ureterálního stentu Resonance® Návod k použití
DA 15	Resonance® ureterstentsæt af metal Brugsanvisning
NL 19	Resonance® metalen ureterale stentset Gebruiksaanwijzing
ET 23	Kusejuha metallstendi komplekt Resonance® Kasutusjuhised
FR 27	Set d'endoprothèse métallique urétérale Resonance® Mode d'emploi
DE 31	Resonance® Metall-Ureterstentset Gebrauchsanweisung
EL 35	Σετ μεταλλικής ουρητηρικής ενδοπρόσθεσης Resonance® Οδηγίες χρήσης
HU 39	Resonance® fém ureterális sztentkészlet Használati utasítás
IT 43	Set con stent ureterale metallico Resonance® Istruzioni per l'uso
LV 47	Resonance® metāla ureterālā stenta komplekts Lietošanas instrukcija
LT 51	Resonance® metalinio šlapimtakio stento rinkinys Naudojimo nurodymai
NO 55	Resonance® metallisk ureteralt stentsett Bruksanvisning
PL 59	Zestaw metalowego stentu moczowodowego Resonance® Instrukcja użycia
PT 63	Conjunto de stent ureteral metálico Resonance® Instruções de utilização
RO 67	Set al stentului ureteral metalic Resonance® Instrucțiuni de utilizare
SK 71	Súprava s kovovým ureterálnym stentom Resonance® Návod na použitie
ES 75	Equipo de stent ureteral metálico Resonance® Instrucciones de uso
SV 79	Resonance® metalluretärstentset Bruksanvisning

RESONANCE® METALLIC URETERAL STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Resonance Metallic Ureteral Stent Set contains (set components may vary):

- Resonance® metal stent with pigtail straightener
- Radiopaque introduction catheter and sheath
- HiWire® Nitinol Core wire guide (not included in the United States; optional component outside the United States)

INTENDED USE

Used for temporary stenting of the ureter in adult patients with extrinsic ureteral obstruction. Intended for one-time use.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- These stents are not intended as permanent indwelling devices.
- The stent must not remain indwelling more than twelve (12) months. If the patient's status permits, the stent may be replaced with a new stent.
- Patients should be checked at regular intervals utilizing techniques such as abdominal X-ray (KUB film). Patients using calcium supplements must be more closely monitored for possible stent encrustation. The stent must be removed if encrustation hampers drainage.
- Improper handling of the stent prior to insertion into the ureter may harm the functionality of the stent. Bending, stretching or any other type of improper handling may deform the stent. It is important that the stent is handled with care.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- Change in urine viscosity may hamper drainage
- Hematuria and Incontinence may indicate fistula formation

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated that the Resonance® Metallic Ureteral Stent is MR Conditional. It can be scanned safely under the following conditions:

At 1.5 Tesla Static Magnetic Field	At 3 Tesla Static Magnetic Field
• Spatial gradient field of 450 Gauss/cm or less; and	• Spatial gradient field of 720 Gauss/cm or less; and
• Maximum whole-body-averaged specific absorption rate (SAR) of 1.5 W/kg for 20 min. of scanning	• Maximum whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 min. of scanning
In non-clinical testing the stent produced a temperature rise during scanning of less than 0.80 °C at an SAR of 1.5 W/kg for 20 min.	In non-clinical testing the stent produced a temperature rise during scanning of less than 1.3 °C at an SAR of 3 W/kg for 15 min.

Note: MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the stent. The image artifact extends approximately 16 mm from the device when scanned in non clinical testing using: GRE sequence in a 3.0 Tesla Siemens Magnetom Trio, A Tim System (Software Numaris/4) with the body coil. Therefore it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with indwelling ureteral stents include:

- Diminished urine drainage / stent occlusion
- Fistula formation including Arteriuoreteral Fistula
- Hemorrhage
- Hydronephrosis

- Infection
- Loss of renal function
- Pain/discomfort
- Perforation of kidney, renal pelvis, ureter and /or bladder
- Peritonitis
- Pyuria
- Stent degradation/fracture
- Stent dislodgement/migration
- Stent encrustation
- Tissue ingrowth
- Ureteral reflux
- Urinary symptoms (frequency, urgency, incontinence, dysuria, hematuria)
- Urinary tract tissue erosion

INSTRUCTIONS FOR USE

Illustrations

Note: Use with wire guide with 0.038" diameter.

Note: Traditional plastic stents drain through a finite number of side ports. Resonance is designed with an infinite number of side ports (the coils) whereby urine may pass through or over any of the coils.

Note: Stent may be utilized with either an antegrade or retrograde placement technique. The figures shown (1 and 2) are for retrograde placement.

1. Over a previously placed wire guide, pass introduction catheter and sheath into appropriate position utilizing fluoroscopic guidance or direct vision (**Fig. 1**).

2. With introduction catheter and sheath in proper position, remove introduction catheter to allow passage of the stent (**Fig. 1**).

3. Use the pigtail straightener to introduce the stent to the sheath.

Note: Once stent is in the sheath pull the pigtail straightener away from the hub.

4. Use the introduction catheter to advance the stent through the sheath (**Fig. 2**) until the numbered marker corresponding to the stent length being deployed reaches the hub on the sheath (**Fig. 3**). At this point the first pigtail will have fully deployed from the sheath.

Note: The sheath tip is radiopaque to aid placement.

5. To prevent further deployment of the stent in the kidney hold the introduction catheter in place and retract the sheath.

6. When the hub of the sheath aligns with the proximal ink mark on the introduction catheter the second pigtail is about to deploy. At this point retract both the sheath and introduction catheter completely. (**Fig. 2**)

7. The stent may be removed using conventional cystoscopic techniques utilizing forceps or graspers. **Note:** Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered

HOW SUPPLIED

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization. Do not use this device for any purpose other than the stated intended use. Store in a dark, dry location away from temperature extremes. Inventory rotation of sterile products is essential. Verify the expiration date on the package label prior to using the product. If the expiration date has lapsed, do not use or resterilize the device. If the product package is open or damaged when received, do not use this device.

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© COPYRIGHT COOK 2012
COOK IRELAND LTD.