- EN
3 **Resonance® Metallic Ureteral Stent Set**
Instructions for Use
- BG
6 **Комплект метален уретерален стент Resonance®**
Инструкции за употреба
- HR
9 **Komplet metalnog ureteralnog stenta Resonance®**
Upute za uporabu
- CS
12 **Souprava kovového ureterálního stentu Resonance®**
Návod k použití
- DA
15 **Resonance® ureterstentsæt af metal**
Brugsanvisning
- NL
18 **Set met Resonance® metalen ureterale stent**
Gebruiksaanwijzing
- ET
22 **Kusejuha metallstendi komplekt Resonance®**
Kasutusjuhised
- FR
25 **Set d'endoprothèse métallique urétérale Resonance®**
Mode d'emploi
- DE
28 **Resonance® Metall-Ureterstentset**
Gebrauchsanweisung
- EL
31 **Σετ μεταλλικής ουρητηρικής ενδοπρόσθεσης Resonance®**
Οδηγίες χρήσης
- HU
35 **Resonance® fém ureterális sztentkészlet**
Használati utasítás
- IT
38 **Set con stent ureterale metallico Resonance®**
Istruzioni per l'uso
- LV
41 **Resonance® metāla ureterālā stenta komplekts**
Lietošanas instrukcija
- LT
44 **„Resonance®“ metalinio šlapimtakio stento rinkinys**
Naudojimo nurodymai
- NO
47 **Resonance® metallisk ureteralt stentsett**
Bruksanvisning
- PL
50 **Zestaw metalowego stentu moczowodowego Resonance®**
Instrukcja użycia
- PT
54 **Conjunto de stent ureteral metálico Resonance®**
Instruções de utilização
- RO
57 **Set al stentului ureteral metalic Resonance®**
Instrucțiuni de utilizare
- SK
60 **Súprava s kovovým ureterálnym stentom Resonance®**
Návod na použitie
- ES
63 **Equipo de stent ureteral metálico Resonance®**
Instrucciones de uso
- SV
66 **Resonance® metalluretärstentset**
Bruksanvisning

Fig. 1


Fig. 2



Fig. 3

RESONANCE® METALLIC URETERAL STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Resonance Metallic Ureteral Stent Set contains (set components may vary):

- Resonance® metal stent with pigtail straightener
- Radiopaque introduction catheter and sheath

Note: The Resonance Metallic Ureteral Stent is constructed of coiled wire allowing urine to pass through or over any of the coils.

INTENDED USE

Used for temporary stenting of the ureter in adult patients with extrinsic ureteral obstruction. Intended for one-time use.

CONTRAINDICATIONS

There are no known contraindications.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with indwelling ureteral stents include: allergic reaction to nickel • bladder spasm • diminished urine drainage/stent occlusion • fever • fistula formation including arterioureteral fistula • hemorrhage • hydronephrosis • infection • insufficient urine drainage • loss of renal function • pain/discomfort • perforation of kidney, renal pelvis, ureter and/or bladder • peritonitis • pyuria • stent degradation/fracture • stent dislodgement/migration • stent encrustation • stent failure • tissue ingrowth • ureteral reflux • urinary symptoms (frequency, urgency, incontinence, dysuria, hematuria) • urinary tract tissue erosion.

PRECAUTIONS

This product is intended for use by physicians trained and experienced in urology techniques. Standard techniques for urology procedures should be employed.

Do not use this device for any purpose other than the stated intended use.

The Resonance stent must only be used with the positioning system provided and vice versa.

Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.

Improper handling of the stent prior to insertion into the ureter may harm the functionality of the stent. Bending, stretching or any other type of improper handling may deform the stent. It is important that the stent is handled with care.

WARNINGS

- These stents are not intended as permanent indwelling devices.
- The stent must not remain indwelling more than twelve (12) months. If the patient's status permits, the stent may be replaced with a new stent.
- Patients should be checked at regular intervals utilizing techniques such as abdominal X-ray (KUB film). Patients using calcium supplements must be more closely monitored for possible stent encrustation. The stent must be removed if encrustation hampers drainage.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Change in urine viscosity may hamper drainage.
- Hematuria and Incontinence may indicate fistula formation.
- If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- The stent contains nickel which may cause an allergic reaction in individuals with nickel sensitivity.

MR INFORMATION

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated that the Resonance® Metallic Ureteral Stent is MR Conditional. It can be scanned safely under the following conditions:

At 1.5 Tesla Static Magnetic Field	At 3 Tesla Static Magnetic Field
<ul style="list-style-type: none">• Spatial gradient field of 450 Gauss/cm or less; and	<ul style="list-style-type: none">• Spatial gradient field of 720 Gauss/cm or less; and
<ul style="list-style-type: none">• Maximum whole-body-averaged specific absorption rate (SAR) of 1.5 W/kg for 20 min. of scanning	<ul style="list-style-type: none">• Maximum whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 min. of scanning
In non-clinical testing the stent produced a temperature rise during scanning of less than 0.80 °C at an SAR of 1.5 W/kg for 20 min.	In non-clinical testing the stent produced a temperature rise during scanning of less than 1.3 °C at an SAR of 3 W/kg for 15 min.

Note: MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the stent. The image artifact extends approximately 16 mm from the device when scanned in non-clinical testing using: GRE sequence in a 3.0 Tesla Siemens Magnetom Trio, A Tim System (Software Numaris/4) with the body coil. Therefore it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

INSTRUCTIONS FOR USE

Illustrations

Note: Use with wire guide with 0.038" diameter.

Note: Stent may be utilized with either an antegrade or retrograde placement technique. The figures shown (**Figs. 1 and 2**) are for retrograde placement.

1. Using a baseline pyelogram, estimate the proper stent length. Accurate measurement enhances drainage efficiency and patient comfort. To facilitate this, the introduction catheter may be utilized for the injection of contrast media using a syringe.
2. Over a previously placed wire guide, pass introduction catheter and sheath into appropriate position utilizing fluoroscopic guidance or direct vision (**Fig. 1**). Rigid cystoscopes with a sheath size of 21Fr or larger allow for insertion of an introduction catheter and sheath with a maximum outside diameter as specified on the product label.
3. With introduction catheter and sheath in proper position, remove introduction catheter and wire guide to allow passage of the stent (**Fig. 1**).
4. Use the pigtail straightener to introduce the stent to the sheath.
Note: Once stent is in the sheath pull the pigtail straightener away from the hub.
5. Use the introduction catheter to advance the stent through the sheath (**Fig. 2**) until the numbered marker corresponding to the stent length being deployed reaches the hub on the sheath (**Fig. 3**). At this point the first pigtail will have fully deployed from the sheath.
Note: The sheath tip is radiopaque to aid placement.
6. To prevent further deployment of the stent in the kidney hold the introduction catheter in place and retract the sheath.
7. When the hub of the sheath aligns with the proximal ink mark on the introduction catheter the second pigtail is about to deploy. At this point retract both the sheath and introduction catheter completely. (**Fig. 2**)
8. The stent may be removed using conventional cystoscopic techniques utilizing forceps or graspers.

Upon completion of the procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dark, dry location away from temperature extremes. Inventory rotation of sterile products is essential. Verify the expiration date on the package label prior to using the product.

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© COPYRIGHT COOK 2019 ©
COOK IRELAND LTD.