COOK
MEDICAL

CE
0123

EN
3
BG
5
CS
7
DA
9
NL
11
ET
13
FR
15
DE
17
EL
19
HU
21
IT
23
LV
25
LT
27
NO
29
PL
31
PT
33
RO
35
SK
37
ES
39
SV
41

Zenker's Diverticulum Overtube

Instructions for Use

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ZENKER'S DIVERTICULUM OVERTUBE

INTENDED USE

The Zenker's Diverticulum Overtube is intended to act as a diverticuloscope for the purpose of creating a diverticulotomy in patients with Zenker's Diverticulum.

DEVICE DESCRIPTION

The Zenker's Diverticulum Overtube is a flexible silicone diverticuloscope with two distal flaps that straddle the diverticulum septum to protect the anterior esophageal and posterior diverticular walls and provide effective exposure of the area to be cut during endoscopic diverticulotomy. The device is supplied in only one size to facilitate insertion into the esophagus over a flexible endoscope.

CONTRAINDICATIONS

Those specific to upper GI endoscopic diverticulotomy.

Diverticulum less than 2cm in depth. In these cases, the shorter flap of the overtube will not fit into the diverticulum and therefore prevents proper isolation of the septum.

WARNINGS

This device is designed for single use only. Attempts to reprocess, re-sterilize, and/or reuse may lead to device failure and/or transmission of disease.

Visually inspect the packaging and device. If the packaging is opened or damaged when received, do not use. If a device abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

PRECAUTIONS

The product is intended for use by physicians trained and experienced in upper GI endoscopy and endoscopic diverticulotomy techniques.

Standard techniques for upper GI endoscopic diverticulotomy should be employed.

POTENTIAL ADVERSE EVENTS

Potential adverse events include but are not limited to: bleeding, fever, infection, pain/discomfort, perforation, those specific to upper GI endoscopy.

INSTRUCTIONS FOR USE

Illustrations

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device confirming it is free from damage prior to use.
3. Lubricate the outer surface of the endoscope with water-based lubricant approved for medical use.
4. Place the device over the endoscope with the flaps oriented towards the distal end of the endoscope.
5. Advance the endoscope into the opening of the patient's mouth and locate the diverticulum using the endoscopic view.
6. Once the diverticulum is located, advance the device over the endoscope into the patient with the shorter flap oriented towards the diverticulum and the longer flap towards the esophagus.
7. Continue to advance the device into the patient until it straddles the diverticulum, i.e. the longer flap is in the esophagus and the shorter flap is in the pharyngeal pouch.

NOTE: It may be necessary to manipulate the distal end of the endoscope towards the pharyngeal pouch in order to direct the shorter flap of the device into the correct position.

8. If the device is positioned incorrectly or resistance is encountered, retract and/or realign as necessary under direct endoscopic observation.
9. Once the device has straddled the diverticulum septum, an appropriate cutting tool can be inserted under endoscopic observation.
10. Retract the device from the patient leaving the endoscope in position.
11. Retract the endoscope from the patient.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Store in a dry location, away from temperature extremes.

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© COPYRIGHT COOK 2023
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.