EN
3**Sonnet® Polypectomy Snare**

Instructions for Use

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Sonnet® Polypectomy Snare

INTENDED USE

The Sonnet® Polypectomy Snare is a monopolar electrosurgical device used with an electrosurgical unit for endoscopic polypectomy. This device is supplied sterile and intended for single use only.

NOTES

Do not use this device for any purpose other than the stated intended use.

Inventory rotation of sterile products is essential. Verify the expiration date on the package label prior to using the product. If the expiration date has lapsed, do not use.

If the product package is open or damaged when received, do not use this device.

Cook devices must be stored in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the polypectomy site.

Contraindications to polypectomy include, but are not limited to: coagulopathy and insufficiently prepped bowel.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to medication, respiratory depression or arrest, cardiac arrhythmia or arrest. Complications which can occur with endoscopic polypectomy include, but are not limited to: transmural burns, thermal injury to the patient, explosion.

PRECAUTION

The coordination of endoscope accessory channel size with compatible devices is essential in order to obtain optimal results during a procedure. Please refer to the product package label for the minimum channel size required for this device.

CAUTIONS

U.S. Federal Law restricts this device to sale by or on the order of a physician.

A thorough understanding of the technical principles, clinical applications and risks associated with GI endoscopy and polypectomy is necessary before using this device. Polypectomy snares should only be used by, or under the supervision of, physicians thoroughly trained in therapeutic endoscopy and polypectomy.

This device meets the recognized standard for high frequency electrosurgical leakage current (ANSI/AAMI HF 18-2001). Safe and effective electrosurgery is dependent not only on equipment design but, to a large extent, on factors under the control of the operator. Any electrosurgical accessory constitutes a potential electrical hazard to

the patient and the operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

In order to enhance safety and effectiveness, it is important that these cautions, as well as the instructions supplied with the electrosurgical unit, be read, understood and followed when using this device.

Before using this device, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure. Examine all components to be used during the procedure. Do not use a device that has been cut, burned or damaged. Damaged device insulation may cause unsafe currents in either the patient or operator. Do not use a handle that is cracked, broken or has a damaged metal active cord post.

Polypectomy should be performed under direct endoscopic visualization.

Switch the electrosurgical unit to the "off" position when it is not in use.

INSTRUCTIONS FOR USE

1. Upon removing the device from the package, uncoil the device and visually inspect with particular attention to kinks, bends, breaks or fraying. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
2. Inspect the active cord. The cord must be free of kinks, bends, breaks and exposed wires to allow for accurate transfer of current. If an abnormality is noted, do not use the active cord.
3. With the electrosurgical unit off, prepare the equipment. Securely connect the active cord to the device handle and electrosurgical unit. The active cord fittings should fit snugly into both the device handle and the electrosurgical unit. Following the instructions from the electrosurgical unit manufacturer, position the patient return electrode and connect it to the electrosurgical unit.
4. Remove the foam insert from the handle of the device. Fully retract and extend the snare to confirm smooth operation of the device.
5. When the polyp is in endoscopic view, introduce the sheath and retracted snare into the endoscope accessory channel. **Caution: To ensure patient safety, the power to the electrosurgical unit should remain off until the snare is properly positioned around the polyp.**
6. Advance the device, in small increments, until it is endoscopically viewed exiting the endoscope.
7. Advance the snare wire out of the sheath and position it around the polyp to be removed. **Warning: When applying current, the tissue must be isolated from the surrounding mucosa. Failure to isolate the tissue may cause fulguration of normal mucosa and/or perforation. Contact of the snare wire with the endoscope during electrosurgery may cause grounding, which could result in injury to the patient and/or operator as well as damage to the endoscope and/or snare wire.**
8. Following the electrosurgical unit manufacturer's instructions for settings, verify the desired settings and activate the electrosurgical unit. **Note:** The maximum rated input voltage for this device is 2 kV peak-to-peak for cut mode and 5 kV peak-to-peak for coagulation mode.
9. Proceed with polypectomy.

10. Upon completion of the polypectomy, turn the electrosurgical unit off. Retract the snare into the sheath and remove the device from the endoscope.
11. Retrieve the polyp and prepare the specimen per institutional guidelines.
12. Upon completion of the procedure, disconnect the active cord from the device handle and dispose of this device per institutional guidelines for biohazardous medical waste.
13. Disconnect the active cord from the electrosurgical unit. Wipe the active cord with a damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping the active cord tightly may damage the device.

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Sonnet is a trademark of Cook Endoscopy.
Cook is a registered trademark of Cook Incorporated.
© COPYRIGHT COOK 2013 ©
COOK IRELAND LTD.