- EN
4 **ShortShot® Saeed Hemorrhoidal Multi-Band Ligator with TriView® Anoscope**
Instructions for Use
- BG
6 **ShortShot® Saeed хемороидален многолентов лигатор с TriView® аноскоп**
Инструкции за употреба
- CS
9 **Vícegumičkový ligátor hemoroidů ShortShot® Saeed s anoskopem TriView®**
Návod k použití
- DA
11 **ShortShot® Saeed ligator til hæmorider med flere bånd og med TriView® anoskop**
Brugsanvisning
- NL
13 **ShortShot® Saeed hemorrhoidale multiband ligator met TriView® anoscoop**
Gebruiksaanwijzing
- ET
16 **Hemorroidide mitme paelaga ligeerimisvahend ShortShot® Saeed koos anoskoobiga TriView®**
Kasutusjuhised
- FR
18 **Ligateur multiple pour hémorroïdes ShortShot® de Saeed avec anoscope TriView®**
Mode d'emploi
- DE
21 **ShortShot® Saeed Multiband-Ligator für Hämorrhoiden mit TriView®-Anoskop**
Gebrauchsanweisung
- EL
24 **Απολινωτής αιμορροΐδων ShortShot® Saeed πολλαπλών δακτυλίων με πρωκτοσκόπιο TriView®**
Οδηγίες χρήσης
- HU
27 **ShortShot® Saeed többgyűrűs aranyérligátor TriView® anoszkóppal**
Használati utasítás
- IT
30 **Dispositivo di legatura multi-elastico per emorroidi ShortShot® Saeed con anoscopio TriView®**
Istruzioni per l'uso

LV
33

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Linia zębata

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INTENDED USE

The ShortShot® Saeed Hemorrhoidal MultiBand Ligator with TriView® Anoscope is used to ligate internal hemorrhoids facilitated by an anoscope, for use in adult patients only.

NOTES

This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than the intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to joints, cracks or breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization. Store in a dry location, at room temperature.

CAUTION

U.S. Federal Law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

Contraindications include those specific to primary anoscope procedure to be performed in gaining access to desired banding site.

Contraindications specific to hemorrhoidal banding include, but are not limited to: coagulopathy • anal stricture • active inflammatory bowel disease in the distal rectum.

Use of ligation bands is contraindicated in patients with a known hypersensitivity to latex.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with the procedural use of this device include, but are not limited to: allergic reaction to latex • allergic reaction to medication • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • nausea • pain/discomfort • pelvic sepsis • perforation • perineal sepsis • perirectal sepsis • respiratory depression or arrest • stricture formation and obstruction • thrombosis • urinary retention.

PRECAUTIONS

This device may be used with an alternative anoscope if preferred. However, inner lumen diameter of anoscope must be compatible with ShortShot® to obtain optimal results during a procedure. Ensure ShortShot® maneuvers within lumen of alternative anoscope with adequate clearance to avoid dislodging bands.

Banding should begin at the most proximal location from anal sphincter and proceed distally because passing ShortShot® over a previously-placed band may dislodge band.

Prior to using this device a routine anoscopic examination is recommended to confirm the diagnosis requiring treatment of hemorrhoids. This exam may be performed with the accompanying TriView® anoscope in accordance with the instructions for use, or an alternative anoscope may be used at preference of physician. If alternative anoscope is used, please consult relevant packaging, label and Instructions for Use.

Hemorrhoidal banding may result in severe pain if the procedure is performed below the dentate line.

OTHER EQUIPMENT NEEDED

- Medical Grade Suction Equipment
- Source of Illumination
- Water-Soluble lubricant

PREPARATION

1. Anoscope - Prior to procedure, test manoeuvrability of ShortShot® barrel in lumen of anoscope.
2. Ligation Bander - Connect suction tube from suction unit to graduated connector at bottom of ligator handle. Place finger on suction port and verify suction is present at tip of ligator barrel. Ensure settings on suction unit are appropriate for suction of hemorrhoids into tip of ligator barrel. Device is now ready for ligation of hemorrhoids.

INSTRUCTIONS FOR USE

Illustrations

1. Anoscope with obturator comes assembled in package. Remove obturator from anoscope and reinsert before use to ensure proper function. Lubricate anoscope and obturator and insert them into the patient's anus. Rotate anoscope in anus so handle is aligned upwards towards patient's spine (see **Fig. 1**).
2. Remove obturator from anoscope and inspect lumen. Up to three internal hemorrhoids will be aligned in the three slots of anoscope. Slots are designed to coincide with anatomical location of hemorrhoids in most adults. It may be necessary to rotate anoscope a small amount in either direction to encourage any swollen hemorrhoids to fall into place. **Note:** This anoscope may be used in patients with fewer than three hemorrhoids.
3. Insert device into lumen of the anoscope (**Fig. 2**). **Caution: Do not cover suction port prior to achieving desired ligation position.**
4. After full contact is made between barrel tip of ligator and hemorrhoid, suction hemorrhoid into ligator by covering suction port.
5. Maintain suction of hemorrhoid by keeping suction port covered. Deploy band by slowly rotating spool downward until tension is released (**Fig. 3**).
6. Uncover suction port of ligator to relieve suction on hemorrhoid. This will allow ligated hemorrhoid to be released from tip of device. Withdraw device.
7. Repeat ligation process on remaining hemorrhoids, if any. If more than 4 bands are required to complete treatment, use an additional device.
8. Upon completion of ligation procedure, remove ligator and anoscope from patient. It is not necessary to reinsert obturator into anoscope to remove.
9. Dispose of devices as per institutional guidelines for biohazardous medical waste.

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Cook is a registered trademark of Cook Incorporated.
© COPYRIGHT COOK 2019 ©
COOK IRELAND LTD.