- EN 5 **EchoTip® Ultra Celiac Plexus Neurolysis Needle**  
Instructions for Use
- BG 8 **Игла за невролиза в слънчевия сплит EchoTip® Ultra**  
Инструкции за употреба
- CS 11 **Jehla EchoTip® Ultra pro provádění neurolyzy plexus coeliacus**  
Návod k použití
- DA 14 **EchoTip® Ultra nål til neurolyse af plexus coeliacus**  
Brugsanvisning
- NL 17 **Echotip® Ultra naald voor neurolyse van de plexus coeliacus**  
Gebruiksaanwijzing
- ET 20 **Κόηυόόηεπόίμικυ νευρολύσινόελ EchoTip® Ultra**  
Kasutusjuhised
- FR 23 **Aiguille pour neurolyse du plexus coeliaque EchoTip® Ultra**  
Mode d'emploi
- DE 26 **EchoTip® Ultra Nadel zur Neurolyse des Plexus coeliacus**  
Gebrauchsanweisung
- EL 29 **Βελόνα νευρόλυσης κοιλιακού πλέγματος EchoTip® Ultra**  
Οδηγίες χρήσης
- HU 32 **EchoTip® Ultra plexus coeliacus neurolysis tű**  
Használati utasítás
- IT 35 **Ago per la neurolisi del plesso celiaco EchoTip® Ultra**  
Istruzioni per l'uso
- LV 38 **EchoTip® Ultra plexus coeliacus neirolyzes adata**  
Lietošanas instrukcija
- LT 41 **„EchoTip® Ultra“ saulės rezginio neurolyzės adata**  
Naudojimo nurodymai
- NO 44 **EchoTip® Ultra-nål for solar pleksus-nevrolyse**  
Bruksanvisning
- PL 47 **Igła EchoTip® Ultra do neurolyzy splotu trzewnego**  
Instrukcja użycia
- PT 50 **Agulha para neurólise do plexo celiaco EchoTip® Ultra**  
Instruções de utilização
- RO 53 **Ac EchoTip® Ultra de neuroliză a plexului cardiac**  
Instrucțiuni de utilizare
- SK 56 **Ihla na neurolyzu solárneho plexu EchoTip® Ultra**  
Návod na použitie
- ES 59 **Aguja de neurólisis del plexo celiaco EchoTip® Ultra**  
Instrucciones de uso
- SV 62 **EchoTip® Ultra nål för celiac plexus-neurolysis**  
Bruksanvisning


Fig. 1

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

**Fig. 2**

## ECHOTIP® ULTRA CELIAC PLEXUS NEUROLYSIS NEEDLE

**CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### INTENDED USE

This device is used to deliver neurolysing agents to the celiac plexus under guidance by endoscopic ultrasound.

### DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 6 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an ultrasound endoscope. The device has a maximum needle extension length of 8 cm and sheath adjustment length of 5 cm. The device has a cone shaped needle tip design with side holes to administer a uniform radial spray of neurolysing agent. The device is available in one needle size: 20 ga.

### CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired injection site.

Contraindications to celiac plexus neurolysis include, but are not limited to: coagulopathy • use of anticoagulants • aspirin or nonsteroidal anti-inflammatory drugs (NSAIDs).

### POTENTIAL ADVERSE EVENTS

Potential adverse events associated with gastrointestinal endoscopy and celiac plexus neurolysis include, but are not limited to: airway obstruction • allergic reaction to medication • allergic reaction to nickel • aspiration • cardiac arrhythmia or arrest • discomfort • fever • gastroparesis • hemorrhage • hypotension • infection • inflammation • lower extremity weakness • nausea • paraplegia • paresthesias • perforation • renal puncture • respiratory depression or arrest • sepsis • transient diarrhea • transient increasing pain • vomiting.

### WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the endoscope.

Ensure the stylet is fully inserted when advancing the needle into the target site.

A thorough understanding of the technical principles, clinical applications and risks associated with endoscopic ultrasound and celiac plexus neurolysis is necessary before using this device. The use of prophylactic antibiotics is recommended for endoscopic ultrasound-guided celiac plexus neurolysis. The use of general anesthesia may be needed in some cases.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

## INSTRUCTIONS FOR USE

## Illustrations

1. Identify the desired injection site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting on the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 1**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 2**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, to the desired injection site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional neurolysis injections are required later.

9. If required, inadvertent puncture of a blood vessel may be determined by attempting to aspirate blood. A Luer lock syringe may be attached to the Luer fitting on the needle handle for this purpose and then removed.

10. Attach the Luer lock fitting of a previously prepared syringe, containing local anaesthetic and neurolysing agents, securely to the Luer fitting on the needle handle. Instill volume required while observing the ultrasound view.

11. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.

12. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

13. For additional neurolysing injections at alternative sites on the celiac plexus, disconnect the syringe from the Luer fitting on the needle handle and gently reinsert the stylet. Repeat Steps 2 through 12 of the "Instructions for Use". **Note:** Prior to reinserting stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.

**Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.**

## **HOW SUPPLIED**

Store in a dry location, away from temperature extremes.

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© COPYRIGHT COOK 2019 ©
**COOK IRELAND LTD.**