EN
4FR
7ES
10**EchoTip Ultra® Celiac Plexus Neurolysis Needle**

Instructions for Use

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

Fig. 2

ECHOTIP ULTRA® CELIAC PLEXUS NEUROLYSIS NEEDLE

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 6 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an ultrasound endoscope. The device has a maximum needle extension length of 8 cm and sheath adjustment length of 5 cm. The device has a cone shaped needle tip design with side holes to administer a uniform radial spray of neurolysing agent. The device is available in one needle size: 20 ga.

Device compatibility

The EchoTip Ultra Celiac Plexus Neurolysis Needle is compatible with:

- Minimum endoscope accessory channel of 2.8mm.
- Syringe with a male Luer connection fitting.

INTENDED USE

This device is used to deliver neurolysing agents to the celiac plexus under guidance by endoscopic ultrasound.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired injection site.

Contraindications to celiac plexus neurolysis: coagulopathy • use of anticoagulants, aspirin, or nonsteroidal anti-inflammatory drugs (NSAIDs).

WARNINGS

- Not for use in the heart or vascular system.
- The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.
- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

- Refer to package label for minimum channel size required for this device.
- The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement, or withdrawal of the device. Failure to retract the needle may result in damage to the endoscope.
- Ensure the stylet is fully inserted when advancing the needle into the target site.
- A thorough understanding of the technical principles, clinical applications, and risks associated with endoscopic ultrasound and celiac plexus neurolysis is necessary before using this device.
- The use of prophylactic antibiotics is recommended for endoscopic ultrasound-guided celiac plexus neurolysis.
- The use of general anesthesia may be needed in some cases.
- Do not use this device for any purpose other than the stated intended use.
- Use of this device is restricted to a trained healthcare professional (Rx only).
- The product is intended for use by physicians trained and experienced in endoscopic ultrasound and celiac plexus neurolysis techniques.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with gastrointestinal endoscopy and celiac plexus neurolysis: airway obstruction • allergic reaction to medication • allergic reaction to nickel • aspiration • cardiac arrhythmia or arrest • discomfort • fever • gastroparesis • hemorrhage • hypotension • infection • inflammation • ischemia • lower extremity weakness • nausea • organ infarction • paraplegia • paresthesias • perforation • renal puncture • respiratory depression or arrest • sepsis • transient diarrhea • transient increasing pain • vomiting.

HOW SUPPLIED

This device is supplied ethylene oxide (EO) sterilized in a tray with a Tyvek lid.

STORAGE

Store in a dry location, away from temperature extremes.

INSTRUCTIONS FOR USE

Illustrations

1. Identify the desired injection site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting on the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 1**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 2**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, to the desired injection site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope, and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage, or malfunction.
8. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional neurolysing injections are required later.
9. If required, inadvertent puncture of a blood vessel may be determined by attempting to aspirate blood. A Luer lock syringe may be attached to the Luer fitting on the needle handle for this purpose and then removed.
10. Attach the Luer lock fitting of a previously prepared syringe, containing local anaesthetic and neurolysing agents, securely to the Luer fitting on the needle handle. Instill volume required while observing the ultrasound view.
11. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
12. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

13. For additional neurolysing injections at alternative sites on the celiac plexus, disconnect the syringe from the Luer fitting on the needle handle and gently reinsert the stylet. Repeat Steps 2 through 12 of the “Instructions for Use”. **Note:** Prior to reinserting stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.

DISPOSAL OF DEVICE

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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© COOK 2024
COOK IRELAND LTD.