- EN
2
Zilver® 518™ & 635™ Biliary Stent
Instructions for Use
- BG
10
Билиарни стентове Zilver® 518™ и 635™
Инструкции за употреба
- CS
19
Biliární stenty Zilver® 518™ a 635™
Návod k použití
- DA
28
Zilver® 518™ og 635™ galdestent
Brugsanvisning
- DE
37
Zilver® 518™ und 635™ Gallengangstents
Gebrauchsanweisung
- EL
46
Ενδοπρόσθεση χοληφόρων 518™ και 635™ Zilver®
Οδηγίες χρήσης
- ES
55
Stents biliares Zilver® 518™ y 635™
Instrucciones de uso
- ET
64
Biliaarsed stendid Zilver® 518™ ja 635™
Kasutusjuhend
- FR
73
Endoprothèse biliaire Zilver® 518™ et 635™
Mode d'emploi
- HU
82
Zilver® 518™ és 635™ epeúti sztent
Használati utasítás
- IT
91
Stent biliari Zilver® 518™ e 635™
Istruzioni per l'uso
- LT
100
„Zilver®“ 518™ ir 635™ tulžies lataku stentas
Naudojimo nurodymai
- LV
109
Zilver® 518™ un 635™ žultsvadu stents
Lietošanas instrukcija
- NL
118
Zilver® 518™ en 635™ biliaire stents
Gebruiksaanwijzing
- NO
127
Zilver® 518™ og 635™ gallestenter
Bruksanvisning
- PL
135
Stent Zilver® 518™ i 635™ do dróg żółciowych
Instrukcja użycia
- PT
144
Stents biliares Zilver® 518™ e 635™
Instruções de utilização
- RO
153
Stent biliar Zilver® 518™ și 635™
Instrucțiuni de utilizare
- SV
162
Zilver® 518™ och 635™ biliär stent
Bruksanvisning
- SK
170
Biliární stent Zilver® 518™ a 635™
Návod na použitie


ZILVER® 518™ & 635™ BILIARY STENT

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Do not re-sterilize.

DEVICE DESCRIPTION

- a. Handle
- b. Hub
- c. Safety Lock
- d. Introducer Catheter
- e. Introducer Tip
- f. Side-arm Flushing Port
- g. Metal Cannula
- h⁽¹⁾. Radiopaque Marker on Delivery System (on 518 delivery system only)
- h⁽²⁾. Radiopaque Marker on Delivery System
- i. Gold Radiopaque Markers on the Stent

The Zilver 518 & 635 Biliary Stents are self-expandable stents made of nitinol. They are flexible, slotted tubes that provide strength and flexibility in the biliary duct upon deployment. Post deployment, the stents impart an outwardly strong radial force upon the inner lumen of the duct, establishing patency in the designated stent region.

The Zilver 518 Biliary Stent is preloaded in a 5.0 French (1.67mm) delivery catheter and the Zilver 635 Biliary Stent is preloaded in a 6.0 French (2.0mm) delivery catheter. Stent deployment is controlled by retracting the handle while holding the metal cannula stationary.

Note: Not all products listed in this IFU are available in all regions.

INTENDED USE

The Zilver 518 and 635 Biliary Stent has been designed for use in palliation of malignant neoplasms in the biliary tree.

MR CONDITIONS

This symbol means the device is MR conditional.

Non-clinical testing has demonstrated that the Zilver Biliary Stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Whole body averaged specific absorption rate (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes of scanning (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla), respectively.

In non-clinical testing, the Zilver Biliary Stent produced maximum temperature rises of 0.1, 3.8, 0.8, and 0.1 degrees C (for a single stent at 1.5 Tesla, a pair of overlapping stents at 1.5 Tesla, a single stent at 3 Tesla, and a pair of overlapping stents at 3 Tesla, respectively) at whole body averaged specific absorption rates (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) of MR scanning in a 1.5 Tesla/64 MHz General Electric MR scanner and a 3 Tesla Excite General Electric MR scanner.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Zilver Biliary Stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant. Heating in the MR environment for stents with fractured struts is not known.

CONTRAINDICATIONS

- Total biliary obstructions and postoperative strictures that cannot be dilated to permit passage of the introducer catheter.
- Stenting of a perforated duct where leakage from the duct could be exacerbated by placement of a stent.
- Patients with bleeding disorders.
- Severe ascites.

WARNINGS

- This stent is **not intended to be removed**. Attempts to remove the stent after placement may cause damage to the surrounding mucosa. This stent is considered to be a permanent implant.
- This device is **not** intended to be deployed through the wall of a previously placed or existing metal stent. Doing so could result in difficulty or inability to remove the introducer.

- Possible allergic reactions to nitinol should be considered.
- The safety and effectiveness of this device for use in the vascular system have not been established.
- After stent placement, alternate methods of treatment such as chemotherapy and irradiation may increase the risk of: a) stent migration due to tumor shrinkage, b) stent erosion of the tissue, and/or c) mucosal bleeding.

PRECAUTIONS

- Manipulation of the product requires high-resolution fluoroscopic control.
- The stent has not been designed to inhibit tumor ingrowth. Tissue may grow through stents.
- Ensure that the safety lock is not inadvertently removed until final stent release.
- Deploy the stent over an extra-stiff or ultra-stiff wire guide.
- Repositioning of the device after deployment is not possible since the introducer catheter cannot be re-advanced over the stent once deployment begins.
- This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.
- Do not attempt to remove the stent from the introducer system before use.
- Do not attempt to push the handle away from the hub during deployment.
- Do not rotate any part of the system during deployment.
- Do not expose the delivery system to organic solvents (e.g., alcohol).
- Do not use power injection systems with the delivery system.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur: Allergic reaction to contrast or medication • Allergic reaction to Nitinol • Bile duct ulceration • Bile leakage • Biloma • Cardiac arrhythmia or arrest • Cholangitis • Cholecystitis • Cholestasis • Death (other than due to normal disease progression) • Fever • Hemothorax • Hematoma • Hemorrhage • Infection/abscess at access site • Inflammation • Liver abscess • Nausea/Vomiting • Obstruction of the pancreatic duct • Pain / Discomfort • Pancreatitis • Perforation • Peritonitis • Pleural effusions • Pleuritis • Pneumonia • Pneumothorax • Recurrent obstructive jaundice • Respiratory depression or arrest • Sepsis • Stent fracture • Stent migration • Stent misplacement • Stent occlusion • Trauma to the biliary duct or duodenum • Tumor overgrowth • Tumor seeding • Vascular injury (including portal or hepatic vein or artery).

PRODUCT RECOMMENDATIONS

Placement of this biliary stent requires advanced skills in interventional biliary procedures. The following instructions give technical guidance, but do not obviate formal training in the use of the device.

Wire Guide Selection

The use of an extra-support wire guide is recommended.

- 0.018 inch (0.46 mm) for the Zilver 518
- 0.035 inch (0.89 mm) for the Zilver 635

Selection of Stent

Measure the length of the target stricture to determine the length of the stent required.

Allow for the area proximal and distal to the tumor to be covered with the stent to protect against interference from further tumor ingrowth.

Measure the diameter of the reference lumen (proximal and distal to the stricture) and refer to Stent Size Selection Table below.

Note: Not all product sizes are available in all markets.

Stent Size Selection Table

6 French (2.0 mm)						
40 cm Delivery System						
Stent length (mm)	20	30	40	60	70	80
Stent Outer Diameter (mm)	6		x	x		x
	8		x	x		x
	9		x	x		x
	10		x	x		x
	12		x	x		x
	14		x	x		x


INSTRUCTIONS FOR USE

This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous transhepatic biliary drainage catheters and metal stents should be employed.

1. Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control.
2. Introduce the extra support wire guide (.018 inch [0.46 mm] for the Zilver 518 or .035 inch [0.89 mm] for the Zilver 635) through the access catheter across the distal segment of the target lesion.
3. Remove the access catheter, leaving the wire guide in place.
4. Use the 1ml syringe included in the inner package to flush the introducer catheter with saline through the side-arm flushing port.

- Use the 1ml syringe to flush the wire guide lumen with saline through the hub.
- Using fluoroscopic visualization, insert the introducer catheter over the wire guide. Position the catheter so that the distal and proximal radiopaque markers on the stent (i) are at the desired position. The stent is now ready to be deployed. (**Fig. 1**)


Fig. 1

Deployment of the Stent

- Before deployment it is important to straighten the proximal part of the introducer catheter as much as possible and to keep the handle in a stable position.
- The stent expansion must be performed under fluoroscopic control.

Fig. 2

3. Hold the hub (**b**) on the metal cannula (**g**) steady. Prior to deploying the stent, remove the red safety lock (**c**). (**Fig. 2**). Hold the hub end stationary. The stent will be deployed as you pull the handle (**a**) toward the hub (**b**). (**Figs. 3 and 4**)

The introducer catheter cannot be re-advanced over the stent during deployment.

Fig. 3




Fig. 4

4. The stent is fully deployed when the handle (a) reaches the hub (b) (Fig. 5).

5. Perform a final cholangiogram to ensure proper placement of the device.

NOTE: Stents bridging the papilla should extend approximately 5 mm beyond the papilla and into the duodenum after deployment.

6. Remove the delivery system.

NOTE: If resistance is met during the withdrawal of the delivery system, re-advance the outer sheath to its pre-deployment position. Withdraw the system as one unit.

7. After stent placement, a temporary external drainage catheter may be placed over the pre-positioned wire guide. This provides access to the stent if an additional treatment is needed.

Multiple Stent Placement

If placement of more than one stent is required in a patient, to cover the length of the lesion, the following recommendations should be considered:

- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem should slightly overlap to allow for complete coverage of the lesion.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2022
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.