
EN 6	<b>Zilver® Stent Vascular Design</b> Instructions for Use
BG 11	<b>Васкуларен стент Zilver®</b> Инструкции за употреба
CS 17	<b>Cévní stent Zilver®</b> Návod k použití
DA 22	<b>Zilver® stent vaskulært design</b> Brugsanvisning
NL 27	<b>Zilver® vasculaire stent</b> Gebruiksaanwijzing
ET 33	<b>Vaskulaarne stent Zilver®</b> Kasutusjuhised
FR 38	<b>Endoprothèse Zilver® à usage vasculaire</b> Mode d'emploi
DE 44	<b>Zilver® Gefäßstent</b> Gebrauchsanweisung
EL 50	<b>Αγγειακό σχέδιο ενδοπρόσθεσης Zilver®</b> Οδηγίες χρήσης
HU 56	<b>Vaszkuláris kialakítású Zilver® sztent</b> Használati utasítás
IT 62	<b>Stent Zilver® per applicazioni vascolari</b> Istruzioni per l'uso
LV 68	<b>Zilver® vaskulāras konstrukcijas stents</b> Lietošanas instrukcija
LT 73	<b>„Zilver®“ kraujagyslėms skirtos konstrukcijos stentas</b> Naudojimo nurodymai
NO 78	<b>Zilver® stent vaskulær design</b> Bruksanvisning
PL 83	<b>Stent Zilver® skonstruowany do użytku naczyniowego</b> Instrukcja użycia
PT 89	<b>Stent de design vascular Zilver®</b> Instruções de utilização
RO 95	<b>Stent Zilver® cu aplicație vasculară</b> Instrucțiuni de utilizare
SK 100	<b>Vaskulárny stent Zilver®</b> Návod na použitie
ES 105	<b>Stent Zilver® de diseño vascular</b> Instrucciones de uso
SV 111	<b>Zilver® stent med vaskulär design</b> Bruksanvisning

Fig. 1

- a. Handle
- b. Hub
- c. Safety lock
- d. Introducer catheter
- e. Introducer tip
- f. Flushing port
- g. Metal cannula
- h. Radiopaque markers
- i. Radiopaque marker

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

Fig. 5

**ZILVER® STENT VASCULAR DESIGN**

**CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

**Do not re-sterilize.**

**DEVICE DESCRIPTION**

The Zilver® Stent Vascular Design (ZIV) is a self-expandable stent made of nitinol. It is a flexible, slotted tube that is designed to provide support while maintaining flexibility in the vessel upon deployment. Post-deployment, the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region.

The stent has radiopaque markers at both ends (h).

The stent is supplied preloaded in a 5.0 French (1.67 mm) or 6.0 French (2.0 mm) introducer catheter. A radiopaque marker (i) on the distal tip of the sheath is applied to visualize deployment of the stent. Stent deployment is controlled by means of a hand held device.

Refer to Stent Size Selection Tables below.

5 French (1.67 mm)													
Stent length (mm)		20		30		40		60		70		80	
Delivery System (cm)		80	125	80	125	80	125	80	125	80	125	80	125
Stent Outer Diameter (mm)	4	x	x	x	x	x	x	x	x			x	x
	5	x	x	x	x	x	x	x	x			x	x
	6	x	x	x	x	x	x	x	x			x	x
	7	x	x	x	x	x	x	x	x			x	x
	8	x	x	x	x	x	x	x	x			x	x
	9	x	x	x	x	x	x	x	x			x	x
	10	x	x	x	x	x	x	x	x		x	x	x

6 French (2.0 mm)													
Stent length (mm)		20		30		40		60		70		80	
Delivery System (cm)		80	125	80	125	80	125	80	125	80	125	80	125
Stent Outer Diameter (mm)	4												
	5	x	x	x	x	x	x	x	x	x		x	x
	6	x	x	x	x	x	x	x	x	x		x	x
	7	x	x	x	x	x	x	x	x	x		x	x
	8	x	x	x	x	x	x	x	x	x		x	x
	9	x	x	x	x	x	x	x	x	x		x	x
	10	x	x	x	x	x	x	x	x	x		x	x
	12			x	x	x	x	x	x	x		x	x
	14			x	x	x	x	x	x	x		x	x

### INTENDED USE

The product is intended for use in the iliac arteries for the following treatments:

- Arteriosclerotic stenosis
- Total occlusions that have been recanalized

The product provides mechanical support to maintain constant blood flow of the vessel.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

Placement of the stent requires advanced skills in interventional procedures. The following instructions will give technical guidance but do not obviate formal training in the use of the device.

Standard techniques for placement of an arterial access sheath and wire guide should be employed during use of the vascular stent.

### CONTRAINDICATIONS

Patients with bleeding diathesis or other disorder, e.g. peptic ulceration or recent cerebrovascular accident, limiting the use of antiplatelet and anticoagulant therapy.

Patients with total arterial occlusion.

Coronary and below the knee use is not recommended.

### WARNINGS

- Persons allergic to nitinol (nickel titanium) may suffer an allergic reaction to this implant.

### PRECAUTIONS

- Manipulation of the product requires high resolution fluoroscopic control.
- The lesion and vessel can be predilated before placement of the stent.
- Do not try to remove the stent from the introducer system before use.
- Ensure that the safety lock is not inadvertently removed prior to stent release.

## MR CONDITIONS

Non-clinical testing has demonstrated that the Zilver Vascular Stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Whole body averaged specific absorption rate (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes of scanning (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla), respectively.

In non-clinical testing, the Zilver Vascular Stent produced maximum temperature rises of 0.1, 3.8, 0.8, and 0.1 degrees C (for a single stent at 1.5 Tesla, a pair of overlapping stents at 1.5 Tesla, a single stent at 3 Tesla, and a pair of overlapping stents at 3 Tesla, respectively) at whole body averaged specific absorption rates (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) of MR scanning in a 1.5 Tesla/64 MHz General Electric MR scanner, a 1.5 Tesla Magnetom Siemens Medical Solutions MR Scanner (to evaluate a pair of overlapping stents), and a 3 Tesla Excite General Electric MR scanner.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Zilver Vascular Stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant. Heating in the MRI environment for stents with fractured struts is not known.

## POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to:

- Abrupt stent closure
- Allergic reaction to nitinol
- Amputation
- Angina/coronary ischemia
- Arrhythmia
- Arterial aneurysm
- Arterial rupture
- Arteriovenous fistula
- Atheroembolization (Blue Toe Syndrome)
- Death
- Embolism
- Fever
- Hematoma/hemorrhage
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection/abscess formation at access site
- Intimal injury/dissection
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Myocardial infarction
- Occlusion
- Pseudoaneurysm formation
- Renal failure
- Restenosis of the stented artery
- Septicemia/bacteraemia

- Stent malapposition
- Stent migration
- Stent strut fracture
- Stroke
- Spasm
- Thrombosis
- Tissue necrosis
- Vessel trauma / vessel rupture
- Worsened claudication/rest pain

## **CHEMICAL TOLERANCES**

- The stent must not come in contact with metal other than nitinol. Please pay attention to this if the patient has an arterial stent already.

## **PRODUCT RECOMMENDATIONS**

### **Arterial access**

For arterial access, we recommend the use of an access set that accepts the introducer catheter (5.0 French [1.67 mm] or 6.0 French [2.0 mm]).

### **Wire Guide Selection**

For vascular introduction a .018-inch (0.46 mm) wire guide is recommended for the 5.0 French (1.67 mm), and a .035-inch (0.89 mm) wire guide for the 6.0 French (2.0 mm) introducer catheter.

### **PTA Balloon Selection**

For pre- and postdilation, an appropriately sized balloon catheter is recommended.

### **Stent Selection**

The chosen stent diameter should be at least 1 mm larger than the diameter of the reference vessel. Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control. Measure the length of the target lesion to determine the length of the stent required. Allow for the proximal and distal aspects of the stent to cover the entire target area.

## **PROCEDURE**

### **Access**

1. Upon removal from the package, inspect the product to ensure no damage has occurred.
2. Introduce the recommended access catheter using standard percutaneous technique.
3. Introduce the wire guide through the access catheter across the distal segment of the target lesion.
4. Remove the access catheter, leaving the wire guide in place.
5. Predilate if necessary. Remove the balloon catheter leaving the wire guide in place.
6. Determine proper stent size.
7. Prior to insertion, flush the delivery system with saline through the flushing port.
8. Insert the delivery system over the wire guide.
9. Place the delivery system under fluoroscopy to determine that the radiopaque markers on the delivery system (h) are at the desired position. The stent is now ready to be deployed (Fig. 2).

NOTE: If multiple stents are required to cover the length of the lesion, please refer to the Multiple Stent Placement section of these instructions for use for further recommendations.

### **Deployment of the stent**

1. Before deployment, it is important to straighten the proximal part of the delivery system as much as possible and to keep the handle in a stable position.
2. The stent expansion must be performed under fluoroscopic control.
3. Remove the safety lock (Fig. 3a). Gently slide the handle towards the hub (Fig. 3b) and the stent will deploy. The radiopaque markers on the delivery system indicate the progress of the deployment.
4. Relocation of the stent during deployment is only possible to a more proximal position; carefully pull the system towards you (Fig. 4). When proper location has been achieved, continue deployment.  
When the handle reaches the hub, the stent has been fully deployed (Fig. 5).  
Post stenting dilation can now be performed if required. Remove the wire guide.
5. Perform completion angiogram to ensure proper placement of the device, remove the introducer catheter.
6. Delivery System Removal - Do not advance sheath after the stent has been deployed. Delivery system can be removed without the need to recapture tip. Check for delivery system integrity post removal from the patient.

### **Multiple Stent Placement**

If placement of more than one stent is required, the following recommendations should be considered:

- The stent must not be in contact with metals other than nitinol. Please take this into consideration if the patient already has a metal stent.
- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations, i.e. a second stent should be placed proximal to the previously placed stent.
- Stents placed in tandem should slightly overlap.

Upon completion of the procedure dispose of device per institutional guidelines for biohazardous medical waste.

### **HOW SUPPLIED**

Supplied sterile by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

### **REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COPYRIGHT COOK 2018 ©
**COOK IRELAND LTD.**