EN
5FR
9**Zilver® 518™ and 635™ Biliary Stents**

Instructions for Use

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
Fig. 5

ZILVER® 518™ AND 635™ BILIARY STENTS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Do not re-sterilize.

DEVICE DESCRIPTION

The Zilver 518 and 635 Biliary Stents are self-expandable stents made of nitinol. They are flexible, slotted tubes that provide strength and flexibility in the biliary duct upon deployment. Post deployment, the stents impart an outwardly strong radial force upon the inner lumen of the duct, establishing patency in the designated stent region.

The Zilver 518 Biliary Stent is preloaded in a 5.0 French (1.67 mm) delivery catheter and the Zilver 635 Biliary Stent is preloaded in a 6.0 French (2.0 mm) delivery catheter. Stent deployment is controlled by retracting the handle while holding the metal cannula stationary.

INTENDED USE

The Zilver 518 and 635 Biliary Stents are intended for palliation of malignant neoplasms in the biliary tree.

MR SAFETY INFORMATION

This symbol means the device is MR conditional

Non-clinical testing has demonstrated the Zilver® 518 and Zilver® 635 Biliary Self-Expanding Stents are MR conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3 T or 1.5 T
- Maximum spatial field gradient of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, non-clinical testing results (including the effect of perfusion) indicate the Zilver® 518 and Zilver® 635 Biliary Self-Expanding Stents are expected to produce a maximum temperature rise of less than 5.3°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 3 mm from the Zilver® 518 and Zilver® 635 Biliary Self-Expanding Stents when imaged with a gradient echo pulse sequence and a 3 T MR system.

CONTRAINDICATIONS

- Total biliary obstructions and postoperative strictures that cannot be dilated to permit passage of the introducer catheter
- Stenting of a perforated duct where leakage from the duct could be exacerbated by placement of a stent
- Patients with bleeding disorders
- Severe ascites

WARNINGS

- This stent is **not intended to be removed**. Attempts to remove the stent after placement may cause damage to the surrounding mucosa. This stent is considered to be a permanent implant.
- This device is **not** intended to be deployed through the wall of a previously placed or existing metal stent. Doing so could result in difficulty or inability to remove introducer.
- Possible allergic reactions to nitinol should be considered.
- The safety and effectiveness of this device for use in the vascular system have not been established.

PRECAUTIONS

- Manipulation of the product requires high-resolution fluoroscopic control.
- The stent has not been designed to inhibit tumor ingrowth. Tissue may grow through stents.
- Ensure that the red safety lock is not inadvertently removed until final stent release.
- Deploy the stent over an extra stiff or ultra stiff wire guide.
- Repositioning of the device after deployment is not possible since the introducer catheter cannot be re-advanced over the stent once deployment begins.
- This device is designed for single use only. Attempts to reprocess, re-sterilize, and/or reuse may lead to device failure and/or transmission of disease.
- Do not attempt to remove the stent from the introducer system before use.
- Do not attempt to push the handle away from the hub during deployment.
- Do not rotate any part of the system during deployment.
- Do not expose the delivery system to organic solvents (e.g., alcohol).
- Do not use power injection system with the delivery system.
- Do not advance the sheath after the stent has been deployed.
- Following stent deployment, if resistance is met during the withdrawal of the delivery system, gently re-advance the outer sheath to its pre-deployment position (while continuing to fluoroscopically monitor the stent). Carefully remove the delivery system and wire guide together as a unit. If resistance is still encountered during removal of the delivery system and wire guide as a unit, remove the wire guide, delivery system and access sheath together as a unit.

PRODUCT RECOMMENDATIONS

Placement of this biliary stent requires advanced skills in interventional biliary procedures. The following instructions give technical guidance, but do not obviate formal training in the use of the device.

Wire Guide Selection

The use of an extra support wire guide is recommended.

- .018 inch (0.46 mm) for the Zilver 518
- .035 inch (0.89 mm) for the Zilver 635

Selection of Stent

Measure the length of the target stricture to determine the length of the stent required.

Allow for the area proximal and distal to the tumor to be covered with the stent to protect against interference from further tumor ingrowth.

Measure the diameter of the reference lumen (proximal and distal to the stricture) and refer to Stent Size Selection Table below.

Stent Size Selection Table

Stent Outer Diameter (mm)	5 French (1.67 mm)					6 French (2.0 mm)				
	Delivery System 40, 80 & 125cm					Delivery System 40, 80 & 125cm				
	Stent Length (mm)									
	20	30	40	60	80	20	30	40	60	80
4	×	×	×	×	×	×	×	×	×	×
5	×	×	×	×	×	×	×	×	×	×
6	×	×	×	×	×	×	×	×	×	×
7	×	×	×	×	×	×	×	×	×	×
8	×	×	×	×	×	×	×	×	×	×
9	×	×	×	×	×	×	×	×	×	×
10	×	×	×	×	×	×	×	×	×	×
12						×	×	×	×	×
14						×	×	×	×	×

INSTRUCTIONS FOR USE

Illustrations

This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of percutaneous transhepatic biliary drainage catheters and metal stents should be employed.

1. Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control.
2. Introduce the extra support wire guide (.018 inch [0.46 mm] for the Zilver 518 or .035 inch [0.89 mm] for the Zilver 635) through the access catheter across the distal segment of the target lesion.
3. Remove the access catheter, leaving the wire guide in place.
4. Use the 1 ml syringe included in the inner package to flush the introducer catheter with saline through the side-arm flushing port.
5. Use the 1 ml syringe to flush the wire guide lumen with saline through the hub.
6. Using fluoroscopic visualization, insert the introducer catheter over the wire guide.

Position the catheter so that the distal and proximal radiopaque markers on the stent (i) are at the desired position. The stent is now ready to be deployed. (Fig. 1)

Deployment of the Stent

1. Before deployment it is important to straighten the proximal part of the introducer catheter as much as possible and to keep the handle in a stable position.
2. The stent expansion must be performed under fluoroscopic control.
3. Hold the hub (b) on the metal cannula (g) steady. Prior to deploying the stent, remove the red safety lock (c). (Fig. 2)
4. Hold the hub end stationary. The stent will be deployed as you pull the handle (a) toward the hub (b). (Figs. 3 and 4)

CAUTION: The introducer catheter cannot be re-advanced over the stent during deployment.

5. The stent is fully deployed when the handle (a) reaches the hub (b). (Fig. 5)
6. Perform a final cholangiogram to ensure proper placement of the device.

NOTE: Stents bridging the papilla should extend no more than 5 mm beyond the papilla and into the duodenum after deployment.

7. Remove the delivery system. **NOTE:** Delivery system can be removed without recapturing tip. If recapture of the tip is required during removal of the delivery system, this must be performed by holding the handle (a) steady and pulling the hub on the metal cannula (b) backwards. Ensure the handle is not pushed forward during recapture.
8. After stent placement, a temporary external drainage catheter may be placed over the pre-positioned wire guide. This provides access to the stent if an additional treatment is needed.

Multiple Stent Placement

If placements of multiple stents are required in a patient, to cover the length of the lesion, the following recommendations should be considered:

- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem should slightly overlap to allow for complete coverage of the lesion.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© COOK 2024
COOK IRELAND LTD.