COOK®
MEDICAL

CE
0123

EN 14	Biliary Stents and Sets Instructions for Use
BG 19	Билиарни стентове и комплекти Инструкции за употреба
CS 24	Biliární stenty a soupravy Návod k použití
DA 29	Galdestenter og sæt Brugsanvisning
DE 34	Gallengangstents und -sets Gebrauchsanweisung
EL 40	Ενδοπροσθέσεις χοληφόρων και σέτ Οδηγίες χρήσης
ES 46	Stents y equipos biliares Instrucciones de uso
ET 51	Biliaarsed stendid ja komplektid Kasutusjuhend
FR 56	Sets et endoprothèses biliaires Mode d'emploi
HR 61	Bilijarni stentovi i kompleti Upute za uporabu
HU 66	Epeüti sztentek és készletek Használati utasítás
IT 71	Stent biliari e set di stent biliari Istruzioni per l'uso
LT 77	Tulžies latakų stentai ir rinkiniai Naudojimo instrukcija
LV 82	Žultsvadu stenti un komplekti Lietošanas instrukcija
NL 86	Biliaire stents en sets Gebruiksaanwijzing
NO 92	Gallestenter og stentsett Bruksanvisning
PL 97	Stenty do dróg żółciowych i zestawy Instrukcja użycia
PT 102	Stents biliares e Conjuntos Instruções de utilização
RO 107	Stenturi biliare și seturi de stenturi biliare Instrucțiuni de utilizare
SK 113	Biliárny stent a súpravy Návod na použitie
SL 118	Stenti za žolčevod in sklopi stentov Navodila za uporabo

SR
123

Bilijarni stentovi i kompleti

Uputstvo za upotrebu

SV
128

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
Bruksanvisning

TR
132

Biliyer Stentler ve Setler

Kullanma Talimatı

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SR Slika 6: Zaštita preklopa

SV Figur 6: Flikskydd

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EN The indicated stent length is the distance between the flaps/pigtails.

- BG Посочваната дължина на стента е разстоянието между крилцата/извитите краища с форма на свинска опашка (pigtail).
- CS Uvedená délka stentu představuje vzdálenost mezi křídélky nebo pigtaily.
- DA Den angivne stentlængde er afstanden mellem flapperne/grisehalerne.
- DE Die angegebene Stentlänge bezieht sich auf den Abstand zwischen den Flaps/Pigtails.
- EL Το ενδεικνυόμενο μήκος ενδοπρόσθεσης είναι η απόσταση μεταξύ των γλωττίδων/σπειροειδών άκρων.
- ES La longitud indicada del stent es la distancia entre las lengüetas o los pigtails.
- ET Stendi näidatud pikkus on keelte / rõngasse keeratud otste vahekaugus.
- FR La longueur indiquée de l'endoprothèse est la distance entre les rabats/pigtails.
- HR Indicirana duljina stenta predstavlja udaljenost između zalistaka/zavojnica.
- HU A megadott sztenthossz a rögzítőszárnyak/pigtail elemek közötti távolsággal egyenlő.
- IT La lunghezza dello stent indicata corrisponde alla distanza tra alette/pigtail.
- LT Nurodytas stento ilgis yra atstumas tarp atraminių liežuvelių / riestų galiukų.
- LV Indicētais stenta garums ir attālums starp atlokiem/saritinātajiem galiem.
- NL De aangegeven stentlengte is de afstand tussen de flappen/pigtails.
- NO Den angitte stentlengden er avstanden mellom flikene/pigtailene.
- PL Wskazana długość stentu to odległość pomiędzy wąsami/końcówkami pigtail.
- PT O comprimento de stent indicado é a distância entre as abas/espirais.
- RO Lungimea indicată a stentului reprezintă distanța dintre aripioare/regiunile pigtail.
- SK Uvádzaná dĺžka stentu je vzdialenosť medzi výbežkami/špirálami
- SL Navedena dolžina stenta je razdalja med zavihkoma/zavitima koncema.
- SR Naznačena dužina stenta je rastojanje između preklopa / „pigtail“ delova.
- SV Den indikerade stentlängden är avståndet mellan flikar/pigtails.
- TR Belirtilen stent uzunluğu, flepler/pigtail'ler arasındaki mesafedir.

ZIMMON®, COTTON-LEUNG®, COTTON-LEUNG® SOF-FLEX®, COTTON-HUIBREGTSE®, ST-2 SOEHENDRA TANNENBAUM®, AND PLASTIC BILIARY STENT (PBS) BILIARY STENTS AND STENT SETS

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

This IFU covers Zimmon® (ZEBD, ZSO), Cotton-Leung® (CLBS, CLSO), Cotton-Leung® Sof-Flex® (CLSO-SF), Cotton-Huibregtse® (CHBS, CHBSO), ST-2 Soehendra Tannenbaum® (TTSO), and Plastic Biliary Stent (PBS) biliary stents and stent sets, as represented by Figures 1–4.

These devices are available in two device configurations:

- Stent only [with or without pigtail straightener/flap protector (positioning sleeve*)]
- Stent set, which includes a combination of a stent (with or without pigtail straightener/flap protector) and an introducer (pushing catheter and, where applicable, a guiding catheter).

*Note: The words flap protector and positioning sleeve refer to the same item.

Performance Characteristics

The function and key features of these components are described below.

A variety of biliary stents that provide drainage are available across the device range to accommodate various patient anatomies. They are offered in French sizes of between 5 Fr and 11.5 Fr, in various shape configurations, and in labelled lengths of between 2 cm and 21 cm. Biliary stents have the following features:

- Atraumatic tip – The ends of the stent are designed to be atraumatic to the anatomy
- Sideports – Aid drainage (not applicable to ST-2 Soehendra Tannenbaum)
- Anti-migration features – Aim to help maintain the position of the stent once it is placed by providing mechanical interference with the anatomy
- Radiopacity – The stent materials are radiopaque to facilitate visibility under fluoroscopy.

Flap protector – This is a tube used to collapse the flap(s) as it is introduced into the endoscope.

Pigtail straightener – This is a tube used to straighten the pigtail to aid in the introduction of the wire guide.

Guiding catheter – The stent is loaded onto the guiding catheter. The function of the guiding catheter is to guide the stent to its intended location.

Guiding catheters contain one or more of the following features:

- Atraumatic tip – The end of the guiding catheter is designed to be atraumatic to the anatomy
- Radiopacity – The guiding catheter tubing material is radiopaque. The guiding catheter also contains marker bands to further enhance its fluoroscopic visibility
- Hub – The guiding catheter contains a hub with a connection to facilitate the attachment of a standard Luer syringe.

Pushing catheter – The function of the pushing catheter is to advance the stent over a pre-positioned wire guide and, where applicable, a guiding catheter to its intended location, and to maintain the position of the stent as it's being deployed. It contains the following feature:

- Proximal identifier – An identifier at its proximal end. It is used to identify the orientation in which the device must be used.

See product label for all dimensions.

Device Compatibility

The Zimmon, ST-2 Soehendra Tannenbaum, Cotton-Leung, Cotton-Leung Sof-Flex, Cotton-Huibregtse and Plastic Biliary Stent (PBS) biliary stents are compatible with the following:

- Endoscope with 3.2 mm, 3.7 mm, or 4.2 mm accessory channel (refer to label for selection of the appropriate accessory channel)
- 0.035" wire guide (refer to label for the wire guide size)
- Endoscopic cap or wire guide locking device (10 Fr Sof-Flex products are not compatible with Fusion wire guide locking device)

- Stent retriever or forceps
- Cook Medical pushing catheter
- Cook Medical guiding catheter
- Contrast media
- Sterile water or saline
- Standard Luer syringe
- Water-soluble lubricant

Table 1 lists the compatible introducers for biliary stents supplied as stent only. **Note:** Stent sets CLBS, ZEBD, and CHBS contain an introducer within the sets.

Table 1: Compatible Cook Introducers/Introduction Systems for Cook Biliary Stents

Cook Biliary Stent (Fr refers to stent diameter in French size)	Compatible Cook Introducers/Introduction Systems ^a						
	Fusion® Pushing Catheter	Pushing Catheter	Guiding Catheter	Huibregtse Guiding Catheter	Stent Introduction System	Oasis® Introduction System	Fusion® Oasis® Introduction System
5 Fr CLSO, ZSO	FS-PC-5	PC-5, PC-5E	—	—	—	—	—
6 Fr ZSO	—	PC-6	—	—	—	—	—
7 Fr CHBSO, CLSO, CLSO-SF, PBS ^b , ZSO	FS-PC-7	PC-7, PC-7E	—	—	—	—	—
8 Fr ZSO	—	PC-8	—	—	—	—	—
8.5 Fr CHBSO ^c , CLSO ^c , PBS ^b , TTSO	—	PC-8.5	GC-5	—	SIS-8.5	OA-8.5, OA-8.5E	FS-OA-8.5
10 Fr CHBSO ^c , CLSO ^c , PBS ^b , TTSO	—	PC-10	GC-6	HGC-6	SIS-10	OA-10, OA-10E	FS-OA-10
10 Fr CLSO-SF	—	PC-10	GC-5	—	—	—	FS-OA-10
10 Fr ZSO	—	PC-10	—	—	—	—	—
11.5 Fr CHBSO, CLSO, PBS ^b , TTSO	—	PC-11.5	GC-6	HGC-6	SIS-11.5	OA-11.5	—

^aIf using Fusion® Pushing Catheter, Pushing Catheter, Guiding Catheter, Oasis® Introduction System, or Fusion® Oasis® Introduction System, refer to the introducer/introduction system IFU for the instructions on device usage.

^bPBS devices are not available in all markets, including the US.

^cWhen using the Fusion® Oasis® Introduction System (FS-OA) with the CLSO and CHBSO stent families, they are not compatible with stent lengths greater than 14 cm.

Qualitative and Quantitative Information

The materials for the stent implant are outlined in Table 2.

Table 2: Stent Implant Materials

Product	Qualitative Information Device Material		Quantitative Information Weight (g)
5 Fr Cotton-Leung (CLSO) and Zimmon (ZSO)	Polyethylene	Stent	Up to 0.40
6 Fr Zimmon (ZSO)	Polyethylene	Stent	Up to 0.58
7 Fr Cotton-Leung (CLSO), Plastic Biliary Stent (PBS), Cotton-Huibregtse (CHBSO), and Zimmon (ZSO)	Ethylene-vinyl acetate (EVA) copolymer	Stent	Up to 1.04
7 Fr Cotton-Leung Sof-Flex (CLSO-SF)	Polyurethane	Stent	Up to 0.64
8 Fr Zimmon (ZSO)	Polyethylene	Stent	Up to 1.12
8.5 Fr Cotton-Leung (CLSO) and Cotton-Huibregtse (CHBSO)	Polyethylene	Stent	Up to 0.76
8.5 Fr Plastic Biliary Stent (PBS)	EVA copolymer	Stent	Up to 0.76
8.5 Fr ST-2 Soehendra Tannenbaum (TTSO)	Polytetrafluoroethylene	Stent	Up to 1.39
10 Fr ST-2 Soehendra Tannenbaum (TTSO)	Polytetrafluoroethylene	Stent	Up to 1.42
10 Fr Cotton-Leung (CLSO), Plastic Biliary Stent (PBS), and Cotton-Huibregtse (CHBSO)	EVA copolymer	Stent	Up to 1.19
10 Fr Zimmon (ZSO)	Polyethylene	Stent	Up to 2.05
10 Fr Cotton-Leung Sof-Flex (CLSO-SF)	Polyurethane	Stent	Up to 1.16
11.5 Fr Cotton-Leung (CLSO), Plastic Biliary Stent (PBS), and Cotton-Huibregtse (CHBSO)	EVA copolymer	Stent	Up to 1.59
11.5 Fr ST-2 Soehendra Tannenbaum (TTSO)	Polytetrafluoroethylene	Stent	Up to 1.92

Patient Population

Adult patients requiring biliary stenting for obstruction. The nature of the underlying pathology is prevalent in varying patients, therefore the devices are indicated for patients with biliary duct obstruction caused by common bile duct stones, malignant biliary obstruction, and benign or malignant strictures.

Intended User

This device is intended for use by physicians trained and experienced in endoscopic retrograde cholangiopancreatography (ERCP) techniques.

Standard techniques for the endoscopic placement of stents should be employed.

Contact with Body Tissue

The device is tissue contacting in line with intended use.

Operating Principles

The biliary stent operates by providing a lumen through which biliary fluid can drain. The stent has an anti-migration feature which provides mechanical interference with the anatomy to help maintain the position of the stent. The stent is radiopaque to facilitate visibility under fluoroscopy. Sideport(s) on the stent aid drainage by providing additional openings for fluid flow. The biliary stent is placed over a pre-positioned wire guide using an introducer. An introducer consists of a pushing catheter only or a pushing catheter and a guiding catheter combined. The guiding catheter contains marker bands which are radiopaque to facilitate visibility under fluoroscopy and guides stent placement. The pushing catheter operates by pushing the stent into place along a pre-positioned wire guide. To protect the flaps and straighten the pigtail on some of the stents a flap protector or pigtail straightener is provided. The flap protector or pigtail straightener operates by providing a lumen through which the stent can be inserted. The flap protector collapses the flaps on the stent before inserting into the endoscope. The pigtail straightener straightens the pigtail to aid in the introduction of the wire guide.

INTENDED USE

This device is used to drain obstructed biliary ducts.

INDICATIONS FOR USE

Endoscopic biliary stent placement for biliary drainage of obstructed ducts that could be caused by common bile duct stones, malignant biliary obstruction, benign or malignant strictures or other obstructed biliary conditions requiring drainage.

CLINICAL BENEFITS

Drainage of the biliary duct.

CONTRAINDICATIONS

Those specific to ERCP.

Inability to pass a wire guide, guiding catheter, and/or stent through the obstructed area.

WARNINGS

- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or migration and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
- These biliary stents are intended up to a maximum of 3 months indwell.

PRECAUTIONS

- A complete diagnostic evaluation of the patient prior to use to determine the proper stent size should be carried out.
- Refer to the label for the minimum channel size required for this device.
- Refer to the device label for the selection of the appropriate wire guide.
- The stent must be placed under fluoroscopic monitoring.
- The flap protector (if any) or the pigtail straightener (where supplied) is not intended for use in the accessory channel of the endoscope.

- Care must be exercised when straightening the pigtail curls (if any) in order to avoid kinking or breaking the stent.
- Do not use excessive force to advance the stent.
- Refer to Figures 1 to 4 for the intended placed orientation of the stent in the biliary duct.
- Periodic evaluation of the device is recommended during indwell period.
- Select the Cook stent introducer system (if not included) in the appropriate French size. See Table 1 for guidance.
- Sphincterotomy is not necessary for the device placement.
- Dislodgement of a placed stent is possible when attempting additional procedures.
- Do not use this device for any purpose other than the stated intended use.
- Store the device in a dry location.
- Use of this device restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in ERCP techniques.
- Standard techniques for the endoscopic placement of stents should be employed.
- Where used on a label, these symbols indicate the following:

Not compatible with Fusion wire guide locking device.

Use with FS-OA-10 or PC-10 and GC-5.

MRI SAFETY INFORMATION

This symbol means the stent is MR Safe.

The Zimmon® Biliary Stent, ST-2 Soehendra Tannenbaum® Biliary Stent, Cotton-Huibregtse® Biliary Stent, Cotton-Leung® Biliary Stent (including Sof-Flex®), and Plastic Biliary Stent are MR Safe.

In non-clinical testing, the image artifact caused by the device extends approximately 7 mm from the stent when imaged with a gradient echo pulse sequence and a 3 T MR system.

POTENTIAL ADVERSE EVENTS

Those associated with ERCP: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • liver abscess • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Those associated with biliary stent placement: fever • obstruction of the pancreatic duct • pain/discomfort • stent migration • stent occlusion • trauma to the biliary duct or duodenum.

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a peel-open pouch.

Sof-Flex Stents only: Following sterilization, the pouch is further packaged into a box.

This device is accompanied by an implant card that should be given to the patient after it has been completed by the healthcare professional.

DEVICE PREPARATION

1. Prior to introducing any devices, ensure working channel of endoscope is lubricated with water or a water-based lubricant.
2. Referring to the device package label, ensure that the endoscope has a channel size greater than or equal to the minimum channel size required to operate the device.
3. Prior to using the device(s) visually inspect it/them for abnormalities that could result in improper function.
4. Attach endoscopic cap or wire guide locking device to endoscope.
5. For 10 Fr Cotton-Leung Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant. Retain the flap protector for use when introducing the stent flaps into the accessory channel.

INSTRUCTIONS FOR USE

I. FOR NON-PIGTAIL STENTS:

1. Gently ensure the full extension of all the flaps.
2. Load the flap protector onto the duodenal flap end of the stent.

Illustrations

For 5 Fr to 7 Fr stents:

- 3a. Introduce the stent, ductal end first, and the flap protector onto a pre-positioned wire guide.
- 3b. Advance the pushing catheter over the wire guide to advance the stent into the accessory channel. As the front flap of the stent enters the accessory channel, keep the flap protector over the back flap until the stent is completely in the accessory channel.

For 8.5 Fr and larger stents:

- 3c. Remove the hub from the end of the guiding catheter, then introduce the guiding catheter into the accessory channel over a pre-positioned wire guide. For 10 Fr Cotton-Leung® Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant. Retain the flap protector for use when introducing the stent flaps into the accessory channel.
- 3d. With the elevator open, advance until a sufficient length of the guiding catheter is above the obstruction.
- 3e. Introduce the stent, ductal end first, and the flap protector onto the guiding catheter and the pre-positioned wire guide.
- 3f. Advance the pushing catheter over the wire guide and the guiding catheter to advance the stent into the accessory channel. As the front flap of stent enters the accessory channel, keep the flap protector over the back flap until the stent is completely in the accessory channel.

Then refer to steps 4-9 below**II. FOR PIGTAIL STENTS:**

1. Straighten the ductal pigtail. Use the pigtail straightener where provided.
2. Introduce the stent, ductal end first, and the pigtail straightener (where supplied) onto the pre-positioned wire guide until the straightener reaches the second curl.
3. Advance the pushing catheter over the wire guide to advance the pigtail stent into the accessory channel.

Then refer to steps 4-9 below

4. As the pushing catheter advances the stent completely into the accessory channel, slide the flap protector/pigtail straightener (where supplied) back over the pushing catheter until it reaches the end of the catheter, keeping it clear of the accessory channel.
5. Advance the pushing catheter in small increments until the stent is in the desired position.
6. Fluoroscopically and endoscopically confirm the desired stent position.
Note: For non-pigtail stents of 8.5 Fr or greater only: If desired, replace the hub on the guiding catheter before flushing the contrast media.
7. After confirming the stent position, gently remove the wire guide and the guiding catheter (where supplied) from the endoscope while maintaining the position of the stent with the pushing catheter.
8. Gently remove the pushing catheter from the accessory channel.
9. These stents may be removed under standard endoscopic techniques.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed in accordance with institutional guidelines. The introducer (if supplied) should be coiled for disposal.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

When available, the EUDAMED website (<https://ec.europa.eu/tools/eudamed>), along with the BUDI for this product (0827002CIRL202007013010B6), can be used to locate the Summary of Safety and Clinical Performance (SSCP) for this product.

Patient-facing information may be accessed at cookmedical.eu/patient-implant-information.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2024
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.