
EN 3	Biliary Stent Instructions for Use
CZ 4	Biliární stent Návod na použití
DA 6	Galdestent Brugsanvisning
NL 8	Biliaire stent Gebruiksaanwijzing
FR 10	Endoprothèse biliaire Mode d'emploi
DE 12	Gallenstent Gebrauchsanweisung
EL 13	Ενδοπρόσθεση χοληφόρων Οδηγίες χρήσης
HU 16	Epeúti sztent Használati utasítás
IT 17	Stent biliare Istruzioni per l'uso
NO 19	Gallestent Bruksanvisning
PL 21	Stent do dróg żółciowych Instrukcja użycia
PT 23	Stent biliar Instruções de utilização
RO 25	Stent biliar Instrucțiuni de utilizare
ES 27	Stent biliar Instrucciones de uso
SV 29	Biliär stent Bruksanvisning

Duodenal end 7 mm flaps

7 mm křídélka duodenálního konce
 7 mm flapper ved duodenalende
 7 mm slippen aan duodenale uiteinde
 Rabats de 7 mm à l'extrémité duodénale
 7-mm-Zungen am Duodenum-Ende
 Γλωττίδες 7 mm άκρου δωδεκαδάκτυλου
 Duodenalis vég 7 mm-es rögzítőszárnyai
 Alette da 7 mm sull'estremità duodenale
 7 mm fliker ved duodenalenden
 Wąsy na końcówce dwunastniczej 7 mm
 Projecções de 7 mm da extremidade duodenal
 Aripioare de 7 mm ale capătului duodenal
 Lengüetas de 7 mm del extremo duodenal
 Duodenala 7 mm ändflikar

15 mm flaps

15 mm křídélka
 15 mm flapper
 15 mm slippen
 Rabats de 15 mm
 15-mm-Zungen
 Γλωττίδες 15 mm
 15 mm-es rögzítőszárnyak
 Alette da 15 mm
 15 mm fliker
 Wąsy 15 mm
 Projecções de 15 mm
 Aripioare de 15 mm
 Lengüetas de 15 mm
 15 mm flikar

Positioning sleeve

Polohovací rukávec
 Positioneringshylster
 Positioneringshuls
 Gaine de positionnement
 Positionierhülse
 Χιτώνιο τοποθέτησης
 Pozicionálöhüvely
 Manicotto di posizionamento
 Posisjoneringshylse
 Rękaw umiejscawiający
 Manga de posicionamento
 Manşon de poziționare
 Manguito posicionador
 Placeringshylsa

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
INTENDED USE

This device is used to drain obstructed biliary ducts.

NOTES

This device is intended for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

This device must be stored in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

For 10 FR Sof-Flex® stents, lubricate the guiding catheter with water-soluble lubricant.

Retain the positioning sleeve (if any)* for use when introducing the stent flaps into the accessory channel.

CONTRAINDICATIONS

Those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Inability to pass wire guide or stent through obstructed area.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum, obstruction of the pancreatic duct, stent migration.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

This stent must be placed under fluoroscopic monitoring.

The positioning sleeve* or the pigtail straightener* is not intended for use in the accessory channel of the endoscope.

Care must be exercised when straightening the pigtail curls* in order to avoid kinking or breaking the stent.

The tapered tip end of the stent* or side holes* must be positioned in the common bile duct while the other end remains in the duodenum.

The longer flaps of the ST-2 Biliary Stent should be positioned in the duct while the shorter flaps remain in the duodenum. (See fig. 1)

This device should not be left indwelling for more than three months or as directed by a physician. Periodic evaluation is recommended.

Select the Cook stent introducer system (if not included) in the appropriate French size.

WARNING

For stents with radiopaque bands, MRI compatibility has not been established.

INSTRUCTIONS FOR USE**I. FOR NON-PIGTAIL STENTS:**

1. Gently ensure full extension of all side flaps*.
2. Load positioning sleeve onto duodenal flap end of stent.

Illustrations

For 5 to 7 FR stents:

3.a. Introduce stent, tapered tip first, and positioning sleeve onto a pre-positioned wire guide.

For 8.5 FR and larger stents:

3.b. Remove Tuohy-Borst adapter from end of guiding catheter, then introduce guiding catheter into accessory channel over a pre-positioned wire guide.

3.c. Introduce stent, tapered tip first, and positioning sleeve onto guiding catheter and pre-positioned wire guide.

THEN REFER TO STEPS 4-7 BELOW

II. FOR PIGTAIL STENTS:

1. Advance pigtail straightener along shaft of stent to tapered tip end in order to straighten pigtail curl. (See fig. 2)

2. Introduce stent, tapered tip first, and pigtail straightener onto pre-positioned wire guide until straightener reaches second curl.

3. Advance pushing catheter over wire guide to advance pigtail stent into accessory channel.

THEN REFER TO STEPS 4-7 BELOW

4. Advance guiding catheter* and/or pushing catheter in 1-2 cm increments until stent is in desired position.

5. Fluoroscopically and endoscopically confirm desired stent position.

Note (for non-pigtail stents only): Replace Tuohy-Borst adapter onto end of guiding catheter*. Inject contrast, if desired, to fluoroscopically visualize stent position.

6. After confirming stent position, gently remove wire guide, then guiding catheter* from endoscope while maintaining position of stent with pushing catheter.

7. Gently remove pushing catheter from accessory channel.

Upon completion of procedure, dispose of introducer components per institutional guidelines for biohazardous medical waste.

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Sof-Flex and Cook are registered trademarks of Cook Incorporated.
© COPYRIGHT COOK 2015 ©
COOK IRELAND LTD.