EN 3	<b>Zilver® Vena™ Venous Self-Expanding Stent</b> Instructions for Use
BG 14	<b>Венозен саморазширяващ се стент Zilver® Vena™</b> Инструкции за употреба
CS 25	<b>Žilní samoexpandující stent Zilver® Vena™</b> Návod na použití
DA 36	<b>Zilver® Vena™ selvekspanderende venestent</b> Brugsanvisning
NL 47	<b>Zilver® Vena™ veneuze zelfexpanderende stent</b> Gebruiksaanwijzing
ET 58	<b>Iselaienev veenistent Zilver® Vena™</b> Kasutusjuhised
FR 69	<b>Endoprothèse veineuse auto-expansible Zilver® Vena™</b> Mode d'emploi
DE 80	<b>Zilver® Vena™ selbstepandierender Venenstent</b> Gebrauchsanweisung
EL 91	<b>Αυτοεκτεινόμενη φλεβική ενδοπρόσθεση Zilver® Vena™</b> Οδηγίες χρήσης
HU 103	<b>Zilver® Vena™ öntáguló vénaszent</b> Használati utasítás
IT 114	<b>Stent venoso autoespansibile Zilver® Vena™</b> Istruzioni per l'uso
LV 125	<b>Zilver® Vena™ venozais stents ar pašpaplašināšanās funkciju</b> Lietošanas instrukcija
LT 136	<b>„Zilver® Vena™“ veninis savaimė išsiplečiantis stentas</b> Naudojimo nurodymai
NO 147	<b>Zilver® Vena™ venøs selvekspanderende stent</b> Bruksanvisning
PL 158	<b>Samorozprężający stent żylny Zilver® Vena™</b> Instrukcja użycia
PT 169	<b>Stent venoso auto-expansível Zilver® Vena™</b> Instruções de utilização
RO 180	<b>Stent venos auto-expansiv Zilver® Vena™</b> Instrucțiuni de utilizare
SK 191	<b>Venózny samorozširovací stent Zilver® Vena™</b> Návod na použitie

ES  
202

## **Stent venoso autoexpandible Zilver® Vena™**

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SV  
213

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Bruksanvisning

## ZILVER® VENA™ VENOUS SELF-EXPANDING STENT

**Do not re-sterilize.**

### DEVICE DESCRIPTION

The Zilver® Vena™ Venous Stent is a self-expanding, flexible, slotted-tube nitinol stent. Post-deployment, the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region.

Upon deployment, the stent provides support, while maintaining flexibility in the vessel.

The Zilver® Vena™ Venous Stent is available in the following sizes:

7 French (2.3 mm) Delivery System							
Stent Length (mm)		60		100		140	
Delivery System (cm)		80	120	80	120	80	120
Stent Inner Diameter (mm)	14	x	x	x	x	x	x
	16	x	x	x	x	x	x

The stent comes preloaded in a 7Fr (2.3 mm) delivery catheter. Hand-loading of the stent is not possible. Stent deployment is controlled by retraction of the handle while holding the metal cannula stationary.

### INTENDED USE

The Zilver® Vena™ Venous Stent is intended for use in the iliofemoral veins for the treatment of symptomatic venous outflow obstruction.

### CONTRAINDICATIONS

Patients with total venous occlusion that cannot be dilated to allow passage of the introducer system or wire guide.

### WARNINGS

Nitinol (nickel titanium) may cause allergic reactions in some patients.

### PRECAUTIONS

- This product should only be used by physicians trained and experienced in diagnostic and interventional vascular techniques. Standard techniques for interventional vascular procedures should be employed.

- Manipulation of the Zilver® Vena™ Venous Stent requires high-resolution fluoroscopic control.
- Do not use power injection systems with the delivery system.
- Prior to the procedure, the patient's underlying condition should be assessed for compatibility with anticipated procedural and post-procedural antiplatelet/anticoagulation therapy.
- Use in patients with a history of contrast sensitivity is not recommended unless the patient can be adequately premedicated.
- Arterial use is not recommended.

### **Stent Handling**

- Do not attempt to remove the stent from the delivery system before use.
- Do not expose any part of the delivery system to organic solvents (e.g., alcohol).
- Carefully inspect the sterile package and stent system prior to use to verify that neither has been damaged during shipment.
- Use the stent system prior to the expiration date specified on the package.
- The device is designed for single use only. Attempts to reprocess, re-sterilize, and/or reuse may lead to device failure and/or transmission of disease.

### **Stent Placement**

- Ensure that the safety lock is not inadvertently removed prior to stent release.
- Do not rotate any part of the system during deployment.
- Repositioning of the device after deployment is not possible since the introducer catheter cannot be re-advanced over the stent once deployment begins.
- Once stent deployment has begun, the stent must be fully deployed.

### **Stent/System Removal**

- Do not advance sheath after stent has been deployed. Delivery system can be removed without the need to recapture tip.

### **Post Implant**

- Antiplatelet/anticoagulant therapy should be administered during and after procedure according to institutional standard of care.
- Use caution when re-crossing a stent to avoid stent damage or migration.

## MRI INFORMATION

Non-clinical testing has demonstrated that the Zilver® Vena™ Venous stent is MR Conditional according to ASTM F2503. A patient with this stent may be scanned safely under the following conditions.

### Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less
- Spatial magnetic gradient field of 1580 Gauss/cm or less

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

### MRI-Related Heating

- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of scanning (i.e., per scanning sequence);
  - 1.5 Tesla Temperature Rise  
In non-clinical testing, the Zilver® Vena™ Venous stent produced maximum temperature rises of 2.6 and 3.8 °C (for a single stent and pair of overlapping stents) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens Magnetom) at an MR system reported whole-body-averaged SAR of 3.8 W/kg (associated with a calorimetry measured whole body averaged value of 3.1 W/kg).
  - 3.0 Tesla Temperature Rise  
In non-clinical testing, the Zilver® Vena™ Venous stent produced a maximum temperature rises of 2.2 and 4.6 °C (for a single stent and pair of overlapping stents) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3 Tesla System (General Electric Excite) at an MR system reported whole-body-averaged SAR of 3.0 W/kg (associated with a calorimetry measured whole body averaged value of 2.8 W/kg).

## **Image Artifact**

MR image quality may be compromised if the area of interest is within approximately 7 mm of the position of Zilver® Vena™ Venous stent as found during non-clinical testing using T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

## **For US Patients Only**

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International  
2323 Colorado Avenue  
Turlock, CA 95382, USA

Phone: 888-633-4298 (toll free)  
209-668-3333 from outside the US

Fax: 209-669-2450

Web: [www.medicalert.org](http://www.medicalert.org)

## **POTENTIAL ADVERSE EVENTS**

Potential adverse events that may occur include, but are not limited to, the following:

- Abrupt stent closure
- Allergic reaction to anticoagulant and/or antithrombotic therapy or contrast medium
- Allergic reaction to Nitinol
- Arrhythmia
- Arteriovenous fistula
- Death
- Embolism
- Fever
- Hematoma/hemorrhage at access site
- Hypersensitivity reactions

- Hypertension
- Hypotension, nausea or symptoms of a vasovagal response
- Infection/abscess formation at access site
- Intimal injury/dissection
- Myocardial infarction (MI)
- Pseudoaneurysm formation
- Pulmonary embolism
- Renal failure
- Restenosis or thrombosis of the stented vein
- Retinal, retroperitoneal, gastrointestinal, or genitourinary bleeding associated with anticoagulant therapy
- Septicemia/bacteremia
- Spasm
- Stent malapposition
- Stent migration
- Stent strut fracture
- Stroke
- Tissue necrosis
- Vessel perforation/rupture
- Worsened pain

## **PRODUCT RECOMMENDATIONS**

### **Venous Access**

For venous access, the use of an access set that accepts a 7.0 French (2.3 mm) introducer catheter is recommended.

### **Wire Guide Selection**

The use of a .035-inch (0.89 mm) wire guide is recommended.

### **PTA Balloon Selection**

For pre- and post-dilation, an appropriately sized balloon catheter is recommended.

### **Stent Selection**

The chosen stent diameter should be at least 1 mm larger than the estimated diameter of the target vein.

## HOW SUPPLIED

The Zilver® Vena™ Venous Stent is supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## ZILVER® VENA™ VENOUS STENT (Fig. 1)


Fig. 1

- a. Handle
- b. Hub
- c. Safety Lock
- d. Delivery System: Outer Sheath
- e. Tip of Delivery System Inner Catheter
- f. Side-arm Flushing Port
- g. Metal Cannula
- h. Radiopaque Marker on the Delivery System
- i. Gold Radiopaque Markers

The stent is supplied preloaded in a 7.0 French (2.3 mm) introducer catheter. A radiopaque marker (h) on the distal tip of the sheath is used to visualize deployment of the stent. Stent deployment is controlled by means of a hand-held device.

## Multiple Stent Placement

If placement of more than one stent is required, the following recommendations should be considered:

- When more than one stent is required, resulting in stent-to-stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.

- In relation to the operator, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem must overlap to allow for complete coverage of the lesion.

## INSTRUCTIONS FOR USE

### Stent Sizing

1. Determine the proper stent size after complete diagnostic evaluation. The stent diameter should be at least 1 mm larger than the estimated diameter of the target vein. Measure the total length of the obstructed vein to determine the length of the segment that requires the stent and the length of the stent needed. Stent from the most distal end of the lesion to the most proximal end of the lesion in a continuous fashion.

**NOTE:** If multiple stents are required to cover the length of the lesion, please refer to the Multiple Stent Placement section of these instructions for use for further recommendations.

### Introduction of the Stent

1. Gain access at the appropriate site utilizing an introducer sheath or guiding catheter. Refer to product labeling for required introducer sheath or guiding catheter size.
2. Introduce an .035 inch (0.89 mm) wire guide through the introducer sheath or guiding catheter across the distal segment of the target lesion.
3. Predilatation before stent placement is at the discretion of the physician. Remove the balloon catheter, leaving the wire guide in place.



Fig. 2

4. Immediately before placing the delivery system into the body, use a 1 ml syringe to flush the delivery system with saline through the side-arm flushing port. Flush until a few drops of saline exit the distal tip, between the delivery system outer sheath (d) and inner catheter (e). (**Fig. 2**)



5. Use the 1 ml syringe to flush the wire guide lumen of the stent delivery system with saline through the hub. (**Fig. 3**)

**Fig. 3**

6. Insert the delivery system over the wire guide.

7. Under fluoroscopy, advance the delivery system past the most distal segment of the lesion.

**Fig. 4**

8. Pull back on the stent delivery system under fluoroscopy until the radiopaque markers on the stent (i) are at the desired position. The stent is now ready to be deployed. (**Fig. 4**)

## Deployment of the Stent

1. Before deployment, it is important to straighten the proximal part of the delivery system as much as possible, remove any slack in the delivery system and to keep the handle in a stable position.
2. The stent expansion must be performed under fluoroscopic control.


Fig. 5

3. Hold the hub (b) on the metal cannula (g) steady. To deploy the stent, remove the red safety lock (c). (**Fig. 5**)




Fig. 6

4. Hold the hub end stationary. The stent will be deployed as you pull the handle (a) toward the hub (b). (**Fig. 6**) The radiopaque marker on the delivery system (h) indicates the progress of the deployment.



**NOTE:** Full deployment of the stent length will occur when the distal end of the sheath has been retracted past the proximal part of the stent.

**Fig. 7**

5. As deployment occurs, continue sliding the handle (a) toward the hub (b) in a slow, smooth and consistent fashion. (**Fig. 7**)

**NOTE:** Once stent deployment has begun, the stent must be fully deployed. Repositioning of the Zilver® Vena™ Venous Stent is not possible since the delivery system's outer sheath cannot be re-advanced over the stent once deployment begins. Refer to the Multiple Stent Placement section of these instructions for use for information on missed lesions.

**Fig. 8**

6. The stent is fully deployed when the handle (a) reaches the hub (b). (**Fig. 8**)

## **Post Stent Deployment**

**1.** Delivery System Removal - Do not advance sheath after the stent has been deployed. Delivery system can be removed without the need to recapture tip. Check for delivery system integrity after removal from the patient.

**2.** Perform a venogram to ensure full deployment of the device.

**NOTE:** If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation can be performed at the discretion of the physician. Select an appropriate size balloon catheter and dilate the lesion with conventional technique. The inflation diameter of the balloon used for post dilatation should approximate the diameter of the reference vessel. Remove the balloon from the patient.

**3.** Introducer sheath and wire guide may be removed at this point.

**4.** Close the entry wound as appropriate.

## **REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COPYRIGHT COOK 2017©
**COOK IRELAND LTD.**