
COOK®
MEDICAL

CE
0123

EN 9	EchoTip Ultra® Endobronchial HD Ultrasound Needle for Olympus Scopes Instructions for Use
BG 13	Ендобронхиална ултразвукова игла с висока разделителна способност EchoTip Ultra® за ендоскопи Olympus Инструкции за употреба
CS 18	Endobronchiální jehla EchoTip Ultra® pro ultrazvuk s vysokým rozlišením, pro endoskopy Olympus Návod k použití
DA 22	EchoTip Ultra® endobronkial HD ultralydsnål til Olympus skoper Brugsanvisning
DE 26	EchoTip Ultra® endobronchiale HD-Ultraschallnadel für Olympus-Endoskope Gebrauchsanweisung
EL 31	Ενδοβρογχική βελόνα υπερήχων υψηλής ανάλυσης EchoTip Ultra® για ενδοσκόπια Olympus Οδηγίες χρήσης
ES 36	Aguja de ecografía endobronquial de alta definición EchoTip Ultra® para endoscopios Olympus Instrucciones de uso
ET 41	Endobronhiaalne kõrglahutusega ultrahelinõel EchoTip Ultra® endoskoopide Olympus jaoks Kasutusjuhend
FI 45	Endobronkiaalin EchoTip Ultra® -HD-ultraäänineula Olympus-skooppeja varten Käyttöohjeet
FR 50	Aiguille pour écho-endoscopie haute définition endobronchique EchoTip Ultra® pour les endoscopes Olympus Mode d'emploi
HU 54	EchoTip Ultra® endobronchialis nagy felbontású ultrahangos tű Olympus endoszkópokhoz Használati utasítás
IT 59	Ago ecografico endobronchiale ad alta definizione EchoTip Ultra® per endoscopi Olympus Istruzioni per l'uso
LT 64	„EchoTip Ultra®“ endobronchinio didelės skiriamosios gebos ultragarso adata, skirta naudoti su „Olympus“ endoskopais Naudojimo instrukcija
LV 68	EchoTip Ultra® endobronhiālā augstas precizitātes ultraskaņas adata Olympus endoskopiem Lietošanas instrukcija
NL 72	EchoTip Ultra® endobronchiale HD ultrasone naald voor Olympus-scopen Gebruiksaanwijzing
NO 77	EchoTip Ultra® endobronkial HD-ultralýdnål for Olympus-skop Bruksanvisning
PL 81	Igła EchoTip Ultra® do ultrasonografii wewnątrzoskrzelowej o wysokiej rozdzielczości do endoskopów Olympus Instrukcja użycia
PT 86	Agulha de ultrassons endobrônquica EchoTip Ultra® de alta definição para endoscópios Olympus Instruções de utilização
RO 90	Ac ecografic HD endobronșic EchoTip Ultra® pentru endoscoape Olympus Instrucțiuni de utilizare
SK 95	Endobronchiálna ultrazvuková ihla s vysokým rozlíšením EchoTip Ultra® pre endoskopy Olympus Návod na použitie
SV 99	EchoTip Ultra® endobronkial HD ultraljudsnål för Olympus endoskop Bruksanvisning
TR 103	Olympus Skopları için EchoTip Ultra® Endobronşiyal HD Ultrason İğnesi Kullanma Talimatı

I F U 0 0 5 1 - 1 1

Fig. 1**EN Biopsy Port Seal**

BG	Уплътнение на отвора за биопсия
CS	Těsnění biotického portu
DA	Forsegling på biopsiport
DE	Abdichtung des Biopsieanschlusses
EL	Στεγανωτικό παρέμβυσμα θύρας βιοψίας
ES	Junta del acceso de biopsia
ET	Biopsiapordi tihend
FI	Biopsiaportin tiivistä
FR	Capuchon étanche de l'orifice pour biopsie
HU	Biopsziányílás tömítése
IT	Guarnizione del raccordo per biopsia
LT	Biopsijos angos uždoris
LV	Biopsijas pieslēgvieta blīve
NL	Afsluiting van de biopsiepoort
NO	Biopsiportforsegling
PL	Uszczelka portu biopsji
PT	Vedante do orifício de biopsia
RO	Sigiliul portului de biopsie
SK	Uzáver biotického portu
SV	Biopsiportstätning
TR	Biopsi Port Mührü

EN Luer Fitting

BG	Сглобка тип Luer
CS	Spojka Luer
DA	Luer-fitting
DE	Luer-Anschluss
EL	Εξάρτημα Luer
ES	Conexión luer
ET	Luer-liitmik
FI	Luer-liitin
FR	Raccord Luer
HU	Luer-illeszték
IT	Attacco Luer
LT	Luerio jungtis
LV	Luera savienojums
NL	Luer-fitting
NO	Luer-kobling
PL	Złącze luer
PT	Conector Luer
RO	Fiting Luer
SK	Spojka Luer
SV	Luer-koppling
TR	Luer Bağlantı

EN Locking Clip

BG	Блокиращ клипс
CS	Aretovací svorka
DA	Låseklemme
DE	Arretierclip
EL	Κλιπ ασφάλισης
ES	Abrazadera de fijación
ET	Kinnitusklamber
FI	Lukitusnippistin
FR	Pince de verrouillage
HU	Rögzítőkapocs
IT	Clip di bloccaggio
LT	Fiksuojamasis spaustukas
LV	Noslēgskava
NL	Borgclip
NO	Låseklips
PL	Zacisk blokujący
PT	Grampo de fixação
RO	Clemă blocantă
SK	Poistná objímka
SV	Låsklämna
TR	Kilitleme Klipsi

EN Adapter Top

BG	Горна част на адаптера
CS	Horní strana adaptéru
DA	Øverste del af adapter
DE	Oberteil des Adapters
EL	Κορυφή του προσαρμογέα
ES	Parte superior del adaptador
ET	Adapteri ülemine osa
FI	Sovittimen yläosa
FR	Haut de l'adaptateur
HU	Adapter teteje
IT	Sommità dell'adattatore
LT	Adapterio viršus
LV	Adaptēra augšdaļa
NL	Adaptertop
NO	Adaptertopp
PL	Część górna łącznika
PT	Parte de cima do adaptador
RO	Capacul adaptorului
SK	Vrch adaptéra
SV	Adaptorns övre del
TR	Adaptörün Üstü

EN Adapter Base

BG	Долна част (основа) на адаптера
CS	Základna adaptéru
DA	Nederste del af adapter
DE	Unterteil des Adapters
EL	Βάση του προσαρμογέα
ES	Base del adaptador
ET	Adapteri alumine osa
FI	Sovittimen pohjaosa
FR	Base de l'adaptateur
HU	Adapter alapja
IT	Base dell'adattatore
LT	Adapterio pagrindas
LV	Adaptēra pamatne
NL	Adapterbasis
NO	Adapterbase
PL	Podstawa łącznika
PT	Base do adaptador
RO	Baza adaptorului
SK	Podstavec adaptéra
SV	Adaptorns bas
TR	Adaptörün Tabanı

Fig. 2


Fig. 3



Fig. 4

EN Plunger

BG Бутало
CS Pist
DA Stempel
DE Kolben
EL Έμβολο
ES Émbolo
ET Kolb
FI Mäntä
FR Poussoir
HU Dugattyú
IT Stantuffo
LT Stūmoklis
LV Virzulis
NL Plunjer
NO Stempel
PL Tłoczek
PT Émbolo
RO Piston
SK Piest
SV Kolv
TR Piston

EN Stopcock "open"

BG Кранче „отворено“
CS „Otevřený“ uzavírací kohout
DA Stophane „åben“
DE Absperrhahn „geöffnet“
EL Στρόφιγγα «ανοικτή»
ES Llave de paso «abierta»
ET Sulgemisventiil „avatud“
FI Sulkuhana „auki“
FR Robinet « ouvert »
HU Elzárócsap nyitott „állásban“
IT Rubinetto “aperto”
LT „Atsuktas“ čiaupas
LV Noslēgkrāns „atvērts“
NL Afsluitkraan ‘open’
NO Stoppekran «åpen»
PL Kranik w pozycji „otwartej”
PT Torneira de passagem “aberta”
RO Robinet „deschis”
SK „Otvorený” kohútik
SV Kran “öppen”
TR Stopkok “açık”

Fig. 5



Fig. 6

EN Sliding Sheath Adjuster

BG Плъзгащ се регулатор на защитната обвивка
CS Posuvné nastavovací zařízení sheathu
DA Glidende sheath-justering
DE Verschiebbarer Schleuseneinsteller
EL Συρόμενος προσαρμογέας θηκariού
ES Ajustador deslizante de la vaina
ET Liughülsi regulaator
FI Holkin liukusäädin
FR Ajusteur de gaine coulissant
HU Elcsúsztható hüvelyhossz-beállító
IT Regolatore scorrevole della guaina
LT Slankusis vamzdelio reguliatorius
LV Pārīdāmais ievadapvalka regulētājs
NL Schuifregeling van de sheath
NO Glidende hylsejusteringsmekanisme
PL Przesuwany regulator koszulki
PT Regulador deslizante da bainha
RO Ajustor glisant al tecii
SK Posuvný nastavovač puzdra
SV Glidande hylsjusterare
TR Kayar Klif Ayarlayıcı

EN Needle Handle

BG Дръжка на иглата
CS Rukojeť jehly
DA Nålegreb
DE Nadelgriff
EL Λαβή βελόνας
ES Mango de la aguja
ET Nõela käepide
FI Neulan kahva
FR Poignée de l'aiguille
HU Tű fogantyúja
IT Impugnatura dell'ago
LT Adatos kotelis
LV Adatas rorktis
NL Naaldhandvat
NO Nålhandtak
PL Uchwyt igły
PT Punho da agulha
RO Mâner al acului
SK Rukováť ihly
SV Nålhandtag
TR İğne Sapi

EN Sheath Reference Mark

BG Означение за дължината на защитната обвивка
CS Referenční značka sheathu
DA Referencemærke på sheath
DE Schleusenbezugsmarke
EL Σήμανση αναφοράς θηκariού
ES Marca de referencia de la vaina
ET Hülsi tähis
FI Holkin viitemerkintä
FR Repère de référence sur la gaine
HU Hüvelyhossz referenciajele
IT Indicatore di riferimento della guaina
LT Vamzdelio atskaitos žymė
LV Ievadapvalka atskaites atzime
NL Referentiemarkering van de sheath
NO Referansemerke på hylsen
PL Znacznik referencyjny wysunięcia koszulki
PT Marca de referência da bainha
RO Marcaj de referință al tecii
SK Referenčná značka puzdra
SV Hylsans referensmarkering
TR Klif Referans İşareti

EN Safety Ring

BG Предпазен пръстен
CS Pojistný prstenec
DA Sikkerhedsring
DE Sicherheitsring
EL Δακτύλιος ασφαλείας
ES Anillo de seguridad
ET Kaitserõngas
FI Turvarengas
FR Bague de sécurité
HU Biztonsági gyűrű
IT Anello di sicurezza
LT Apsauginis žiedas
LV Drošības gredzens
NL Veiligheidsring
NO Sikkerhetsring
PL Pierscień bezpieczeństwa
PT Anel de segurança
RO Inel de siguranță
SK Bezpečnostný krúžok
SV Säkerhetslåsring
TR Güvenlik Halkası

EN Stylet Hub

BG Втулка на стилето
CS Ústí styletu
DA Stiletmuffe
DE Mandrinsatz
EL Ομφαλός του στείλεου
ES Conector del estilete
ET Stiletikinnitus
FI Mandriinin kanta
FR Embase du stylet
HU Mandrin végidoma
IT Pomello del mandrino
LT Zondo kaištis
LV Stiletes galviņa
NL Aanzetstuk van het stilet
NO Stiletteste
PL Złączka mandrynu
PT Conector do estilete
RO Acoperitoarea stiletului
SK Hrdlo styletu
SV Mandrängfättning
TR Stile Göbeği

Fig. 7

ECHOTIP ULTRA® ENDOBRONCHIAL HD ULTRASOUND NEEDLE FOR OLYMPUS SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The EchoTip Ultra® Endobronchial HD Ultrasound Needle for Olympus Scopes is available in two different needle sizes (22 and 25 gage) and contains the following components:

- Needle
- Syringe
- Accessory channel adapter

Performance Characteristics

The function and key features of these components are described below:

Endobronchial needle

- Handle – The handle contains adjustable components to allow the user to adjust the extension of the needle and the sheath. The handle safety ring slides and locks the needle at the desired needle length extension (0–5 cm). The '0 cm' reference mark ensures complete needle retraction within the sheath. The sliding sheath adjuster allows the user to set the required sheath length extension (0–3 cm) for their Olympus endobronchial ultrasound (EBUS) scope. The '0 cm' reference mark denotes the minimum sheath length available. The base of the sliding sheath adjuster has a Luer lock fitting, which facilitates the attachment of the device to the scope accessory channel using the accessory channel adapter.
- 4.1 Fr Sheath – The sheath covers the needle during introduction, advancement, and withdrawal of the device from the scope accessory channel. It protects the scope accessory channel from damage by the needle tip.
- Needle – The stainless steel needle has a dimpling pattern to provide visualization of the needle tip under endobronchial ultrasound. The needle has a beveled tip to puncture and facilitate sample collection from the target site.
- Stylet – The stylet is inserted through the length of the needle. The distal end of the stylet is beveled to match the bevel of the needle tip. The stylet provides support to the needle during advancement into the target site, prevents sample collection outside of the target site, and may also be used for sample retrieval following collection if desired.

Syringe – The device is supplied with a 10 mL syringe to utilize the standard vacuum technique for sample collection, if desired. The Luer lock fitting at the tip of the syringe can be attached to the Luer fitting on the handle. The syringe has a plunger and a stopcock to create a vacuum for sample collection. It is supplied with 5 mL and 10 mL clips attached to the plunger to allow the user to select the preferred volume for sample aspiration.

Accessory channel adapter – The accessory channel adapter consists of a Luer fitting, a biopsy port seal, and a locking clip. The accessory channel adapter attaches to the scope accessory channel port of an Olympus EBUS scope to function as a Luer fitting, to which the needle can be connected.

Device Compatibility

The EchoTip Ultra Endobronchial HD Ultrasound Needle for Olympus Scopes is compatible with:

- Olympus EBUS scope with 2 mm minimum accessory channel
- Accessory channel adapter
- 10 mL Standard Luer syringe

Patient Population

Adult patients requiring fine needle aspiration (FNA) of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or the gastrointestinal tract.

Intended User

Use of this device is restricted to a trained healthcare professional.

Contact with Body Tissue

This device is tissue contacting in line with intended use.

Operating Principles

Olympus EBUS scopes do not have a Luer fitting to allow needle attachment directly to the accessory channel port. The EchoTip Ultra Endobronchial HD Ultrasound Needle for Olympus Scopes is supplied with an accessory channel adapter to function as a Luer fitting, to which the needle can be connected. The device is inserted into the scope accessory channel and attached when the fittings meet.

The sheath length can be adjusted depending on the length of the EBUS scope using the sliding sheath adjuster. The needle length extension can be adjusted as required. The device is used to puncture the target site to obtain a sample. If using a standard vacuum syringe technique, the stylet is then removed, and the syringe is attached to the Luer fitting on the needle handle. Turning the syringe stopcock to the "open" position allows negative pressure to aspirate the sample. The handle is gently moved in small increments back and forth within the target site to obtain a sample. If using the syringe, the stopcock is returned to the "closed" position. The needle is retracted into the sheath and the device is disconnected from the scope accessory channel. Upon completion of the procedure the adapter is removed from the accessory channel port.

INTENDED USE

This device is used to sample targeted submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or gastrointestinal tract through the accessory channel of an ultrasound bronchoscope for fine needle aspiration (FNA).

INDICATIONS FOR USE

This device is indicated for endobronchial ultrasound guided tissue sampling of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or the gastrointestinal tract, which is used for subsequent pathological examination to aid in disease diagnosis and patient management.

CLINICAL BENEFITS

This device facilitates successful sampling of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or gastrointestinal tract to aid in disease diagnosis and patient management.

CONTRAINDICATIONS

Those specific to the primary endobronchial procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications: coagulopathy.

WARNINGS

- Not for use in the heart or vascular system.
- The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.
- This single-use device is not designed for re-use. Attempts to reprocess, resterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- Visually inspect the integrity of the sterile packaging. Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Visually inspect the device with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.
- This device contains nickel, which may cause allergic reaction in individuals with nickel sensitivity.

PRECAUTIONS

- The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement, or withdrawal of the device. Failure to retract the needle may result in damage to the scope.
- Ensure the stylet is fully inserted when advancing the needle into the target site.
- Do not use this device for any purpose other than stated intended use.
- Use of this device is restricted to a trained healthcare professional.
- The product is intended for use by physicians trained and experienced in endobronchial ultrasound techniques.
- When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.
- This device is intended for use with an Olympus EBUS scope.
- Store the device in a dry location.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with an endobronchial ultrasound procedure: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • pain/discomfort • perforation • respiratory depression or arrest.

Those associated with the device: allergic reaction to nickel • damage to blood vessels • death (other than due to normal disease progression) • fever • hemorrhage • infection • inflammation • nerve damage • pain/discomfort • perforation • pneumoperitoneum • sepsis • septicemia/bacteremia • tumor seeding.

HOW SUPPLIED

These devices are supplied ethylene oxide (EO) sterilized in a tray with a Tyvek lid.

PLACEMENT OF ACCESSORY CHANNEL ADAPTER

Illustrations

1. Examine the adapter and note features (**Fig. 1**).
2. Place the adapter, in the open position, in line with the scope accessory channel metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel metal hub (**Fig. 2**).
4. Close the adapter ensuring it is locked (**Fig. 3**).
5. The adapter will now function as a Luer fitting to which the needle can be connected. **Note:** If irrigation of the accessory channel is desired/necessary, irrigate through the closed adapter prior to insertion or after removal of the needle.

SYSTEM PREPARATION

1. Visually inspect the packaging confirming it is unopened and free from damage prior to use.
2. Visually inspect the device components confirming they are free from damage prior to use.
3. Examine the syringe. It has two plunger locks that must be depressed to advance the plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in the "open" position, aligned with the syringe (**Fig. 4**).
4. Prepare the syringe as follows:
 - a. With the stopcock in the "open" position, depress the plunger locks and fully advance the plunger into the syringe.
 - b. Turn the stopcock to the "closed" position.
 - c. Pull the plunger back until it is locked into place at the desired setting, creating a vacuum.
 - d. Set the prepared syringe aside until aspiration is required.

INSTRUCTIONS FOR USE

1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory

- channel of the scope (**Fig. 5**). **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
 4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
 5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the accessory channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 6**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
 6. While maintaining position of ultrasound scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 7**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
 7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
 8. Standard vacuum syringe techniques may be applied for sample collection (Steps 9–11) or, if desired, other techniques that may or may not incorporate the use of this stylet may be used.
 9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
 10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
 11. Turn the stopcock to the “open” position aligned with the syringe, allowing the negative pressure in the syringe to aspirate the sample. Move the needle handle back and forth within the target site. **Note:** Do not remove the needle from the target site during sample collection.
 12. Upon completion of sample collection, if using the syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
 13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
 14. Unlock the thumbscrew on the safety ring and extend the needle.
 15. Use a new syringe or stylet to expel the specimen, then prepare per institutional guidelines.
 16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
 17. For additional sample collection from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
 18. Additional samples may be obtained by repeating Steps 2 through 16 of the “Instructions for Use”.
 19. Upon completion of procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
 20. To remove the adapter from the scope, pull back on the locking clip, lift the top and slide off the scope.

DEVICE DISPOSAL

This device may be contaminated with potentially infectious substances of human origin and should be coiled for disposal in accordance with institutional guidelines.

PATIENT COUNSELLING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2025
William A. Cook Australia Pty Ltd
COOK IRELAND LTD.