
EN 9	<b>EchoTip® Ultra Endobronchial HD Ultrasound Needle for Olympus Scopes</b> Instructions for Use
BG 12	<b>Ендобронхиална ултразвукова игла с висока разделителна способност EchoTip® Ultra за ендоскопи Olympus</b> Инструкции за употреба
CS 17	<b>Endobronchiální jehla EchoTip® Ultra pro ultrazvuk s vysokým rozlišením, pro endoskopy Olympus</b> Návod k použití
DA 21	<b>EchoTip® Ultra endobronkial HD ultralydsnål til Olympus skoper</b> Brugsanvisning
DE 25	<b>EchoTip® Ultra endobronchiale HD-Ultraschallnadel für Olympus-Endoskope</b> Gebrauchsanweisung
EL 29	<b>Ενδοβρογχική βελόνα υπερήχων υψηλής ανάλυσης EchoTip® Ultra για ενδοσκόπια Olympus</b> Οδηγίες χρήσης
ES 34	<b>Aguja de ecografía endobronquial de alta definición EchoTip® Ultra para endoscopios Olympus</b> Instrucciones de uso
ET 39	<b>Endobronhiaalne kõrglahutusega ultrahelinõel EchoTip® Ultra endoskoopidele Olympus</b> Kasutusjuhised
FR 42	<b>Aiguille ultrasonore HD endobronchique EchoTip® Ultra pour les endoscopes Olympus</b> Mode d'emploi
HU 47	<b>EchoTip® Ultra endobronchialis, nagy felbontású, echogén tű Olympus endoszkópokhoz</b> Használati utasítás
IT 51	<b>Ago ecografico endobronchiale ad alta definizione EchoTip® Ultra per endoscopi Olympus</b> Istruzioni per l'uso
LT 55	<b>„EchoTip® Ultra“ endobronchinio didelės skiriamosios gebos ultragarso adata, skirta naudoti su „Olympus“ endoskopais</b> Naudojimo nurodymai
LV 59	<b>EchoTip® Ultra endobronhiālā augstas precizitātes ultraskaņas adata Olympus endoskoopiem</b> Lietošanas instrukcija
NL 63	<b>EchoTip® Ultra endobronchiale HD ultrasone naald voor Olympus endoscopen</b> Gebruiksaanwijzing
NO 68	<b>EchoTip® Ultra endobronkial HD-ultraljudsnål for Olympus-skop</b> Bruksanvisning
PL 72	<b>Igła EchoTip® Ultra do ultrasonografii wewnątrzoskrzelowej o wysokiej rozdzielczości do endoskopów Olympus</b> Instrukcja użycia
PT 76	<b>Agulha de ultrassons endobrônquica EchoTip® Ultra de alta definição para endoscópios Olympus</b> Instruções de utilização
RO 80	<b>Ac ecografic HD endobronșic EchoTip® Ultra pentru endoscoape Olympus</b> Instrucțiuni de utilizare
SK 84	<b>Endobronchiálna ultrazvuková ihla s vysokým rozlišením EchoTip® Ultra pre endoskopy Olympus</b> Návod na použitie
SV 88	<b>EchoTip® Ultra endobronkial HD ultraljudsnål för Olympus endoskop</b> Bruksanvisning
TR 92	<b>Olympus Skopları için EchoTip® Ultra Endobronşiyal HD Ultrason İğnesi</b> Kullanma Talimatı

**Fig. 1**

**EN Biopsy Port Seal**

BG Уплътнение на отвора за биопсия  
 CS Těsnění bioptického portu  
 DA Forsegling på biopsiport  
 DE Abdichtung des Biopsieanschlusses  
 EL Στεγανωτικό παρέμβυσμα θύρας βιοψίας  
 ES Junta del acceso de biopsia  
 ET Biopsiapordi tihend  
 FR Capuchon étanche de l'orifice pour biopsie  
 HU Biopsziányílás tömítése  
 IT Guarnizione del raccordo per biopsia  
 LT Biopsijos angos uždoris  
 LV Biopsijas pieslēgvieta blīve  
 NL Afsluiting biopsiepoort  
 NO Biopsiportforsegling  
 PL Uszczelka portu biopsji  
 PT Vedante do orifício de biopsia  
 RO Sigiliul portului de biopsie  
 SK Uzáver bioptického portu  
 SV Biopsiportstättning  
 TR Biyopsi Port Mührü

**EN Adapter Top**

BG Горна част на адаптера  
 CS Horní strana adaptéru  
 DA Øverste del af adapter  
 DE Oberteil des Adapters  
 EL Κορυφή του προσαρμογέα  
 ES Parte superior del adaptador  
 ET Adapteri ülemine osa  
 FR Haut de l'adaptateur  
 HU Adapter teteje  
 IT Sommità dell'adattatore  
 LT Adapterio viršus  
 LV Adapters augšdaļa  
 NL Adaptertop  
 NO Adaptertopp  
 PL Część górna łącznika  
 PT Parte de cima do adaptador  
 RO Capacul adaptorului  
 SK Vrch adaptéra  
 SV Adapters övre del  
 TR Adaptörün Tabanı

**EN Luer Connector**

BG Конектор тип Luer  
 CS Konektor Luer  
 DA Luer-konnektor  
 DE Luer-Anschluss  
 EL Συνδετικό Luer  
 ES Conector Luer  
 ET Luer-tüüpi ühendus  
 FR Connecteur Luer  
 HU Luer-csatlakozó  
 IT Connettore Luer  
 LT Luerio jungtis  
 LV Luer tipa savienotājs  
 NL Luer-connector  
 NO Luer-kobling  
 PL Łącznik Luer  
 PT Conector Luer  
 RO Conector Luer  
 SK Luerov konektor  
 SV Luer-anlutning  
 TR Luer Konektör


**EN Adapter Base**

BG Долна част (основа) на адаптера  
 CS Základna adaptéru  
 DA Nederste del af adapter  
 DE Unterteil des Adapters  
 EL Βάση του προσαρμογέα  
 ES Base del adaptador  
 ET Adapteri alumine osa  
 FR Base de l'adaptateur  
 HU Adapter alapja  
 IT Base dell'adattatore  
 LT Adapterio pagrindas  
 LV Adapters pamatne  
 NL Adapterbasis  
 NO Adapterbase  
 PL Podstawa łącznika  
 PT Base do adaptador  
 RO Baza adaptorului  
 SK Podstavec adaptéra  
 SV Adapters bas  
 TR Adaptörün Tabanı


**EN Locking Clip**

BG Блокиращ клипс  
 CS Aretovací svorka  
 DA Låseklemme  
 DE Arretierclip  
 EL Κλιπ ασφάλισης  
 ES Abrazadera de fijación  
 ET Kinnitusklamber  
 FR Pince de verrouillage  
 HU Rögzítőkapocs  
 IT Clip di bloccaggio  
 LT Fiksuojamasis spaustukas  
 LV Noslēgskava  
 NL Borgclip  
 NO Låseklips  
 PL Zacisk blokujący  
 PT Grampo de fixação  
 RO Clemă blocantă  
 SK Poistná objímka  
 SV Låsklämma  
 TR Kilitleme Klipsi


**Fig. 2**

**Fig. 3**

**Fig. 4**


**EN Plunger**

BG	Буџало
CS	Pist
DA	Stempel
DE	Kolben
EL	Έμβολο
ES	Émbolo
ET	Kolb
FR	Piston
HU	Dugattyú
IT	Stantuffo
LT	Stūmoklis
LV	Virzulis
NL	Plunjer
NO	Stempel
PL	Tłoczek
PT	Émbolo
RO	Piston
SK	Piest
SV	Kolv
TR	Piston


**EN Stopcock "open"**

BG	Кранче "отворено"
CS	„Otevřený“ uzavírací kohout
DA	Stophane "åben"
DE	Absperrhahn geöffnet
EL	Στρόφιγγα ανοικτή
ES	Llave de paso «abierta»
ET	Sulgemisventiil avatud
FR	Robinet ouvert
HU	Elzárócsap nyitott állásban
IT	Rubinetto aperto
LT	Kranelis „atviras“
LV	Slēgkrāns atvērts
NL	Afsluitkraan 'open'
NO	Stoppekran åpen
PL	Kranik w pozycji otwartej
PT	Torneira de passagem "aberta"
RO	Robinet „deschis“
SK	Otvorený kohútik
SV	Kran "öppen"
TR	Stopkok "açık"

**Fig. 5**

**Fig. 6**

**EN Sliding Sheath Adjuster**

- BG Плъзгащ се регулатор на защитната обвивка
- CS Posuvné nastavovací zařizení sheathu
- DA Glidende sheath-justering
- DE Verschiebbarer Schleuseneinsteller
- EL Συρόμενος προσαρμογέας θηκariου
- ES Ajustador deslizante de la vaina
- ET Liughülsi regulaator
- FR Ajusteur de gaine coulissant
- HU Elcsúszatható hüvelyhossz-beállító
- IT Regolatore scorrevole della guaina
- LT Slankusis vamzdėlio reguliatorius
- LV Pārīdāmais ievadapvalka regulētājs
- NL Schuifregeling van de sheath
- NO Glidende hylsejusteringsmekanisme
- PL Przesuwany regulator koszulki
- PT Regulador deslizante da bainha
- RO Ajustor glisant al tecii
- SK Posuvný nastavovač puzdra
- SV Glidande hylsjusterare
- TR Kayar Kilif Ayarlayıcı

**EN Sheath Reference Mark**

- BG Означение за дължината на защитната обвивка
- CS Referenční značka sheathu
- DA Referencemærke på sheath
- DE Schleusenbezugsmarke
- EL Σήμανση αναφοράς θηκariου
- ES Marca de referencia de la vaina
- ET Hülsi tähis
- FR Repère de référence sur la gaine
- HU Hüvelyhossz referencijele
- IT Indicatore di riferimento della guaina
- LT Vamzdėlio atskaitos žymė
- LV Ievadapvalka atskaites atzīme
- NL Referentiemarkering van de sheath
- NO Referansemerke på hylsen
- PL Znacznik referencyjny wysunięcia koszulki
- PT Marca de referência da bainha
- RO Marcaj de referință al tecii
- SK Referenčná značka puzdra
- SV Hylsans referensmarkering
- TR Kilif Referans İşareti

**EN Needle Handle**

- BG Дръжка на иглата
- CS Rukojeť jehly
- DA Nålegreb
- DE Nadelgriff
- EL Λαβή βελόνας
- ES Mango de la aguja
- ET Nõela käepide
- FR Poignée de l'aiguille
- HU Tú fogantyúja
- IT Impugnatura dell'ago
- LT Adatos kotelis
- LV Adatas rokturis
- NL Veiligheidsring
- NO Nålehåndtak
- PL Uchwyt igły
- PT Punho da agulha
- RO Mânerul acului
- SK Rukováť ihly
- SV Nålhandtag
- TR İğne Sapi



**EN Safety Ring**

- BG Предпазен пръстен
- CS Pojistný prstenec
- DA Sikkerhedsring
- DE Sicherheitsring
- EL Δακτύλιος ασφαλείας
- ES Anillo de seguridad
- ET Kaitserõngas
- FR Bague de sécurité
- HU Biztonsági gyűrű
- IT Anello di sicurezza
- LT Apsauginis žiedas
- LV Drošības gredzens
- NL Veiligheidsring
- NO Sikkerhetsring
- PL Pierścień bezpieczeństwa
- PT Anel de segurança
- RO Inel de siguranță
- SK Bezpečnostný krúžok
- SV Säkerhetslåsring
- TR Güvenlik Halkası

**EN Stylet Hub**

- BG Втулка на стилето
- CS Nástavec stiletu
- DA Stiletmuffe
- DE Mandrinsatz
- EL Ομφαλός του στείλεου
- ES Conector del estilete
- ET Stiletikinnitus
- FR Embase du stylet
- HU Mandrin végidoma
- IT Pomello del mandrino
- LT Zondo kaištis
- LV Stiletas galviņa
- NL Aanzetstuk van het stilet
- NO Stilettfeste
- PL Glowica mandrynu
- PT Conector do estilete
- RO Acoperitoarea stiletului
- SK Hrdlo styletu
- SV Mandrängfättning
- TR Stile Göbeği

**Fig. 7**

## **ECHOTIP® ULTRA ENDOBRONCHIAL HD ULTRASOUND NEEDLE FOR OLYMPUS SCOPES**

**CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.**

### **INTENDED USE**

This device is used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration (FNA).

### **DEVICE DESCRIPTION**

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an Olympus endobronchial ultrasound (EBUS) scope using an adapter. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 ml syringe to utilize the standard vacuum technique for sample collection, if desired, and an accessory channel adapter. The device is available in two different needle sizes: 22 and 25 ga.

### **CONTRAINDICATIONS**

Contraindications include those specific to the primary endoscopic ultrasound procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

### **POTENTIAL ADVERSE EVENTS**

Potential adverse events associated with bronchial endoscopy and gastrointestinal endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

### **WARNINGS**

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Ensure the stylet is fully inserted when advancing the needle into the target site.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with an Olympus EBUS scope.

## PLACEMENT OF ACCESSORY CHANNEL ADAPTER

**Illustrations**

1. Examine the adapter and note features (**Fig. 1**).
2. Place the adapter, in the open position, in line with the scope accessory channel metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel metal hub (**Fig. 2**).
4. Close the adapter ensuring it is locked (**Fig. 3**).
5. The adapter will now function as a Luer lock to which the needle can be connected. **Note:** If irrigation of accessory channel is desired/necessary, irrigate through closed adapter prior to insertion or after removal of needle.

## SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in "open" position, aligned with syringe (**Fig. 4**).
2. Prepare syringe as follows:
  - a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
  - b. Turn stopcock to "closed" position.
  - c. Pull plunger back until it is locked into place at desired setting, creating a vacuum.
  - d. Set prepared syringe aside until aspiration is required.

## INSTRUCTIONS FOR USE

1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope (**Fig. 5**).  
**Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (**Fig. 6**). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of ultrasound scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (**Fig. 7**). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for sample collection (Steps 9-11) or, if desired, other techniques that may or may not incorporate the use of this stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.

11. Turn the stopcock to the “open” position aligned with the syringe, allowing the negative pressure in the syringe to aspirate the sample. Move the needle handle back and forth within the target site. **Note:** Do not remove the needle from the target site during sample collection.
12. Upon completion of sample collection, if using the syringe, turn the stopcock to the “closed” position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
14. Unlock the thumbscrew on the safety ring and extend the needle.
15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
17. For additional sample collection from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
18. Additional samples may be obtained by repeating Step 2 through 16 of the “Instructions for Use”.
19. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
20. To remove the adapter from the scope, pull back on the locking clip, lift the top and slide off the scope.

**Upon completion of the procedure, dispose of the device, including the adapter, per institutional guidelines for biohazardous medical waste.**

#### **HOW SUPPLIED**

Store in a dry location, away from temperature extremes.

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© COOK 2021
William A. Cook Australia Pty Ltd
**COOK IRELAND LTD.**